

ADI Back Order Form - Retail

Stealth Products
Solutions For Life®

32 Wedgwood Road, Bicester, Oxon,
OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk

September 2023

Client Information

PO#: ¹	
NHS:	
Date:	
Contact:	
E-mail: ²	
Phone: ²	
Fax:	
Mark for:	

Order Type*: ☐ Quote ☐ Order

* Required Fields

1. Required if this is an order

2. One of the options is required



1 Back Support Style

☐ **ALBR**

Aluminum
Scapula Cut-out
back
support



☐ **AFTBR**

Aluminum
Flat Top
back
support



☐ **CFBR**

Carbon
Fibre
back
support



2 Back Width

- ☐ **10" (25.40cm)** Only available with Low Back Height
- ☐ **11" (27.94cm)** Only available with Low Back Height
- ☐ **12" (30.48cm)** Only available with Low or Medium Back Height
- ☐ **13" (33.02cm)** Only available with Low, Medium or Standard Back Heights

- ☐ **14" (35.56cm)** Only available with Low, Medium or Standard Back Heights
- ☐ **15" (38.10cm)** Only available with Low, Medium or Standard Back Heights
- ☐ **16" (40.64cm)** Available with all back heights
- ☐ **17" (43.18cm)** Available with all back heights

- ☐ **18" (45.72cm)** Available with all back heights
- ☐ **19" (48.26cm)** Available with all back heights
- ☐ **20" (50.80cm)** Available with all back heights

Please note: Carbon Fibre Back only available in widths 12" to 19"

3 Back Height

- ☐ **L Low 10" (25.40cm)**
Not available in Extra Deep 6" Contour

- ☐ **S Standard 16" (40.64cm)**
Not available in Extra Deep 6" Contour if Back Width is less than 16"

- ☐ **X Extra Tall 20" (50.80cm)**
16" back widths and above

- ☐ **M Medium 13" (33.02cm)**
Not available in Extra Deep 6" Contour.

- ☐ **T Tall 18" (45.72cm)**
16" back width and above

Please note: Carbon Fibre Back are only available in Low, Medium & Standard Heights

4 Back Support Contour

☐ **A Active 2.5" (6.35cm)**



☐ **D Deep 4" (10.16cm)**



☐ **6 Extra Deep 6" (15.24cm)**



Please note: Carbon Fibre Backs are only available in Active & Deep contours

PART NUMBER CONFIGURATOR

Please tick one box in each section and then add the coloured code to the relevant box below to generate your ADI part number. For example an Aluminium Flat Top 16" width, with Standard Height, Deep Contour & 2 Point Pro Hardware would be **AFTBR16SD2**


5

Back Mounting Options

- ☐ **W** Width Adjustable Elite ☐ **4** 4 Point Pro ☐ **R** Quick Release
- ☐ **2** 2 Point Pro ☐ **F** Fixed Elite *Only available with Carbon Fibre Backs*

QR Mounting WC20 Locking Pin & Lanyard Kit.....£115.00 (part number BR6005)

This option must be selected with Quick Release hardware for occupied transportation



6

Pricing (Aluminium Backs)

- ☐ **Active ADI Backrest**.....£650.00 ea. ☐ **Deep ADI Backrest**.....£725.00 ea. ☐ **Extra Deep ADI Backrest**.....£725.00 ea.

Prices shown are for any height and width of scapula or flat top aluminum backrest, and include any fitting kit. Please refer to your Stealth Specialist for Carbon Fibre Backrest Pricing

7

ADI Lateral Options

- | | |
|---|--|
| ☐ Left Fixed Lateral Support 3"H x 3"D £145.00 ea.
<i>Part: LATADJSTSM</i> | ☐ Left Swing-Away Lateral Support, 3"H x 3"D £235.00 ea.
<i>Part: LATADJSWSM</i> |
| ☐ Right Fixed Lateral Support, 3"H x 3"D £145.00 ea.
<i>Part: LATADJSTSM</i> | ☐ Right Swing-Away Lateral Support, 3"H x 3"D £235.00 ea.
<i>Part: LATADJSWSM</i> |
| ☐ Left Fixed Lateral Support, 3"H x 4"D £145.00 ea.
<i>Part: LATADJSTMD</i> | ☐ Left Swing-Away Lateral Support, 3"H x 4"D £235.00 ea.
<i>Part: LATADJSWMD</i> |
| ☐ Right Fixed Lateral Support, 3"H x 4"D £145.00 ea.
<i>Part: LATADJSTMD</i> | ☐ Right Swing-Away Lateral Support, 3"H x 4"D £235.00 ea.
<i>Part: LATADJSWMD</i> |
| ☐ Left Fixed Lateral Support, 4.5"H x 4.5"D £145.00 ea.
<i>Part: LATADJSTLG</i> | ☐ Right Swing-Away Lateral Support, 4.5"H x 4.5"D £235.00 ea.
<i>Part: LATADJSWLG</i> |
| ☐ Right Fixed Lateral Support, 4.5"H x 4.5"D £145.00 ea.
<i>Part: LATADJSTLG</i> | ☐ Right Swing-Away Lateral Support, 4.5"H x 4.5"D £235.00 ea.
<i>Part: LATADJSWLG</i> |
| ☐ Left Fixed Lateral Support, 4"H x 6"D £145.00 ea.
<i>Part: LATADJSTXLG</i> | ☐ Left Swing-Away Lateral Support, 4"H x 6"D £235.00 ea.
<i>Part: LATADJSWXLG</i> |
| ☐ Right Fixed Lateral Support, 4"H x 6"D £145.00 ea.
<i>Part: LATADJSTXLG</i> | ☐ Right Swing-Away Lateral Support, 4"H x 6"D £235.00 ea.
<i>Part: LATADJSWXLG</i> |

8

Outer Covers


- ☐ **Coolcore**.....Standard
- ☐ **AirMesh**.....N/C
- ☐ **StealthTek**.....N/C
- ☐ **Spare CoolCore Outer Cover**.....£74.00 ea.
- ☐ **Spare AirMesh Outer Cover**.....£74.00 ea.
- ☐ **Spare StealthTek Outer Cover**.....£74.00 ea.

All backs are supplied with a Coolcore outer cover as standard. If an Air Mesh or Stealth Tek outer cover is required instead, please tick one of the N/C options.

9


Accessories

- ☐ **Privacy Flap to fit cover widths 12"-14"**..... £41.00 ea.
Part: BR6002
- ☐ **Privacy Flap to fit cover widths 15" - 17"**.....£41.00 ea.
Part: BR6003
- ☐ **Privacy Flap to fit cover widths 18" - 20"**.....£41.00 ea.
Part: BR6004
- ☐ **PSIS Pad Lumbar Support, Grey**.....£7.50 ea.
Part: BR6000



Estimated Total Value ^

^ Estimated value based on estimate breakdown above



Back Support Shell Specifications


Contour Width & Depth

Active Contour																						
width	10"		11"		12"		13"		14"		15"		16"		17"		18"		19"		20"	
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
A	7.5	19.1	8.5	21.6	9.5	24.1	10.5	26.7	11.375	28.9	12.25	31.1	13.5	34.3	14.375	36.5	15.375	39.1	16.375	41.6	17.25	43.8
B	7	17.8	8	20.3	9	22.9	10	25.4	10.75	27.3	11.75	29.8	12.75	32.4	13.75	34.9	14.75	37.5	15.75	40.0	16.75	42.5
C	2.1	5.3	2.1	5.3	2.1	5.3	2.1	5.3	2.25	5.7	2.25	5.7	2.25	5.7	2.25	5.7	2.25	5.7	2.25	5.7	2.25	5.7

Deep Contour																						
width	10"		11"		12"		13"		14"		15"		16"		17"		18"		19"		20"	
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
A	7.5	19.1	8.5	21.6	9.5	24.1	10.5	26.7	12.625	32.1	13.75	34.9	14.25	36.2	15.25	38.7	16.75	42.5	17.75	45.1	18.5	47.0
B	6	15.2	7	17.8	8	20.3	9	22.9	11	27.9	12	30.5	13	33.0	14	35.6	15	38.1	16	40.6	17	43.2
C	3.625	9.2	3.625	9.2	3.625	9.2	3.625	9.2	4.25	10.8	4.25	10.8	4.25	10.8	4.25	10.8	4.25	10.8	4.25	10.8	4.25	10.8

Extra Deep Contour										
width	16"		17"		18"		19"		20"	
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
A	15.5	39.4	16.5	41.9	17.5	44.5	18.5	47.0	19.5	49.5
B	14.75	37.5	15.75	40.0	16.75	42.5	16.75	42.5	17.75	45.1
C	6	15.2	6	15.2	6	15.2	6	15.2	6	15.2

- A. Actual Back Width
- B. Trunk Width
- C. Depth





Back Support Shell Specifications

Lateral Wing Height

Low Deep (10" height / 4" contour)																						
Width	10"		11"		12"		13"		14"		15"		16"		17"		18"		19"		20"	
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
A	4.25	10.8	4.5	11.4	5	12.7	5.25	13.3	5.5	14.0	5.5	14.0	5.5	14.0	5.5	14.0	5.5	14.0	5.5	14.0	5.5	14.0
B	5	12.7	5	12.7	6.25	15.9	6.25	15.9	6.25	15.9	6.25	15.9	6.25	15.9	6.25	15.9	6.25	15.9	6.25	15.9	6.25	15.9
C	10	25.4	10	25.4	10	25.4	10	25.4	10	25.4	10	25.4	10	25.4	10	25.4	10	25.4	10	25.4	10	25.4

Medium Deep (13" height / 4" contour)																		
Width	12"		13"		14"		15"		16"		17"		18"		19"		20"	
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
A	6.5	16.5	6.5	16.5	6.75	17.1	6.75	17.1	6.75	17.1	6.75	17.1	6.75	17.1	6.75	17.1	6.75	17.1
B	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8
C	13	33.0	13	33.0	13	33.0	13	33.0	13	33.0	13	33.0	13	33.0	13	33.0	13	33.0

Standard Deep (16" height / 4" contour)														
Width	14"		15"		16"		17"		18"		19"		20"	
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
A	7.25	18.4	7.25	18.4	7.25	18.4	7.25	18.4	7.25	18.4	7.25	18.4	7.25	18.4
B	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8
C	16	40.6	16	40.6	16	40.6	16	40.6	16	40.6	16	40.6	16	40.6

Standard Extra Deep (16" height / 6" contour)										
Width	16"		17"		18"		19"		20"	
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
A	7.25	18.4	7.25	18.4	7.25	18.4	7.25	18.4	7.25	18.4
B	6	15.2	6	15.2	6	15.2	6	15.2	6	15.2
C	16	40.6	16	40.6	16	40.6	16	40.6	16	40.6

Tall Deep (18" height / 4" contour)										
Width	16"		17"		18"		19"		20"	
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
A	9.5	24.1	9.5	24.1	9.5	24.1	9.5	24.1	9.5	24.1
B	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8
C	18	45.7	18	45.7	18	45.7	18	45.7	18	45.7

Tall Extra Deep (18" height / 6" contour)										
Width	16"		17"		18"		19"		20"	
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
A	9.5	24.1	9.5	24.1	9.5	24.1	9.5	24.1	9.5	24.1
B	6	15.2	6	15.2	6	15.2	6	15.2	6	15.2
C	18	45.7	18	45.7	18	45.7	18	45.7	18	45.7

Extra Tall Deep (20" height / 4" contour)										
Width	16"		17"		18"		19"		20"	
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
A	11	27.9	11	27.9	11	27.9	11	27.9	11	27.9
B	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8
C	20	50.8	20	50.8	20	50.8	20	50.8	20	50.8

Extra Tall Extra Deep (20" height / 6" contour)										
Width	16"		17"		18"		19"		20"	
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
A	11	27.9	11	27.9	11	27.9	11	27.9	11	27.9
B	6	15.2	6	15.2	6	15.2	6	15.2	6	15.2
C	20	50.8	20	50.8	20	50.8	20	50.8	20	50.8

