

ASPEN WITH TB FLEX SEATING

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
	Postcode
Phone	Email
PO No	Order Ref

DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Clients Weight	

BASE MODEL

£4,000

- 4 Colour Options
- Height Adjustable Removable Armrests
- 70° Swing-Away Legrests
- Angle Adjustable Fold-Forward Backrest
- 14" Flat Free Drive Tyres
- 3" x 9" Flat Free Front Castors



- ATX Suspension System
- NE+ Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:
21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres.	Standard
3" x 14" Pneumatic Black Tyres.	FOC

FRONT CASTORS

3" x 9" Flat Free Black Tyres.	Standard
3" x 8" Flat Free Black Tyres.	£55

COLOURS

Back in Black	Standard
Electric Blue	FOC
Silver Bullet	FOC
Red Corvette	FOC

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT WIDTH

15" Width	FOC	16" Width	FOC
17" Width	FOC	18" Width	Standard
19" Width	FOC	20" Width	FOC

SEAT DEPTH

15" Depth	FOC	16" Depth	FOC
17" Depth	FOC	18" Depth	Standard
19" Depth	FOC	20" Depth	FOC

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	£100
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£115
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£115
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£135
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i>	£135
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

Manual Settable Tilt. <i>0° to 30° tool adjustable with a telescopic tube.</i>	Standard		
Powered Seat Tilt. <i>0° to 30°. Uses 1 power option.</i>	£300		
Powered Seat Lift & Tilt. <i>10" lift & 50° tilt. Not available with VR2. Uses 2 power options.</i>	£1,895.00		
Manual Tilt Seat to Floor Height. <i>From seat pan.</i>	17.5"	18"	18.5"
Powered Tilt Seat to Floor Height. <i>From seat pan.</i>	17.5"	18" <i>Changes Tilt range to -5° to 25°</i>	18.5" <i>Changes Tilt range to -10° to 20°</i>
Powered Lift & Tilt Seat to Floor Height. <i>From seat pan.</i>	17.5"	18"	18.5"

STANDARD BACK SUPPORT POSITIONING

Manual Angle Settable, Fold-Forward. <i>Adj. range 95° to 125°. Will not fold with 10" lift.</i>	Standard
Attendant Controlled Manual Recline. <i>Adj. range 95° to 125°. Not suitable for occupied transportation.</i>	£150
Powered Back Recline. <i>Adj. range 95° to 125°. Uses 1 power option.</i>	£345

STANDARD BACK SUPPORTS

Tension Adjustable. <i>Air-mesh cover. Adjustment range 19" to 23" from seat board. OBAS** required to mount headrest.</i>	Standard
TRU-Comfort 2. <i>4-way stretch cover. Adjustable range 24.5" to 28" from seat board.</i>	£135
ADI Aluminium Back. <i>Skin protection, 4 way stretch, Coolcore cover.</i>	POA
Omit Back Support. <i>Supplied with back canes and brace bar only.</i>	-£25

Enter ADI Back Support part number here for quotation

STEALTH PRO BACK SUPPORTS

Standard, Coolcore Cover. <i>Solid back, 9 cm contour.</i>	£227			
Active, Coolcore Cover. <i>Solid back, 12 cm contour.</i>	£341			
Extra Deep Active, Coolcore Cover. <i>Solid back, 17 cm contour.</i>	£448			
Active Flex, Coolcore Cover. <i>Tension adjustable straps. 12 cm contour.</i>	£559			
Back Height	50 cm. <i>Actual height.</i>	Standard	57 cm. <i>Actual height.</i>	FOC

The back support width will match the seat width selected in seat dimensions.
Not available with paediatric seat dimensions.

ARM SUPPORTS

Std Flip-up Height Adjustable. <i>Includes contour arm pads. Do not select arm pads below.</i>	Standard			
TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan.</i>	£85			
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>	£279			
Arm Support Pads	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC

LEG SUPPORTS

70° Swing-away Legrests. <i>Adjustment range 12" to 15".</i>	Standard		
60° to 90° Tool Adjustable Swing-away Legrests. <i>Adjustment range 11" to 16.5".</i>	£330		
Manual Elevating Swing-away Legrests. <i>Adjustment range 15" to 20".</i>	£250		
Power Elevating Swing-away Legrests. <i>Adjustment range 14.5" to 19.5". Min. seat width 16". Uses 2 power options.</i>	£650		
Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away legrest hanger for mounting.</i>	£300		
Limb Support	LEFT	RIGHT	
Seat pan to heel measurement	(Leave blank for standard factory)		Adjustment ranges above are from seat pan to footplates. Individual legrests supplied with heel loops. Elevating legrests supplied with calf supports.

STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£219
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	£219
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i>	£327
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£88
Push Handles Head Mount Adapter OBAS**. <i>Required to mount head support with TA Back or Active Flex Pro Back.</i>	£266

Please visit stealthproducts.com for alternative head supports.

TRUNK SOLUTIONS

Thoracic Lateral Supports.	LEFT	£150	Lateral Thoracic Swing-away Hardware.	£150	RIGHT			
		£35	3" H x 5" L Curved Lateral Pad.	£35				
		£35	4" H x 6" L Curved Lateral Pad.	£35				
Stealth 4-Point. Harness	Extra Small	£137	Classic	STD*	Dynamic	STD*	Vinyl Face	FOC
	Small	£137		Contoured		FOC		Structured
	Medium	£137	*STD = Standard option **Classic Only					
	Large	£137						
	Extra Large	£137						
Chest Strap.	Extra Small	£87	Small	£87	Side Release	STD*		
	Medium	£87	Large	£87			Velcro Release	£19
	Extra Large	£87						

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.				Standard (STD)			
Retractable Lap Belt, Push Button Release.				£156			
Stealth Positioning 2-Point Padded Pelvic Belt.	Extra Small	£78	Single Pull	Standard	Side Release	Standard	
	Small	£78	Dual Pull	FOC	Push Button	FOC	
	Medium	£78	Rear Pull	FOC			
	Large	£78					
	Extra Large	£78					
Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull.		Extra Small	£120	Small	£120	Push Button	STD*
		Medium	£120	Large	£120	Side Release	£20

Please visit stealthproducts.com for alternative lateral, pelvic & lower extremity solutions.

*STD = Standard option

NON-EXPANDABLE CONTROLS

NE + Controls. <i>With no power options.</i>		Standard
NE + Controls. <i>2 power option ready.</i>		£125
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	£100
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

EXPANDABLE CONTROLS

Q-Logic 3e Controls. <i>No power options.</i>	£225
Q-Logic 3e Controls. <i>3 power option ready.</i>	£300
Q-Logic 3e Controls. <i>5 power option ready.</i>	£500
Q-Logic 3 EX Controls. <i>No power options.</i>	£275
Q-Logic 3 EX Controls. <i>3 power option ready.</i>	£350
Q-Logic 3 EX Controls. <i>5 power option ready.</i>	£550
Q-Logic 3 Enhanced Display. <i>Includes 2 egg switches & flexible display mount.</i>	£1074
SCIM Module Kit. <i>Allows Q-Logic 3 EX joystick display to be used for speciality controls. Not available with 3e controls.</i>	£316

iAccess OLED 4 Button Switch Module. <i>Mounted to opposite armrest unless specified.</i>	£300	Standard	If iAccess is selected choose either standard or Toggle lever option.
		Toggle Lever	
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard	
	Omit Hand Control. <i>Often selected when ED is used as master remote with special device.</i>	-£100	
	Attendant Control. <i>Hand control mounted at rear.</i>	£100	
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430	

Please visit stealthproducts.com for alternative control options such as mini joysticks, head arrays & switches.

HAND CONTROL MOUNTS

Fixed.		Standard					
Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i>		£75					
Retract 4 Swing-away. <i>Full 180° Movement.</i>		£150					
Hand Control Joystick Handle.	Mounting Position.	Left	FOC	Right	Standard		
		Standard Carrot	STD	T-Bar	£35	Extended I-Handle	£35
		Mushroom	£35	Medium Ball	£35	3" Goal Post	£50
		Chin Cup	£35	Soft Ball	£35	4" Goal Post	£50

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
4 Mph Program with Lights and indicators.	£450
6 Mph Program with Lights and indicators.	£450
Lights & Indicator Surcharge. <i>Required if no power options have been selected.</i>	£260

ACCESSORIES

Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379
Quantum Back Pack.	£110
XLR USB Charger Port.	£17
Magnetic to XLR Charger Adapter. <i>Only with Q-Logic 3e.</i>	£49
Kerb Climber.	£50

NOTES & SPECIAL REQUESTS

PRICE

**ASPEN WITH
TB FLEX SEATING**



BASE PRICE	£4,000
COST OPTIONS	
TOTAL	

5



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
RETAIL UK - 11/09/24 V1

