

ASPEN WITH TB FLEX SEATING

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
	Postcode
Phone	Email
PO No	Order Ref

DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Clients Weight	

BASE MODEL

£4,450

- 4 Colour Options
- Height Adj. Removable Arm Supports
- 70° Swing-Away Leg Supports
- Angle Adj. Fold-Forward Back Support
- 14" Flat Free Drive Tyres
- 3" x 9" Flat Free Front Castors

- ATX Suspension System
- QLNES Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:
21 st 6lbs. (136 kg)



Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres.	Standard
3" x 14" Pneumatic Black Tyres.	FOC

FRONT CASTORS

3" x 9" Flat Free Black Tyres.	Standard
3" x 8" Flat Free Black Tyres.	£55

COLOURS

Back in Black	Standard
Electric Blue	FOC
Silver Bullet	FOC
Red Corvette	FOC

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT WIDTH

15"Width	FOC	16"Width	FOC
17"Width	FOC	18"Width	Standard
19"Width	FOC	20"Width	FOC

SEAT DEPTH

15"Depth	FOC	16"Depth	FOC
17"Depth	FOC	18"Depth	Standard
19"Depth	FOC	20"Depth	FOC

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	£99
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£126
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£126
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£148
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i>	£148
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£209
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£209

Please note that seat cushions may not be available in all width and depth combinations.
Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

Manual Settable Tilt. 0° to 30° tool adjustable with a telescopic tube.				Standard
Powered Seat Tilt. 0° to 30°. Uses 1 power option.				£500
Powered Seat Lift & Tilt. 10" lift & 50° tilt. Not available with VR2. Uses 2 power options.				£1950
Manual Tilt Seat to Floor Height. From seat pan.	17.5"	18"	18.5"	
Powered Tilt Seat to Floor Height. From seat pan.	17.5"	18" Changes Tilt range to -5° to 25°	18.5" Changes Tilt range to -10° to 20°	
Powered Lift & Tilt Seat to Floor Height. From seat pan.	17.5"	18"	18.5"	

STANDARD BACK SUPPORT POSITIONING

Manual Angle Settable, Fold-Forward. <i>Adj. range 95° to 125°. Will not fold with 10" lift.</i>	Standard
Attendant Controlled Manual Recline. <i>Adj. range 95° to 125°. Not suitable for occupied transportation.</i>	£150
Powered Back Recline. <i>Adj. range 95° to 125°. Uses 1 power option.</i>	£400

STANDARD BACK SUPPORTS

Tension Adjustable. <i>Air-mesh cover. Adjustment range 19" to 23" from seat board. OBAS** required to mount headrest.</i>	Standard
TRU-Comfort 2. <i>4-way stretch cover. Adjustable range 24.5" to 28" from seat board.</i>	£135
ADI Aluminium Back. <i>Skin protection, 4 way stretch, Coolcore cover.</i>	POA
Enter ADI Back Support part number here for quotation	

STEALTH PRO BACK SUPPORTS

Standard, Coolcore Cover. <i>Solid back, 9 cm contour.</i>				£227	
Active, Coolcore Cover. <i>Solid back, 12 cm contour.</i>				£341	
Extra Deep Active, Coolcore Cover. <i>Solid back, 17 cm contour.</i>				£448	
Active Flex, Coolcore Cover. <i>Tension adjustable straps. 12 cm contour.</i>				£559	
●●●●●	Back Height	50 cm. <i>Actual height.</i>	Standard	57 cm. <i>Actual height.</i>	FOC

The back support width will match the seat width selected in seat dimensions.
Not available with paediatric seat dimensions.

ARM SUPPORTS

TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan.</i>				Standard
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>				£279
Arm Support Pads	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC

LEG SUPPORTS

70° Swing-Away. <i>Extension range 12" to 15". Minimum seat width 15".</i>	Standard	
60° to 90° Tool Adjustable Swing-Away. <i>6"W x 8"D foot supports. Extension range 11" to 19.5". Min. seat width 15".</i>	£330	
Manual Elevating Swing-Away. <i>Extension range 15" to 20". Min. seat width 15"</i>	£425	
Power Elevating Articulating Swing-Away. <i>5"W x 7"D foot supports. Ext. range 9" to 21". Min. seat width 15". Uses 2 power options</i>	£850	
Powered Leg Support Operation Method.	Independent. <i>Uses 2 power options.</i>	Combined. <i>Uses 1 power option.</i>

❶ Enter seat pan to heel measurement.

If NO measurement is supplied, standard factory length option will be fitted.

Stealth Residual Limb Support (each). *Includes a 60-90° Swing-away legrest hanger for mounting.*

LEFT

£300

Limb Support

£300

RIGHT

Extension ranges above are from seat pan to foot supports. Individual leg supports supplied with heel loops.
Elevating leg supports supplied with standard calf supports.

STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£230
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	£230
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i>	£342
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£92
Push Handles Head Mount Adapter OBAS**. <i>Required to mount head support with TA Back or Active Flex Pro Back.</i>	£266

Please visit stealthproducts.com for alternative head supports.

TRUNK SOLUTIONS

Thoracic Lateral Supports.	LEFT	£175	Lateral Thoracic Swing-away Hardware.				£175	RIGHT			
		£42	3" H x 5" L Curved Lateral Pad.				£42				
		£42	4" H x 6" L Curved Lateral Pad.				£42				
Stealth 4-Point. Harness	Extra Small	£143	}	Classic	STD*	}	Dynamic	STD*	}	Vinyl Face	FOC
	Small	£143		Contoured	FOC		Structured	FOC		Zipper Face**	FOC
	Medium	£143		*STD = Standard option **Classic Only							
	Large	£143									
	Extra Large	£143									
Chest Strap.	Extra Small	£91	Small	£91	}	Side Release				STD*	
	Medium	£91	Large	£91		Velcro Release				£20	
	Extra Large	£91									

*STD = Standard option **Classic Only

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.						Standard						
Retractable Lap Belt, Push Button Release.						£169						
Stealth Positioning 2-Point Padded Pelvic Belt.	Extra Small	£84	}	Single Pull	Standard	}	Side Release	Standard				
	Small	£84		Dual Pull	FOC		Push Button	FOC				
	Medium	£84		Rear Pull	FOC							
	Large	£84										
	Extra Large	£84										
Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull.						Extra Small	£143	Small	£143	}	Push Button	STD*
						Medium	£143	Large	£143		Side Release	£21

Please visit stealthproducts.com for alternative lateral, pelvic & lower extremity solutions.

*STD = Standard option

NON-EXPANDABLE CONTROLS

QLNES Controls. <i>With no power options.</i>		Standard
QLNES Controls. <i>2 power option ready.</i>		£125
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	£100
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

EXPANDABLE CONTROLS

Q-Logic 3e Controls. <i>No power options.</i>		£300
Q-Logic 3e Controls. <i>3 power option ready.</i>		£350
Q-Logic 3e Controls. <i>5 power option ready.</i>		£500
Q-Logic 3 EX Controls. <i>No power options.</i>		£350
Q-Logic 3 EX Controls. <i>3 power option ready.</i>		£375
Q-Logic 3 EX Controls. <i>5 power option ready.</i>		£550
Q-Logic 3 Enhanced Display. <i>Includes 2 egg switches & flexible display mount. Not available with 3e controls.</i>		£1074
SCIM Module Kit. <i>Allows Q-Logic 3 EX joystick display to be used for speciality controls. Not available with 3e controls.</i>		£316
iAccess OLED 4 Button Switch Module. <i>Mounted to opposite armrest unless specified.</i>		£300
		Standard Toggle Lever
		If iAccess is selected choose either standard or Toggle leaver option.
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Omit Hand Control. <i>Often selected when ED is used as master remote with special device.</i>	-£100
	Attendant Control. <i>Hand control mounted at rear.</i>	£100
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

Please visit stealthproducts.com for alternative control options such as mini joysticks, head arrays & switches.

HAND CONTROL MOUNTS

Fixed.						Standard
Swing-away.						£75
Advanced Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i>						£100
Mounting Position.	Left	FOC	Right	Standard		
Hand Control Joystick Handle.	Standard Carrot	STD	T-Bar	£35	Extended I-Handle	£35
	Mushroom	£35	Medium Ball	£35	3" Goal Post	£50
	Chin Cup	£35	Soft Ball	£35	4" Goal Post	£50

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
4 Mph Program with Lights and indicators.	£450
6 Mph Program with Lights and indicators.	£450
Lights & Indicator Surcharge. <i>Required if no power options have been selected.</i>	£260

ACCESSORIES

Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379
Quantum Back Pack.	£110
XLR USB Charger Port.	£17
Magnetic to XLR Charger Adapter. <i>Only with Q-Logic 3e.</i>	£49
Kerb Climber.	£50

NOTES & SPECIAL REQUESTS

PRICE

**ASPEN WITH
TB FLEX SEATING**



BASE PRICE	£4,450
COST OPTIONS	
TOTAL	

5



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
RETAIL UK - 19/03/25 V1

