

## EDGE 3 STRETTO WITH CONTOURED SEATING

### DEALERS DETAILS

Dealers Name \_\_\_\_\_

Order Date \_\_\_\_\_

Address \_\_\_\_\_

Postcode \_\_\_\_\_

Phone \_\_\_\_\_ Email \_\_\_\_\_

PO No \_\_\_\_\_ Order Ref \_\_\_\_\_

### DELIVERY DETAILS

Dealers Address \_\_\_\_\_

Postcode \_\_\_\_\_

Quantum Specialist \_\_\_\_\_

Ⓢ Clients Weight (required for suspension settings) \_\_\_\_\_

### BASE MODEL

**£3,670**

- SRS Suspension System
- 12.5" Flat Free Drive Tyres
- Flip-up Armrests
- 70° Swing-Away Legrests
- Contoured Seat with Fold-Forward Backrest
- 4 Mph (6 Km/h) Motors



- 14 Colour Options
- Max User Weight: 21 st 6lbs. (136 kg)
- 2 x 62 ah Batteries
- NE+ Control System
- 2" x 5" Flat Free Front Castors
- 2" x 6" Flat Free Rear Castors

Standard = Standard Factory Feature  
 FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 12.5" Flat Free Black Tyres.	Standard
3" x 12.5" Pneumatic Black Tyres.	FOC
3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	£50
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	£50
4" x 14" Pneumatic Black Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£220
4" x 14" Flat Free Black Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£320

### COLOURS

Back in Black	Standard	Lemon Crush*	FOC
Green Machine	FOC	Night Train*	FOC
Midnight Blue	FOC	Evolution Orange*	FOC
Electric Blue	FOC	Military Green*	FOC
Blue Jean	FOC	Grey Street*	FOC
Deep Purple	FOC	Raspberry Beret*	FOC
Red Corvette	FOC	Custom Hydro Dip*	£600
Pretty in Pink	FOC		

\* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



## CONTOURED SEAT PACKAGES

16"W x 16"D	FOC	18"W x 18"D	Standard	20"W x 20"D	FOC
-------------	-----	-------------	----------	-------------	-----

Seat package includes standard headrest and auto buckle lap belt.

## SEAT CUSHIONS

No Seat Cushion Required	FOC
Contoured Seat Cushion. <i>136 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	£70
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	£130
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	£150
Solution SPP: Skin Protection, Positioning, 4-way stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£170

Seat cushion dimensions will match seat package selected above.

## SEAT POSITIONING

Manual Settable Tilt. <i>0° to 30° tool adjustable with a telescopic tube.</i>	Standard		
iLevel Powered Seat Lift. <i>12" lift. Uses 1 power option. Not available with VR2.</i>	£1,995		
Manual Seat to Floor Height. <i>From seat pan.</i>	16.5"	17.5"	18.5"
iLevel Powered Lift Seat to Floor Height. <i>From seat pan.</i>	18"	19"	20"

## ARM SUPPORTS

Std Flip-up Height Adjustable. <i>Includes contour arm pads. Do not select arm pads below.</i>	Standard			
TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan.</i>	£100			
Single Post Height Adjustable. <i>Adjustment range 11"-17" from seat pan.</i>	£150			
Arm Support Pads	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC

## LEG SUPPORTS

70° Swing-away Legrests. <i>Adjustment range 12" to 15".</i>	Standard		
60° to 90° Tool Adjustable Swing-away Legrests. <i>Adjustment range 11" to 16.5".</i>	£350		
Manual Elevating Swing-Away Legrests. <i>Inc. calf pads. Adjustment range 15" to 20". Min. seat width 15".</i>	£250		
Manually Adjustable Centre Mount Foot Platform. <i>Inc. calf pads. Adj. range 10" to 16". 5"W x 11"L Footplates. Min. seat width 15".</i>	£350		
Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away legrest hanger for mounting.</i>	£300		
Limb Support	<b>LEFT</b>	<b>RIGHT</b>	Adjustment ranges above are from seat pan to footplates. Individual legrests are supplied with heel loops.
Seat pan to heel measurement	(Leave blank for standard factory)		

## HEAD SUPPORTS

Contoured Seating Headrest.	Standard
-----------------------------	----------

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative head supports.

## PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.	Standard
--	----------

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative lap belts.

## NON-EXPANDABLE CONTROLS

NE + Controls. <i>With no power options.</i>	Standard
NE + Controls. <i>2 power option ready.</i>	£125
VR2 Controls. <i>With no power options.</i>	FOC
VR2 Controls. <i>2 power option ready.</i>	£125

## HAND CONTROL MOUNTS

Fixed.	Standard
Swing-away. <i>Swings in &amp; out with height, tension and pivot adjustment.</i>	£75
Mounting Position.	Left FOC Right Standard

## SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
4 Mph Program with Lights and indicators.	£450
6 Mph Program with Lights and indicators.	£450
Lights & Indicator Surcharge. <i>Required if no power options have been selected.</i>	£260
LED Fender Lights. <i>Includes single switch button mounted to side rail. Not road legal.</i>	£195

## ACCESSORIES

Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379
Quantum Back Pack.	£110
Additional XLR USB Charger Port.	£17

## NOTES & SPECIAL REQUESTS

## PRICE

**EDGE 3 STRETTO WITH  
CONTOURED SEATING**



BASE PRICE	<b>£3,670</b>
COST OPTIONS	
TOTAL	

3



**Quantum®**  
32 Wedgwood Road, Bicester,  
Oxon, OX26 4UL  
Phone: 01869 324600  
Fax: 01869 323070  
Email: sales@quantumrehab.co.uk  
Web: www.quantumrehab.co.uk  
**RETAIL UK - 17/03/23 V1**



**Crash Tested**



V1: 14/06/2022