

EDGE 3 STRETTO WITH CONTOURED SEATING

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Ⓢ Clients Weight (required for suspension settings) _____

BASE MODEL

£3,850

- SRS Suspension System
- 12.5" Flat Free Drive Tyres
- Flip-up Armrests
- 70° Swing-Away Legrests
- Contoured Seat with Fold-Forward Backrest
- 4 Mph (6 Km/h) Motors



- 14 Colour Options
- Max User Weight: 21 st 6lbs. (136 kg)
- 2 x 62 ah Batteries
- QLNES Control System
- 2" x 5" Flat Free Front Castors
- 2" x 6" Flat Free Rear Castors

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 12.5" Flat Free Black Tyres.	Standard
3" x 12.5" Pneumatic Black Tyres.	FOC
3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	£50
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	£50
4" x 14" Pneumatic Black Tyres, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£220
4" x 14" Flat Free Black Tyres, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£320

COLOURS

Back in Black	Standard	Lemon Crush*	FOC
Green Machine	FOC	Night Train*	FOC
Midnight Blue	FOC	Evolution Orange*	FOC
Electric Blue	FOC	Military Green*	FOC
Blue Jean	FOC	Grey Street*	FOC
Deep Purple	FOC	Raspberry Beret*	FOC
Red Corvette	FOC	Custom Hydro Dip*	£600
Pretty in Pink	FOC		

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



CONTOURED SEAT PACKAGES

16"W x 16"D	FOC	18"W x 18"D	Standard	20"W x 20"D	FOC
-------------	-----	-------------	----------	-------------	-----

Seat package includes standard headrest and auto buckle lap belt.

SEAT CUSHIONS

No Seat Cushion Required	FOC
Contoured Seat Cushion. <i>136 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	£70
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	£130
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	£150
Solution SPP: Skin Protection, Positioning, 4-way stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£170

Seat cushion dimensions will match seat package selected above.

SEAT POSITIONING

Static Seat.	Standard
iLevel Powered Seat Lift. <i>12" lift. Uses 1 power option. Not available with VR2. Not available with lights and indicators.</i>	£1,995
Static Seat to Floor Height. <i>From seat pan.</i>	16.5" 17.5" 18.5"
iLevel Powered Lift Seat to Floor Height. <i>From seat pan.</i>	18" 19" 20"

ARM SUPPORTS

Std Flip-up Height Adjustable. <i>Includes contour arm pads. Do not select arm pads below.</i>	Standard
TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan.</i>	£100
Single Post Height Adjustable. <i>Adjustment range 11"-17" from seat pan.</i>	£150
Arm Support Pads	10" Standard Standard 10" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC
	14" Standard FOC 14" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC

LEG SUPPORTS

70° Swing-away Legrests. <i>Adjustment range 12" to 15".</i>	Standard
60° to 90° Tool Adjustable Swing-away Legrests. <i>Adjustment range 11" to 16.5".</i>	£350
Manual Elevating Swing-Away Legrests. <i>Inc. calf pads. Adjustment range 15" to 20". Min. seat width 15".</i>	£250
Manually Adjustable Centre Mount Foot Platform. <i>Inc. calf pads. Adj. range 10" to 16". 5"W x 11"L Footplates. Min. seat width 15".</i>	£350
Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away legrest hanger for mounting.</i>	£300
Limb Support	LEFT RIGHT
Seat pan to heel measurement	(Leave blank for standard factory)

Adjustment ranges above are from seat pan to footplates. Individual legrests are supplied with heel loops.

HEAD SUPPORTS

Contoured Seating Headrest.	Standard
-----------------------------	----------

Please visit stealthproducts.com for alternative head supports.

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.	Standard
--	----------

Please visit stealthproducts.com for alternative lap belts.

NON-EXPANDABLE CONTROLS

QLNES Controls. <i>With no power options.</i>	Standard
QLNES Controls. <i>2 power option ready.</i>	£125
VR2 Controls. <i>With no power options. Not available with fender lights.</i>	FOC

HAND CONTROL MOUNTS

Fixed.	Standard
Basic swing-away.	£75
Advanced swing-away. <i>Chain driven. Swings in & out with height, tension and pivot adjustment.</i>	£100
Mounting Position.	Left FOC Right Standard

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
4 Mph Program with Lights and indicators.	£450
6 Mph Program with Lights and indicators. <i>Not available with iLevel.</i>	£450
Lights & Indicator Surcharge. <i>Required if no power options have been selected.</i>	£260
LED Fender Lights. <i>Includes single switch button mounted to seat rail. Not road legal. Not available with VR2 controls.</i>	£195

ACCESSORIES

Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379
Quantum Back Pack.	£110
Additional XLR USB Charger Port.	£17

NOTES & SPECIAL REQUESTS

PRICE

**EDGE 3 STRETTO WITH
CONTOURED SEATING**



BASE PRICE

£3,850

COST OPTIONS

TOTAL

3



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
RETAIL UK - 28/02/25 V1.1

