

## EDGE 3 STRETTO WITH CONTOURED SEATING

### DEALERS DETAILS

|              |           |
|--------------|-----------|
| Dealers Name |           |
| Order Date   |           |
| Address      |           |
|              | Postcode  |
| Phone        | Email     |
| PO No        | Order Ref |

### DELIVERY DETAILS

|   |          |
|---|----------|
| Dealers Address                                     |          |
|   |          |
|   | Postcode |
| Quantum Specialist                                  |          |
| Ⓢ Clients Weight (required for suspension settings) |          |

### BASE MODEL

£4,000

SRS Suspension System  
12.5" Flat Free Drive Tyres  
Flip-up Armrests  
70° Swing-Away Legrests  
Contoured Seat with Fold-Forward Backrest  
4 Mph (6 Km/h) Motors

14 Colour Options  
Max User Weight: 21 st 6lbs. (136 kg)  
2 x 62 ah Batteries  
QLNES Control System  
2" x 5" Flat Free Front Castors  
2" x 6" Flat Free Rear Castors



Standard = Standard Factory Feature  
FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

|   |          |
|---|----------|
| 3" x 12.5" Flat Free Black Tyres.   | Standard |
| 3" x 12.5" Pneumatic Black Tyres.   | FOC      |
| 3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>                              | £50      |
| 3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>                              | £50      |
| 4" x 14" Pneumatic Black Tyres, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i> | £220     |
| 4" x 14" Flat Free Black Tyres, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i> | £320     |

### COLOURS

|                |          |                   |      |
|----------------|----------|-------------------|------|
| Back in Black  | Standard | Lemon Crush*      | FOC  |
| Green Machine  | FOC      | Night Train*      | FOC  |
| Midnight Blue  | FOC      | Evolution Orange* | FOC  |
| Electric Blue  | FOC      | Military Green*   | FOC  |
| Blue Jean      | FOC      | Grey Street*      | FOC  |
| Deep Purple    | FOC      | Raspberry Beret*  | FOC  |
| Red Corvette   | FOC      | Custom Hydro Dip* | £600 |
| Pretty in Pink | FOC      |                   |      |

\* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



## CONTOURED SEAT PACKAGES

| 16"W x 16"D   | FOC | 18"W x 18"D | Standard | 20"W x 20"D | FOC |
|---|-----|-------------|----------|-------------|-----|
| Seat package includes standard headrest and auto buckle lap belt. |     |             |          |             |     |

## SEAT CUSHIONS

|  |          |
|--|----------|
| No Seat Cushion Required   | FOC      |
| Contoured Seat Cushion. <i>136 Kg weight limit.</i>  | Standard |
| Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>                              | £83      |
| Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>                                      | £165     |
| TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>             | £182     |
| Solution SPP: Skin Protection, Positioning, 4-way stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i> | £209     |

Seat cushion dimensions will match seat package selected above.

## SEAT POSITIONING

|   |       |       |          |
|---|-------|-------|----------|
| Static Seat.  |       |       | Standard |
| iLevel Powered Seat Lift. <i>12" lift. Uses 1 power option. Not available with VR2. Not available with lights and indicators.</i> |       |       | £1,995   |
| Static Seat to Floor Height. <i>From seat pan.</i>  | 16.5" | 17.5" | 18.5"    |
| iLevel Powered Lift Seat to Floor Height. <i>From seat pan.</i>   | 18"   | 19"   | 20"      |

## ARM SUPPORTS

|  |              |          |  |          |
|--|--------------|----------|--|----------|
| Std Flip-up Height Adjustable. <i>Includes contour arm pads. Do not select arm pads below.</i> |              |          |  | Standard |
| TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan.</i>             |              |          |  | £100     |
| Single Post Height Adjustable. <i>Adjustment range 11"-17" from seat pan.</i>                  |              |          |  | £150     |
| Arm Support Pads   | 10" Standard | Standard | 10" Waterfall. <i>Not with Std flip-up or TB Flex.</i> | FOC      |
|  | 14" Standard | FOC      | 14" Waterfall. <i>Not with Std flip-up or TB Flex.</i> | FOC      |

## LEG SUPPORTS

|   |                                    |       |  |
|---|------------------------------------|-------|--|
| 70° Swing-away Legrests. <i>Adjustment range 12" to 15".</i>  |                                    |       | Standard   |
| 60° to 90° Tool Adjustable Swing-away Legrests. <i>Adjustment range 11" to 16.5".</i>   |                                    |       | £350   |
| Manual Elevating Swing-Away Legrests. <i>Inc. calf pads. Adjustment range 15" to 20". Min. seat width 15".</i>                            |                                    |       | £425   |
| Manually Adjustable Centre Mount Foot Platform. <i>Inc. calf pads. Adj. range 10" to 16". 5"W x 11"L Footplates. Min. seat width 15".</i> |                                    |       | £350   |
| Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away legrest hanger for mounting.</i>                                    |                                    |       | £300   |
| Limb Support  | LEFT                               | RIGHT | Adjustment ranges above are from seat pan to footplates. Individual legrests are supplied with heel loops. |
| Seat pan to heel measurement  | (Leave blank for standard factory) |       |  |

## HEAD SUPPORTS

|                             |          |
|-----------------------------|----------|
| Contoured Seating Headrest. | Standard |
|-----------------------------|----------|

## PELVIC SOLUTIONS

|  |          |
|--|----------|
| Auto Buckle Lap Belt, Push Button Release. | Standard |
|--|----------|

Please visit [stealthproducts.com](https://stealthproducts.com) for alternative lap belts.

## NON-EXPANDABLE CONTROLS

|   |          |
|---|----------|
| QLNES Controls. <i>With no power options.</i>                                 | Standard |
| QLNES Controls. <i>2 power option ready.</i>                                  | £125     |
| VR2 Controls. <i>With no power options. Not available with fender lights.</i> | FOC      |

## HAND CONTROL MOUNTS

|  |                         |
|--|-------------------------|
| Fixed.   | Standard                |
| Basic swing-away.  | £75                     |
| Advanced swing-away. <i>Chain driven. Swings in &amp; out with height, tension and pivot adjustment.</i> | £100                    |
| Mounting Position.   | Left FOC Right Standard |

## SPEED PROGRAM & LIGHTS PACKAGE

|  |          |
|--|----------|
| 4 Mph Program. <i>Lights and indicators not included.</i>  | Standard |
| 4 Mph Program with Lights and indicators.  | £450     |
| 6 Mph Program with Lights and indicators. <i>Not available with iLevel.</i>  | £450     |
| Lights & Indicator Surcharge. <i>Required if no power options have been selected.</i>  | £260     |
| LED Fender Lights. <i>Includes single switch button mounted to seat rail. Not road legal. Not available with VR2 controls.</i> | £195     |

## ACCESSORIES

|  |      |
|--|------|
| Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i> | £379 |
| Quantum Back Pack.   | £110 |
| Additional XLR USB Charger Port.   | £17  |

## NOTES & SPECIAL REQUESTS

## PRICE

### EDGE 3 STRETTO WITH CONTOURED SEATING



BASE PRICE

**£4,000**

COST OPTIONS

TOTAL

3



**Quantum®**  
32 Wedgwood Road, Bicester,  
Oxon, OX26 4UL  
Phone: 01869 324600  
Fax: 01869 323070  
Email: sales@quantumrehab.co.uk  
Web: www.quantumrehab.co.uk  
**RETAIL UK - 01/04/25 V1.1**

