

## EDGE 3 STRETTO WITH TB FLEX SEATING

### DEALERS DETAILS

Dealers Name	
Order Date	
Address	
Postcode	
Phone	Email
PO No	Order Ref

### DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Ⓢ Clients Weight (required for suspension settings)	

### BASE MODEL

**£5,750**

SRS Suspension  
12.5" Flat Free Drive Tyres  
Height Adjustable Removable Armrests  
70° Swing-Away Legrests  
Angle Adjustable Flip Forward Backrest  
4 Mph (6 Km/h) Motors



14 Colour Options  
Max User Weight: 21 st 6lbs. (136 kg)  
2 x 40ah batteries  
QLNES Control System  
2x5" Front Castors, Flat Free  
2x6" Rear Castors, Flat Free

Standard = Standard Factory Feature  
FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 12.5" Flat Free Black Tyres.	Standard
3" x 12.5" Pneumatic Black Tyres.	FOC
3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	£50
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	£50
4" x 14" Pneumatic Black Tyres, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£220
4" x 14" Flat Free Black Tyres, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£320

### COLOURS

Back in Black	Standard	Lemon Crush*	FOC
Green Machine	FOC	Night Train*	FOC
Midnight Blue	FOC	Evolution Orange*	FOC
Electric Blue	FOC	Military Green*	FOC
Blue Jean	FOC	Grey Street*	FOC
Deep Purple	FOC	Raspberry Beret*	FOC
Red Corvette	FOC	Custom Hydro Dip*	£600
Pretty in Pink	FOC		

\* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



## PAEDIATRIC SEAT DIMENSIONS

10"Width	FOC	10"Depth	FOC
11"Width	FOC	11"Depth	FOC
12"Width	FOC	12"Depth	FOC
13"Width	FOC	13"Depth	FOC
14"Width	FOC	14"Depth	FOC

## STANDARD SEAT DIMENSIONS

15"Width	FOC	15"Depth	FOC
16"Width	FOC	16"Depth	FOC
17"Width	FOC	17"Depth	FOC
18"Width	Standard	18"Depth	Standard
19"Width	FOC	19"Depth	FOC
20"Width	FOC	20"Depth	FOC

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	£99
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£126
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£126
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£148
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i>	£148
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£209
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£209

Please note that seat cushions may not be available in all width and depth combinations.  
Refer to our Stealth Cushion Order Forms for size availability.

## SEAT POSITIONING

Manual Settable Tilt. <i>0° to 30° tool adjustable with a telescopic tube.</i>	Standard		
Powered Seat Tilt. <i>0° to 30°. Uses 1 power option.</i>	£500		
Powered Seat Lift & Tilt. <i>10" lift &amp; 50° tilt. 15" to 20" seat widths only. Uses 2 power options.</i>	£1950		
Manual Seat to Floor Height. <i>From seat pan.</i>	16.5"	17.5"	18.5"
Powered Tilt Seat to Floor Height. <i>From seat pan.</i>	17.5"	18.5"	19.5"
Powered Lift & Tilt Seat to Floor Height. <i>From seat pan.</i>	18"	19"	20"

## PAEDIATRIC BACK SUPPORTS & POSTIONING *Compatible with paediatric seat dimensions.*

Paediatric Fixed Tool Adj. with Synergy Back Support. <i>Adj. range 95°-125°. Not available with 15"-20" seat widths.</i>				Standard
Paediatric Fold-forward with Synergy Back Support. <i>Release levers. Not suitable for occupied transport. Not available with 15"-20" seat widths.</i>				£100
ADI Aluminium Back. <i>Skin protection, 4 way stretch, Coolcore cover.</i>				POA
Back Height. From seat pan.	15"H. <i>10" to 12" &amp; 14" seat widths only.</i>	16"H. <i>10" to 12" &amp; 14" seat widths only.</i>	17"H. <i>13" &amp; 14" seat widths only.</i>	
	18"H. <i>13" &amp; 14" seat widths only.</i>	19"H. <i>14" seat width only.</i>	20"H. <i>14" seat width only.</i>	
	21"H. <i>14" seat width only.</i>	22"H. <i>14" seat width only.</i>	23"H. <i>14" seat width only.</i>	
	24"H. <i>14" seat width only.</i>	25"H. <i>14" seat width only.</i>		

## STANDARD BACK SUPPORT POSITIONING *Compatible with standard seat dimensions.*

Manual Angle Settable, Fold-Forward. <i>Adj. range 95° to 125°. Will not fold with 10" lift. Not available with 10"-14" seat widths.</i>	Standard
Attendant Controlled Manual Recline. <i>Adjustment range 95° to 125°. Not available with 10" lift or 10"-14" seat widths. Not suitable for occupied transportation.</i>	£150
Powered Back Recline. <i>Adj. range 95° to 125°. Will not fold with 10" lift. Not available with 10"-14" seat widths. Uses 1 power option.</i>	£400

Not available with paediatric seat dimensions.

## STANDARD BACK SUPPORTS *Compatible with standard seat dimensions.*

Tension Adjustable. <i>Air-mesh cover. Adjustment range 19" to 23".</i>	Standard
TRU-Comfort 2. <i>4-way stretch cover. Adjustable range 24.5" to 28" from seat pan.</i>	£170
ADI Aluminium Back. <i>Skin protection, 4 way stretch, Coolcore cover.</i>	POA

Not available with paediatric seat dimensions.

## STEALTH PRO BACK SUPPORTS *Compatible with standard seat dimensions.*

Standard, Coolcore Cover. <i>Solid back, 9 cm contour.</i>				£227
Active, Coolcore Cover. <i>Solid back, 12 cm contour.</i>				£341
Extra Deep Active, Coolcore Cover. <i>Solid back, 17 cm contour.</i>				£448
Active Flex, Coolcore Cover. <i>Tension adjustable straps. 12 cm contour.</i>				£559
Back Height	50 cm. <i>Actual height.</i>	Standard	57 cm. <i>Actual height.</i>	FOC

The back support width will match the seat width selected in seat dimensions. Not available with paediatric seat dimensions.

## ARM SUPPORTS

TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan. Minimum seat depth 15".</i>				Standard
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>				£279
Single Post Height Adjustable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i>				£100
Paediatric Single Post Height Adjustable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i>				£100
Arm Support Pads	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC

## LEG SUPPORTS

70° Swing-Away. <i>Extension range 12" to 15". Minimum seat width 15".</i>	Standard	
60° to 90° Tool Adjustable Swing-Away. <i>6"W x 8"D foot supports. Extension range 11" to 19.5". Min. seat width 15".</i>	£330	
Manually Adjustable Centre Mount Foot Platform. <i>5"W x11"L Footplates. Extension range 10" to 19". Min. seat width 15".</i>	£340	
Manual Elevating Swing-Away. <i>Extension range 15" to 20". Min. seat width 15"</i>	£425	
Power Elevating Articulating Swing-Away. <i>5"W x 7"D foot supports. Ext. range 9" to 21". Min. seat width 15". Uses 2 power options</i>	£850	
Powered Leg Support Operation Method.	Independent. <i>Uses 2 power options.</i>	Combined. <i>Uses 1 power option.</i>
Power Elevating Articulating Centre Mount Foot Platform. <i>5"W x11"L Footplates. Extension range 10" to 19". Min. seat width 15"</i>	£950	

① Enter seat pan to heel measurement.

*If NO measurement is supplied, standard factory length option will be fitted.*

Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away legrest hanger for mounting.</i>	LEFT	£300	Limb Support	£300	RIGHT
--	------	------	--------------	------	-------

Extension ranges above are from seat pan to foot supports. Individual leg supports supplied with heel loops. Elevating leg supports supplied with standard calf supports.

## STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£230
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	£230
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i>	£342
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£92
Push Handles Head Mount Adapter OBAS**. <i>Required to mount head support with TA Back or Active Flex Pro Back.</i>	£266

Please visit [stealthproducts.com](https://stealthproducts.com) for alternative head supports.

## TRUNK SOLUTIONS

Thoracic Lateral Supports.	LEFT	£175	Lateral Thoracic Swing-away Hardware.				£175	RIGHT
		£42	3" H x 5" L Curved Lateral Pad.				£42	
		£42	4" H x 6" L Curved Lateral Pad.				£42	
4-Point Harness.	Extra Small	£143	Classic	STD*	Dynamic	STD*	Vinyl Face	FOC
	Small	£143	Contoured	FOC	Structured	FOC	Zipper Face**	FOC
	Medium	£143						
	Large	£143						
	Extra Large	£143						
Chest Strap.	Extra Small	£91	Small	£91	Side Release		Standard	
	Medium	£91	Large	£91	Velcro Release		£20	
	Extra Large	£91						

\*\*Classic Only.

\*STD = Standard option.

## PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.								Standard (STD)
Retractable Lap Belt, Push Button Release.								£169
Stealth Positioning 2-Point Padded Pelvic Belt.	Extra Small	£84	Small	£84	Single Pull	STD*	Side Release	STD*
	Medium	£84	Large	£84	Dual Pull	FOC	Push Button	FOC
	Extra Large	£84			Rear Pull	FOC		
Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull.	Extra Small		£143	Small	£143	Push Button	STD*	
	Medium		£143	Large	£143	Side Release	£20	

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative lateral, pelvic & lower extremity solutions

\*STD = Standard option

## NON-EXPANDABLE CONTROLS

QLNES Controls. <i>With no power options.</i>		Standard
QLNES Controls. <i>2 power option ready.</i>		£125
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	£100
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

## EXPANDABLE CONTROLS

Q-Logic 3e Controls. <i>No power options.</i>		£300
Q-Logic 3e Controls. <i>3 power option ready.</i>		£350
Q-Logic 3e Controls. <i>5 power option ready.</i>		£500
Q-Logic 3 EX Controls. <i>No power options.</i>		£250
Q-Logic 3 EX Controls. <i>3 power option ready.</i>		£375
Q-Logic 3 EX Controls. <i>5 power option ready.</i>		£550
Q-Logic 3 Enhanced Display. <i>Includes 2 egg switches &amp; flexible display mount. Not available with 3e controls.</i>		£1074
SCIM Module Kit. <i>Allows Q-Logic 3 EX joystick display to be used for speciality controls. Not available with 3e controls.</i>		£316
iAccess OLED 4 Button Switch Module. <i>Mounted to opposite armrest unless specified.</i>		£300
Control Method.	Standard Toggle Lever	If iAccess is selected choose either standard or Toggle lever option.
	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Omit Hand Control. <i>Often selected when ED is used as master remote with special device.</i>	-£100
	Attendant Control. <i>Hand control mounted at rear.</i>	£100
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative control options such as mini joysticks, head arrays & switches.

## HAND CONTROL MOUNTS

Fixed.	Standard				
Basic swing-away.	£75				
Advanced swing-away. <i>Chain driven. Swings in &amp; out with height, tension, and pivot adjustment.</i>	£100				
Mounting Position.	Left	FOC	Right	Standard	
	Standard Carrot	STD	T-Bar	£35	Extended I-Handle £35
Hand Control Joystick Handle.	Mushroom	£35	Medium Ball	£35	3" Goal Post £50
	Chin Cup	£35	Soft Ball	£35	4" Goal Post £50

## SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program.	Standard
4 Mph Program with Lights and indicators.	£450
6 Mph Program with Lights and indicators.	£450
Lights & Indicator Surcharge. <i>Required if no power options have been selected.</i>	£260
LED Fender Lights. <i>Includes single switch button mounted to side rail. Not road legal.</i>	Standard
Omit LED Fender Lights.	FOC

## ACCESSORIES

Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379
Quantum Back Pack.	£110
Additional XLR USB Charger Port.	£17
Magnetic to XLR Charger Adapter. <i>Only with Q-Logic 3e.</i>	£49

## NOTES & SPECIAL REQUESTS

## PRICE

### EDGE 3 STRETTO WITH TB FLEX SEATING



BASE PRICE

**£5,750**

COST OPTIONS

TOTAL

5



**Quantum®**  
32 Wedgwood Road, Bicester,  
Oxon, OX26 4UL  
Phone: 01869 324600  
Fax: 01869 323070  
Email: sales@quantumrehab.co.uk  
Web: www.quantumrehab.co.uk  
**RETAIL UK - 01/04/25 V1.1**

