

## EDGE 3 STRETTO WITH TB FLEX SEATING

### DEALERS DETAILS

Dealers Name \_\_\_\_\_

Order Date \_\_\_\_\_

Address \_\_\_\_\_

Postcode \_\_\_\_\_

Phone \_\_\_\_\_ Email \_\_\_\_\_

PO No \_\_\_\_\_ Order Ref \_\_\_\_\_

### DELIVERY DETAILS

Dealers Address \_\_\_\_\_

Postcode \_\_\_\_\_

Quantum Specialist \_\_\_\_\_

Ⓢ Clients Weight (required for suspension settings) \_\_\_\_\_

### BASE MODEL

**£5,675**

- SRS Suspension
- 12.5" Flat Free Drive Tyres
- Height Adjustable Removable Armrests
- 70° Swing-Away Legrests
- Angle Adjustable Flip Forward Backrest
- 4 Mph (6 Km/h) Motors

- 14 Colour Options
- Max User Weight: 21 st 6lbs. (136 kg)
- 2 x 40ah batteries
- NE+ Control System
- 2x5" Front Castors, Flat Free
- 2x6" Rear Castors, Flat Free



Standard = Standard Factory Feature  
 FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 12.5" Flat Free Black Tyres.	Standard
3" x 12.5" Pneumatic Black Tyres.	FOC
3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	£50
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	£50
4" x 14" Pneumatic Black Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£220
4" x 14" Flat Free Black Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£320

### COLOURS

Back in Black	Standard	Lemon Crush*	FOC
Green Machine	FOC	Night Train*	FOC
Midnight Blue	FOC	Evolution Orange*	FOC
Electric Blue	FOC	Military Green*	FOC
Blue Jean	FOC	Grey Street*	FOC
Deep Purple	FOC	Raspberry Beret*	FOC
Red Corvette	FOC	Custom Hydro Dip*	£600
Pretty in Pink	FOC		

\* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



## PAEDIATRIC SEAT DIMENSIONS

10" Width	FOC	10" Depth	FOC
11" Width	FOC	11" Depth	FOC
12" Width	FOC	12" Depth	FOC
13" Width	FOC	13" Depth	FOC
14" Width	FOC	14" Depth	FOC

## STANDARD SEAT DIMENSIONS

15" Width	FOC	15" Depth	FOC
16" Width	FOC	16" Depth	FOC
17" Width	FOC	17" Depth	FOC
18" Width	Standard	18" Depth	Standard
19" Width	FOC	19" Depth	FOC
20" Width	FOC	20" Depth	FOC

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	£100
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£115
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£115
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£135
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i>	£135
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

## SEAT POSITIONING

Manual Settable Tilt. <i>0° to 30° tool adjustable with a telescopic tube.</i>	Standard		
Powered Seat Tilt. <i>0° to 30°. Uses 1 power option.</i>	£300		
Powered Seat Lift & Tilt. <i>10" lift &amp; 50° tilt. 15" to 20" seat widths only. Uses 2 power options.</i>	£1,895.00		
Manual Seat to Floor Height. <i>From seat pan.</i>	16.5"	17.5"	18.5"
Powered Tilt Seat to Floor Height. <i>From seat pan.</i>	17.5"	18.5"	19.5"
Powered Lift & Tilt Seat to Floor Height. <i>From seat pan.</i>	18"	19"	20"

## PAEDIATRIC BACK SUPPORTS & POSTIONING

Paediatric Fixed Tool Adj. with Synergy Back Support. <i>Adj. range 95°-125°. Not available with 15"-20" seat widths.</i>	Standard		
Paediatric Fold-forward with Synergy Back Support. <i>Release levers. Not suitable for occupied transport. Not available with 15"-20" seat widths.</i>	£100		
ADI Aluminium Back. <i>Skin protection, 4 way stretch, Coolcore cover.</i>	POA		
Omit Back Support. <i>Supplied with back canes and brace bar only.</i>	-£25		
Back Height. From seat pan.	15"H. <i>10" to 12" &amp; 14" seat widths only.</i>	16"H. <i>10" to 12" &amp; 14" seat widths only.</i>	17"H. <i>13" &amp; 14" seat widths only.</i>
	18"H. <i>13" &amp; 14" seat widths only.</i>	19"H. <i>14" seat width only.</i>	20"H. <i>14" seat width only.</i>
	21"H. <i>14" seat width only.</i>	22"H. <i>14" seat width only.</i>	23"H. <i>14" seat width only.</i>
	24"H. <i>14" seat width only.</i>	25"H. <i>14" seat width only.</i>	

## STANDARD BACK SUPPORT POSITIONING

Manual Angle Settable, Fold-Forward. <i>Adj. range 95° to 125°. Will not fold with 10" lift. Not available with 10"-14" seat widths.</i>	Standard
Attendant Controlled Manual Recline. <i>Adjustment range 95° to 125°. Not available with 10" lift or 10"-14" seat widths. Not suitable for occupied transportation.</i>	£150
Powered Back Recline. <i>Adj. range 95° to 125°. Will not fold with 10" lift. Not available with 10"-14" seat widths. Uses 1 power option.</i>	£200

Not available with paediatric seat dimensions.

## STANDARD BACK SUPPORTS

Tension Adjustable. <i>Air-mesh cover. Adjustment range 19" to 23".</i>	Standard
TRU-Comfort 2. <i>4-way stretch cover. Adjustable range 24.5" to 28" from seat pan.</i>	£170
ADI Aluminium Back. <i>Skin protection, 4 way stretch, Coolcore cover.</i>	POA
Omit Back Support. <i>Supplied with back canes and brace bar only.</i>	-£25

Not available with paediatric seat dimensions.

## STEALTH PRO BACK SUPPORTS

Standard, Coolcore Cover. <i>Solid back, 9 cm contour.</i>	£227
Active, Coolcore Cover. <i>Solid back, 12 cm contour.</i>	£341
Extra Deep Active, Coolcore Cover. <i>Solid back, 17 cm contour.</i>	£448
Active Flex, Coolcore Cover. <i>Tension adjustable straps. 12 cm contour.</i>	£559

Back Height	50 cm. <i>Actual height.</i>	Standard	57 cm. <i>Actual height.</i>	FOC
-------------	------------------------------	----------	------------------------------	-----

The back support width will match the seat width selected in seat dimensions. Not available with paediatric seat dimensions.

## ARM SUPPORTS

Std Flip-up Height Adjustable. <i>Includes contour arm pads. Do not select arm pads below. Minimum seat depth 15".</i>	Standard
TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan. Minimum seat depth 15".</i>	FOC
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>	£279
Single Post Height Adjustable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i>	£100
Paediatric Single Post Height Adjustable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i>	£100

Arm Support Pads	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC

## LEG SUPPORTS

70° Swing-away Legrests. <i>Adjustment range 12" to 15". Minimum seat width 15".</i>	Standard
90° Paediatric Swing-away Legrests. <i>Adjustment range 7" to 11".</i>	£150
60° to 90° Tool Adjustable Swing-away Legrests. <i>Adjustment range 11" to 16.5". Min. seat width 15".</i>	£330
Manual Elevating Swing-Away Legrests. <i>Inc. calf pads. Adjustment range 15" to 20". Min. seat width 15".</i>	£250
Power Elevating Swing-away Legrests. <i>Inc. calf pads. Adjustment range 14.5" to 19.5". Min. seat width 16". Uses 2 power options.</i>	£650
Power Elevating Articulating Swing-Away Legrests. <i>Inc. calf pads. Adj. range 15" to 20". Min. seat width 15". Uses 2 power options.</i>	£1300
Manually Adjustable Centre Mount Foot Platform. <i>Inc. calf pads. Adj. range 10" to 16". 5"W x 11"L Footplates. Min. seat width 15".</i>	£340
Paediatric Adj. Centre Mount Foot Platform. <i>Inc. calf pads. Adj. range 8"-12". 5"W x 8"L Footplates. 12"-14" Seat widths only.</i>	£340
Power Elevating Articulating Centre Mount Foot Platform. <i>Inc. calf pads. Adj. range 10" to 16". Min. seat width 16".</i>	£642
Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away legrest hanger for mounting.</i>	£300

Limb Support	<b>LEFT</b>	<b>RIGHT</b>
Seat pan to heel measurement		(Leave blank for standard factory)

Adjustment ranges above are from seat pan to footplates. Individual legrests are supplied with heel loops.

## STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£219
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	£219
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i>	£327
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£88
Push Handles Head Mount Adapter OBAS**. <i>Required to mount head support with TA Back or Active Flex Pro Back.</i>	£266

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative head supports.

## TRUNK SOLUTIONS

Thoracic Lateral Supports.	LEFT	£150	Lateral Thoracic Swing-away Hardware.		£150	RIGHT		
		£35	3" H x 5" L Curved Lateral Pad.		£35			
		£35	4" H x 6" L Curved Lateral Pad.		£35			
4-Point Harness.	Extra Small	£137	Classic	STD*	Dynamic	STD*	Vinyl Face	FOC
	Small	£137		Contoured		FOC		Structured
	Medium	£137						
	Large	£137						
	Extra Large	£137						
Chest Strap.	Extra Small	£87	Small	£87	Side Release		Standard	
	Medium	£87	Large	£87	Velcro Release		£19	
	Extra Large	£87						

\*\*Classic Only.

\*STD = Standard option.

## PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.				Standard (STD)					
Retractable Lap Belt, Push Button Release.				£156					
Stealth Positioning 2-Point Padded Pelvic Belt.	Extra Small	£78	Small	£78	Single Pull	STD*	Side Release	STD*	
	Medium	£78	Large	£78		Dual Pull		FOC	Push Button
	Extra Large	£78			Rear Pull	FOC			
Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull.	Extra Small		£120	Small		£120	Push Button		STD*
	Medium		£120	Large		£120	Side Release		£20

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative lateral, pelvic & lower extremity solutions

\*STD = Standard option

## NON-EXPANDABLE CONTROLS

NE + Controls. <i>With no power options.</i>		Standard
NE + Controls. <i>2 power option ready.</i>		£125
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	£100
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

## EXPANDABLE CONTROLS

Q-Logic 3e Controls. <i>No power options.</i>		£225
Q-Logic 3e Controls. <i>3 power option ready.</i>		£300
Q-Logic 3e Controls. <i>5 power option ready.</i>		£500
Q-Logic 3 EX Controls. <i>No power options.</i>		£275
Q-Logic 3 EX Controls. <i>3 power option ready.</i>		£350
Q-Logic 3 EX Controls. <i>5 power option ready.</i>		£550
Q-Logic 3 Enhanced Display. <i>Includes 2 egg switches &amp; flexible display mount.</i>		£1074
SCIM Module Kit. <i>Allows Q-Logic 3 EX joystick display to be used for speciality controls. Not available with 3e controls.</i>		£316
iAccess OLED 4 Button Switch Module. <i>Mounted to opposite armrest unless specified.</i>		£300
		Standard Toggle Lever
If iAccess is selected choose either standard or Toggle lever option.		
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Omit Hand Control. <i>Often selected when ED is used as master remote with special device.</i>	-£100
	Attendant Control. <i>Hand control mounted at rear.</i>	£100
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative control options such as mini joysticks, head arrays & switches.

## HAND CONTROL MOUNTS

Fixed.						Standard
Swing-away. <i>Swings in &amp; out with height, tension and pivot adjustment.</i>						£75
Retract 4 Swing-away. <i>Full 180° Movement.</i>						£150
Mounting Position.	Left	FOC	Right	Standard		
Hand Control Joystick Handle.	Standard Carrot	STD	T-Bar	£35	Extended I-Handle	£35
	Mushroom	£35	Medium Ball	£35	3" Goal Post	£50
	Chin Cup	£35	Soft Ball	£35	4" Goal Post	£50

## SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
4 Mph Program with Lights and indicators.	£450
6 Mph Program with Lights and indicators.	£450
Lights & Indicator Surcharge. <i>Required if no power options have been selected.</i>	£260
LED Fender Lights. <i>Includes single switch button mounted to side rail. Not road legal.</i>	Standard
Omit LED Fender Lights.	FOC

## ACCESSORIES

Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379
Quantum Back Pack.	£110
Additional XLR USB Charger Port.	£17
Magnetic to XLR Charger Adapter. <i>Only with Q-Logic 3e.</i>	£49

## NOTES & SPECIAL REQUESTS

## PRICE

### EDGE 3 STRETTO WITH TB FLEX SEATING



BASE PRICE **£5,675**

COST OPTIONS

TOTAL

5



**Quantum®**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: sales@quantumrehab.co.uk  
 Web: www.quantumrehab.co.uk  
**RETAIL UK - 27/03/24 V1.1**



**Crash Tested**



V1: 14/06/2022

