

EDGE 3 STRETTO WITH TB4 SEATING

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
	Postcode
Phone	Email
PO No	Order Ref

DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
ⓘ Clients Weight (required for suspension settings)	

BASE MODEL

£14,950

14 Colour Options
TB3 Flip-back Arm Supports
Centre Mount Foot Platform
TB4 Seating with Anterior Tilt
12.5" Flat Free Drive Tyres
2" x 5" Flat Free Front Castors

2" x 6" Flat Free Rear Castors
SRS Suspension System
Q-Logic 3 Control System
2 x 62 ah Batteries
4 Mph (6 Km/h) Motors
Max User Weight: 21 st 6lbs. (136 kg)



Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 12.5" Flat Free Black Tyres.	Standard
3" x 12.5" Pneumatic Black Tyres.	FOC
3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	£50
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	£50

COLOURS

Back in Black	Standard	Lemon Crush*	FOC
Green Machine	FOC	Night Train*	FOC
Midnight Blue	FOC	Evolution Orange*	FOC
Electric Blue	FOC	Military Green*	FOC
Blue Jean	FOC	Grey Street*	FOC
Deep Purple	FOC	Raspberry Beret*	FOC
Red Corvette	FOC	Custom Hydro Dip*	£600
Pretty in Pink	FOC		

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

16"-20"W.

20"-22"W.

Set width at:

SEAT PACKAGES - SEAT DEPTH

14"-16"D.

16"-22"D.

Set depth at:

STEALTH SEAT CUSHIONS

No Seat Cushion Required.	FOC
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	Standard
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	£82
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	£99
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£126
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£126
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£148
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i>	£148
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£209
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£209

Please note that seat cushions may not be available in all width and depth combinations.
Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

TRU-Balance 4 iLevel Powered Seat Lift and Tilt. <i>12" lift & 50° tilt. Up to 3.5 Mph drive speed while elevated. Contains positional feedback.</i>	Standard
Seat to Floor Height. <i>From seat pan.</i>	18" 18.5" 19"

ANTERIOR TILT OPTIONS

10° Anterior Tilt Angle.	Standard
20° Anterior Tilt Angle. <i>Requires a chest bar and knee block assembly.</i>	£600
30° Anterior Tilt Angle. <i>Requires a chest bar and knee block assembly.</i>	£800
Turn-off Anterior Tilt Angle. <i>Only available with 10° anterior tilt through programming.</i>	FOC
Knee Block Assembly. <i>Required with 20° and 30° anterior tilt. Optional with 10° anterior tilt.</i>	£300
Small Chest Bar. <i>14"-20" back widths. Required with 20° and 30° anterior tilt. Optional with 10° anterior tilt.</i>	£250
Medium Chest Bar. <i>17"-22" back widths. Required with 20° and 30° anterior tilt. Optional with 10° anterior tilt.</i>	£250

BACKREST SPECIFICATIONS

TRU-Balance 4 Back Support System. <i>Select TB4 back options below.</i>	Standard
--------------------------------------------------------------------------	----------

Includes a metal active contour back support shell, TRU-Comfort 2 back cushion, and a pair of plastic lateral wedge inserts.

TB4 BACK SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Lateral Wedges. <i>Select size in Back Support Dimensions</i>	Standard
Privacy Flap. <i>12.75"W. Closes gap between back support and seat pan.</i>	£50
Stealth Adapter Interface Plate. <i>SSAPT4. For mounting third party back supports to TB4 back system. Not required for Q-Back 2.</i>	£358
ADI Q-Back 2. <i>Skin protection, 4 way stretch, Coolcore cover. Mounts directly to TB4 system.</i>	POA

Enter ADI Q-Back 2 part number here for quotation.

TB4 BACK SUPPORT POSTIONING

TRU-Balance 4 Powered Recline. *165° with sheer reduction. Uses 1 power option.*

Standard

BACK SUPPORT DIMENSIONS

Back Support can be +/- 4" of seat width with TB4 Back Support System.

Back Support Width Package.		14"-16"W.	16"-20"W.	18"-22"W.
Back Support Width.	14"W.	16"W.	17"W.	18"W.
	19"W.	20"W.	21"W.	22"W.
Back Support Height. <i>From seat pan.</i>	19"H.	20"H.	21"H.	22"H.
	23"H.	24"H.	25"H.	26"H.
	27"H.	28"H.	29"H.	30"H.
	31"H.	32"H.		

ARM SUPPORTS

TB3 Flip-back. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>				Standard
Down Posts for TB3 Flip-back.				£100
Arm Support Pads.	10" Standard.	Standard	10" Falcon Gel.	£40
	14" Standard.	FOC	14" Falcon Gel.	£40
	17" Standard. <i>Minimum seat depth 16"</i>	£50		

CENTRE MOUNTED LEG SUPPORTS

Power Elevating Articulating Dual Actuator Centre Mount Foot Platform. <i>Uses 2 power options.</i>				Standard
Lower Leg Extensions.	10"-16".			Standard
<i>From seat pan.</i>				
Footplates.	10"W x 6"D Split Footplates.	FOC	10"W x 8"D Split Footplates.	FOC
	10"W x 11"D Split Footplates.			Standard
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>			Standard
	Reverse Dartex Calf Pad Covers.			£42
Footplate Accessories.	4" Width Extension Kit. <i>Not available with 2" length kit or 10"W x 6"D footplates.</i>			£16
	2" Length Extension Kit. <i>Not available with 4" width kit or 10"W x 6"D footplates.</i>			£16
	Footplate Connector.	£10	Outer Edge Supports.	£10
	Wheel Bumpers.	£10	Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.	

STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>	£230
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£230
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	£230
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M/ACP150M.</i>	£342
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£92

Please visit stealthproducts.com for alternative head supports.

STEALTH TRUNK SOLUTIONS

UniLink Thoracic Lateral Supports. <i>Supplied as a pair.</i>	Lateral Thoracic Swing-away Hardware. £444		Additional 1.5" Link (each).		£54
	3" H x 5" L Curved Lateral Pad. £84		Additional 2" Link (each).		£54
	4" H x 6" L Curved Lateral Pad. £84		Additional 3" Link (each).		£54
4-Point Harness.	Extra Small. £143	Classic. STD*	Dynamic. STD*	Vinyl Face. FOC	
	Small. £143	Contoured. FOC	Structured. FOC	Zipper Face.** FOC	
	Medium. £143				**Classic Only.
	Large. £143				
	Extra Large. £143				
Chest Strap.	Extra Small. £91	Small. £91	Side Release.	Standard	
	Medium. £91	Large. £91	Velcro Release.	£20	
	Extra Large. £91				

*STD = Standard option.

STEALTH PELVIC & THIGH SOLUTIONS

UniLink Thigh Supports.	LEFT	£185	UniLink Hip Guide Long Removable.	£185	RIGHT
		£185	Unilink Hip Guide Short Removable.	£185	
		£51	4x4 Pelvic Thigh Pad.	£51	
		£51	4x8 Pelvic Thigh Pad.	£51	
		£51	4x12 Pelvic Thigh Pad.	£51	
Knee Adductors.	LEFT	£179	Swing-away Track Mounted.	£179	RIGHT
		£38	3" Round Knee Pad.	£38	
		£38	4" Round Knee Pad.	£38	
Pelvic Belts.	Auto Buckle Lap Belt, Push Button Release.				Standard
	Retractable Lap Belt, Push Button Release.				£169
Stealth Positioning 2-Point Padded Pelvic Belt.	Extra Small. £84	Small. £84	Single Pull. STD*	Side Release. STD*	
	Medium. £84	Large. £84	Dual Pull. FOC	Push Button. FOC	
	Extra Large. £84		Rear Pull. FOC		
Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull.	Extra Small. £143	Small. £143	Push Button. STD*		
	Medium. £143	Large. £143	Side Release. £21		
Medial Supports (Pommels).	Fixed Flip-down Hardware. <i>MTHW-660.</i>		£189	Small Wedge. <i>MTMD-12M.</i>	£108
	Medium Wedge. <i>MTMD-13M.</i>		£108	Large Wedge. <i>MTMD-14M.</i>	£108
	Adjustable Wedge. <i>MTMD-42M.</i>		£209		

Please visit stealthproducts.com for alternative lateral, pelvic & lower extremity solutions. *STD = Standard option

EXPANDABLE CONTROLS

Q-Logic 3e Controls. <i>5 power option ready.</i>	Standard
Q-Logic 3 EX Controls. <i>5 power option ready.</i>	£100
Q-Logic 3 Enhanced Display. <i>Includes 2 egg switches & flexible display mount. Not available with Q-Logic 3e.</i>	£1,074
IR Blaster. <i>Infrared module for Q-Logic 3 EX hand control.</i>	£50
SCIM Module Kit. <i>Allows Q-Logic 3 EX joystick display to be used for speciality controls. Not available with 3e controls.</i>	£316
iAccess OLED 4 Button Switch Module. <i>Mounted to opposite armrest unless specified.</i>	Standard
Push Button iAccess	<i>If iAccess is selected choose either Push Button or Toggle lever option.</i>
Toggle Lever iAccess	
Omit iAccess OLED 4 Button Switch Module.	FOC
Control Method.	
Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
Omit Hand Control. <i>Often selected when ED is used as master remote with special device.</i>	-£100
Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

Please visit stealthproducts.com for alternative control options such as mini joysticks, head arrays & switches.

HAND CONTROL MOUNTS

Fixed.						Standard
Basic swing-away.						£75
Advanced swing-away. <i>Chain driven. Swings in & out with height, tension, and pivot adjustment.</i>						£100
Mounting Position.	Left	FOC	Right	Standard		
Hand Control Joystick Handle.	Standard Carrot	STD	T-Bar	£35	Extended I-Handle £35	
	Mushroom	£35	Medium Ball	£35	3" Goal Post £50	
	Chin Cup	£35	Soft Ball	£35	4" Goal Post £50	

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
4 Mph Program with Lights and indicators.	£450
6 Mph Program with Lights and indicators.	£450
Accu-Trac Motors. <i>Requires Q-Logic 3 EX. Delivers enhanced traction and tracking, especially with switched drive controls. Only available with 6 mph speed program and lights package.</i>	£450
LED Fender Lights. <i>Includes single switch button mounted to side rail. Not road legal.</i>	Standard
Omit LED Fender Lights.	FOC

ACCESSORIES

USB Mobile Device Charger. <i>Please specify mounting position. i.e arm support, backrest, seat rail.</i>	£110
Cup Holder.	£20
Personal Item Hooks.	£15
Cane/Crutch Holder.	£50
Mobile Phone Holder.	£40
Back Pack Holders.	£30
Clothes Guards.	£100
Glove Box.	£30

ACCESSORIES CONTINUED

Hydration System.	£125
Oxygen Holder.	£50
Quick Release Transfer Bars.	£125
Stealth Tablet Holder.	£342
Push Handles.	£50
Accessory Bar.	£50
Back-up Camera.	£300
Quantum Back Pack.	£110
Magnetic to XLR Charger Adapter. <i>Only with Q-Logic 3e.</i>	£49
Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379
Wooden Trim Kit for Tray. <i>Only with swing-away tray table.</i>	£100

NOTES & SPECIAL REQUESTS

PRICE

EDGE 3 STRETTO WITH TB4 SEATING



BASE PRICE

£14,950

COST OPTIONS

TOTAL



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
RETAIL UK - 01/04/25 V1

