

EDGE 3 WITH TB3 SEATING

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
	Postcode
Phone	Email
PO No	Order Ref

DELIVERY DETAILS

Dealers Address	
	Postcode
Quantum Specialist	
Ⓢ Clients Weight	

BASE MODEL

£6,750

- 14 Colour Options
- Height Adjustable Removable Armrests
- Centre Mount Foot Platform
- TB3 Seating
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- SRS Suspension System
- QLNES Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	Standard
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC
4" x 14" Pneumatic Black Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£220
4" x 14" Flat Free Black Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£320

COLOURS

Back in Black	Standard	Lemon Crush*	FOC
Green Machine	FOC	Night Train*	FOC
Midnight Blue	FOC	Evolution Orange*	FOC
Electric Blue	FOC	Military Green*	FOC
Blue Jean	FOC	Grey Street*	FOC
Deep Purple	FOC	Raspberry Beret*	FOC
Red Corvette	FOC	Custom Hydro Dip*	£600
Pretty in Pink	FOC		

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

12"-16"W.
14"-16"W.
16"-20"W
20"-22"W

Set width at:

SEAT PACKAGES - SEAT DEPTH

12"-16"D.
14"-16"D.
16"-20"D
20"-22"D

Set depth at:

Minimum seat size of 14"W x 14"D required for manually adjustable and powered centre mount foot platform when iLevel is selected.

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	Standard
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	£82
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	£99
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£126
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£126
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£148
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i>	£148
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£209
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£209

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

Static Seat.				Standard
TRU-Balance 3 Powered Seat Tilt. <i>0° to 50° COG tilt. Uses 1 power option.</i>				£950
iLevel Powered Seat Lift & Tilt. <i>12" lift & 50° tilt. Up to 4.5 Mph drive speed while elevated. Uses 2 power options.</i>				£2250
Static Seat to Floor Height. <i>From seat pan.</i>	16.5"	17.5"	18.5"	
Powered Tilt Seat to Floor Height. <i>From seat pan.</i>	17.5"	18.5"	19.5"	
Powered iLevel Lift & Tilt Seat to Floor Height. <i>From seat pan.</i>	17.5"	18.5"	19.5"	

BACKREST SPECIFICATIONS

TRU-Balance 3 Back Support System. <i>Select TB3 back options below. Ignore all 'Back Cane' options.</i>	Standard
Back Canes Only System. <i>Not available with Power recline or TB3 flip-back arm supports. Move directly to 'Back Cane Only Support Options' if selected.</i>	FOC

Select one option above. Both options above include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a set of wedge inserts.

TB3 BACK SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i>	Standard
Privacy Flap. <i>12.75"W. Closes gap between back support and seat pan.</i>	£50
Stealth Adapter Interface Plate. <i>SSAPT4. For mounting third party back supports to TB3 back system. Not required for Q-Back 2.</i>	£358
ADI Q-Back 2. <i>Skin protection, 4 way stretch, Coolcore cover. Mounts directly to TB3 system.</i>	POA

Enter ADI Q-Back 2 part number here for quotation.

TB3 BACK SUPPORT POSTIONING

TRU-Balance 3 Static Back. <i>Angle can be preset manually.</i>	Standard
TRU-Balance 3 Flip-forward Back. <i>Not available with TB3 flip-back arm supports or powered recline.</i>	£150
TRU-Balance 3 Powered Recline. <i>175° with sheer reduction. Uses 1 power option.</i>	£700

BACK CANE ONLY SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions.</i>	Standard
Privacy Flap. <i>12.75"W. Closes gap between back support and seat pan.</i>	£50
ADI Back Support. <i>Skin protection, 4 way stretch, Coolcore cover. Mounts to canes with cane clamps.</i>	POA

Enter ADI Back Support part number here for quotation

Back Cane Only Support Positioning.	Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i>					Standard
	Flip-forward Back Canes. <i>Not suitable for occupied transport or with Endo mounts.</i>					£200
	2" Endo Cane Mounts. <i>1" each side. Reduces back width by 2". Minimum seat width 12".</i>					£40
	4" Endo Cane Mounts. <i>2" each side. Reduces back width by 4". Minimum seat width 14".</i>					£40
Back Cane Height. <i>From seat pan.</i>	16"H.	17"H.	18"H.	19"H.	20"H.	
	21"H.	22"H.	23"H.	24"H.	25"H.	

BACK SUPPORT DIMENSIONS

Compatible with TB3 and Back Canes options above. Ignore Back Support Dimensions if an ADI back support is selected in TB3 or Back Canes sections. Back Support can be +/- 4" of seat width with TB3 Back Support System, but must be same as seat width with Back Canes Only System unless -2" or -4" Endo Mounts have been selected.

Backrest Support Width Package.	12"-14"W.	14"-16"W.	16"-20"W.
	18"-22"W.	20"-24"W.	
Back Support Width.	12"W.	14"W.	16"W.
	18"W.	19"W.	20"W.
	22"W.	24"W.	21"W.
Back Support Height. <i>From seat pan.</i>	14"H. <i>12"W back only.</i>	15"H. <i>12"W back only.</i>	16"H.
	18"H.	19"H.	20"H.
	22"H.	23"H.	24"H.
	26"H.	27"H.	28"H.
	30"H.	31"H.	32"H.

19" minimum backrest height when powered recline is selected.

ARM SUPPORTS

2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>			Standard	
Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i>			FOC	
Single Post Drop-in Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i>			FOC	
TB3 Flip-back. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>			£350	
Down Posts for TB3 Flip-back. <i>Only with TB3 Flip-back arm supports. Minimum seat/back width 16".</i>			£100	
Arm Support Pads.	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	17" Standard. <i>TB3 only. Minimum seat depth 16".</i>	£50	Falcon Gel Mounting Plates. <i>Required for 2-post, Single Post & Paediatric arm supports</i>	£40
	10" Falcon Gel	£40		
	14" Falcon Gel	£40		

CENTRE MOUNTED LEG SUPPORTS

Base Mounted Flip-up Foot Platform. <i>Mounted to chassis. Will not move with powered lift or tilt.</i>					FOC		
Manually Adjustable Centre Mount Foot Platform. <i>Minimum seat width 14".</i>					Standard		
Paediatric Adjustable Centre Mount Foot Platform. <i>12"-14" seat widths only.</i>					FOC		
Power Elevating Articulating Centre Mount Foot Platform. <i>Minimum seat width 14". Uses 1 power option.</i>					£950		
Power Elevating Articulating & Lower To Floor Centre Mount Foot Platform. <i>Minimum seat width 14". Lowers to the floor to aid standing transfers. Uses 2 power options.</i>					£1300		
Lower Leg Extensions. <i>From seat pan.</i>		8"-12". <i>Paediatric centre mount only.</i>	FOC	10"-16". *	FOC	13"-19". *	STD
Footplates.	10"W x 6"D Split Footplates. <i>Paediatric centre mount only.</i>					FOC	
	10"W x 8"D Split Footplates.		FOC	10"W x 11"D Split Footplates.		Standard	
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>		Standard	}	Calf Panel. <i>18"W x 8"H.</i>		£42
	Reverse Dartex Calf Pad Covers.		£42		Reverse Dartex Calf Panel Cover.		£62
Footplate Accessories.	4" Width Extension Kit. <i>Not available with 2" length kit or 10"W x 6"D footplates.</i>					£16	
	2" Length Extension Kit. <i>Not available with 4" width kit or 10"W x 6"D footplates.</i>					£16	
	Footplate Connector.		£10	Outer Edge Supports.		£10	
	Wheel Bumpers.		£10				

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.

* Manually adjustable only.

INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

60° to 90° Tool Adjustable Swing-away Leg Supports. <i>Select lowers and footplates below. Minimum seat width 15".</i>							£110
60° TO 90°	Lower Leg Extensions. <i>From seat pan.</i>	11"-13"	FOC	13"-16.5"	Standard	15.5"-19.5"	FOC
	Angle & Depth Adj. Foot Support Size.	4"W x 6"D	FOC	5"W x 6"D	FOC	5"W x 8"D	FOC
		6"W x 8"D	STD	6"W x 11"D	FOC	8"W x 11"D	FOC
	Accessories.	Wheel Bumpers					
Manual Elevating Swing-Away Leg Supports. <i>Flip-up footplates. Adj. range 15" to 20". Min. seat width 15".</i>							£700
Power Elevating Articulating Swing-Away Leg Supports. <i>Includes calf pads.</i>							£850
ELEVATING ARTICULATING	Lower Leg Extensions. <i>From seat pan.</i>	9"-15"		FOC	13"-21"		Standard
	Footplate Type.	Angle Adjustable Flip-up. <i>Select size below.</i>					Standard
	Foot Support Size. <i>Includes heel loops.</i>	4"W x 6"D. <i>14" seat widths and above.</i>	FOC	5"W x 7"D. <i>16" seat widths and above.</i>	FOC	6"W x 8"D. <i>18" seat widths and above.</i>	FOC
	Accessories.	Wheel Bumpers		£10	Gel Padded Knee Adductors		£60
LIMB SUPPORT	Powered Leg Support Operation Method.	Independent. <i>Uses 2 power options.</i>			Combined. <i>Uses 1 power option.</i>		
	Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away legrest hanger for mounting.</i>	LEFT	£300	Limb Support	£300	RIGHT	

Certain configurations may result in reduced ground clearance below foot supports. Contact Quantum Technical for information.

STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>	£230
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£230
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	£230
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i>	£342
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£92

Please visit stealthproducts.com for alternative head supports.

STEALTH TRUNK SOLUTIONS

UniLink Thoracic Lateral Supports. <i>Supplied as a pair.</i>	Lateral Thoracic Swing-away Hardware.		£444	Additional 1.5" Link (each).		£52		
	3" H x 5" L Curved Lateral Pad.		£84	Additional 2" Link (each).		£52		
	4" H x 6" L Curved Lateral Pad.		£84	Additional 3" Link (each).		£52		
4-Point Harness.	Extra Small	£143	Classic	STD*	Dynamic	STD*	Vinyl Face	FOC
	Small	£143			Structured	FOC	Zipper Face**	FOC
	Medium	£143						
	Large	£143						
	Extra Large	£143						
Chest Strap.	Extra Small	£91	Small	£91	Side Release		Standard	
	Medium	£91	Large	£91	Velcro Release		£20	
	Extra Large	£91						

*STD = Standard option **Classic Only

STEALTH PELVIC & THIGH SOLUTIONS

UniLink Thigh Supports.	LEFT	£185	UniLink Hip Guide Long Removable.	£185	RIGHT		
		£185	Unilink Hip Guide Short Removable.	£185			
		£51	4x4 Pelvic Thigh Pad.	£51			
		£51	4x8 Pelvic Thigh Pad.	£51			
		£51	4x12 Pelvic Thigh Pad.	£51			
Knee Adductors.	LEFT	£163	Fixed Cane Mounted.	£163	RIGHT		
		£179	Swing-Away Track Mounted.	£179			
		£38	3" Round Knee Pad.	£38			
		£38	4" Round Knee Pad.	£38			
Pelvic Belts.	Auto Buckle Lap Belt, Push Button Release.				Standard		
	Retractable Lap Belt, Push Button Release.				£169		
Positioning 2-Point Padded Pelvic Belt.	Extra Small	£84	}Single Pull	Standard	}Side Release	Standard	
	Small	£84		FOC		Push Button	FOC
	Medium	£84		FOC			
	Large	£84					
	Extra Large	£84					
Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull.	Extra Small	£143	Small	£143	}Push Button	Standard	
	Medium	£143	Large	£143		Side Release	£21
Medial Supports (Pommels).	Fixed Flip-down Hardware. <i>MTHW-660.</i>		£189	Small Wedge. <i>MTMD-12M.</i>		£108	
	Medium Wedge. <i>MTMD-13M.</i>		£108	Large Wedge. <i>MTMD-14M.</i>		£108	
	Adjustable Wedge. <i>MTMD-42M.</i>		£209	Please visit stealthproducts.com for alternative lateral walking & lower extremity solutions.			

Please visit stealthproducts.com for alternative lateral, pelvic & lower extremity solutions.

NON-EXPANDABLE CONTROLS

QLNES Controls. <i>With no power options.</i>		Standard
QLNES Controls. <i>2 power option ready.</i>		£125
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	£100
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

EXPANDABLE CONTROLS

Q-Logic 3e and 3 EX hand controls are supplied with built-in Bluetooth. Q-Logic 3 Enhanced Display is supplied with built-in Bluetooth and Infrared.

Q-Logic 3e Controls. <i>No power options.</i>	£300
Q-Logic 3e Controls. <i>2 power option ready.</i>	£300
Q-Logic 3e Controls. <i>5 power option ready.</i>	£500
Q-Logic 3 EX Controls. <i>No power options.</i>	£350
Q-Logic 3 EX Controls. <i>2 power option ready.</i>	£375
Q-Logic 3 EX Controls. <i>5 power option ready.</i>	£550
Q-Logic 3 Enhanced Display. <i>Includes 2 egg switches & flexible display mount.</i>	£1074
IR Blaster. <i>Infrared module for Q-Logic 3 EX hand control.</i>	£50
SCIM Module Kit. <i>Allows Q-Logic 3 EX joystick display to be used for speciality controls. Not available with 3e controls.</i>	£316

iAccess OLED 4 Button Switch Module.
Mounted to opposite armrest unless specified.

£300

Push Button iAccess
Toggle Lever iAccess

If iAccess is selected choose either
Push Button or Toggle lever option.

Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Omit Hand Control. <i>Often selected when ED is used as master remote with special device.</i>	-£100
	Attendant Control. <i>Hand control mounted at rear.</i>	£125
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

Please visit stealthproducts.com for alternative control options such as mini joysticks, head arrays & switches.

HAND CONTROL MOUNTS

Fixed.						Standard
Swing-away.						£75
Advanced Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i>						£100
Mounting Position.	Left	FOC	Right	Standard		
Hand Control Joystick Handle.	Standard Carrot	Standard	T-Bar	£35	Extended I-Handle £35	
	Mushroom	£35	Medium Ball	£35	3" Goal Post £50	
	Chin Cup	£35	Soft Ball	£35	4" Goal Post £50	

SPEED PROGRAM & LIGHTS PACKAGE

2-Pole Motors & 4 Mph Program. <i>Lights and indicators not included.</i>	Standard
2-Pole Motors & 4 Mph Program with Lights and indicators.	£450
2-Pole Motors & 6 Mph Program with Lights and indicators.	£450
4-Pole Motor Package & 7 Mph Program with Lights and indicators.	£900
Lights & Indicator Surcharge. <i>Required if no power options have been selected. Not available with 3e controls.</i>	£260
Accu-Trac Motors. <i>Requires Q-Logic 3 EX. Delivers enhanced traction and tracking, especially with switched drive controls. Only available with 6 mph speed program and lights package.</i>	£450
LED Fender Lights. <i>Includes single switch button mounted to side rail. Not road legal.</i>	Standard

ACCESSORIES

USB Mobile Device Charger. <i>Please specify mounting position. i.e arm support, backrest, seat rail. Not available with VR2 or R-Net.</i>	£110
Cup Holder. <i>Only with TB3 flip-back arm supports.</i>	£20
Personal Item Hooks. <i>Not with back cane option.</i>	£15
Mobile Phone Holder. <i>Only with TB3 arm supports.</i>	£40
Back Pack Holders. <i>Not with back cane option.</i>	£30
Cane/Crutch Holder. <i>Requires Accessory Bar.</i>	£50
Clothes Guards. <i>Only with TB3 arm supports.</i>	£100
Glove Box.	£30
Hydration System.	£125
Oxygen Holder.	£50
Quick Release Transfer Bars.	£125
Stealth Tablet Holder.	£342
Push Handles. <i>Not with back cane option.</i>	£50
Accessory Bar.	£50
Back-up Camera.	£300
Quantum Back Pack.	£110
Additional XLR USB Charger Port.	£17
Magnetic to XLR Charger Adapter. <i>Only with Q-Logic 3e.</i>	£49
Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379
Prepare Chassis for Dahl Docking System. <i>Dahl system not included.</i>	£100

NOTES & SPECIAL REQUESTS

PRICE

EDGE 3 WITH TB3 SEATING



BASE PRICE

£6,750

COST OPTIONS

TOTAL

7



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
RETAIL UK - 01/04/25 V1

