

**Quantum®**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: sales@quantumrehab.co.uk  
 Web: www.quantumrehab.co.uk



**EDGE 3**  
**Stretto**



April 2022 <sup>(v1.0)</sup>



**Crash Tested**



### Edge 3 Stretto with TB Flex Seating

Dealer Name:	Order Date:
Address:	
Post Code:	
Phone	Email:
PO No:	Order Ref:
Delivery Address:	
Post Code:	
Quantum Specialist:	

### Standard Features

- 17 Colour Options
- Max User Weight: 21 st 6lbs. (136 kg)
- 2 x 40ah batteries
- NE+ Control System
- 2" x 5" Front Castors, Flat Free
- 2" x 6" Rear Castors, Flat Free
- SRS Suspension
- 12.5" Flat Free Drive Tyres
- Height Adjustable Removable Armrests
- 70° Swing-Away Legrests
- Angle Adjustable Flip Forward Backrest
- 4 Mph (6 Km/h) Motors

**Standard = Standard Factory Feature | FOC = Select instead of Standard Factory Feature**

### Base Model

**Edge 3 Stretto with TB Flex Seating - £5,140.00**

### Drive Wheels

- |  |          |
|--|----------|
| <input type="checkbox"/> 2" x 12.5" " Flat Free, Black Tyres | Standard |
| <input type="checkbox"/> 2" x 12.5" Pneumatic, Black Tyres   | FOC      |
| <input type="checkbox"/> 3" x 14" Flat Free, Black Tyres     | £50.00   |
| <input type="checkbox"/> 3" x 14" Pneumatic, Black Tyres     | £50.00   |

### Castor Wheels

- |  |          |
|--|----------|
| <input checked="" type="checkbox"/> 2" x 5" Flat Free, Black Front Tyres | Standard |
| <input checked="" type="checkbox"/> 2" x 6" Flat Free, Black, Rear Tyres | Standard |

### Colours



*\* This colour choice may increase lead time – Please check with Q-rep before ordering  
 Matching back shroud and joystick shroud included where applicable.*

## Seat Dimensions

### Seat Width

- 15"       16"       17"       18"  
 19"       20"

### Seat Depth

- 15"       16"       17"       18"  
 19"       20"

Please note that seat cushions may not be available in all width and depth combinations. Please refer to our Stealth Cushion Order Forms for alternative models and full size range.

## Seat Cushions

Please refer to our Stealth Cushion Order form for available sizes and enter your preferred size here.....W x .....D

- |   |          |
|---|----------|
| <input type="checkbox"/> No Seat Cushion Required   | FOC      |
| <input type="checkbox"/> Stealth Simplicity G: General use, Fluid resistant Stealth-Tek cover - 25 stone weight limit                                       | Standard |
| <input type="checkbox"/> Stealth TRU-Comfort 2 SPP: Skin protection, positioning, 4-way stretch Coolcore cover - 32 stone weight limit                      | £70.00   |
| <input type="checkbox"/> Stealth Glacial SP: Skin protection, 4-way stretch, Coolcore cover - 79.3kg to 181.4kg weight limit                                | £100.00  |
| <input type="checkbox"/> Stealth Solution® SPP: Skin protection, positioning, 4-way stretch, Coolcore cover and urethane liner - 25 stone weight limit      | £100.00  |
| <input type="checkbox"/> Stealth Essence® SPP: Skin protection, positioning, 4-way stretch and Coolcore cover - 35 stone 7lbs weight limit                  | £120.00  |
| <input type="checkbox"/> Stealth Spectrum® Gel SPP: Skin protection, positioning, 4-way stretch, Coolcore cover and urethane liner - 25 stone weight limit  | £170.00  |
| <input type="checkbox"/> Stealth Spectrum® Foam SPP: Skin protection, positioning, 4-way stretch, Coolcore cover and urethane liner - 42 stone weight limit | £170.00  |
| <input type="checkbox"/> Stealth Premiere P: Positioning, 4-way stretch, Coolcore cover - 158.7kg weight limit  | £120.00  |
| <input type="checkbox"/> Stealth Zen SP: Skin Protection, 4-way stretch, Coolcore cover - 79.3 to 181.4 kg weight limit                                     | £60.00   |

## Seat Positioning Adjustments

- |   |          |
|---|----------|
| <input type="checkbox"/> Manual Tilt 0° to 30° Tool adjustable with a telescopic tube   | Standard |
| <input type="checkbox"/> TB-Flex® Power Tilt 0° to 30° Full Drive Lock at 20° when combined with Power Recline. Uses 1 power option | £300.00  |

### Seat to Floor Height

- 16.5"                                       17.5"                                       18.25"

## Backrest Options

- |  |          |
|--|----------|
| <input type="checkbox"/> Tension Adjustable Backrest Air-mesh cover. Range of adjustment 19" to 23" from seat board          | Standard |
| <input type="checkbox"/> TRU-Comfort 2 Backrest 4-way stretch cover. Range of height adjustable 24.5" to 28" from seat board | £170.00  |
| <input type="checkbox"/> ADI Backrest Skin protection, 4 way stretch, Coolcore cover   | POA      |
| <input type="checkbox"/> Omit Backrest Supplied with back canes and brace bar only   | - £25.00 |

## Backrest Positioning Adjustments

- |  |          |
|--|----------|
| <input type="checkbox"/> Manual Settable / Flip Forward Range of adjustments 95° to 125°   | Standard |
| <input type="checkbox"/> Attendant Controlled Manual Recline Range of adjustments 95° to 125°. Not suitable for occupied transportation. | £150.00  |
| <input type="checkbox"/> Powered Recline Range of adjustments 95° to 125°. Uses 1 power option.  | £200.00  |

## Armrests

- |  |          |
|--|----------|
| <input type="checkbox"/> TB Flex Removable Armrest with Quick Release Range of adjustments 9" to 14" | Standard |
| <input type="checkbox"/> TB Flex Removable Armrest with Armrest Lever Range of adjustments 9" to 14" | £85.00   |
| <input type="checkbox"/> 2-Post Flip Up Removable Armrest Range of adjustments 9" to 13"             | £280.00  |

### Armpads

- |                                       |                                       |   |   |
|---------------------------------------|---------------------------------------|---|---|
| <input type="checkbox"/> 10" Standard | <input type="checkbox"/> 14" Standard | <input type="checkbox"/> 10" Waterfall<br><small>Only available with 2-post flip up</small> | <input type="checkbox"/> 14" Waterfall<br><small>Only available with 2-post flip up</small> |
|---------------------------------------|---------------------------------------|---|---|

## Non-Expandable Controls

- |   |          |
|---|----------|
| <input type="checkbox"/> NE + Controls 2 power option ready | Standard |
| <input type="checkbox"/> VR2 Controls 2 power option ready  | FOC      |

## Expandable Controls

<input type="checkbox"/> Q-Logic 3 Controls 2 power option ready	£350.00
<input type="checkbox"/> Q-Logic 3 Controls 5 power option ready	£550.00
<input type="checkbox"/> R-Net Controls 2 power option ready	£350.00
<input type="checkbox"/> R-Net Controls 5 power option ready	£550.00
<input type="checkbox"/> Q-Logic 3 Enhanced Display Includes x2 egg switches and flexible display mount. Does not include Q-Logic 3 joystick.	£1,074.00
<input type="checkbox"/> R-Net Omni 2 Display Includes x2 egg switches and display mount. Does not include R-Net joystick.	£1,074.00
<input type="checkbox"/> SCIM Module Kit Allows Q-Logic 3 joystick display to be used for speciality controls.	£316.00
<input type="checkbox"/> iAccess OLED 4 Button Switch Module Mounted to opposite armrest unless specified	£300.00
<input type="checkbox"/> Push Button iAccess Module    or <input type="checkbox"/> Toggle Lever iAccess Module	

## Joystick Mounting Options

<input type="checkbox"/> Fixed	Standard
<input type="checkbox"/> Retract 4 Swing-away Joystick Arm Kit (Full 180° Movement)	£150.00
<input type="checkbox"/> Swing-Away	£75.00
<input type="checkbox"/> Dual Control	£430.00
<input type="checkbox"/> Attendant Control Only	£100.00

### Mounting Position

Left                      or                       Right

## Legrests

<input type="checkbox"/> 70° Swing-away Legrests Range of adjustments 12" to 15" from seat pan to heel	Standard
<input type="checkbox"/> Manually Adjustable Centre Mount Foot Platform Mounted to Seat. Range of adjustments 10" to 16" from seat pan to heel	£340.00
<input type="checkbox"/> Quantum 60° to 90° Adjustable Swing-away Legrests <input type="checkbox"/> 11"-13"      OR <input type="checkbox"/> 13.5"-16.5" from seat pan to heel	£330.00
<input type="checkbox"/> Manual Elevating Swing-Away Legrests Range of adjustments 15" to 20" from seat pan to heel	£250.00
<input type="checkbox"/> Powered Elevating Swing-Away Legrests Range of adjustments 15" to 20" from seat pan to heel. Uses 2 power options	£650.00

## Headrests

<input type="checkbox"/> Stealth 10" Comfort Plus CP250	£140.00
<input type="checkbox"/> Stealth 14" Comfort Plus CP150	£140.00
<input type="checkbox"/> Stealth 14" Adjustable Comfort Plus ACP150	£250.00

Please visit [www.stealthproducts.com](http://www.stealthproducts.com) for more headrest options

## Lateral and Trunk Supports

### Stealth Swing Away Thoracic Lateral Supports

#### Left Lateral Package

<input type="checkbox"/> Lateral Thoracic Thin Hardware with Mount	£150.00
<b>Left Lateral Pads</b> Choose one option only	
<input type="checkbox"/> 3" H x 5" L Curved Lateral Pad	£35.00
<input type="checkbox"/> 4" H x 6" L Curved Lateral Pad	£35.00

#### Right Lateral Package

<input type="checkbox"/> Lateral Thoracic Thin Hardware with Mount	£150.00
<b>Right Lateral Pads</b> Choose one option only	
<input type="checkbox"/> 3" H x 5" L Curved Lateral Pad	£35.00
<input type="checkbox"/> 4" H x 6" L Curved Lateral Pad	£35.00

## Pelvic Supports

<input type="checkbox"/> Auto Buckle Lap Belt, Push Button Release	Standard
<b>i-Fit Pelvic Belts</b> Supplied un-installed	
<input type="checkbox"/> Padded Offers a large padded surface area	£70.00

## Pelvic Supports Continued

Non-Padded *Narrowed section for smaller support surfaces. Not available with 4-point* £65.00

### Size

Extra Small *Not available with non-padded*       Small       Medium       Large       Extra Large *Not available with non-padded or 4-point.*

### Attachment Method

2-Point       4-Point *Single pull only*       Push Button       Side Release

### Adjustment

Single Pull       Dual Pull

## Speed Program and Lights Package

4 Mph program. *Lights and indicators not included* Standard

4 Mph program with Lights and indicators £450.00

6 Mph program with Lights and indicators £450.00

Lights and Indicator Surcharge *If no power options have been selected, you must tick this box in addition to the lights and indicator package to receive the correct lighting module* £260.00

LED Fender Lights *Includes single switch button for operation, mounted to side rail. Not road legal. Does not include indicators* Standard

Omit LED Fender Lights FOC

## Accessories

Swing-Away Tray Table *Supplied with relevant joystick cut-out, or no cut out if attendant control only* £350.00

## Price

<b>Sub Total:</b>	£5,140.00
<b>Options:</b>	
<b>Total:</b>	

### Notes

1. All measurements given are approximate.
2. By faxing or mailing this order you accept all responsibility for the configuration ordered.
3. Errors made in filling out this order, may result in further charges or delays.
4. If this order has a deadline you must advise us ASAP. We will try to deliver for this date, but cannot accept liability if we are unable to fulfil this order by your deadline.
5. Information relating to crash testing is correct at time of printing.  
Pride Mobility constantly reviews and updates this information, if you require further information about crash testing, please contact your Quantum Specialist