



RETAIL

ESSENCE SPP CUSHION

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
	Postcode
Email	
Phone	Fax

DELIVERY DETAILS

Dealers Address		
Postcode		
PO No		
Order Type	Quote	Order

DETAILS

This cushion offers skin protection and positioning with a Coolcore, 4-way stretch, breathable cover, a urethane liner, and an ABS flat rigidiser.

Low profile Tru-Comfort 2, lateral and pre-ischial shelf. Contoured, HD molded foam base, coccyx and sacral relief.



Weight capacities
Up to 20" - 158.75 kg
22" and Above - 204.11 kg
Overall weight: 1.5 kg - 2.2 kg
Overall height 3.5"
15x15" up to 22x24"
Odd width sizes in 15" to 21"
Custom sizes available

CUSHION

		ADD WIDTH		ADD DEPTH	Width											
					10	12	14	15	16	17	18	19	20	21	22	
Essence SPP	No Wedge. 15"-20"W.	ST-ESSSPP			10											
	No Wedge. 22"W.	ST-ESSSPP			12											
	Small Wedge. 15"-20"W.	ST-ESSSPP			14											
	Small Wedge. 22"W.	ST-ESSSPP			15											
	Med Wedge. 15"-20"W.	ST-ESSSPP			16											
	Medium Wedge. 22"W.	ST-ESSSPP			17											
	Large Wedge. 15"-20"W.	ST-ESSSPP			18											
	Large Wedge. 22"W.	ST-ESSSPP			19											
Extra Covers And Liners.	Coolcore Cover.	ESSSPP-NBLK			20											
	Mesh Cover.	ESSSPP-MBLK			21											
	Polyurethane Liner.	PUL1709080			22											
					24											

Select your size, then complete the Part Number by adding the Width first then Depth.

ACM CONTOURED RIGIDISER FOR USE WITH SLING SEATS ONLY.

15"W x 16"D. CUS142467-1616. £72	17"W x 18"D. CUS142466-1818. £72	19"W x 18"D. CUS142465-2018. £72
15"W x 17"D. CUS142467-1616. £72	17"W x 19"D. CUS142466-1818. £72	19"W x 19"D. CUS142465-2018. £72
15"W x 18"D. CUS142467. £72	17"W x 20"D. CUS142466. £72	19"W x 20"D. CUS142466. £72
15"W x 19"D. CUS142467. £72	17"W x 21"D. CUS142466. £72	19"W x 21"D. CUS142465. £72
15"W x 20"D. CUS142467. £72	17"W x 22"D. CUS142466. £72	19"W x 22"D. CUS142465. £72
16"W x 16"D. CUS142467. £72	18"W x 16"D. CUS142466-1816. £72	20"W x 18"D. CUS142465-2018. £72
16"W x 17"D. CUS142467-1616. £72	18"W x 17"D. CUS142466-1816. £72	20"W x 19"D. CUS142465-2018. £72
16"W x 18"D. CUS142467. £72	18"W x 18"D. CUS142466-1818. £72	20"W x 20"D. CUS142466. £72
16"W x 19"D. CUS142467. £72	18"W x 19"D. CUS142466-1818. £72	20"W x 21"D. CUS142465. £72
16"W x 20"D. CUS142467. £72	18"W x 20"D. CUS142466. £72	20"W x 22"D. CUS142465. £72
17"W x 16"D. CUS142467-1616. £72	18"W x 21"D. CUS142466. £72	
17"W x 17"D. CUS142467-1616. £72	18"W x 22"D. CUS142466. £72	

NOTES & SPECIAL REQUESTS

Special Order Essence SPP Cushion.

£450

Please detail your special order cushion request.

All special order cushion requests must be approved with the Quantum Custom Department.

Please contact Quantum Inside Sales or Custom Department with questions. Extended lead time will apply.

Cushion part numbers involved will be the cushion part number listed in previous sections modified to meet special order request stated in details below. Drawings are available for sketch details.


FRONT VIEW



RIGHT SIDE VIEW



TOP VIEW



PRICE

TOTAL

2



Stealth Products
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
ESSENCE SPP - 01/04/25 V1

