



RETAIL

GLACIAL SP CUSHION

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Email _____

Phone _____ Fax _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____


PO No _____

Order Type _____ Quote _____ Order _____

DETAILS

This cushion offers skin protection and a Coolcore 4-way stretch, breathable cover.

New SIM G base, fluid resistant overlay, top and sides.
TPU heat-sealed, low profile lightweight cooling gel topper



Weight capacities
158.75 kg: 15"-21" width,
181.43 kg: 22" width+
Overall weight: 1.5kg
Overall height 2.5"
10x10" up to 24x24"
Odd width sizes in 15" up to 21"
Custom sizing not available

CUSHION

		ADD WIDTH	ADD DEPTH	
Glacial SP. 10"-20" W.	ST-GLSP	<input type="text"/>	<input type="text"/>	£209
Glacial SP. 22"-24" W	ST-GLSP	<input type="text"/>	<input type="text"/>	£231
Extra Covers And Liners.	Coolcore Cover.	GLSP-NBLK	<input type="text"/>	£77
	Mesh Cover.	GLSP-MBLK	<input type="text"/>	£77
	Stealth Tek Cover.	GLSP-TBLK	<input type="text"/>	£77
	Urethane Liner.	URL1709070	<input type="text"/>	£44
	Polyurethane Liner.	PUL1709030	<input type="text"/>	£72

	Width												
	10	12	14	15	16	17	18	19	20	21	22	24	
10													
12													
14													
15													
16													
17													
18													
19													
20													
21													
22													
23													
24													

Select your size, then complete the Part Number by adding the Width first then Depth.

ADD SOLID BASE

ABS FLAT RIGIDISER

12"W x 12"D. <i>DWR1010E189C.</i>	£39	17"W x 20"D. <i>DWR1010E234H.</i>	£39	22"W x 18"D. <i>DWR1010E194I.</i>	£39
12"W x 14"D. <i>DWR1010E189E.</i>	£39	18"W x 16"D. <i>DWR1010E192G.</i>	£39	22"W x 20"D. <i>DWR1010E194K.</i>	£39
14"W x 12"D. <i>DWR1010E190C.</i>	£39	18"W x 18"D. <i>DWR1010E192I.</i>	£39	22"W x 22"D. <i>DWR1010E194M.</i>	£39
14"W x 14"D. <i>DWR1010E190E.</i>	£39	18"W x 20"D. <i>DWR1010E192K.</i>	£39	22"W x 24"D. <i>DWR1010E193O.</i>	£39
14"W x 16"D. <i>DWR1010E190G.</i>	£39	18"W x 22"D. <i>DWR1010E192M.</i>	£39	24"W x 18"D. <i>DWR1010E236I.</i>	£39
15"W x 15"D. <i>MEC1804152-1515.</i>	£39	19"W x 17"D. <i>MEC1804152-1917.</i>	£39	24"W x 20"D. <i>DWR1010E236K.</i>	£39
15"W x 16"D. <i>MEC1804152-1515.</i>	£39	19"W x 18"D. <i>MEC1804152-1917.</i>	£39	24"W x 21"D. <i>DWR1010E236L.</i>	£39
15"W x 17"D. <i>DWR1010E232H.</i>	£39	19"W x 19"D. <i>MEC1804152-1919.</i>	£39	24"W x 22"D. <i>DWR1010E236M.</i>	£39
16"W x 16"D. <i>DWR1010E191G.</i>	£39	19"W x 20"D. <i>MEC1804152-1919.</i>	£39	24"W x 23"D. <i>DWR1010E236N.</i>	£39
16"W x 18"D. <i>DWR1010E191I.</i>	£39	20"W x 16"D. <i>DWR1010E193G.</i>	£39	24"W x 24"D. <i>DWR1010E236O.</i>	£39
16"W x 20"D. <i>DWR1010E191K.</i>	£39	20"W x 18"D. <i>DWR1010E193I.</i>	£39		
17"W x 17"D. <i>DWR1010E234H.</i>	£39	20"W x 20"D. <i>DWR1010E193K.</i>	£39		
17"W x 18"D. <i>DWR1010E234H.</i>	£39	20"W x 22"D. <i>DWR1010E193M.</i>	£39		
17"W x 19"D. <i>DWR1010E234J.</i>	£39	22"W x 16"D. <i>DWR1010E193G.</i>	£39		

ACM CONTOURED RIGIDISER FOR USE WITH SLING SEATS ONLY.

16"W x 16"D. <i>CUS142467-1616.</i>	£72	17"W x 20"D. <i>CUS142466.</i>	£72	19"W x 19"D. <i>CUS142465-2018.</i>	£72
16"W x 18"D. <i>CUS142467.</i>	£72	18"W x 16"D. <i>CUS142466-1816.</i>	£72	19"W x 20"D. <i>CUS142465.</i>	£72
16"W x 20"D. <i>CUS142467.</i>	£72	18"W x 18"D. <i>CUS142466-1818.</i>	£72	20"W x 18"D. <i>CUS142465-2018.</i>	£72
17"W x 17"D. <i>CUS142466-1816.</i>	£72	18"W x 20"D. <i>CUS142466.</i>	£72	20"W x 20"D. <i>CUS142465.</i>	£72
17"W x 18"D. <i>CUS142466-1818.</i>	£72	18"W x 22"D. <i>CUS142466.</i>	£72	20"W x 22"D. <i>CUS142465.</i>	£72
17"W x 19"D. <i>CUS142466-1818.</i>	£72	19"W x 18"D. <i>CUS142465-2018.</i>	£72		

NOTES & SPECIAL REQUESTS

PRICE

TOTAL

2



Stealth Products
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
GLACIAL SP - 01/04/25 V1

