

## LIGHTNING WITH SYNERGY SEATING

### DEALERS DETAILS

Dealers Name \_\_\_\_\_

Order Date \_\_\_\_\_

Address \_\_\_\_\_

Postcode \_\_\_\_\_

Phone \_\_\_\_\_ Email \_\_\_\_\_

PO No \_\_\_\_\_ Order Ref \_\_\_\_\_

### DELIVERY DETAILS

Dealers Address \_\_\_\_\_

Postcode \_\_\_\_\_

Quantum Specialist \_\_\_\_\_

Ⓞ Clients Weight \_\_\_\_\_

### BASE MODEL

**£4,995**

- 10 Colour Options
- Height Adjustable Removable Armrests
- 70° Swing-away Leg Supports
- Angle Adjustable Back Support
- Simplicity G Seat Cushion
- 14" Pneumatic Drive Tyres



- 2" x 9" Pneumatic Front Castors
- ATX Suspension System
- NE+ Control System
- 2 x 70 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature  
 FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 14" Pneumatic Grey Tyres.	Standard
3" x 14" Flat Free Grey Tyres.	£150
4" x 14" Pneumatic Grey Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£220
4" x 14" Flat Free Grey Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£320

### FRONT CASTORS

2" x 9" Pneumatic Grey Tyres.	Standard
2" x 9" Flat Free Grey Tyres.	£60

### COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Night Train	FOC	Evolution Orange	FOC
Electric Blue	FOC	Grey Street	FOC
Midnight Blue	FOC	Green Machine	FOC
Deep Purple	FOC	Custom Hydro Dip*	£600
Red Corvette	FOC		

\* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



## SEAT PACKAGES - SEAT WIDTH

14"-18"W.

16"-20"W

20"-22"W

Set width at:

## SEAT PACKAGES - SEAT DEPTH

14"-18"D.

15"-20"D

20"-22"D

Set depth at:

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	Standard
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	£75
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	£100
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£115
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£115
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£135
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i>	£135
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

## SEAT POSITIONING

Static Seat.				Standard
Static Seat to Floor Height. <i>From seat pan.</i>	15.875"	16.875"	17.875"	

## BACKREST SPECIFICATIONS

Back Canes System.						Standard
Back Cane Height.	16"H.	17"H.	18"H.	19"H.	20"H.	
<i>From seat pan.</i>	21"H.	22"H.	23"H.	24"H.	25"H.	

## BACK SUPPORT

Active Contour Back Support with TRU-Comfort 4 Back Cushion and Wedge Inserts.				Standard	
Tension Adjustable Back Upholstery (TAB). <i>May increase lead-time.</i>				£120	
Back Support Height.	16"H.	17"H.	18"H.		
<i>From seat pan. Back support height must be similar to back cane height. Ignore section if using an ADI back support.</i>	19"H.	20"H.	21"H.		
	22"H.	23"H.	24"H.		
	25"H.				
ADI Back Support. <i>Skin protection, 4 way stretch, Coolcore cover. Mounts to back canes with cane clamps.</i>				POA	
Enter ADI Q-Back 2 part number here for quotation.					
Privacy Flap. <i>12.75"W. Closes gap between back support and seat pan. Not available with a TAB.</i>				£50	

## BACK SUPPORT POSITIONING

Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i>	Standard
Flip-forward Back Canes. <i>Not suitable for occupied transport or with Endo mounts.</i>	£200
2" Endo Cane Mounts. <i>Not suitable for occupied transport. 1" each side. Reduces back width by 2". Minimum seat width 12".</i>	£40
4" Endo Cane Mounts. <i>Not suitable for occupied transport 2" each side. Reduces back width by 4". Minimum seat width 14".</i>	£40

## ARM SUPPORTS

Single Post Drop-in Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i>	Standard
Singe Post HD Tool Adjustable Removable. <i>Adjustment range 10"-14" from seat pan. Minimum seat depth 14".</i>	£220
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>	£145
Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i>	FOC
Cantilever Flip-Up. <i>Not available with a TAB. Increases overall width by 7". Height adj. 11"-17". Not available with Flip Forward Backrest.</i>	£280
Arm Support Pads.	10" Standard Standard 14" Standard FOC
	10" Waterfall. FOC 14" Waterfall. FOC
	10" Falcon Gel. <i>Includes mounting plate.</i> £80 14" Falcon Gel. <i>Includes mounting plate.</i> £80

## INDIVIDUALLY MOUNTED LEG SUPPORTS

70° Swing-Away Leg Supports. <i>Flip-up footplates. Adj. range 12" to 15". Minimum seat width 15".</i>	FOC
90° Paediatric Swing-Away Leg Supports. <i>4"W x 6"D angle &amp; depth adj. footplates. Adjustment range 7" to 11".</i>	£110
Manual Elevating Swing-Away Leg Supports. <i>Flip-up footplates. Adj. range 15" to 20". Min. seat width 15".</i>	£130
Power Elevating Swing-Away Leg Supports. <i>Flip-up f/plates. Adj. range 14.5" to 19.5". Min. seat width 16". Uses 2 power options.</i>	£650

## INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

60° to 90°	60° to 90° Tool Adjustable Swing-Away Leg Supports. <i>Select lowers and footplates below. Minimum seat width 15".</i>	£110	
60° TO 90°	Lower Leg Extensions. <i>From seat pan.</i> 11"-13" FOC 13.5"-16.5" Standard 15.5"-19.5" FOC		
	Footplate Type. Flip-up Extruded. <i>8.5"W x 6"D.</i> STD Angle & Depth Adjustable. <i>Select size below.</i> £115 Multi-axis. <i>Select size below.</i> £263		
	Footplate Size. <i>Angle adj. &amp; Multi-axis only.</i> 4"W x 6"D FOC 5"W x 6"D FOC 5"W x 8"D FOC	6"W x 8"D STD 6"W x 11"D FOC 8"W x 11"D FOC	
		Accessories. Wheel Bumpers £10	
	Power Elevating Articulating Swing-Away Leg Supports. <i>Inc. calf pads. Adj. range 15" to 20". Min. seat width 15". Uses 2 power options.</i>	£1250	
POWER ELEVATING ARTICULATING	Lower Leg Extensions. <i>From seat pan.</i> 11"-13" FOC 14.5"-18" Standard 17.5"-21" FOC		
	Footplate Type. Angle & Depth Adjustable. <i>Select size below.</i> Standard		
	Footplate Size. 4"W x 6"D FOC 5"W x 6"D FOC 5"W x 8"D FOC	6"W x 8"D Standard 6"W x 11"D FOC 8"W x 11"D FOC	
		Accessories. Small Calf Pads Standard Wheel Bumpers £10	
	Large Calf Pads FOC Gel Padded Knee Adductors £60		

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information. Adjustment ranges above are from seat pan to footplates. Individual legrests are supplied with heel loops.

Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away leg support hanger for mounting.</i>	LEFT	£300	Limb Support	£300	RIGHT
--	------	------	--------------	------	-------

## STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>	£219
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£219
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	£219
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i>	£327
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£88

## STEALTH TRUNK SOLUTIONS

Thoracic Lateral Supports.	<b>LEFT</b>	£425	Lateral Thoracic Swing-away Hardware.		£425	<b>RIGHT</b>			
		£72	3" H x 5" L Curved Lateral Pad.		£72				
		£72	4" H x 6" L Curved Lateral Pad.		£72				
Thoracic Lateral Support Extras.	Additional 1.5" Link (each).					£52			
	Additional 2" Link (each).					£52			
	Additional 3" Link (each).					£52			
4-Point Harness.	Extra Small	£137	Classic	STD*	Dynamic	STD*	Vinyl Face	FOC	
	Small	£137		Contoured		FOC		Structured	FOC
	Medium	£137							
	Large	£137							
	Extra Large	£137							
Chest Strap.	Extra Small	£87	Small	£87	Side Release		Standard		
	Medium	£87	Large	£87	Velcro Release		£19		
	Extra Large								£87

\*STD = Standard option \*\*Classic Only

## STEALTH PELVIC & THIGH SOLUTIONS

UniLink Thigh Supports.	<b>LEFT</b>	£176	UniLink Hip Guide Long Removable.		£176	<b>RIGHT</b>		
		£176	Unilink Hip Guide Short Removable.		£176			
		£48	4x4 Pelvic Thigh Pad.		£48			
		£48	4x8 Pelvic Thigh Pad.		£48			
		£48	4x12 Pelvic Thigh Pad.		£48			
Knee Adductors.	<b>LEFT</b>	£156	Fixed Cane Mounted.		£156	<b>RIGHT</b>		
		£170	Swing-Away Track Mounted.		£170			
		£36	3" Round Knee Pad.		£36			
		£36	4" Round Knee Pad.		£36			
Pelvic Belts.	Auto Buckle Lap Belt, Push Button Release.					Standard		
	Retractable Lap Belt, Push Button Release.					£156		
Positioning 2-Point Padded Pelvic Belt.	Extra Small	£78	Single Pull	Standard	Side Release	Standard		
	Small	£78		Rear Pull		FOC	Push Button	FOC
	Medium	£78	Dual Pull		FOC			
	Large	£78						
	Extra Large	£78						
Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull.	Extra Small	£120	Small	£120	Push Button		Standard	
	Medium	£120	Large	£120	Side Release		£20	
Medial Supports (Pommels).	Fixed Flip-down Hardware. <i>MTHW-660.</i>		£180	Small Wedge. <i>MTMD-12.</i>		£103		
	Medium Wedge. <i>MTMD-13.</i>		£103	Large Wedge. <i>MTMD-14.</i>		£103		
	Adjustable Wedge. <i>MTMD-42.</i>		£199					

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative lateral, pelvic & lower extremity solutions.

## NON-EXPANDABLE CONTROLS

NE + Controls. <i>With no power options.</i>	Standard	
NE + Controls. <i>2 power option ready.</i>	£125	
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	£100
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

## EXPANDABLE CONTROLS

Q-Logic 3e and 3 EX hand controls are supplied with built-in Bluetooth. Q-Logic 3 Enhanced Display is supplied with built-in Bluetooth and Infrared.

Q-Logic 3e Controls. <i>No power options.</i>	£225
Q-Logic 3e Controls. <i>2 power option ready.</i>	£300
Q-Logic 3 EX Controls. <i>No power options.</i>	£275
Q-Logic 3 EX Controls. <i>2 power option ready.</i>	£350
Q-Logic 3 Enhanced Display. <i>Includes 2 egg switches &amp; flexible display mount.</i>	£1074
IR Blaster. <i>Infrared module for Q-Logic 3 EX hand control.</i>	£50
SCIM Module Kit. <i>Allows Q-Logic 3 EX joystick display to be used for speciality controls. Not available with 3e controls.</i>	£316

iAccess OLED 4 Button Switch Module.  
*Mounted to opposite armrest unless specified.*

£300

Push Button iAccess  
Toggle Lever iAccess

If iAccess is selected choose either  
Push Button or Toggle lever option.

Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Omit Hand Control. <i>Often selected when ED is used as master remote with special device.</i>	-£100
	Attendant Control. <i>Hand control mounted at rear.</i>	£125
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

## HAND CONTROL MOUNTS

Fixed.	Standard					
Swing-away. <i>Swings in &amp; out with height, tension and pivot adjustment.</i>	£75					
Retract 4 Swing-away. <i>Full 180° Movement.</i>	£150					
Mounting Position.	Left	FOC	Right	Standard		
Hand Control Joystick Handle.	Standard Carrot	Standard	T-Bar	£35	Extended I-Handle	£35
	Mushroom	£35	Medium Ball	£35	3" Goal Post	£50
	Chin Cup	£35	Soft Ball	£35	4" Goal Post	£50

## SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
4 Mph Program with Lights and indicators.	£450
6 Mph Program with Lights and indicators.	£450
Lights & Indicator Surcharge. <i>Required if no power options have been selected.</i>	£260
Accu-Trac Motors. <i>Requires Q-Logic 3 EX. Delivers enhanced traction and tracking, especially with switched drive controls. Only available with 6 mph speed program and lights package.</i>	£450
LED Fender Lights. <i>Includes single switch button mounted to side rail. Not road legal.</i>	£190

## ACCESSORIES

USB Mobile Device Charger. <i>Please specify mounting position. i.e arm support, backrest, seat rail. Not available with VR2 or R-Net.</i>	£110
Glove Box.	£30
Hydration System.	£125
Oxygen Holder.	£50
Stealth Tablet Holder.	£342
Back-up Camera .	£300
Magnetic to XLR Charger Adapter. <i>Only with Q-Logic 3e.</i>	£49
Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379
Wooden Trim Kit for Tray. <i>Only with swing-away tray table.</i>	£100

## NOTES & SPECIAL REQUESTS

## PRICE

**LIGHTNING WITH  
SYNERGY SEATING**



**BASE PRICE** **£4,995**

**COST OPTIONS**

**TOTAL**

6



**Quantum®**  
32 Wedgwood Road, Bicester,  
Oxon, OX26 4UL  
Phone: 01869 324600  
Fax: 01869 323070  
Email: sales@quantumrehab.co.uk  
Web: www.quantumrehab.co.uk  
**RETAIL UK - 30/08/23 V1**



**Crash  
Tested**



V1: 14/06/2022

