



## RETAIL

### PREMIERE P CUSHION

#### DEALERS DETAILS

Dealers Name \_\_\_\_\_

Order Date \_\_\_\_\_

Address \_\_\_\_\_

Postcode \_\_\_\_\_

Email \_\_\_\_\_

Phone \_\_\_\_\_ Fax \_\_\_\_\_

#### DELIVERY DETAILS

Dealers Address \_\_\_\_\_

Postcode \_\_\_\_\_

PO No \_\_\_\_\_

Order Type \_\_\_\_\_ Quote \_\_\_\_\_ Order \_\_\_\_\_

### DETAILS

This cushion offers positioning with a Coolcore 4-way stretch, breathable cover.

Gray X-Soft Sunmate FRG foam topper, high density positioning wedges built in for offloading



Weight capacities  
 79.37 kg: 10"-14" width  
 158.75 kg: 15"-21" width,  
 226.79 kg: 22" width+  
 Overall weight: 1.5kg  
 Overall height 3.5"  
 14x14" up to 22x22"  
 Odd sizes available  
 Custom sizes available

### CUSHION

		ADD WIDTH	ADD DEPTH	
Premiere P. 14"-20"W.		ST-PREMP	<input type="checkbox"/>	£210
Premiere P. 22"W.		ST-PREMP	<input type="checkbox"/>	£230
Cushion Cover.	Coolcore.	PREMP-NBLK	<input type="checkbox"/>	STD
	Mesh.	PREMP-MBLK	<input type="checkbox"/>	FOC
Fluid Proof Liner.	Urethane.	URL1709070	<input type="checkbox"/>	£30
	Polyurethane.	PUL1709030	<input type="checkbox"/>	£48
Extra Covers And Liners.	Coolcore Cover.	PREMP-NBLK	<input type="checkbox"/>	£70
	Mesh Cover.	PREMP-MBLK	<input type="checkbox"/>	£44
	Urethane Liner.	URL1709070	<input type="checkbox"/>	£30
	Polyurethane Liner.	PUL1709030	<input type="checkbox"/>	£48

		Width							
		14	15	16	17	18	19	20	22
Depth	13								
	14								
	15								
	16								
	17								
	18								
	19								
	20								
22									

Select your size, then complete the Part Number by adding the Width first then Depth.

## ADD SOLID BASE

### ABS FLAT RIGIDISER

14"W x 13"D. <i>DWR1010E190C.</i> £34	16"W x 19"D. <i>DWR1010E191I.</i> £34	19"W x 17"D. <i>MEC1804152-1917.</i> £34
14"W x 14"D. <i>DWR1010E190E.</i> £34	16"W x 20"D. <i>DWR1010E191K.</i> £34	19"W x 18"D. <i>MEC1804152-1917.</i> £34
14"W x 15"D. <i>DWR1010E190E.</i> £34	17"W x 16"D. <i>MEC1804152-1715.</i> £34	19"W x 19"D. <i>MEC1804152-1917.</i> £34
14"W x 16"D. <i>DWR1010E190G.</i> £34	17"W x 17"D. <i>DWR1010E234H.</i> £34	20"W x 16"D. <i>DWR1010E193G.</i> £34
15"W x 14"D. <i>DWR1010E190G.</i> £34	17"W x 18"D. <i>DWR1010E234H.</i> £34	20"W x 17"D. <i>DWR1010E193G.</i> £34
15"W x 15"D. <i>MEC1804152-1515.</i> £34	17"W x 19"D. <i>DWR1010E234J.</i> £34	20"W x 18"D. <i>DWR1010E193G.</i> £34
15"W x 16"D. <i>MEC1804152-1515.</i> £34	18"W x 16"D. <i>DWR1010E192G.</i> £34	20"W x 19"D. <i>DWR1010E193I.</i> £34
15"W x 17"D. <i>DWR1010E232H.</i> £34	18"W x 17"D. <i>DWR1010E192G.</i> £34	20"W x 20"D. <i>DWR1010E193K.</i> £34
16"W x 16"D. <i>DWR1010E191G.</i> £34	18"W x 18"D. <i>DWR1010E192G.</i> £34	22"W x 18"D. <i>DWR1010E194I.</i> £34
16"W x 17"D. <i>DWR1010E191G.</i> £34	18"W x 19"D. <i>DWR1010E192I.</i> £34	22"W x 20"D. <i>DWR1010E194K.</i> £34
16"W x 18"D. <i>DWR1010E191I.</i> £34	18"W x 20"D. <i>DWR1010E192K.</i> £34	22"W x 22"D. <i>DWR1010E194M.</i> £34

Available with optional ABS Flat Rigidizer or ACM Rigidizer.

### ACM CONTOURED RIGIDISER SUITABLE FOR SLING SEATS

16"W x 16"D. <i>CUS142467-1616.</i> £56	18"W x 16"D. <i>CUS142466-1816.</i> £56	20"W x 18"D. <i>CUS142465-2018.</i> £56
16"W x 18"D. <i>CUS142467.</i> £56	18"W x 18"D. <i>CUS142466-1818.</i> £56	20"W x 20"D. <i>CUS142465.</i> £56
16"W x 20"D. <i>CUS142467.</i> £56	18"W x 20"D. <i>CUS142466.</i> £56	

For use on manual wheelchairs only.

## NOTES & SPECIAL REQUESTS

Special Order Stealth Premiere P Cushion.

£380

Please detail your special order cushion request.

All special order cushion requests must be approved with the Quantum Custom Department.

Please contact Quantum Inside Sales or Custom Department with questions. Extended lead time will apply.

Cushion part numbers involved will be the cushion part number listed in previous sections modified to meet special order request stated in details below. Drawings are available below for sketch details.


FRONT VIEW



RIGHT SIDE VIEW



TOP VIEW



## PRICE

TOTAL

2



### Stealth Products

32 Wedgwood Road, Bicester,  
Oxon, OX26 4UL

Phone: 01869 324600

Fax: 01869 323070

Email: sales@quantumrehab.co.uk

Web: www.quantumrehab.co.uk

PREMIERE P - 02/11/23 V1



V1: 14/06/2022

