

Q1450 WITH CONTOURED SEATING

DEALERS DETAILS

Dealers Name _____
 Order Date _____
 Address _____
 _____ Postcode _____
 Phone _____ Email _____
 PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____
 _____ Postcode _____
 Quantum Specialist _____
 Ⓢ Clients Weight _____

BASE MODEL

£8,995

- 5 Colour Options
- Height Adjustable Removable Armrests
- Flip-up Centre Mount Foot Platform
- Angle Adjustable Back Support
- Spectrum Foam SPP Seat Cushion
- 3" x 14" Flat Free Drive Tyres
- 2" x 9" Flat Free Rear Castors



- 1.5" Anti-tip Wheels
- Swing-arm Suspension System
- VR2 Control System
- 2 x 70 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:
42 st 8lbs. (272 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Grey Tyres. Standard

REAR CASTORS

2" x 9" Flat Free Grey Tyres. Standard

COLOURS

Back in Black	Standard
Electric Blue	FOC
Red Corvette	FOC
Evolution Orange	FOC
Pretty in Pink	FOC

Please download our free PQ-365 app to view colours and other product features in more detail.



SEATING OPTIONS

Pinchless Hinge High Back Captain's Seat. <i>Includes a head support and contoured seat cushion.</i>	22"W-20"D.	FOC	22"W-22"D.	FOC	22"W-24"D.	Standard
	24"W-20"D.	FOC	24"W-22"D.	FOC	24"W-24"D.	FOC
HD High Back Captain's Seat. <i>Includes a head support, lumbar cut-out, static hinges, contoured seat cushion.</i>			26"W-20"D.	£770	26"W-22"D	£770
Pinchless Hinge High Back Contoured Seat. <i>Includes a head support and solid seat pan.</i>					22"W-20"D	FOC
Special Order Contoured Seat.	Seat type		Width		Depth	POA

STEALTH SEAT CUSHIONS

Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	STD
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may only ordered with the Contoured Pinchless Hinge High Back seat and the solid seat pan option.

SEAT POSITIONING

Static Seat.				Standard
Static Seat to Floor Height. <i>From seat pan.</i>	17.5"	18.5"	19.5"	

ARM SUPPORTS

Bariatric Single Post Tool Height Adjustable Removable. <i>Includes Bariatric full length arm pads.</i>				Standard
Single Post HD Height Adjustable Removable.				£130
HD & Bariatric Height Package. <i>From seat pan.</i>	8"-12"	FOC	10"-14"	Standard
Single Post HD Arm Support Pads.	10" Standard. <i>Not available with Bariatric.</i>			FOC
	14" Standard. <i>Not available with Bariatric.</i>			Standard
	10" Waterfall. <i>Not available with Bariatric.</i>			FOC
	14" Waterfall. <i>Not available with Bariatric.</i>			FOC
Flip-up Arm Supports <i>Includes arm pads.</i>	Full Length, Flip-up, Height Adjustable, Removable. <i>Required with pinchless hinge high back. Not available with HD seat</i>			£130
	Desk Length, Flip-up, Height Adjustable. <i>Removable Required with pinchless hinge high back. Not available with HD seat.</i>			£130

CENTRE MOUNTED LEG SUPPORTS

17"W x 13"D Base Mounted Flip-up Foot Platform. <i>Mounted to chassis. Will not move with power tilt.</i>	Standard
17"W x 16"D Base Mounted Flip-up Foot Platform. <i>Mounted to chassis. Will not move with power tilt.</i>	FOC

INDIVIDUALLY MOUNTED LEG SUPPORTS

70° HD Drop-in Swing-away Legrests.						£340
70° HD Lower Extensions. <i>From seat pan.</i>	8.5"-12".	FOC	13"-18".	Standard	16"-21".	FOC
70° HD Footplates.	6"W x 8"D.	Standard	6"W x 11"D.	FOC	8"W x 11"D.	FOC

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information. Adjustment ranges above are from seat pan to footplates. Individual legrests are supplied with heel loops.

Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away leg support hanger for mounting.</i>	LEFT	£300	Limb Support	£300	RIGHT
--	-------------	------	--------------	------	--------------

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.						Standard
Auto Buckle Lap Belt Length.	60".	Standard	70".	FOC	80".	FOC
	90".	FOC				
Retractable Lap Belt, Push Button Release.						£156
Stealth Positioning 2-Point Padded Pelvic Belt.	Medium	£78	} Single Pull	Standard	} Push Button	Standard
	Large	£78		FOC		
	Extra Large	£78		FOC		

NON-EXPANDABLE CONTROLS

VR2 Controls. <i>With no power options.</i>		Standard
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard

EXPANDABLE CONTROLS

Q-Logic 3e and 3 EX hand controls are supplied with built-in Bluetooth. Q-Logic 3 Enhanced Display is supplied with built-in Bluetooth and Infrared.

Q-Logic 3e Controls. <i>No power options.</i>		£225			
Q-Logic 3 EX Controls. <i>No power options.</i>		£275			
Q-Logic 3 Enhanced Display. <i>Includes 2 egg switches & flexible display mount.</i>		£1074			
IR Blaster. <i>Infrared module for Q-Logic 3 EX hand control.</i>		£50			
SCIM Module Kit. <i>Allows Q-Logic 3 EX joystick display to be used for speciality controls. Not available with 3e controls.</i>		£316			
iAccess OLED 4 Button Switch Module. <i>Mounted to opposite armrest unless specified.</i>	£300	<table border="0"> <tr> <td rowspan="2">} Push Button iAccess</td> <td rowspan="2">} Toggle Lever iAccess</td> <td rowspan="2">} If iAccess is selected choose either Push Button or Toggle lever option.</td> </tr> <tr> </tr> </table>	} Push Button iAccess	} Toggle Lever iAccess	} If iAccess is selected choose either Push Button or Toggle lever option.
} Push Button iAccess	} Toggle Lever iAccess	} If iAccess is selected choose either Push Button or Toggle lever option.			
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard			

HAND CONTROL MOUNTS

Fixed.						Standard
Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i>						£75
Mounting Position.	Left	FOC	Right	Standard		
	Standard Carrot	Standard	T-Bar	£35	Extended I-Handle	£35
Hand Control Joystick Handle.	Mushroom	£35	Medium Ball	£35	3" Goal Post	£50
	Chin Cup	£35	Soft Ball	£35	4" Goal Post	£50

SPEED PROGRAM

4 Mph Program.

Standard

NOTES & SPECIAL REQUESTS

PRICE

**Q1450 WITH
CONTOURED SEATING**



BASE PRICE

£8,995

COST OPTIONS

TOTAL

4



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
RETAIL UK - 13/08/24 V1

