

## Q1450 WITH SYNERGY SEATING

### DEALERS DETAILS

Dealers Name \_\_\_\_\_

Order Date \_\_\_\_\_

Address \_\_\_\_\_

Postcode \_\_\_\_\_

Phone \_\_\_\_\_ Email \_\_\_\_\_

PO No \_\_\_\_\_ Order Ref \_\_\_\_\_

### DELIVERY DETAILS

Dealers Address \_\_\_\_\_

Postcode \_\_\_\_\_

Quantum Specialist \_\_\_\_\_

Ⓞ Clients Weight \_\_\_\_\_

### BASE MODEL

**£10,195**

- 5 Colour Options
- Height Adjustable Removable Armrests
- Flip-up Centre Mount Foot Platform
- Angle Adjustable Back Support
- Spectrum Foam SPP Seat Cushion
- 3" x 14" Flat Free Drive Tyres
- 2" x 9" Flat Free Rear Castors



- 1.5" Anti-tip Wheels
- Swing-arm Suspension System
- VR2 Control System
- 2 x 70 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:  
42 st 8lbs. (272 kg)

Standard = Standard Factory Feature  
FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 14" Flat Free Grey Tyres.

Standard

### REAR CASTORS

2" x 9" Flat Free Grey Tyres.

Standard

### COLOURS

Back in Black	Standard
Electric Blue	FOC
Red Corvette	FOC
Evolution Orange	FOC
Pretty in Pink	FOC

Please download our free PQ-365 app to view colours and other product features in more detail.



## SEAT PACKAGES - SEAT WIDTH

20"-24"W.
24"-26"W
26"-28"W
28"-30"W
30"-32"W
Set width at:

## SEAT PACKAGES - SEAT DEPTH

20"-24"D.
24"-28"D
Set depth at:

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	STD
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

## SEAT POSITIONING

Static Seat.	Standard
HD Powered Seat Tilt. <i>0° to 45°. Uses 1 power option.</i>	£1600
Static Seat to Floor Height. <i>From seat pan.</i>	16.5" 17.5" 18.5"
HD Tilt Seat to Floor Height. <i>From seat pan.</i>	19.5" 20.5"

## SYNERGY BACK SUPPORTS

Synergy Visco Padded Back. <i>Height range 17" to 25" tall from seat pan.</i>	Standard
HD Sling Back. <i>Height range 14" to 25" tall from seat pan.</i>	£150
Enter height preference here.	

## BACK SUPPORT POSITIONING

Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i>	Standard
1" Endo Cane Mounts. <i>Reduces backrest width by 2" compared to seat width.</i>	£40
2" Endo Cane Mounts. <i>Reduces backrest width by 4" compared to seat width.</i>	£40
Back Cane Height.	16"H. 17"H. 18"H. 19"H. 20"H.
<i>From seat pan.</i>	21"H. 22"H. 23"H. 24"H. 25"H.

## ARM SUPPORTS

HD Single Post Tool Height Adjustable Removable.	Standard	
Bariatric Single Post Tool Height Adjustable Removable.	£150	
HD & Bariatric Height Package. <i>From seat pan.</i>	8"-12" FOC 10"-14" Standard	
Cantilever Flip-Up Armrests. <i>Adj. range 11"-17". Not available with HD sling back support. Increases overall width by 7".</i>	£150	
Arm Support Pads.	10" Standard. <i>Not available with Bariatric.</i>	FOC
	14" Standard. <i>Not available with Bariatric.</i>	Standard
	10" Waterfall. <i>Not available with Bariatric.</i>	FOC
	14" Waterfall. <i>Not available with Bariatric.</i>	FOC
	Bariatric. <i>Bariatric arm supports only.</i>	FOC

## CENTRE MOUNTED LEG SUPPORTS

17"W x 13"D Base Mounted Flip-up Foot Platform. <i>Mounted to chassis. Will not move with power tilt. Ignore other centre mount options if selected.</i>		FOC
17"W x 16"D Base Mounted Flip-up Foot Platform. <i>Mounted to chassis. Will not move with power tilt. Ignore other centre mount options if selected.</i>		Standard
Adjustable Centre Mount Foot Platform. <i>Mounted to seat.</i>		£244
Footplates.	17"W x 13"D One Piece Footplate.	FOC
	17"W x 16"D One Piece Footplate.	FOC
Calf Support.	Individual Calf Pads. <i>8"W x 8"H tapering to 6"W</i>	Standard } Calf Panel. <i>18"W x 8"H.</i> £42
	Reverse Dartex Calf Pad Covers. £42	Reverse Dartex Calf Panel Cover. £62

## INDIVIDUALLY MOUNTED LEG SUPPORTS

70° HD Drop-in Swing-away Legrests.		£340				
70° HD Lower Extensions. <i>From seat pan.</i>	8.5"-12".	FOC	13"-18".	Standard	16"-21".	FOC
	70° HD Footplates.	6"W x 8"D.	Standard	6"W x 11"D.	FOC	8"W x 11"D.

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information. Adjustment ranges above are from seat pan to footplates. Individual legrests are supplied with heel loops.

Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away leg support hanger for mounting.</i>	<b>LEFT</b>	£300	Limb Support	£300	<b>RIGHT</b>
--	-------------	------	--------------	------	--------------

## STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>	£219
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£219
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	£219
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i>	£327
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£88

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative head supports.

## PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.		Standard				
Auto Buckle Lap Belt Length.	60".	Standard	70".	FOC	80".	FOC
	90".	FOC				
Retractable Lap Belt, Push Button Release.		£156				
Stealth Positioning 2-Point Padded Pelvic Belt.	Medium	£78	} Single Pull	Standard	} Push Button	Standard
	Large	£78		FOC		
	Extra Large	£78		FOC		

## NON-EXPANDABLE CONTROLS

VR2 Controls. <i>With no power options.</i>	Standard
VR2 Controls. <i>2 power option ready.</i>	£125
Control Method. Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard

## EXPANDABLE CONTROLS

Q-Logic 3e and 3 EX hand controls are supplied with built-in Bluetooth. Q-Logic 3 Enhanced Display is supplied with built-in Bluetooth and Infrared.

Q-Logic 3e Controls. <i>No power options.</i>	£225
Q-Logic 3e Controls. <i>3 power option ready.</i>	£300
Q-Logic 3 EX Controls. <i>No power options.</i>	£275
Q-Logic 3 EX Controls. <i>3 power option ready.</i>	£350
Q-Logic 3 Enhanced Display. <i>Includes 2 egg switches &amp; flexible display mount.</i>	£1074
IR Blaster. <i>Infrared module for Q-Logic 3 EX hand control.</i>	£50
SCIM Module Kit. <i>Allows Q-Logic 3 EX joystick display to be used for speciality controls. Not available with 3e controls.</i>	£316

iAccess OLED 4 Button Switch Module.  
*Mounted to opposite armrest unless specified.*

£300

Push Button iAccess  
Toggle Lever iAccess

If iAccess is selected choose either  
Push Button or Toggle lever option.

Control Method. Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
Omit Hand Control. <i>Often selected when ED is used as master remote with special device.</i>	-£100

## HAND CONTROL MOUNTS

Fixed.	Standard					
Swing-away. <i>Swings in &amp; out with height, tension and pivot adjustment.</i>	£75					
Mounting Position.	Left	FOC	Right	Standard		
Hand Control Joystick Handle.	Standard Carrot	Standard	T-Bar	£35	Extended I-Handle	£35
	Mushroom	£35	Medium Ball	£35	3" Goal Post	£50
	Chin Cup	£35	Soft Ball	£35	4" Goal Post	£50

## SPEED PROGRAM

4 Mph Program.	Standard
----------------	----------

## NOTES & SPECIAL REQUESTS

---

### PRICE

**Q1450 WITH  
SYNERGY SEATING**



<b>BASE PRICE</b>	<b>£10,195</b>
<b>COST OPTIONS</b>	
<b>TOTAL</b>	



**Quantum®**  
32 Wedgwood Road, Bicester,  
Oxon, OX26 4UL  
Phone: 01869 324600  
Fax: 01869 323070  
Email: [sales@quantumrehab.co.uk](mailto:sales@quantumrehab.co.uk)  
Web: [www.quantumrehab.co.uk](http://www.quantumrehab.co.uk)  
**RETAIL UK - 13/08/24 V1**



V1: 14/06/2022

