

Q4 WITH TB FLEX SEATING

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
Postcode	
Phone	Email
PO No	Order Ref

DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Clients Weight	

BASE MODEL

£3,940

7 Colour Options
Height Adjustable Removable Armrests
70° Swing-Away Legrests
Angle Adjustable Fold-Forward Backrest
14" Flat Free Drive Tyres
2" x 6" Flat Free Front Castors



2" x 6" Flat Free Rear Castors
ATX Suspension System
NE+ Control System
2 x 62 ah Batteries
4 Mph (6 Km/h) Motors
Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres.	Standard
3" x 14" Pneumatic Black Tyres.	FOC

COLOURS

Back in Black	Standard	Red Corvette	FOC
Midnight Blue	FOC	Custom Hydro Dip*	£600
Electric Blue	FOC	* Extended lead-time	
Deep Purple	FOC	Please download our free PQ-365 app to view colours and other product features in more detail.	
Pretty in Pink	FOC		
Evolution Orange	FOC		



SEAT WIDTH

15" Width	FOC	16" Width	FOC
17" Width	FOC	18" Width	Standard
19" Width	FOC	20" Width	FOC

SEAT DEPTH

15" Depth	FOC	16" Depth	FOC
17" Depth	FOC	18" Depth	Standard
19" Depth	FOC	20" Depth	FOC

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	£100
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£115
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£115
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£135
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i>	£135
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability. The nearest size will be provided unless otherwise specified

SEAT POSITIONING

Manual Settable Tilt. <i>0° to 30° tool adjustable with a telescopic tube.</i>				Standard
Powered Seat Tilt. <i>0° to 30°. Uses 1 power option.</i>				£300
Manual Seat to Floor Height. <i>From seat pan.</i>	16"	17"	18"	
Powered Tilt Seat to Floor Height. <i>From seat pan.</i>	16"	17" <i>Changes Tilt range to -5° to 25°</i>	18" <i>Changes Tilt range to -10° to 20°</i>	

STANDARD BACK SUPPORT POSITIONING

Manual Angle Settable, Fold-Forward. <i>Adjustment range 95° to 125°.</i>	Standard
Attendant Controlled Manual Recline. <i>Adjustment range 95° to 125°. Not suitable for occupied transportation.</i>	£150
Powered Back Recline. <i>Adjustment range 95° to 125°. Uses 1 power option.</i>	£200

STANDARD BACK SUPPORTS

TRU-Comfort 2. <i>4-way stretch cover. Adjustable range 24.5" to 28" from seat board.</i>	Standard
Tension Adjustable. <i>Air-mesh cover. Adjustment range 19" to 23" from seat board. OBAS** required to mount headrest.</i>	FOC
ADI Aluminium Back. <i>Skin protection, 4 way stretch, Coolcore cover.</i>	POA
Omit Back Support. <i>Supplied with back canes and brace bar only.</i>	-£25
Enter ADI Back Support part number here for quotation.	

STEALTH PRO BACK SUPPORTS

Standard, Coolcore Cover. <i>Solid back, 9 cm contour.</i>				£227
Active, Coolcore Cover. <i>Solid back, 12 cm contour.</i>				£341
Extra Deep Active, Coolcore Cover. <i>Solid back, 17 cm contour.</i>				£448
Active Flex, Coolcore Cover. <i>Tension adjustable straps. 12 cm contour.</i>				£559
Back Height	50 cm.	Actual height.	Standard	57 cm.
				Actual height.
				FOC

The back support width will match the seat width selected in seat dimensions.

ARM SUPPORTS

Std Flip-up Height Adjustable. <i>Includes contour arm pads. Do not select arm pads below.</i>				Standard
TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan.</i>				£85
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>				£279
Arm Support Pads	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC

LEG SUPPORTS

70° Swing-away Legrests. <i>Adjustment range 12" to 15".</i>	Standard
Base Mounted Foot Platform. <i>Mounted to chassis. Will not move with power tilt.</i>	FOC
60° to 90° Tool Adjustable Swing-away Legrests. <i>Adjustment range 11" to 16.5".</i>	£330
Manual Elevating Swing-Away Legrests. <i>Adjustment range 15" to 20".</i>	£250
Power Elevating Swing-away Legrests. <i>Adjustment range 14.5" to 19.5". Min. seat width 16". Uses 2 power options.</i>	£650
Manually Adjustable Centre Mount Foot Platform. <i>Adjustment range range 10" to 16". 5"W x 11"L Footplates.</i>	£340
Power Elevating Articulating Centre Mount Foot Platform. <i>Adj. range 10" to 16". Min. seat width 16". Uses 2 power options.</i>	£642
Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away legrest hanger for mounting.</i>	£300

Limb Support

LEFT

RIGHT

Seat pan to heel measurement

(Leave blank for standard factory)

Adjustment ranges above are from seat pan to footplates. Individual legrests supplied with heel loops. Elevating legrests supplied with calf supports.

STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£219
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	£219
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i>	£327
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£88
Push Handles Head Mount Adapter OBAS**. <i>Required to mount head support with TA Back or Active Flex Pro Back.</i>	£266

Please visit stealthproducts.com for alternative head supports.

TRUNK SOLUTIONS

Thoracic Lateral Supports.	LEFT	£150	Lateral Thoracic Swing-away Hardware.		£150	RIGHT			
		£35	3" H x 5" L Curved Lateral Pad.		£35				
		£35	4" H x 6" L Curved Lateral Pad.		£35				
Stealth 4-Point. Harness	Extra Small	£137	Small	£137	}	Classic	STD*	Contoured	FOC
	Medium	£137	Large	£137		Dynamic	STD*	Structured	FOC
	Extra Large	£137				Vinyl Face	FOC	Zipper Face**	FOC
Chest Strap.	Extra Small	£87	Small	£87	}	Side Release			STD*
	Medium	£87	Large	£87		Velcro Release			£19
	Extra Large	£87							

*STD = Standard option **Classic Only

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.						Standard (STD)			
Retractable Lap Belt, Push Button Release.						£156			
Stealth Positioning 2-Point Padded Pelvic Belt.	Extra Small	£78	Small	£78	}	Single Pull	STD*	Dual Pull	FOC
	Medium	£78	Large	£78		Rear Pull	FOC	Side Release	STD*
	Extra Large	£78				Push Button	FOC		
Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull.			Extra Small	£120	}	Small	£120	Push Button	STD*
			Medium	£120		Large	£120	Side Release	£20

Please visit stealthproducts.com for alternative lateral, pelvic & lower extremity solutions *STD = Standard option

NON-EXPANDABLE CONTROLS

NE + Controls.	FOC
NE + Controls. <i>2 power option ready.</i>	£125
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i> Standard
	Attendant Control. <i>Hand control mounted at rear.</i> £100
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i> £430

HAND CONTROL MOUNTS

Fixed.					Standard	
Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i>					£75	
Retract 4 Swing-away. <i>Full 180° Movement.</i>					£150	
Mounting Position.	Left	FOC	Right	Standard		
Hand Control Joystick Handle.	Standard Carrot	STD	T-Bar	£35	Extended I-Handle	£35
	Mushroom	£35	Medium Ball	£35	3" Goal Post	£50
	Chin Cup	£35	Soft Ball	£35	4" Goal Post	£50

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
4 Mph Program with Lights and indicators.	£450
6 Mph Program with Lights and indicators.	£450
Lights & Indicator Surcharge. <i>Required if no power options have been selected.</i>	£260

ACCESSORIES

Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379
Quantum Back Pack.	£110
XLR USB Charger Port.	£17

NOTES & SPECIAL REQUESTS

PRICE

**Q4 WITH TB
FLEX SEATING**



BASE PRICE **£3,940**

COST OPTIONS

TOTAL

4



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
RETAIL UK - 09/08/23 V1



V1: 14/06/2022



**Crash
Tested**

