

Q4 WITH TB FLEX SEATING

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Clients Weight _____

BASE MODEL

£4,450

- 7 Colour Options
- Height Adjustable Removable Armrests
- 70° Swing-Away Legrests
- Angle Adjustable Fold-Forward Backrest
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- ATX Suspension System
- NE+ Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

| | |
|---------------------------------|----------|
| 3" x 14" Flat Free Black Tyres. | Standard |
| 3" x 14" Pneumatic Black Tyres. | FOC |

COLOURS

| | | | |
|------------------|----------|--|------|
| Back in Black | Standard | Red Corvette | FOC |
| Midnight Blue | FOC | Custom Hydro Dip* | £600 |
| Electric Blue | FOC | * Extended lead-time | |
| Deep Purple | FOC | Please download our free PQ-365 app to view colours and other product features in more detail. | |
| Pretty in Pink | FOC | | |
| Evolution Orange | FOC | | |



SEAT WIDTH

| | | | |
|-----------|-----|-----------|----------|
| 15" Width | FOC | 16" Width | FOC |
| 17" Width | FOC | 18" Width | Standard |
| 19" Width | FOC | 20" Width | FOC |

SEAT DEPTH

| | | | |
|-----------|-----|-----------|----------|
| 15" Depth | FOC | 16" Depth | FOC |
| 17" Depth | FOC | 18" Depth | Standard |
| 19" Depth | FOC | 20" Depth | FOC |

STEALTH SEAT CUSHIONS

| | |
|---|----------|
| No Seat Cushion Required | FOC |
| Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i> | Standard |
| Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i> | FOC |
| TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i> | £100 |
| Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i> | £115 |
| Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i> | £115 |
| Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i> | £135 |
| Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i> | £135 |
| Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i> | £190 |
| Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i> | £190 |

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability. The nearest size will be provided unless otherwise specified

SEAT POSITIONING

| | | | |
|--|----------|--------------------------------------|---------------------------------------|
| Manual Settable Tilt. <i>0° to 30° tool adjustable with a telescopic tube.</i> | Standard | | |
| Powered Seat Tilt. <i>0° to 30°. Uses 1 power option.</i> | £300 | | |
| Manual Seat to Floor Height. <i>From seat pan.</i> | 16" | 17" | 18" |
| Powered Tilt Seat to Floor Height. <i>From seat pan.</i> | 16" | 17" Changes Tilt range to -5° to 25° | 18" Changes Tilt range to -10° to 20° |

STANDARD BACK SUPPORT POSITIONING

| | |
|---|----------|
| Manual Angle Settable, Fold-Forward. <i>Adjustment range 95° to 125°.</i> | Standard |
| Attendant Controlled Manual Recline. <i>Adjustment range 95° to 125°. Not suitable for occupied transportation.</i> | £150 |
| Powered Back Recline. <i>Adjustment range 95° to 125°. Uses 1 power option.</i> | £200 |

STANDARD BACK SUPPORTS

| | |
|--|----------|
| TRU-Comfort 2. <i>4-way stretch cover. Adjustable range 24.5" to 28" from seat board.</i> | Standard |
| Tension Adjustable. <i>Air-mesh cover. Adjustment range 19" to 23" from seat board. OBAS** required to mount headrest.</i> | FOC |
| ADI Aluminium Back. <i>Skin protection, 4 way stretch, Coolcore cover.</i> | POA |
| Omit Back Support. <i>Supplied with back canes and brace bar only.</i> | -£25 |
| Enter ADI Back Support part number here for quotation. | |

STEALTH PRO BACK SUPPORTS

| | | | | |
|---|------------------------------|----------|------------------------------|-----|
| Standard, Coolcore Cover. <i>Solid back, 9 cm contour.</i> | £227 | | | |
| Active, Coolcore Cover. <i>Solid back, 12 cm contour.</i> | £341 | | | |
| Extra Deep Active, Coolcore Cover. <i>Solid back, 17 cm contour.</i> | £448 | | | |
| Active Flex, Coolcore Cover. <i>Tension adjustable straps. 12 cm contour.</i> | £559 | | | |
| Back Height | 50 cm. <i>Actual height.</i> | Standard | 57 cm. <i>Actual height.</i> | FOC |

The back support width will match the seat width selected in seat dimensions.

ARM SUPPORTS

| | | | | |
|--|--------------|----------|--|-----|
| Std Flip-up Height Adjustable. <i>Includes contour arm pads. Do not select arm pads below.</i> | Standard | | | |
| TB Flex Removable, Quick Release. <i>Adjustment range 9" to 14" from seat pan.</i> | £85 | | | |
| 2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i> | £279 | | | |
| Arm Support Pads | 10" Standard | Standard | 10" Waterfall. <i>Not with Std flip-up or TB Flex.</i> | FOC |
| | 14" Standard | FOC | 14" Waterfall. <i>Not with Std flip-up or TB Flex.</i> | FOC |

LEG SUPPORTS

| | |
|---|----------|
| 70° Swing-away Legrests. <i>Adjustment range 12" to 15".</i> | Standard |
| Base Mounted Foot Platform. <i>Mounted to chassis. Will not move with power tilt.</i> | FOC |
| 60° to 90° Tool Adjustable Swing-away Legrests. <i>Adjustment range 11" to 16.5".</i> | £330 |
| Manual Elevating Swing-Away Legrests. <i>Adjustment range 15" to 20".</i> | £250 |
| Power Elevating Swing-away Legrests. <i>Adjustment range 14.5" to 19.5". Min. seat width 16". Uses 2 power options.</i> | £650 |
| Manually Adjustable Centre Mount Foot Platform. <i>Adjustment range range 10" to 16". 5"W x 11"L Footplates.</i> | £340 |
| Power Elevating Articulating Centre Mount Foot Platform. <i>Adj. range 10" to 16". Min. seat width 16". Uses 2 power options.</i> | £642 |
| Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away legrest hanger for mounting.</i> | £300 |

| | | |
|------------------------------|-------------|------------------------------------|
| Limb Support | LEFT | RIGHT |
| Seat pan to heel measurement | | (Leave blank for standard factory) |

Adjustment ranges above are from seat pan to footplates. Individual legrests supplied with heel loops. Elevating legrests supplied with calf supports.

STEALTH HEAD SUPPORTS

| | |
|---|------|
| UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i> | £219 |
| UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i> | £219 |
| UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i> | £327 |
| Upgrade Mounting Hardware from HMO475-17M to UHB3D33200. | £88 |
| Push Handles Head Mount Adapter OBAS**. <i>Required to mount head support with TA Back or Active Flex Pro Back.</i> | £266 |

Please visit stealthproducts.com for alternative head supports.

TRUNK SOLUTIONS

| | | | | | | | | |
|----------------------------|-------------|------|---------------------------------------|------|----------------|------|---------------|------|
| Thoracic Lateral Supports. | LEFT | £150 | Lateral Thoracic Swing-away Hardware. | £150 | RIGHT | | | |
| | | £35 | 3" H x 5" L Curved Lateral Pad. | £35 | | | | |
| | | £35 | 4" H x 6" L Curved Lateral Pad. | £35 | | | | |
| Stealth 4-Point. Harness | Extra Small | £137 | Small | £137 | Classic | STD* | Contoured | FOC |
| | Medium | £137 | Large | £137 | Dynamic | STD* | Structured | FOC |
| | Extra Large | £137 | | | Vinyl Face | FOC | Zipper Face** | FOC |
| Chest Strap. | Extra Small | £87 | Small | £87 | Side Release | | | STD* |
| | Medium | £87 | Large | £87 | Velcro Release | | | £19 |
| | Extra Large | £87 | | | | | | |

*STD = Standard option **Classic Only

PELVIC SOLUTIONS

| | | | | | | | | |
|---|----------------|------|-------|------|--------------|------|--------------|------|
| Auto Buckle Lap Belt, Push Button Release. | Standard (STD) | | | | | | | |
| Retractable Lap Belt, Push Button Release. | £156 | | | | | | | |
| Stealth Positioning 2-Point Padded Pelvic Belt. | Extra Small | £78 | Small | £78 | Single Pull | STD* | Dual Pull | FOC |
| | Medium | £78 | Large | £78 | Rear Pull | FOC | Side Release | STD* |
| | Extra Large | £78 | | | Push Button | FOC | | |
| Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull. | Extra Small | £120 | Small | £120 | Push Button | STD* | | |
| | Medium | £120 | Large | £120 | Side Release | £20 | | |

Please visit stealthproducts.com for alternative lateral, pelvic & lower extremity solutions *STD = Standard option

NON-EXPANDABLE CONTROLS

| | | |
|---|---|----------|
| NE + Controls. | FOC | |
| NE + Controls. <i>2 power option ready.</i> | £125 | |
| Control Method. | Hand Control. <i>Select position in Hand Control Mounts below.</i> | Standard |
| | Attendant Control. <i>Hand control mounted at rear.</i> | £100 |
| | Dual Controls. <i>Hand control and secondary attendant control at rear.</i> | £430 |

HAND CONTROL MOUNTS

| | | | | | | |
|---|-----------------|-----|-------------|----------|-------------------|-----|
| Fixed. | Standard | | | | | |
| Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i> | £75 | | | | | |
| Retract 4 Swing-away. <i>Full 180° Movement.</i> | £150 | | | | | |
| Mounting Position. | Left | FOC | Right | Standard | | |
| Hand Control Joystick Handle. | Standard Carrot | STD | T-Bar | £35 | Extended I-Handle | £35 |
| | Mushroom | £35 | Medium Ball | £35 | 3" Goal Post | £50 |
| | Chin Cup | £35 | Soft Ball | £35 | 4" Goal Post | £50 |

SPEED PROGRAM & LIGHTS PACKAGE

| | |
|---|----------|
| 4 Mph Program. <i>Lights and indicators not included.</i> | Standard |
| 4 Mph Program with Lights and indicators. | £450 |
| 6 Mph Program with Lights and indicators. | £450 |
| Lights & Indicator Surcharge. <i>Required if no power options have been selected.</i> | £260 |

ACCESSORIES

| | |
|--|------|
| Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i> | £379 |
| Quantum Back Pack. | £110 |
| XLR USB Charger Port. | £17 |

NOTES & SPECIAL REQUESTS

PRICE

**Q4 WITH TB
FLEX SEATING**



BASE PRICE

£4,450

COST OPTIONS

TOTAL

4



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
RETAIL UK - 11/12/23 V1



Crash Tested



V1: 14/06/2022