

Q6 EDGE HD WITH SYNERGY SEATING

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
	Postcode
Phone	Email
PO No	Order Ref

DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Clients Weight	

BASE MODEL

£5,995

- 10 Colour Options
- Height Adjustable Removable Armrests
- Flip-up Centre Mount Foot Platform
- Angle Adjustable Back Support
- Spectrum Foam SPP Seat Cushion
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- ATX Suspension System
- NE+ Control System
- 2 x 70 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight: 35 st 7lbs. (226 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	Standard
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC
4" x 14" Pneumatic Black Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£220
4" x 14" Flat Free Black Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£320

COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Night Train	FOC	Evolution Orange	FOC
Electric Blue	FOC	Grey Street	FOC
Midnight Blue	FOC	Green Machine	FOC
Deep Purple	FOC	Custom Hydro Dip*	£600
Red Corvette	FOC		

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

20"-24"W
24"-28"W
26"-28"W
Set width at:

SEAT PACKAGES - SEAT DEPTH

15"-20"D
20"-24"D
24"-28"D
Set depth at:

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>226.7 Kg weight limit.</i>	STD
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>204 Kg weight limit.</i>	FOC
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>340 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

Static Seat.	Standard
HD Powered Seat Tilt. <i>0° to 45°. Uses 1 power option.</i>	£1600
Static Seat to Floor Height. <i>From seat pan.</i>	16.75" 17.75" 18.75"
HD Tilt Seat to Floor Height. <i>From seat pan.</i>	18.5" 19.5" 20.5"

SYNERGY BACK SUPPORTS

Synergy Visco Padded Back. <i>Height range 17" to 25" tall from seat pan.</i>	Standard
HD Sling Back. <i>Height range 14" to 25" tall from seat pan.</i>	£200
Enter height preference here.	
Back Cane Height. <i>From seat pan.</i>	14"H. 15"H. 16"H. 17"H. 18"H. 19"H.
	20"H. 21"H. 22"H. 23"H. 24"H. 25"H.

SYNERGY BACK SUPPORT POSTIONING

Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i>	Standard
1" Endo Cane Mounts. <i>Reduces backrest width by 2" compared to seat width.</i>	£40
2" Endo Cane Mounts. <i>Reduces backrest width by 4" compared to seat width.</i>	£40

ARM SUPPORTS

HD Single Post Tool Height Adjustable Removable.	Standard
Bariatric Single Post Tool Height Adjustable Removable.	£220
HD & Bariatric Height Package. <i>From seat pan.</i>	8"-12" FOC 10"-14" Standard
Cantilever Flip-Up Armrests. <i>Adj. range 11"-17". Not available with HD sling back support. Increases overall width by 7".</i>	£280
Arm Support Pads	10" Standard. <i>Not available with Bariatric.</i> FOC 10" Waterfall. <i>Not available with Bariatric.</i> FOC
	14" Standard. <i>Not available with Bariatric.</i> Standard 14" Waterfall. <i>Not available with Bariatric.</i> FOC
	Bariatric. <i>Bariatric arm supports only.</i> FOC

CENTRE MOUNT LEG SUPPORTS

Base Mounted Flip-up Foot Platform. <i>Mounted to chassis. Will not move with power tilt. Ignore other centre mount options if selected.</i>					Standard
Manually Adjustable Centre Mount Foot Platform.					£326
Manual Adjustable Lower Leg Extensions. <i>From seat pan.</i>	10"-16".	Standard	13.5"-19.5".	FOC	
Power Elevating Articulating Centre Mount Foot Platform. <i>Uses 1 power option.</i>					£1320
Power Elevating Articulating & Lower To Floor Centre Mount Foot Platform. <i>Minimum seat width 14". Only with Q-Logic 3. Lowers to the floor to aid standing transfers.</i>					£1790
Powered Elevating Lower Leg Extensions. <i>From seat pan.</i>	11"-17".	Standard	13.5"-19.5".	FOC	
Manual Adjustable & Powered Elevating Footplates.	10"W x 8"D Split Footplates.	FOC	10"W x 11"D Split Footplates.	Standard	
Calf Support.	Individual Small Calf Pads. <i>6"W x 6"H.</i>	STD	Small Calf Panel. <i>18"W x 8"H.</i>		£84
	Individual Large Calf Pads. <i>6"W x 8"H.</i>	FOC	Large Calf Panel. <i>24"W x 8"H.</i>		£84
	Reverse Dartex Calf Pad Covers.	£42	Reverse Dartex Calf Panel Cover.		£62
Footplate Accessories.	4" Width Extension Kit. <i>Not available with 2" length kit or 10"W x 6"D footplates.</i>				£16
	2" Length Extension Kit. <i>Not available with 4" width kit or 10"W x 6"D footplates.</i>				£16
	Footplate Connector.				£10
	Outer Edge Supports.				£10
	Wheel Bumpers.				£10

INDIVIDUALLY MOUNTED LEG SUPPORTS

70° HD Drop-in Swing-away Legrests.						£340
70° HD Lower Extensions. <i>From seat pan.</i>	8.5"-12".	FOC	13"-18".	Standard	16"-21".	FOC
70° HD Footplates.	6"W x 8"D.	Standard	6"W x11"D.	FOC	8"W x 11"D.	FOC
HD Manual Elevating Swing-Away Legrests. <i>Adjustment range 15" to 20". Min. seat width 15"</i>						£380
HD Manual Elevating Lower Extensions. <i>From seat pan.</i>	8.5"-12".	FOC	13"-18".	Standard	16"-21".	FOC
HD Manual Elevating Footplates.	Flip-up Extruded. <i>7.25"W x 6"D</i>					FOC
	Flip-up Extruded. <i>8.5"W x 6"D</i>					Standard
	6"W x 8"D.	£60	6"W x11"D.	£60	8"W x 11"D.	£60

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.
Adjustment ranges above are from seat pan to footplates. Individual legrests are supplied with heel loops.

Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away legrest hanger for mounting.</i>	LEFT	£300	Limb Support	£300	RIGHT
--	------	------	--------------	------	-------

STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£219
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	£219
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i>	£327
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£88

Please visit stealthproducts.com for alternative head supports.

TRUNK SOLUTIONS

Thoracic Lateral Supports.	LEFT	£175	Lateral Thoracic Swing-away Hardware.				£175	RIGHT			
		£36	4" H x 6" L Curved Lateral Pad.				£36				
Stealth 4-Point Harness.	Medium	£136	}	Classic	STD*	}	Dynamic	STD*	}	Vinyl Face	FOC
	Large	£136		Contoured	FOC		Structured	FOC		Zipper Face**	FOC
	Extra Large	£136									
Chest Strap.	Medium		}		£87	}	Side Release		}	Standard	
	Large				£87		Velcro Release			£19	
	Extra Large				£87						

*STD = Standard option

**Classic Only

*STD = Standard option **Classic Only

STEALTH PELVIC & THIGH SOLUTIONS

UniLink Thigh Supports.	LEFT	£176	UniLink Hip Guide Long Removable.		£176	RIGHT					
		£176	Unilink Hip Guide Short Removable.		£176						
		£48	4x4 Pelvic Thigh Pad.		£48						
		£48	4x8 Pelvic Thigh Pad.		£48						
		£48	4x12 Pelvic Thigh Pad.		£48						
Knee Adductors.	LEFT	£156	Fixed Cane Mounted.		£156	RIGHT					
		£170	Swing-Away Track Mounted.		£170						
		£36	3" Round Knee Pad.		£36						
		£36	4" Round Knee Pad.		£36						
Stealth 4-Point Harness.	Medium	£137	} Classic	STD*	} Dynamic	STD*	} Vinyl Face	FOC			
	Large	£137		Contoured		FOC		Structured	FOC	Zipper Face**	FOC
	Extra Large	£137									
Chest Strap.	Medium		£87	} Side Release	Standard						
	Large		£87		} Velcro Release	£19					
	Extra Large		£87								
*STD = Standard option **Classic Only											

*STD = Standard option **Classic Only

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.						Standard
Auto Buckle Lap Belt Length.	38".	FOC	50".	FOC	60".	Standard
	70".	FOC	80".	FOC	90".	FOC
Retractable Lap Belt, Push Button Release.						£156
Stealth Positioning 2-Point Padded Pelvic Belt.	Large	£78	Single Pull	Standard	Push Button	Standard
	Extra Large	£78	Dual Pull	FOC		
			Rear Pull	FOC		
Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull.	Medium		£120	Push Button		Standard
	Large		£120			
Medial Supports (Pommels).	Fixed Flip-down Hardware. <i>MTHW-660.</i>		£180	Small Wedge. <i>MTMD-12.</i>	£103	
	Medium Wedge. <i>MTMD-13.</i>		£103	Large Wedge. <i>MTMD-14.</i>	£103	
	Adjustable Wedge. <i>MTMD-42.</i>		£199	Please visit stealthproducts.com for alternative lateral pelvic & lower extremity solutions.		

Please visit stealthproducts.com for alternative lateral, pelvic & lower extremity solutions.

NON-EXPANDABLE CONTROLS

NE + Controls. <i>With no power options.</i>		Standard
NE + Controls. <i>2 power option ready.</i>		£125
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	£100
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

EXPANDABLE CONTROLS

Q-Logic 3e Controls. <i>No power options.</i>		£225
Q-Logic 3e Controls. <i>2 power option ready.</i>		£300
Q-Logic 3e Controls. <i>5 power option ready.</i>		£500
Q-Logic 3 EX Controls. <i>No power options.</i>		£275
Q-Logic 3 EX Controls. <i>2 power option ready.</i>		£350
Q-Logic 3 EX Controls. <i>5 power option ready.</i>		£550
Q-Logic 3 Enhanced Display. <i>Includes 2 egg switches & flexible display mount.</i>		£1074
SCIM Module Kit. <i>Allows Q-Logic 3 EX joystick display to be used for speciality controls. Not available with 3e controls.</i>		£316
iAccess OLED 4 Button Switch Module. <i>Mounted to opposite armrest unless specified.</i>	£300	<div> <div>Push Button iAccess</div> <div>Toggle Lever iAccess</div> </div> If iAccess is selected choose either Push Button or Toggle lever option.
R-Net Controls. <i>With no power options.</i>		£275
R-Net Controls. <i>1 power option ready.</i>		£350
R-Net Controls. <i>5 power option ready.</i>		£550
Omni 2 Displays. <i>With built-in Infrared. Includes 2 egg switches & flip-down display mount.</i>		£1074
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Omit Hand Control. <i>Often selected when ED is used as master remote with special device.</i>	-£100
	Attendant Control. <i>Hand control mounted at rear.</i>	£125
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

Please visit stealthproducts.com for alternative control options such as mini joysticks, head arrays & switches.

HAND CONTROL MOUNTS

Fixed.					Standard
Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i>					£75
Retract 4 Swing-away. <i>Full 180° Movement.</i>					£150
Mounting Position.	Left	FOC	Right	Standard	
Hand Control Joystick Handle.	Standard Carrot	Standard	T-Bar	£35	Extended I-Handle £35
	Mushroom	£35	Medium Ball	£35	3" Goal Post £50
	Chin Cup	£35	Soft Ball	£35	4" Goal Post £50

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
4 Mph Program with Lights and indicators.	£450
Lights & Indicator Surcharge. <i>Required if no power options have been selected.</i>	£260

ACCESSORIES

USB Mobile Device Charger. <i>Q-Logic only. Mounting to seat rail.</i>	£110
Glove Box.	£30
Hydration System.	£125
Oxygen Holder.	£50
Stealth Tablet Holder.	£342
Back-up Camera.	£300
Magnetic to XLR Charger Adapter. <i>Only with Q-Logic 3e.</i>	£49
Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379
Wooden Trim Kit for Tray. <i>Only with swing-away tray table.</i>	£100

NOTES & SPECIAL REQUESTS

PRICE

Q6 EDGE HD WITH SYNERGY SEATING



BASE PRICE	£5,995
COST OPTIONS	
TOTAL	



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
RETAIL UK - 05/06/23 V1



**Crash
Tested**



V1: 14/06/2022

