

## Q6 EDGE HD WITH TB3 SEATING

### DEALERS DETAILS

Dealers Name \_\_\_\_\_

Order Date \_\_\_\_\_

Address \_\_\_\_\_

Postcode \_\_\_\_\_

Phone \_\_\_\_\_ Email \_\_\_\_\_

PO No \_\_\_\_\_ Order Ref \_\_\_\_\_

### DELIVERY DETAILS

Dealers Address \_\_\_\_\_

Postcode \_\_\_\_\_

Quantum Specialist \_\_\_\_\_

Clients Weight \_\_\_\_\_

### BASE MODEL

**£7,350**

- 10 Colour Options
- Height Adj. Removable Arm Supports
- Flip-up Centre Mount Foot Platform
- Angle Adjustable Back Support
- Spectrum Foam SPP Seat Cushion
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- ATX Suspension System
- NE+ Control System
- 2 x 70 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:  
32 st 2lbs. (204 kg)

Standard = Standard Factory Feature  
FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 14" Flat Free Charcoal Tyres.	Standard
3" x 14" Pneumatic Charcoal Tyres.	FOC
4" x 14" Pneumatic Charcoal Tyres, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£220
4" x 14" Flat Free Charcoal Tyres, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£320

### COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Night Train	FOC	Evolution Orange	FOC
Electric Blue	FOC	Grey Street	FOC
Midnight Blue	FOC	Green Machine	FOC
Deep Purple	FOC	Custom Hydro Dip*	£600
Red Corvette	FOC		

\* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



## SEAT PACKAGES - SEAT WIDTH

20"-24"W

24"-28"W

28"-32"W

Set width at:

## SEAT PACKAGES - SEAT DEPTH

20"-25"D

Set depth at:

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>226.7 Kg weight limit.</i>	STD
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>204 Kg weight limit.</i>	FOC
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>340 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

## SEAT POSITIONING

Static Seat.	Standard		
HD Powered Seat Tilt. <i>0° to 50°. Uses 1 power option.</i>	£1600		
HD Powered iLevel 10" Lift & Tilt. <i>10" lift &amp; 50° tilt. Up to 3.5 mph drive speed while elevated. 28.5 Stone weight capacity. Maximum seat size 24"W x 24"D. Uses 2 power options.</i>	£2250		
Static Seat to Floor Height. <i>From seat pan.</i>	16.5"	17.5"	18.5"
HD Tilt Seat to Floor Height. <i>From seat pan.</i>	18.5"	19.5"	20.5"
HD Lift & Tilt Seat to Floor Height. <i>From seat pan.</i>	18.5"	19.5"	20.5"

## BACKREST SPECIFICATIONS

TRU-Balance 3 Back Support System. <i>Select TB3 back options below. Ignore all 'Back Cane' options.</i>	Standard
Back Canes Only System. <i>Not available with Power recline or TB3 flip-back arm supports. Move directly to 'Back Cane Only Support Options' if selected.</i>	FOC

Select one option above. Both options above include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a set of wedge inserts.

## TB3 BACK SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i>	Standard
Privacy Flap. <i>12.75"W. Closes gap between back support and seat pan.</i>	£50
Stealth Adapter Interface Plate SSAPT4. For mounting third party back supports to TB3 back system. <i>Not required for Q-Back 2.</i>	£358
ADI Q-Back 2. <i>Skin protection, 4 way stretch, Coolcore cover. Mounts directly to TB3 system.</i>	POA

Enter ADI Q-Back 2 part number here for quotation.

## TB3 BACK SUPPORT POSITIONING

TRU-Balance 3 Static Back. <i>Angle can be preset manually.</i>	Standard
TRU-Balance 3 Flip-forward Back. <i>Not available with TB3 flip-back arm supports or powered recline.</i>	£150
TRU-Balance 3 Powered Recline. <i>175° with sheer reduction. Uses 1 power option.</i>	£850

## BACK CANE ONLY SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i>		Standard			
Privacy Flap. <i>12.75"W. Closes gap between back support and seat pan.</i>		£50			
ADI Back Support. <i>Skin protection, 4 way stretch, Coolcore cover. Mounts to canes with cane clamps.</i>		POA			
Enter ADI Back Support part number here for quotation					
Back Cane Only Support Positioning.	Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i>	Standard			
	Flip-forward Back Canes. <i>Not suitable for occupied transport or with Endo mounts.</i>	£200			
	2" Endo Cane Mounts. <i>1" each side. Reduces back width by 2".</i>	£40			
	4" Endo Cane Mounts. <i>2" each side. Reduces back width by 4".</i>	£40			
Back Cane Height. <i>From seat pan.</i>	16"H.	17"H.	18"H.	19"H.	20"H.
	21"H.	22"H.	23"H.	24"H.	25"H.

## BACK SUPPORT DIMENSIONS

Compatible with TB3 and Back Canes options above. Back Support can be +/- 4" of seat width with TB3 Back Support System, but must be same as seat width with Back Canes Only System unless -2" or -4" Endo Mounts have been selected.

Backrest Support Width Package.	16"-20"W.		20"-24"W.	
	14"-28"W.		28"-32"W.	
Back Support Width.	18"W.	19"W.	20"W.	21"W.
	22"W.	24"W.	26"W.	28"W.
	30"W.	32"W.		
Back Support Height. <i>From seat pan.</i>	16"H.	17"H.	18"H.	19"H.
	20"H.	21"H.	22"H.	23"H.
	24"H.	25"H.	26"H.	27"H.
	28"H.			

19" minimum backrest height when powered recline is selected.

## ARM SUPPORTS

Height Adjustable, Removable, Heavy Duty.		Standard
HD Removable Height Package. <i>From seat pan.</i>	8"-12".	FOC 10"-14".
TB3 Flip-back with Down Posts. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>		£450
Arm Support Pads.	10" Standard	Standard
	14" Standard	FOC
	17" Standard. <i>TB3 only. Minimum seat depth 16".</i>	£50
	10" Waterfall.	FOC
	14" Waterfall.	FOC
	10" Falcon Gel	£40
	14" Falcon Gel	£40
	Falcon Gel Mounting Plates. <i>Required for Height Adjustable HD arm supports.</i>	£40

## CENTRE MOUNTED LEG SUPPORTS

Base Mounted Flip-up Foot Platform. <i>Mounted to chassis. Will not move with powered lift or tilt.</i>		FOC
Manually Adjustable Centre Mount Foot Platform. <i>Mounted to seat.</i>		Standard
Power Elevating Articulating Centre Mount Foot Platform. <i>Mounted to seat. Uses 1 power option.</i>		£950
Power Elevating Articulating & Lower To Floor Centre Mount Foot Platform. <i>Mounted to seat. Minimum seat width 14". Only with Q-Logic 3. Lowers to the floor to aid standing transfers. Uses 2 power options.</i>		£1300
Lower Leg Extensions. <i>From seat pan.</i>	10"-16"	Standard 13"-19" FOC
Footplates.	10"W x 8"D Split Footplates. FOC	10"W x 11"D Split Footplates. Standard
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i> Standard	Calf Panel. <i>18"W x 8"H.</i> £42
	Individual Large Calf Pads. <i>8"W x 8"H tapering to 6"W.</i> FOC	Reverse Dartex Calf Panel Cover. £62
	Reverse Dartex Calf Pad Covers. £42	
Footplate Accessories.	4" Width Extension Kit. <i>Not available with 2" length kit or 10"W x 6"D footplates.</i>	£16
	2" Length Extension Kit. <i>Not available with 4" width kit or 10"W x 6"D footplates.</i>	£16
	Footplate Connector. £10	Outer Edge Supports. £10
	Wheel Bumpers.	£10

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.

## INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

70° Heavy Duty, Drop-in, Swing-away Leg Supports. <i>Select lowers and footplates below.</i>		£360
70° HD	Lower Leg Extensions. <i>From seat pan.</i>	8.5"-12" FOC 13"-18" Standard 16"-21" FOC
	Footplate Type.	Angle & Depth Adjustable. <i>Select size below.</i> Standard
	Footplate Size.	6"W x 8"D Standard 6"W x 11"D FOC 8"W x 11"D FOC
	Accessories.	Wheel Bumpers £10
Heavy Duty, Manual Elevating, Swing-away Leg Supports.		£360
ELEVATING HD	Lower Leg Extensions. <i>From seat pan.</i>	11"-14" FOC 14"-17" Standard 12.5"-16.5" FOC
	Footplate Type.	15.5"-19.5" FOC
		Flip-up Extruded. <i>8.5"W x 6"D.</i> Standard
	Footplate Size. <i>Angle adj. only.</i>	Angle & Depth Adjustable. <i>Select size below.</i> £115
		6"W x 8"D Standard
		6"W x 11"D FOC 8"W x 11"D FOC
Accessories.	Wheel Bumpers £10	
LIMB SUPPORT	Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away leg support hanger for mounting.</i>	LEFT £300 Limb Support £300 RIGHT

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information. Adjustment ranges above are from seat pan to footplates. Individual leg supports supplied with heel loops. Elevating leg supports supplied with calf supports.

## STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£230
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	£230
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i>	£342
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£92

## STEALTH TRUNK SOLUTIONS

UniLink Thoracic Lateral Supports. <i>Supplied as a pair.</i>	Lateral Thoracic Swing-away Hardware. £444		Additional 1.5" Link (each). £54	
	3" H x 5" L Curved Lateral Pad. £84		Additional 2" Link (each). £54	
	4" H x 6" L Curved Lateral Pad. £84		Additional 3" Link (each). £54	
4-Point Harness.	Medium £143	Classic STD* } Contoured FOC }	Dynamic STD* } Structured FOC }	Vinyl Face FOC
	Large £143			Zipper Face** FOC
	Extra Large £143			
Chest Strap.	Medium	£91	Side Release	Standard
	Large	£91	Velcro Release	£20
	Extra Large	£91		

\*STD = Standard option \*\*Classic Only

## STEALTH PELVIC & THIGH SOLUTIONS

UniLink Thigh Supports.	<b>LEFT</b>	£185	UniLink Hip Guide Long Removable.	£185	<b>RIGHT</b>
		£185	Unilink Hip Guide Short Removable.	£185	
		£51	4x4 Pelvic Thigh Pad.	£51	
		£51	4x8 Pelvic Thigh Pad.	£51	
		£51	4x12 Pelvic Thigh Pad.	£51	
Knee Adductors.	<b>LEFT</b>	£163	Fixed Cane Mounted.	£163	<b>RIGHT</b>
		£179	Swing-away Track Mounted.	£179	
		£38	3" Round Knee Pad.	£38	
		£38	4" Round Knee Pad.	£38	
Pelvic Belts.	Auto Buckle Lap Belt, Push Button Release.			Standard	
	Retractable Lap Belt, Push Button Release.			£169	
Stealth Positioning 2-Point Padded Pelvic Belt.	Large £84	Single Pull } Rear Pull }	Standard } FOC }	Push Button } Standard }	
	Extra Large £84		Dual Pull } FOC }		
	Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull.	Medium	£143	Push Button	Standard
Large		£143			
Medial Supports (Pommels).	Fixed Flip-down Hardware. <i>MTHW-660.</i> £189		Small Wedge. <i>MTMD-12M.</i>	£108	
	Medium Wedge. <i>MTMD-13M.</i> £108		Large Wedge. <i>MTMD-14M.</i>	£108	
	Adjustable Wedge. <i>MTMD-42M.</i> £209				

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative lateral, pelvic & lower extremity solutions

## NON-EXPANDABLE CONTROLS

NE+ Controls. <i>With no power options. Not available with lights &amp; indicators.</i>	Standard	
NE+ Controls. <i>2 power option ready. Not available with lights &amp; indicators.</i>	£125	
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	£100

## EXPANDABLE CONTROLS

Q-Logic 3e Controls. <i>No power options.</i>	£300
Q-Logic 3e Controls. <i>3 power option ready.</i>	£350
Q-Logic 3e Controls. <i>5 power option ready.</i>	£500
Q-Logic 3 EX Controls. <i>No power options.</i>	£350
Q-Logic 3 EX Controls. <i>3 power option ready.</i>	£375
Q-Logic 3 EX Controls. <i>5 power option ready.</i>	£550
Q-Logic 3 Enhanced Display. <i>Includes 2 egg switches &amp; flexible display mount. Not available with 3e controls.</i>	£1074
SCIM Module Kit. <i>Allows Q-Logic 3 EX joystick display to be used for speciality controls. Not available with 3e controls.</i>	£316

iAccess OLED 4 Button Switch Module.  
*Mounted to opposite armrest unless specified.  
 Supplied FOC when iLevel is selected.*

£300

Push Button iAccess  
 Toggle Lever iAccess

If iAccess is selected choose either  
 Push Button or Toggle lever option.

Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Omit Hand Control. <i>Often selected when ED is used as master remote with special device.</i>	−£100
	Attendant Control. <i>Hand control mounted at rear.</i>	£125
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative control options such as mini joysticks, head arrays & switches.

## HAND CONTROL MOUNTS

Fixed.							Standard
Basic swing-away.							£75
Advanced swing-away. <i>Chain driven. Swings in &amp; out with height, tension, and pivot adjustment.</i>							£100
Mounting Position.	Left	FOC	Right	Standard			
Hand Control Joystick Handle.	Standard Carrot	Standard	T-Bar	£35	Extended I-Handle	£35	
	Mushroom	£35	Medium Ball	£35	3" Goal Post	£50	
	Chin Cup	£35	Soft Ball	£35	4" Goal Post	£50	

## SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
4 Mph Program with Lights and indicators. <i>Not available with NE+ controls.</i>	£450
Lights & Indicator Surcharge. <i>Required if no power options have been selected.</i>	£260

## ACCESSORIES

USB Mobile Device Charger. <i>Please specify mounting position. i.e arm support, backrest, seat rail.</i>	£110
Cup Holder. <i>Only with TB3 flip-back arm supports.</i>	£20
Personal Item Hooks (each). <i>Not with back cane option.</i>	£15
Mobile Phone Holder. <i>Only with TB3 arm supports.</i>	£40
Back Pack Holders. <i>Not with back cane option.</i>	£30
Clothes Guards. <i>Only with TB3 arm supports.</i>	£100
Glove Box.	£30
Hydration System.	£125
Oxygen Holder.	£50

## ACCESSORIES CONTINUED

Quick Release Transfer Bars.	£125
Stealth Tablet Holder.	£342
Push Handles. <i>Not with back cane option.</i>	£50
Accessory Bar.	£50
Back-up Camera.	£300
Quantum Back Pack.	£110
Additional XLR USB Charger Port.	£17
Magnetic to XLR Charger Adapter. <i>Only with Q-Logic 3e.</i>	£49
Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379

## NOTES & SPECIAL REQUESTS

## PRICE

**Q6 EDGE HD WITH  
TB3 SEATING**



<b>BASE PRICE</b>	<b>£7,350</b>
<b>COST OPTIONS</b>	
<b>TOTAL</b>	

7



**Quantum®**  
32 Wedgwood Road, Bicester,  
Oxon, OX26 4UL  
Phone: 01869 324600  
Fax: 01869 323070  
Email: sales@quantumrehab.co.uk  
Web: www.quantumrehab.co.uk  
**RETAIL UK - 01/04/25 V1**

