

## Q6 EDGE Z WITH TB3 SEATING

### DEALERS DETAILS

Dealers Name \_\_\_\_\_

Order Date \_\_\_\_\_

Address \_\_\_\_\_

Postcode \_\_\_\_\_

Phone \_\_\_\_\_ Email \_\_\_\_\_

PO No \_\_\_\_\_ Order Ref \_\_\_\_\_

### DELIVERY DETAILS

Dealers Address \_\_\_\_\_

Postcode \_\_\_\_\_

Quantum Specialist \_\_\_\_\_

Ⓞ Clients Weight \_\_\_\_\_

### BASE MODEL

**£7,695**

- 14 Colour Options
- Height Adjustable Removable Armrests
- Centre Mount Foot Platform
- TB3 Seating
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- ATX Suspension System
- NE+ Control System
- 2 x 70 ah Batteries
- 6 Mph ( 9.65 Km/h) Motors
- Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature  
 FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	Standard
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC
4" x 14" Pneumatic Black Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£220
4" x 14" Flat Free Black Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£320

### COLOURS

Back in Black	Standard	Pretty in Pink	FOC
Green Machine	FOC	Night Train*	FOC
Midnight Blue	FOC	Evolution Orange*	FOC
Electric Blue	FOC	Grey Street*	FOC
Deep Purple	FOC	Custom Hydro Dip*	£600
Red Corvette	FOC		

\* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



## SEAT PACKAGES - SEAT WIDTH

12"-16"W.

14"-16"W.

16"-20"W

20"-24"W

Set width at:

## SEAT PACKAGES - SEAT DEPTH

12"-16"D.

14"-16"D.

16"-20"D

20"-24"D

Set depth at:

Minimum seat size of 14"W x 14"D required for manually adjustable and powered centre mount foot platform when iLevel is selected.

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	Standard
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	£75
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	£100
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£115
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£115
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£135
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i>	£135
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

## SEAT POSITIONING

Static Seat.	Standard		
TRU-Balance 3 Powered Seat Tilt. <i>0° to 50° COG tilt. Uses 1 power option.</i>	£872		
iLevel Powered Seat Lift & Tilt. <i>12" lift &amp; 50° tilt. Up to 3.5 Mph drive speed while elevated. Uses 2 power options.</i>	£2170		
Static Seat to Floor Height. <i>From seat pan.</i>	16.5"	17.5"	18.5"
Powered Tilt Seat to Floor Height. <i>From seat pan.</i>	18.5"	19.5"	20.5"
Powered iLevel Lift & Tilt Seat to Floor Height. <i>From seat pan.</i>	18.5"	19.5"	20.5"

## BACKREST SPECIFICATIONS

TRU-Balance 3 Back Support System. <i>Select TB3 back options below. Ignore all 'Back Cane' options.</i>	Standard
Back Canes Only System. <i>Not available with Power recline or TB3 flip-back arm supports. Move directly to 'Back Cane Only Support Options' if selected.</i>	FOC

Select one option above. Both options above include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a set of wedge inserts.

## TB3 BACK SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i>	Standard
Privacy Flap. <i>12.75"W. Closes gap between back support and seat pan.</i>	£50
Stealth Adapter Interface Plate. <i>SSAPT4. For mounting third party back supports to TB3 back system. Not required for Q-Back 2.</i>	£341
ADI Q-Back 2. <i>Skin protection, 4 way stretch, Coolcore cover. Mounts directly to TB3 system.</i>	POA

Enter ADI Q-Back 2 part number here for quotation.

## TB3 BACK SUPPORT POSTIONING

TRU-Balance 3 Static Back. <i>Angle can be preset manually.</i>	Standard
TRU-Balance 3 Flip-forward Back. <i>Not available with TB3 flip-back arm supports, powered recline.</i>	£150
TRU-Balance 3 Powered Recline. <i>175° with sheer reduction. Uses 1 power option.</i>	£600

## BACK CANE ONLY SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions.</i>	Standard
Privacy Flap. <i>12.75"W. Closes gap between back support and seat pan.</i>	£50
ADI Back Support. <i>Skin protection, 4 way stretch, Coolcore cover. Mounts to canes with cane clamps.</i>	POA

Enter ADI Back Support part number here for quotation

Back Cane Only Support Positioning.	Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i>	Standard
	Flip-forward Back Canes. <i>Not suitable for occupied transport or with Endo mounts.</i>	£200
	2" Endo Cane Mounts. <i>1" each side. Reduces back width by 2". Minimum seat width 12".</i>	£40
	4" Endo Cane Mounts. <i>2" each side. Reduces back width by 4". Minimum seat width 14".</i>	£40
Back Cane Height. <i>From seat pan.</i>	16"H.      17"H.      18"H.      19"H.      20"H.	
	21"H.      22"H.      23"H.      24"H.      25"H.	

## BACK SUPPORT DIMENSIONS

Compatible with TB3 and Back Canes options above. Ignore Back Support Dimensions if an ADI back support is selected in TB3 or Back Canes sections. Back Support can be +/- 4" of seat width with TB3 Back Support System, but must be same as seat width with Back Canes Only System unless -2" or -4" Endo Mounts have been selected.

Backrest Support Width Package.	12"-14"W.	14"-16"W.	16"-20"W.
	18"-22"W.	20"-24"W.	
Back Support Width.	12"W.	14"W.	16"W.
	18"W.	19"W.	20"W.
	22"W.	24"W.	
Back Support Height. <i>From seat pan.</i>	14"H. <i>12"W back only.</i>	15"H. <i>12"W back only.</i>	16"H.
	18"H.	19"H.	20"H.
	22"H.	23"H.	24"H.
	26"H.	27"H.	28"H.

19" minimum backrest height when powered recline is selected.

## ARM SUPPORTS

2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>	Standard
Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i>	FOC
Single Post Drop-in Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i>	FOC
Single Post HD Drop-in Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i>	£150
TB3 Flip-back. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>	£350
Down Posts for TB3 Flip-back. <i>Only with TB3 Flip-back arm supports. Minimum seat/back width 16".</i>	£100

Arm Support Pads.	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	17" Standard. <i>TB3 only. Minimum seat depth 16".</i>	£50	Falcon Gel Mounting Plates. <i>Required for 2-post, Single Post &amp; Paediatric arm supports</i>	£40
	10" Falcon Gel	£40	14" Falcon Gel	£40

## CENTRE MOUNTED LEG SUPPORTS

Base Mounted Flip-up Foot Platform. <i>Mounted to chassis. Will not move with powered lift or tilt.</i>		FOC	
Manually Adjustable Centre Mount Foot Platform. <i>Minimum seat width 14"</i>		Standard	
Paediatric Adjustable Centre Mount Foot Platform. <i>12"-14" seat widths only.</i>		FOC	
Power Elevating Articulating Centre Mount Foot Platform. <i>Minimum seat width 14". Uses 1 power option.</i>		£700	
Power Elevating Articulating & Lower To Floor Centre Mount Foot Platform. <i>Minimum seat width 14". Lowers to the floor to aid standing transfers. Uses 2 power options.</i>		£1112	
Lower Leg Extensions. <i>From seat pan.</i>	8"-12" <i>Paediatric centre mount only.</i> FOC	10"-16" <i>Manually adj. &amp; powered elev. only.</i> STD	13"-19" <i>Manually adj. &amp; powered elev. only.</i> STD
Footplates.	10"W x 6"D Split Footplates. <i>Paediatric centre mount only.</i>		FOC
	10"W x 8"D Split Footplates.	FOC	10"W x 11"D Split Footplates. Standard
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>	Standard	Calf Panel. <i>18"W x 8"H.</i> £42
	Reverse Dartex Calf Pad Covers. £42		Reverse Dartex Calf Panel Cover. £62
Footplate Accessories.	4" Width Extension Kit. <i>Not available with 2" length kit or 10"W x 6"D footplates.</i>		£16
	2" Length Extension Kit. <i>Not available with 4" width kit or 10"W x 6"D footplates.</i>		£16
	Footplate Connector.	£10	Outer Edge Supports. £10
	Wheel Bumpers.	£10	Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.

## INDIVIDUALLY MOUNTED LEG SUPPORTS

70° Swing-Away Leg Supports. <i>Flip-up footplates. Adj. range 12" to 15". Minimum seat width 15"</i>	FOC
90° Paediatric Swing-Away Leg Supports. <i>4"W x 6"D angle &amp; depth adj. footplates. Adjustment range 7" to 11"</i>	£110
Manual Elevating Swing-Away Leg Supports. <i>Flip-up footplates. Adj. range 15" to 20". Min. seat width 15"</i>	£130
Power Elevating Swing-Away Leg Supports. <i>Flip-up f/plates. Adj. range 14.5" to 19.5". Min. seat width 16". Uses 2 power options.</i>	£650

## INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

60° to 90°		60° to 90° Tool Adjustable Swing-Away Leg Supports. <i>Select lowers and footplates below. Minimum seat width 15"</i>				£110	
60° TO 90°	Lower Leg Extensions. <i>From seat pan.</i>	11"-13"	FOC	13.5"-16.5"	Standard	15.5"-19.5"	FOC
	Footplate Type.	Flip-up Extruded. <i>8.5"W x 6"D.</i>	STD	Angle & Depth Adjustable. <i>Select size below.</i>	£115	Multi-axis. <i>Select size below.</i>	£263
	Footplate Size. <i>Angle adj. &amp; Multi-axis only.</i>	4"W x 6"D	FOC	5"W x 6"D	FOC	5"W x 8"D	FOC
		6"W x 8"D	STD	6"W x 11"D	FOC	8"W x 11"D	FOC
Accessories.	Wheel Bumpers					£10	
POWER ELEVATING ARTICULATING		Power Elevating Articulating Swing-Away Leg Supports. <i>Inc. calf pads. Adj. range 15" to 20". Min. seat width 15". Uses 2 power options.</i>				£1250	
POWER ELEVATING ARTICULATING	Lower Leg Extensions. <i>From seat pan.</i>	11"-13"	FOC	14.5"-18"	Standard	17.5"-21"	FOC
	Footplate Type.	Angle & Depth Adjustable. <i>Select size below.</i>					Standard
	Footplate Size.	4"W x 6"D	FOC	5"W x 6"D	FOC	5"W x 8"D	FOC
		6"W x 8"D	Standard	6"W x 11"D	FOC	8"W x 11"D	FOC
Accessories.	Small Calf Pads		Standard	Wheel Bumpers		£10	
	Large Calf Pads		FOC	Gel Padded Knee Adductors		£60	
Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information. Adjustment ranges above are from seat pan to footplates. Individual legrests are supplied with heel loops.							

Stealth Residual Limb Support (each). *Includes a 60-90° Swing-away leg support hanger for mounting.*

LEFT

£300

Limb Support

£300

RIGHT

## STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>	£219
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£219
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	£219
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i>	£327
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£88

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative head supports.

## STEALTH TRUNK SOLUTIONS

UniLink Thoracic Lateral Supports. <i>Supplied as a pair.</i>	Lateral Thoracic Swing-away Hardware. £425		Additional 1.5" Link (each).		£52	
	3" H x 5" L Curved Lateral Pad. £72		Additional 2" Link (each).		£52	
	4" H x 6" L Curved Lateral Pad. £72		Additional 3" Link (each).		£52	
4-Point Harness.	Extra Small	£137	Classic STD*	Dynamic STD*	Vinyl Face	FOC
	Small	£137			Contoured FOC	Structured FOC
	Medium	£137				
	Large	£137				
	Extra Large	£137				
Chest Strap.	Extra Small	£87	Small	£87	Side Release	Standard
	Medium	£87	Large	£87	Velcro Release	£19
	Extra Large	£87				

\*STD = Standard option \*\*Classic Only

## STEALTH PELVIC & THIGH SOLUTIONS

UniLink Thigh Supports.	<b>LEFT</b>	£176	UniLink Hip Guide Long Removable.	£176	<b>RIGHT</b>	
		£176	Unilink Hip Guide Short Removable.	£176		
		£48	4x4 Pelvic Thigh Pad.	£48		
		£48	4x8 Pelvic Thigh Pad.	£48		
		£48	4x12 Pelvic Thigh Pad.	£48		
Knee Adductors.	<b>LEFT</b>	£156	Fixed Cane Mounted.	£156	<b>RIGHT</b>	
		£170	Swing-Away Track Mounted.	£170		
		£36	3" Round Knee Pad.	£36		
		£36	4" Round Knee Pad.	£36		
Pelvic Belts.	Auto Buckle Lap Belt, Push Button Release.				Standard	
	Retractable Lap Belt, Push Button Release.				£156	
Positioning 2-Point Padded Pelvic Belt.	Extra Small	£78	Single Pull	Standard	Side Release	Standard
	Small	£78	Rear Pull	FOC	Push Button	FOC
	Medium	£78	Dual Pull	FOC		
	Large	£78				
	Extra Large	£78				
Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull.	Extra Small	£120	Small	£120	Push Button	Standard
	Medium	£120	Large	£120	Side Release	£20
Medial Supports (Pommels).	Fixed Flip-down Hardware. <i>MTHW-660.</i>		£180	Small Wedge. <i>MTMD-12.</i>		£103
	Medium Wedge. <i>MTMD-13.</i>		£103	Large Wedge. <i>MTMD-14.</i>		£103
	Adjustable Wedge. <i>MTMD-42.</i>		£199			

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative lateral, pelvic & lower extremity solutions.

## NON-EXPANDABLE CONTROLS

NE + Controls. <i>With no power options.</i>	Standard	
NE + Controls. <i>2 power option ready.</i>	£125	
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	£100
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

## EXPANDABLE CONTROLS

Q-Logic 3e and 3 EX hand controls are supplied with built-in Bluetooth. Q-Logic 3 Enhanced Display is supplied with built-in Bluetooth and Infrared.

Q-Logic 3e Controls. <i>No power options.</i>	£225
Q-Logic 3e Controls. <i>2 power option ready.</i>	£300
Q-Logic 3e Controls. <i>5 power option ready.</i>	£500
Q-Logic 3 EX Controls. <i>No power options.</i>	£275
Q-Logic 3 EX Controls. <i>2 power option ready.</i>	£350
Q-Logic 3 EX Controls. <i>5 power option ready.</i>	£550
Q-Logic 3 Enhanced Display. <i>Includes 2 egg switches &amp; flexible display mount.</i>	£1074
IR Blaster. <i>Infrared module for Q-Logic 3 EX hand control.</i>	£50
SCIM Module Kit. <i>Allows Q-Logic 3 EX joystick display to be used for speciality controls. Not available with 3e controls.</i>	£316

iAccess OLED 4 Button Switch Module.  
*Mounted to opposite armrest unless specified.*

£300

Push Button iAccess  
Toggle Lever iAccess

If iAccess is selected choose either  
Push Button or Toggle lever option.

Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Omit Hand Control. <i>Often selected when ED is used as master remote with special device.</i>	-£100
	Attendant Control. <i>Hand control mounted at rear.</i>	£125
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative control options such as mini joysticks, head arrays & switches.

## HAND CONTROL MOUNTS

Fixed.	Standard					
Swing-away. <i>Swings in &amp; out with height, tension and pivot adjustment.</i>	£75					
Retract 4 Swing-away. <i>Full 180° Movement.</i>	£150					
Mounting Position.	Left	FOC	Right	Standard		
Hand Control Joystick Handle.	Standard Carrot	Standard	T-Bar	£35	Extended I-Handle	£35
	Mushroom	£35	Medium Ball	£35	3" Goal Post	£50
	Chin Cup	£35	Soft Ball	£35	4" Goal Post	£50

## SPEED PROGRAM & LIGHTS PACKAGE

6 Mph Program with Lights and indicators.	Standard
4 Mph Program. <i>Lights and indicators not included.</i>	FOC
4 Mph Program with Lights and indicators.	FOC
Accu-Trac Motors. <i>Requires Q-Logic 3 EX. Delivers enhanced traction and tracking, especially with switched drive controls. Only available with 6 mph speed program and lights package.</i>	£450
LED Fender Lights. <i>Includes single switch button mounted to side rail. Not road legal.</i>	£190

## ACCESSORIES

USB Mobile Device Charger. <i>Please specify mounting position. i.e arm support, backrest, seat rail. Not available with VR2 or R-Net.</i>	£110
Cup Holder. <i>Only with TB3 flip-back arm supports.</i>	£20
Personal Item Hooks. <i>Not with back cane option.</i>	£15
Cane/Crutch Holder. <i>Accessory bar required.</i>	£50
Mobile Phone Holder. <i>Only with TB3 arm supports.</i>	£40
Back Pack Holders. <i>Not with back cane option.</i>	£30
Clothes Guards. <i>Only with TB3 arm supports.</i>	£100
Glove Box.	£30
Hydration System.	£125
Oxygen Holder.	£50
Quick Release Transfer Bars.	£125
Stealth Tablet Holder.	£342
Push Handles. <i>Not with back cane option.</i>	£50
Accessory Bar.	£50
Back-up Camera.	£300
Quantum Back Pack.	£110
Additional XLR USB Charger Port.	£17
Magnetic to XLR Charger Adapter. <i>Only with Q-Logic 3e.</i>	£49
Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379
Wooden Trim Kit for Tray. <i>Only with swing-away tray table.</i>	£100

## NOTES & SPECIAL REQUESTS

## PRICE

**Q6 EDGE Z WITH  
TB3 SEATING**



**BASE PRICE** **£7,695**

**COST OPTIONS**

**TOTAL**

7



**Quantum®**  
32 Wedgwood Road, Bicester,  
Oxon, OX26 4UL  
Phone: 01869 324600  
Fax: 01869 323070  
Email: sales@quantumrehab.co.uk  
Web: www.quantumrehab.co.uk  
**RETAIL UK - 04/06/23 V1**



**Crash  
Tested**



V1: 14/06/2022