

Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk



QUANTUM



April 2022



Crash Tested



Q6 Edge HD with Synergy® Seating

Dealer Name: _____ Order Date: _____

Address: _____

Post Code: _____

Phone _____ Email: _____

PO No: _____ Order Ref: _____

Delivery Address: _____

Post Code: _____

Quantum Specialist: _____

Standard Features

- Synergy Seating System
- 4" of Seat Width Adjustment
- 4" of Seat Depth Adjustment
- 2 x 70ah batteries
- NE+ Control System
- 2" x 6" Front & Rear Castors, Flat Free
- 3" x 14" Flat Free Drive Tyres
- 12 Colour Options
- 4 Mph (6 Km/h) Motors
- ATX Suspension System
- Single Post Height Adj. Removable Armrests
- Base Mounted Foot Platform
- Maximum User Weight 35 Stone 10lbs
- Synergy Spectrum Foam Seat Cushion
- Tie-Down System for Transportation

Standard = Standard Factory Feature | FOC = Select instead of Standard Factory Feature

Base Model

Q6 Edge HD with Synergy® Seating - £5,720.00

Drive Wheels

- 3" x 14" Flat Free, Dark Grey Tyres Standard
- 3" x 14" Pneumatic, Black Tyres FOC

Castor Wheels

- 2" x 6" Flat Free, Dark Grey Tyres, Front and Rear Standard

Colours



Red Corvette



Back in Black



Pretty in Pink



Lime Light



Silver Bullet



Orange Crush



Electric Blue



White Satin



Deep Purple



Grey Street



Midnight Blue



Yellow Blaze

Seat Dimensions

Select width and depth package and enter desired seat width and depth in that package range. Available in 1" increments.

Seat Width Packages

- 20"W-24"W
- 24"W-26"W

Set Width at: _____

Seat Depth Packages

- 20"D-24"D Range of adjustment 18" - 24"
- 24"D-26"D

Set Depth at: _____

Seat Cushions

Please refer to our Stealth Cushion Order form for available sizes and enter your preferred size here.....W xD

- No Seat Cushion Required FOC
- Stealth Spectrum® Foam SPP: *Skin protection, positioning, 4-way stretch, Coolcore cover and urethane liner - 42 stone weight limit* Standard
- Stealth Simplicity G: *General use, Fluid resistant Stealth-Tek cover - 25 stone weight limit* FOC
- Stealth TRU-Comfort 2 SPP: *Skin protection, positioning, 4-way stretch Coolcore cover - 32 stone weight limit* FOC
- Stealth Solution® SPP: *Skin protection, positioning, 4-way stretch, Coolcore cover and urethane liner - 25 stone weight limit* FOC
- Stealth Spectrum® Gel SPP: *Skin protection, positioning, 4-way stretch, Coolcore cover and urethane liner - 25 stone weight limit* FOC

Synergy® Backrest Specifications

- Synergy® Backrest Assembly *Including Visco padded backrest* Standard
- Tension Adjustable Backrest Upholstery *This option may increase the delivery time* £120.00
- Back Canes Only *Suitable for third party backrests. A backrest will not be supplied* FOC

Backrest/Back Cane Height:

Set Backrest/Back Cane Height at: _____ Available in 1" increments. Please refer to the Synergy Backrest Matrix below for seat width and backrest height compatibility. Backrest Width must match Seat Width selected unless using 1" or 2" Endo brackets in Backrest Positioning Adjustments section, in which case Backrest Width can be 2" or 4" narrower than Seat Width.

Back Height = Seat Pan to Top of Back | Odd heights = 3" gap at bottom, Even heights = 4" gap at bottom.

		Synergy Static Back Width										
		16"	17"	18"	19"	20"	21"	22"	23"	24"	25"	26"
Synergy Back Height	12"											
	13"											
	14"											
	15"											
	16"											
	17"	✓	✓	✓	✓	✓	✓	✓	✓	✓		
	18"	✓	✓	✓	✓	✓	✓	✓	✓	✓		
	19"	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	20"	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	21"	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	22"	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	23"	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
24"	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
25"			✓	✓	✓	✓	✓	✓	✓	✓	✓	

TRU-Comfort Plus® Backrest Specifications

- TRU-Comfort Plus® Backrest Assembly *4-way stretch fabric (Positioning). Only available with seat widths 20", 22" and 24"* £280.00

Backrest Height:

- 27" Tall Standard
- 24" Tall FOC

Back Height = Seat Pan to Top of Back | 3" gap at bottom.

Seat Positioning Adjustments

- Static Seat Standard
- TRU-Balance®2 HD Tilt *0-45 degree centre of gravity (CG) power tilt seating system. Uses 1 power option* £1600.00

Seat to Floor Height

- | | | | |
|-------------------|---|---|---|
| Static Seat | <input type="checkbox"/> 16.75" Seat to Floor | <input type="checkbox"/> 17.75" Seat to Floor | <input type="checkbox"/> 18.75" Seat to Floor |
| TB2 HD Power Tilt | <input type="checkbox"/> 18.5" Seat to Floor | <input type="checkbox"/> 19.5" Seat to Floor | <input type="checkbox"/> 20.5" Seat to Floor |

Backrest Positioning Adjustments

- Static Backrest *Tool angle adjustable from 85° to 120° in 5° increments* Standard
- Additional Static Back Options:**
- Flip Forward Backrest *Not suitable for occupied transportation* £200.00
- 1" Endo Brackets *Reduces backrest width by 2". Not available with TRU-Comfort Plus backrest on seat widths 14"-17"* £40.00
- 2" Endo Brackets *Reduces backrest width by 4". Not available with TRU-Comfort Plus backrest on seat widths 14"-19"* £40.00

Armrests

- Single Post, Quick Height Adjustable Removable Armrests, Pair *Height adjustment from 11" - 17" in 1/2" increments* Standard
- Single Post, Tool Height Adjustable Armrests, Heavy Duty, Removable, Pair *Height adjustment from 10" - 14"* £220.00
- Cantilever Flip-Up Armrests, Pair *Not available with a tension adjustable backrest. Increases overall width by 7". Height adj. 11"-17" depending on backrest height. Not available with Flip Forward Backrest.* £280.00
- 2-Post, Flip Back, Height Adjustable, Removable *Height range 9"-13". Not available with a seat depth less than 16"* £145.00

Armrests Pad Options

<input type="checkbox"/> 14"L x 2"W Straight, Pair	or	<input type="checkbox"/> 10"L x 2"W Straight, Pair	Standard
<input type="checkbox"/> 14"L x 3"W Waterfall, Pair	or	<input type="checkbox"/> 10"L x 3"W Waterfall, Pair	FOC
<input type="checkbox"/> 14"L x 3.5"W Gel Armpad, Pair			£160.00
<input type="checkbox"/> 10"L x 3.5"W Gel Armpad, Pair			£130.00

Headrests

<input type="checkbox"/> Stealth 10" Comfort Plus CP250	£140.00
<input type="checkbox"/> Stealth 14" Comfort Plus CP150	£140.00
<input type="checkbox"/> Stealth 14" Adjustable Comfort Plus ACP150	£250.00
<input type="checkbox"/> Stealth Small Contoured STL-OTSM250	£200.00
<input type="checkbox"/> Stealth Large Contoured STL-OTLG150	£200.00

Lateral and Trunk Supports

Stealth Thoracic Lateral Supports

Left Lateral Package

<input type="checkbox"/> Thin Hardware, Short Link TWBTL-L	£175.00
<input type="checkbox"/> Thin Hardware, Standard Link TWBSTL-L	£175.00

Left Lateral Pads

<input type="checkbox"/> 3" H x 5" L Curved Lateral Pad TWBL-3X5	£35.00
<input type="checkbox"/> 4" H x 6" L Curved Lateral Pad TWBL-4X6	£35.00

Right Lateral Package

<input type="checkbox"/> Thin Hardware, Short Link TWBTL-R	£175.00
<input type="checkbox"/> Thin Hardware, Standard Link TWBSTL-R	£175.00

Right Lateral Pads

<input type="checkbox"/> 3" H x 5" L Curved Lateral Pad TWBL-3X5	£35.00
<input type="checkbox"/> 4" H x 6" L Curved Lateral Pad TWBL-4X6	£35.00

i-Fit Chest Harness Supplied uninstalled

<input type="checkbox"/> Classic Style <small>Offers a large padded surface for the thoracic area</small>	£100.00
<input type="checkbox"/> Contoured Style <small>Narrowed sectioned over the sternum that allows a better fit for females, G-tube users, or those requiring a smaller support surface</small>	£100.00

Size

<input type="checkbox"/> Extra Small	<input type="checkbox"/> Small	<input type="checkbox"/> Medium	<input type="checkbox"/> Large	<input type="checkbox"/> Extra Large
--------------------------------------	--------------------------------	---------------------------------	--------------------------------	--------------------------------------

Construction

<input type="checkbox"/> Dynamic <small>Allows for a degree of trunk movement and softer fit</small>	<input type="checkbox"/> Structured <small>More resistant to high tone, aiding maintenance of correct posture</small>
--	---

Pelvic Supports

Stealth Pelvic Supports may conflict with armrest option. Contact Quantum for more information

Left Pelvic Package

<input type="checkbox"/> Quick Release Removable Hardware <small>INDPART2805</small>	£100.00
<input type="checkbox"/> Flip-Down Hardware TWBL-L	£150.00
<input type="checkbox"/> Additional link for Flip-down Hardware LHW-129	£35.00
<input type="checkbox"/> Full Surface Contact Feature for Flip-Down Hardware <small>LTHW511</small>	£50.00

Left Pelvic Pads

<input type="checkbox"/> 4" H x 4" L Pelvic Pad TWBADD-4x4	£35.00
<input type="checkbox"/> 4" H x 8" L Pelvic Pad TWBADD-4x8	£35.00
<input type="checkbox"/> 4" H x 12" L Pelvic Pad TWBADD-4x12	£35.00

Right Pelvic Package

<input type="checkbox"/> Quick Release Removable Hardware <small>INDPART2805</small>	£100.00
<input type="checkbox"/> Flip-Down Hardware TWBL-R	£150.00
<input type="checkbox"/> Additional link for Flip-down Hardware LHW-129	£35.00
<input type="checkbox"/> Full Surface Contact Feature for Flip-Down Hardware <small>LTHW511</small>	£50.00

Right Pelvic Pads

<input type="checkbox"/> 4" H x 4" L Pelvic Pad TWBADD-4x4	£35.00
<input type="checkbox"/> 4" H x 8" L Pelvic Pad TWBADD-4x8	£35.00
<input type="checkbox"/> 4" H x 12" L Pelvic Pad TWBADD-4x12	£35.00

Stealth Knee Adductors

Left Knee Support Package

<input type="checkbox"/> Cane Mounted Hardware <small>Only available with individual legrests</small> <small>KN1ADD</small>	£145.00
<input type="checkbox"/> Swing-Away Track Mount Hardware <small>Mounts to seat rail</small> <small>KN1ADD-TRK</small> <small>May conflict with tie-down bracket</small>	£160.00

Left Knee Pads With 1/2" medium sunmate

<input type="checkbox"/> 3" Round Knee Adductor Pad KN1ADD3	£17.00
<input type="checkbox"/> 4" Round Knee Adductor Pad KN1ADD4	£17.00

Right Knee Support Package

<input type="checkbox"/> Cane Mounted Hardware <small>Only available with individual legrests</small> <small>KN1ADD</small>	£145.00
<input type="checkbox"/> Swing-Away Track Mount Hardware <small>Mounts to seat rail</small> <small>KN1ADD-TRK</small> <small>May conflict with tie-down bracket</small>	£160.00

Right Knee Pads With 1/2" medium sunmate

<input type="checkbox"/> 3" Round Knee Adductor Pad KN1ADD3	£17.00
<input type="checkbox"/> 4" Round Knee Adductor Pad KN1ADD4	£17.00

Pelvic Supports (continued)

i-Fit Pelvic Belts

- | | |
|---|--------|
| <input type="checkbox"/> Padded <i>Offers a large padded surface area</i> | £70.00 |
| <input type="checkbox"/> Non-Padded <i>Narrowed sectioned for smaller support surface</i> | £65.00 |
- Size**
- | | | | | |
|--|--------------------------------|---------------------------------|--------------------------------|--|
| <input type="checkbox"/> Extra Small
<i>Not available with non-padded</i> | <input type="checkbox"/> Small | <input type="checkbox"/> Medium | <input type="checkbox"/> Large | <input type="checkbox"/> Extra Large
<i>Not available with non-padded or 4-point.</i> |
|--|--------------------------------|---------------------------------|--------------------------------|--|
- Attachment Method**
- | | | | |
|----------------------------------|--|--------------------------------------|---------------------------------------|
| <input type="checkbox"/> 2-Point | <input type="checkbox"/> 4-Point <i>Single pull only</i> | <input type="checkbox"/> Push Button | <input type="checkbox"/> Side Release |
|----------------------------------|--|--------------------------------------|---------------------------------------|
- Adjustment**
- | | |
|--------------------------------------|------------------------------------|
| <input type="checkbox"/> Single Pull | <input type="checkbox"/> Dual Pull |
|--------------------------------------|------------------------------------|

Medial (Pommel) Supports

Stealth Pommel Mounting Hardware *May conflict with centre mount foot platform. Please contact Quantum for more information*

- | | |
|---|---------|
| <input type="checkbox"/> Fixed Modular Removable Medial Hardware <i>Push button release. MTHW-650</i> | £225.00 |
| <input type="checkbox"/> Flip-Down Link Style Medial Hardware. <i>MTHW-620</i> | £170.00 |
| <input type="checkbox"/> Fixed Flip-Down Medial Hardware. <i>MTHW-640</i> | £135.00 |

Stealth Pommel Pads

- | | |
|---|---------|
| <input type="checkbox"/> Small Wedge Pad 3" W x 4" L x 3" H x 2" W. <i>MTMD-12</i> | £95.00 |
| <input type="checkbox"/> Medium Wedge Pad 4.5" W x 5" L x 3" H x 2" W. <i>MTMD-13</i> | £95.00 |
| <input type="checkbox"/> Large Wedge Pad 5" W x 5" L x 4" H x 3" W. <i>MTMD-14</i> | £95.00 |
| <input type="checkbox"/> Adjustable Pommel with Dual 3" x 5" pads <i>Includes slip cover. MTMD-42</i> | £185.00 |

Centre Mount Foot Platforms

- | | |
|--|-----------|
| <input type="checkbox"/> Base Mounted Foot Platform <i>Mounted to chassis, will not move with Power Tilt or Seat Lift</i> | Standard |
| <input type="checkbox"/> Adjustable Centre Mount Foot Platform <i>Mounted to Seat.</i> | £326.00 |
| <input type="checkbox"/> Power Articulating Centre Mount Foot Platform <i>Mounted to Seat. Uses 1 power option.</i> | £1,320.00 |
| <input type="checkbox"/> Power Articulating Centre Mount Foot Platform with independent actuators
<i>Mounted to Seat. Enables powered foot platform to be lowered to the floor for standing transfers. Only available with QLogic 3. Uses 2 power options</i> | £1,790.00 |

Legrest Options

Lower Extensions *Seat pan to Heel*

- | | | | |
|--|----|--|----------|
| <input type="checkbox"/> 10-16"L <i>Seat pan to heel</i> | or | <input type="checkbox"/> 13-19"L <i>Seat pan to heel</i> | Standard |
|--|----|--|----------|

Footplates

- | | | | |
|---|----|--|----------|
| <input type="checkbox"/> 10"W x 11"D Split Footplates | or | <input type="checkbox"/> 10"W x 8"D Split Footplates | Standard |
|---|----|--|----------|

Calf Pads

- | | | | |
|--|----|--|--------|
| <input type="checkbox"/> Small 6"H x 6"W | or | <input type="checkbox"/> Large 8"H x 8"W <i>tapered to 6"W</i> | £40.00 |
|--|----|--|--------|

Calf Pad Covers (Pair) *1/2" Soft Visco Foam*

- | | | |
|---|--|--------|
| <input type="checkbox"/> Air Mesh | <input type="checkbox"/> Reversed Dartex | £42.00 |
| <input type="checkbox"/> Calf Panel 8"H x 18"W. <i>Not available with Calf Pads</i> | | £84.00 |

Calf Panel Extension (each)

£30.00

Calf Panel Cover *1/2" Soft Visco Foam*

- | | | |
|-----------------------------------|--|--------|
| <input type="checkbox"/> Air Mesh | <input type="checkbox"/> Reversed Dartex | £72.00 |
|-----------------------------------|--|--------|

Footplate Extension Kits



- | | |
|---|--------|
| <input type="checkbox"/> 4" Width Extension Kit. <i>Not available with 2" Bell Shaped, Length Extension Kit</i> | £16.00 |
|---|--------|



- | | |
|---|--------|
| <input type="checkbox"/> 2" Bell Shaped, Length Extension Kit. <i>Not available with 4" Width Extension Kit</i> | £16.00 |
|---|--------|

Additional Footplate Accessories:

- | | |
|---|--------|
| <input type="checkbox"/> Footplate Outside Edge, Pair | £10.00 |
| <input type="checkbox"/> Wheel Bumpers, Pair | £10.00 |
| <input type="checkbox"/> Footplate Connector | £10.00 |

Individual Legrests

Choose instead of centre mount foot platform options

- | | | |
|--------------------------|--|----------|
| <input type="checkbox"/> | 60-90° Swing-Away Legrests Tool adjustable in 5° Increments. Factory set at 70°
Footplate Selection <i>Includes heel loops</i> | £310.00 |
| <input type="checkbox"/> | Flip-up, Extruded Footplates, 8.5"W x 6"D | Standard |
| <input type="checkbox"/> | Angle Adjustable Footplates, 6"W x 8"D, Pair <i>Other sizes available</i> | £115.00 |

Please provide the seat pan to heel measurement required, minimum 11", maximum 19.5": _____

- | | | |
|--------------------------|---|---------|
| <input type="checkbox"/> | 70° Heavy duty, drop in, swing-away legrests with angle adjustable footplates.
Footplate Selection <i>Includes heel loops</i> | £340.00 |
| <input type="checkbox"/> | 6"W x 8"D | |
| <input type="checkbox"/> | 8"W x 11"D <i>Other sizes available</i> | |

Please provide the seat pan to heel measurement required, minimum 12", maximum 18": _____

Elevating Legrests

Choose instead of centre mount foot platform options

- | | | |
|--------------------------|---|----------|
| <input type="checkbox"/> | Manual Elevating Swing-Away Legrests <i>Not suitable for occupied transportation</i> | £130.00 |
| <input type="checkbox"/> | Supply base mounted foot platform to allow power chair to remain crash tested <i>Manual elevating legrests must be removed during transport</i> | FOC |
| | Footplate Selection <i>Includes heel loops</i> | |
| <input type="checkbox"/> | Flip-up, Extruded Footplates, 8.5"W x 6"D & 14"-17" Lower leg extensions | Standard |
| <input type="checkbox"/> | Angle Adjustable Footplates, 6" W x 8"D & 14.5"-18" Lower leg extensions | £115.00 |

- | | | |
|--------------------------|---|----------|
| <input type="checkbox"/> | 70° Power Articulating Elevating Swing-Away Legrests 4" articulation range. Uses 2 power options
Footplate Selection <i>Includes heel loops</i> | £1320.00 |
| <input type="checkbox"/> | Angle Adjustable Footplates, 6"W x 8"D <i>Other sizes available</i> | Standard |

Please provide the seat pan to heel measurement required, minimum 11", maximum 18": _____

Accessories:

- | | | | | | |
|--------------------------|------------------------------------|----|--------------------------|------------------------------------|--------|
| <input type="checkbox"/> | Small Calf Pads, Pair 5.35"W x 6"H | or | <input type="checkbox"/> | Large Calf Pads, Pair 5.86"W x 8"H | £88.00 |
| <input type="checkbox"/> | Gel Padded Adductor Buttons, Pair | | | | £30.00 |

Operation Method: Independent Combined

Additional Legrest Options

- | | | | |
|--------------------------|---|--|---------|
| <input type="checkbox"/> | Stump Board support (each) <i>Includes a 60-90° Swing-Away Legrest hanger</i> | <input type="checkbox"/> Left <input type="checkbox"/> Right | £300.00 |
| <input type="checkbox"/> | Residual Limb support (each) | <input type="checkbox"/> Left <input type="checkbox"/> Right | £180.00 |
- If both a left and right support are required, please amend price accordingly*

NE+ Non-Expandable Controls

- | | | |
|--------------------------|--|----------|
| <input type="checkbox"/> | NE+ - No Power Options | Standard |
| <input type="checkbox"/> | NE+ - 2 Power Option Ready <i>Required when 1 or 2 power options are selected.</i> | £125.00 |
| <input type="checkbox"/> | Mount Joystick as Attendant Control Only <i>Please select NE+ power option above</i> | £100.00 |
| <input type="checkbox"/> | Dual Control <i>Please select NE+ power option above</i> | £380.00 |

Q-Logic 3 Expandable Controls (includes Bluetooth)

- | | | |
|--------------------------|---|----------|
| <input type="checkbox"/> | Q-Logic 3 - 2 Power Option Ready <i>Required when 0, 1 or 2 power options are selected</i> | £350.00 |
| <input type="checkbox"/> | Q-Logic 3 - 5 Power Option Ready <i>Required when more than 2 power options are selected</i> | £550.00 |
| <input type="checkbox"/> | User Joystick not required <i>A Q-Logic 3 hand control will not be supplied. An Enhanced Display must be selected.</i> | -£100.00 |
| <input type="checkbox"/> | Bluetooth enabled <i>Please tick this box if you would like the Bluetooth function turned on.</i> | Standard |
| <input type="checkbox"/> | IR Blaster <i>Infrared module for Q-Logic 3 hand control</i> | £50.00 |
| <input type="checkbox"/> | Mount Joystick as Attendant Control Only <i>Please select Q-Logic 3 power option above</i> | £100.00 |
| <input type="checkbox"/> | Dual Control <i>Please select Q-Logic 3 power option above</i> | £430.00 |
| <input type="checkbox"/> | SCIM Module Kit <i>Allows use of joystick display for specialty controls. Q-Logic 3 Joystick must be selected above</i> | £316.00 |
| <input type="checkbox"/> | Q-Logic 3 Enhanced Display <i>Includes x2 egg switches and flexible display mount. Does not include Q-Logic 3 Joystick</i> | £1074.00 |
| <input type="checkbox"/> | iAccess OLED 4 button switch module <i>Supplied as standard with iLevel seat lift. Mounted to opposite armrest unless specified</i> | £300.00 |

Push Button iAccess Module



Toggle lever iAccess Module



Joystick Handles

- | | | |
|---|---|--|
| <input type="checkbox"/> Standard Carrot Handle | | Standard |
| <input type="checkbox"/> Mushroom Handle | <input type="checkbox"/> Extended I-Handle | £35.00 (each) |
| <input type="checkbox"/> T-Bar | <input type="checkbox"/> Soft Ball Handle | <input type="checkbox"/> Bodypoint 3 Inch Goal Post Handle |
| <input type="checkbox"/> Chin Cup | <input type="checkbox"/> Medium Ball Handle | <input type="checkbox"/> Bodypoint 4 Inch Goal Post Handle |

Joystick Mounting Brackets

- | | |
|--|----------|
| <input type="checkbox"/> Fixed Inline | Standard |
| <input type="checkbox"/> Swing-Away | £75.00 |
| <input type="checkbox"/> Multi-Axis | £125.00 |
| <input type="checkbox"/> Retract 4 Swing-Away Joystick Arm Kit <i>Full 180° Movement</i> | £150.00 |

Mounting Position




- | | | |
|---------------------------------|----|--|
| <input type="checkbox"/> Left | or | <input type="checkbox"/> Right |
| <input type="checkbox"/> Inline | or | <input type="checkbox"/> Offset <i>Not available with Fixed Joystick Bracket</i> |

Speed Program and Lights Package

- | | |
|---|----------|
| <input type="checkbox"/> 4 Mph program <i>Lights and Indicators not Included</i> | Standard |
| <input type="checkbox"/> 4 Mph program with Lights and indicators | £450.00 |
| <input type="checkbox"/> Lights and indicator surcharge <i>If no power options have been selected, you must tick this box in addition to the lights and indicator package to receive the correct lighting module.</i> | £260.00 |

Quantum Accessories

Some accessories will be shipped with the unit uninstalled. Multiple accessories can be mounted along the various tracks, some interference may occur.

- | | | | |
|--------------------------|---|---|---------|
| <input type="checkbox"/> |  | Glove Box <i>Supplied uninstalled</i> | £30.00 |
| <input type="checkbox"/> |  | Clothing Guard (pair) <i>Only available with Flip Up armrests</i> | £100.00 |
| <input type="checkbox"/> |  | Quantum Backpack | £40.00 |

Price

Sub Total:	£5,720.00
Options:	
Total:	

Special Instructions

For additional seating options and alternative drive controls please refer to your Quantum Specialist.

Quantum Notes:

1. This product is suitable for occupied transportation in accordance with ISO 7176-19. We strongly recommend that all heavy accessories such as manual ELRs are removed and secured during occupied transportation. We also advise that the use of third party accessories including life Support Systems, Chin Controls and Trays be risk assessed prior to the wheelchair being transported. Full details of suitable restraints for this wheelchair can be obtained from Quantum Technical on 01869 324600 or Email us on: technical@quantumrehab.co.uk
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3 %. Overall measurements will vary based on seating and accessory selections.