

Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk



April 2022



Crash Tested



Quantum Lightning with Synergy® Seating and R-Net Controls

Dealer Name: _____ Order Date: _____

Address: _____

_____ Post Code: _____

Phone _____ Email: _____

PO No: _____ Order Ref: _____

Delivery Address: _____

_____ Post Code: _____

Quantum Specialist: _____

Standard Features

- Synergy Seating System
- 4" of Seat Width Adjustment
- 4" of Seat Depth Adjustment
- 2 x 70ah batteries
- R-Net Joystick
- 3" x 9" Front Castors, Pneumatic
- 3" x 14" Pneumatic Drive Tyres
- 5 Colour Options
- 4 Mph (6 Km/h) Motors
- ATX Suspension System
- Single Post Height Adj. Removable Armrests
- 60°-90° Swing-Away Legrests
- Maximum User Weight 21 Stone 6lbs
- Synergy Simplicity Seat Cushion
- Tie-Down System for Transport

Standard = Standard Factory Feature | FOC = Select instead of Standard Factory Feature

Base Model

Lightning with Synergy® Seating and R-Net Controls - £4,410.00

Drive Wheels

- 3" x 14" Pneumatic, Grey Tyres Standard
- 3" x 14" Flat Free, Grey Tyres £150.00

Castor Wheels

- 3" x 9" Pneumatic, Grey Tyres Standard
- 3" x 9" Flat Free, Grey Tyres £60.00

Colours



Red



Black



Blue



Silver



Orange

Seat Dimensions

Select width and depth package and enter desired seat width and depth in that package range, available in 1" increments. Maximum seat size is 22"W x 22"D.

Seat Width Packages

- 14"W-18"W 16"W-20"W
- Set Width at: _____ (Max seat width 20")

Seat Depth Packages

- 14"D-18"D 16"D-20"D
- Set Depth at: _____ (Max seat depth 20")

Seat Cushions

- | | |
|---|----------|
| <input type="checkbox"/> No Seat Cushion Required | FOC |
| <input type="checkbox"/> Simplicity Cushion: <i>General Use</i> | Standard |
| <input type="checkbox"/> TRU-Comfort™2 4-Way Stretch: <i>Skin Protection</i> | £45.00 |
| <input type="checkbox"/> Solution® Cushion: <i>Skin Protection & Positioning</i> | £100.00 |
| <input type="checkbox"/> Solution® 1 Cushion: <i>Skin Protection</i> | £100.00 |
| <input type="checkbox"/> Spectrum® Gel Cushion: <i>Positioning</i> | £160.00 |
| <input type="checkbox"/> Spectrum® Foam Cushion: <i>Skin Protection & Positioning</i> | £160.00 |

Refer to cushion matrix for size availability and enter your selection here: W..... X D.....

Synergy® Backrest Specifications

- | | |
|---|----------|
| <input type="checkbox"/> Synergy® Backrest Assembly <i>Including Visco padded backrest</i> | Standard |
| <input type="checkbox"/> Tension Adjustable Backrest Upholstery <i>This option may increase the delivery time</i> | £120.00 |
| <input type="checkbox"/> Back Canes Only <i>Suitable for third party backrests. A backrest will not be supplied</i> | FOC |

Backrest/Back Cane Height:

Set Backrest/Back Cane Height at: _____ (Available in 1" increments. Min 13", Max 25") *Please refer to Synergy Backrest Matrix on page * for seat width and backrest height compatibility*

Back Height = Seat Pan to Top of Back | Odd heights = 4" gap at bottom, Even heights = 3" gap at bottom.

TRU-Comfort Plus® Backrest Specifications

- | | |
|---|----------|
| <input type="checkbox"/> TRU-Comfort Plus® Backrest Assembly <i>4-way stretch fabric (Positioning). Only available with seat widths 16", 18", 20" and 22"</i> | £220.00 |
| Backrest Height: | |
| <input type="checkbox"/> 27" Tall | Standard |
| <input type="checkbox"/> 24" Tall | FOC |

Back Height = Seat Pan to Top of Back | 3" gap at bottom.

Seat Positioning Adjustments

- | | |
|--|----------|
| <input type="checkbox"/> Static Seat | Standard |
| <input type="checkbox"/> TRU-Balance®2 Power Tilt <i>0-55 degree centre of gravity (CG) power tilt seating system. Uses 1 power option</i> | £1100.00 |

Seat to Floor Height

- | | | | |
|----------------|--|---|--|
| Static Seat | <input type="checkbox"/> 16.5" Seat to Floor | <input type="checkbox"/> 17.25" Seat to Floor | <input type="checkbox"/> 18" Seat to Floor |
| TB2 Power Tilt | <input type="checkbox"/> 18" Seat to Floor | <input type="checkbox"/> 18.25" Seat to Floor | <input type="checkbox"/> 19.5" Seat to Floor |

Backrest Positioning Adjustments

- | | |
|---|----------|
| <input type="checkbox"/> Static Backrest <i>Tool angle adjustable from 85° to 120° in 5° increments</i> | Standard |
| Additional Static Back Options: | |
| <input type="checkbox"/> Flip Forward Backrest <i>Not suitable for occupied transportation</i> | £200.00 |
| <input type="checkbox"/> 1" Endo Brackets <i>Reduces backrest width by 2". Not available with TRU-Comfort Plus backrest on seat widths 14"-17"</i> | £40.00 |
| <input type="checkbox"/> 2" Endo Brackets <i>Reduces backrest width by 4". Not available with TRU-Comfort Plus backrest on seat widths 14"-19"</i> | £40.00 |
| <input type="checkbox"/> TRU-Balance®2 Power Recline <i>Only available with 16"-22" s/widths. 168° of power recline with full drive lockout at 20°. Uses 1 power option. Limited backrest heights available. See reference chart on page*</i> | £1790.00 |

Armrests

- | | |
|---|----------|
| <input type="checkbox"/> Single Post, Quick Height Adjustable Removable Armrests, Pair <i>Height adjustment from 11" - 17" in 1/2" increments</i> | Standard |
| <input type="checkbox"/> Single Post, Tool Height Adjustable Armrests, Heavy Duty, Removable, Pair <i>Height adjustment from 10"-14"</i> | £220.00 |
| <input type="checkbox"/> Paediatric, Removable, Height Adjustable <i>Height adjustment from 7.75" to 9.75" in 1/2" increments. Not available with full length arm pads</i> | FOC |
| <input type="checkbox"/> Cantilever Flip-Up Armrests, Pair <i>Not available with a tension adjustable backrest. Increases overall width by 7". Height adj. 11"-17" depending on backrest height. Not available with Flip Forward Backrest. Not available with TB2 Power Recline</i> | £280.00 |
| <input type="checkbox"/> 2-Post, Flip Back, Height Adjustable, Removable <i>Height range 9"-13". Not available with a seat depth less than 16"</i> | £145.00 |

Armrests Pad Options

- | | | |
|--|---|----------|
| <input type="checkbox"/> 14"L x 2"W Straight, Pair | <input type="checkbox"/> 10"L x 2"W Straight, Pair | Standard |
| <input type="checkbox"/> 14"L x 3"W Waterfall, Pair | <input type="checkbox"/> 10"L x 3"W Waterfall, Pair | FOC |
| <input type="checkbox"/> 14"L x 3.5"W Gel Armpad, Pair | | £160.00 |
| <input type="checkbox"/> 10"L x 3.5"W Gel Armpad, Pair | | £130.00 |

Offotobock Arm Channel (Each)

- | | |
|---|---------|
| <input type="checkbox"/> Medium (4.75" x 11.5") <i>Desk length armrest required</i> | £140.00 |
| <input type="checkbox"/> Large (4.75" x 14") <i>Full length armrest required</i> | £140.00 |
| <input type="checkbox"/> One Piece (4.75" x 18") <i>Full length armrest required. Not compatible with Hand Pad.</i> | £140.00 |
| <input type="checkbox"/> Offotobock Hand Pad <i>Not available on same side as controller</i> | £50.00 |
| <input type="checkbox"/> Medium <input type="checkbox"/> Large | |

Offotobock Arm Channel Mounting Position

- Left Right Both *Please amend price for pair*

Headrests

- | | |
|---|---------|
| <input type="checkbox"/> Quantum 10" W Headrest | £140.00 |
| <input type="checkbox"/> Stealth Small Contoured 9.75"W x 3"H | £200.00 |
| <input type="checkbox"/> Stealth Large Contoured 11"W x 6"H | £200.00 |
| <input type="checkbox"/> Stealth 10" Comfort Plus | £200.00 |
| <input type="checkbox"/> Stealth 14" Comfort Plus | £200.00 |

Please visit www.stealthproducts.com for more headrest options

Individual Legrests

- | | |
|---|----------|
| <input type="checkbox"/> 60-90° Swing-Away Legrests <i>Tool adjustable in 5° increments. Factory set at 70°</i> | Standard |
| Footplate Selection <i>Includes heel loops</i> | |
| <input type="checkbox"/> Flip-up, Extruded Footplates, 7.25"W x 6"D, Pair | Standard |
| <input type="checkbox"/> Angle Adjustable Footplates, 6"W x 8"D, Pair <i>Other sizes available</i> | £115.00 |

Please provide the seat pan to heel measurement required, minimum 11", maximum 19.5": _____

- | | |
|--|---------|
| <input type="checkbox"/> Heavy duty, drop in, swing-away legrests with angle adjustable footplates | £130.00 |
| Legrest Angle Footplate Selection <i>Includes heel loops</i> | |
| <input type="checkbox"/> 70° <input type="checkbox"/> 80° <input type="checkbox"/> 6"W x 8"D <input type="checkbox"/> 8"W x 11"D <i>Other sizes available</i> | |

Please provide the seat pan to heel measurement required, minimum 12", maximum 18": _____

- | | |
|---|-----|
| <input type="checkbox"/> 90° Paediatric Swing-Away Legrests <i>Seat Pan to Heel range of adjustment 7"-10"L. Includes 4" x 6"D angle adjustable footplates and heel loops</i> | FOC |
|---|-----|

- | | |
|---|----------|
| <input type="checkbox"/> 70° Tapered Legrests <i>Not available with seat widths 17" or less</i> | FOC |
| Footplate Selection <i>Includes heel loops</i> | |
| <input type="checkbox"/> Flip-up, Extruded Footplates, 7.25"W x 6"D, Pair | Standard |
| <input type="checkbox"/> Angle Adjustable Footplates, 6"W x 8"D, Pair | £115.00 |

Please provide the seat pan to heel measurement required, minimum 11", maximum 19.5": _____

Elevating Legrests

Choose instead of centre mount foot platform options

- | | |
|---|----------|
| <input type="checkbox"/> Manual Elevating Swing-Away Legrests <i>Not suitable for occupied transportation</i> | £130.00 |
| Footplate Selection <i>Includes heel loops</i> | |
| <input type="checkbox"/> Flip-up, Extruded Footplates, 7.25" W x 6"D, Pair | Standard |
| <input type="checkbox"/> Angle Adjustable Footplates, 6" W x 8"D, Includes Wheel Bumpers, Pair <i>Other sizes available</i> | £115.00 |

Please provide the seat pan to heel measurement required, minimum 11", maximum 18": _____

Elevating Legrests (Continued)

70° Euro Power Elevating Drop-in Legrests. *Includes calf pads. Uses 2 power options.* £720.00

Footplate Selection *Includes heel loops*

Flip-up, Extruded Footplates, 7.25"W x 6"D, Pair Standard

Angle Adjustable Footplates, 6"W x 8"D, Pair *Other sizes available* £115.00

Please provide the seat pan to heel measurement required, minimum 11", maximum 18": _____

70° Power Articulating Elevating Swing-Away Legrests 4" articulation range. *Uses 2 power options.* £1320.00

Footplate Selection *Includes heel loops*

Angle Adjustable Footplates, 6"W x 8"D, Includes Wheel Bumpers, Pair *Other sizes available* Standard

Please provide the seat pan to heel measurement required, minimum 11", maximum 18": _____

Accessories:

Small Calf Pads, Pair 5.35"W-6"H or Large Calf Pads, Pair 5.86"Wx8"H £88.00

Gel Padded Adductor Buttons, Pair £30.00

Operation Method: Independent Combined

Additional Legrest Options

Stump Board support (each) *Includes a 60-90° Swing-Away Legrest hanger* Left Right £300.00

Residual Limb support (each) Left Right £180.00

If both a left and right support are required, please amend price accordingly

R-Net Expandable Controls

R-Net Joystick Standard

Mount Joystick as Attendant Control Only £100.00

Dual Control £430.00

R-Net Omni Enhanced Display *Includes x2 egg switches and swing-away display mount. Does not include R-Net Joystick* £763.00

Joystick Handles

Standard Carrot Handle Standard

Mushroom Handle Extended I-Handle Bodypoint 3 Inch Goal Post Handle £35.00 (each)

T-Bar Soft Ball Handle Bodypoint 4 Inch Goal Post Handle

Chin Cup Medium Ball Handle

Joystick Mounting Brackets

Fixed, Inline Standard

Swing-Away £75.00

Multi-Axis £125.00

Retract 4 Swing-Away Joystick Arm Kit *Full 180° Movement* £150.00

Mounting Position

Left or Right

Inline or Offset *Not available with Fixed Joystick Bracket*

Speed Program and Lights Package

4 Mph program *Available with VR2 or R-Net. Lights and indicators not included* Standard

4 Mph program, with Lights and indicators £400.00

Quantum Accessories

Some accessories will be shipped with the unit, not installed. Multiple accessories can be mounted along the various tracks, some interference may occur.

 Glove Box *Supplied uninstalled* £30.00

 Clothing Guard (pair) *Only available with Flip Up armrests* £100.00

