



## RETAIL

### SOLUTION SPP CUSHION

#### DEALERS DETAILS

|              |          |
|--------------|----------|
| Dealers Name |          |
| Order Date   |          |
| Address      |          |
|              | Postcode |
| Email        |          |
| Phone        | Fax      |


#### DELIVERY DETAILS

|                 |       |       |
|-----------------|-------|-------|
| Dealers Address |       |       |
|                 |       |       |
| Postcode        |       |       |
| PO No           |       |       |
| Order Type      | Quote | Order |

#### DETAILS

This cushion offers skin protection and positioning with a Coolcore, 4-way stretch, breathable cover and an urthane liner.

Standard viscoelastic foam top with twin cell aqueous viscous gel (heavier gel) insert in ishchial well, mid-level contour.



Weight capacities  
79.37 kg: 10"-14" width  
158.75 kg: 15"-21" width,  
226.79 kg: 22" width+  
Overall weight: 2.4kg  
Overall height 3"  
10x10" up to 24x24"  
Odd sizes in 17" and 19"  
Custom sizes available

#### CUSHION

|                                   |                     | ADD<br>WIDTH | ADD<br>DEPTH         |      |
|-----------------------------------|---------------------|--------------|----------------------|------|
| Solution SPP. 10"-20" W.          |                     | ST-SOLSPP    | <input type="text"/> | £209 |
| Solution SPP. 22"-24" W           |                     | ST-SOLSPP    | <input type="text"/> | £231 |
| Extra<br>Covers<br>And<br>Liners. | Coolcore Cover.     | SOLSPP-NBLK  | <input type="text"/> | £88  |
|                                   | Mesh Cover.         | SOLSPP-MBLK  | <input type="text"/> | £88  |
|                                   | Stealth Tek Cover.  | SOLSPP-TBLK  | <input type="text"/> | £44  |
|                                   | Urethane Liner.     | URL1709070   | <input type="text"/> | £44  |
|                                   | Polyurethane Liner. | PUL1709030   | <input type="text"/> | £72  |

|       |    | Width |    |    |    |    |    |    |    |    |    |    |  |  |
|-------|----|-------|----|----|----|----|----|----|----|----|----|----|--|--|
|       |    | 10    | 12 | 14 | 15 | 16 | 17 | 18 | 19 | 20 | 22 | 24 |  |  |
| Depth | 10 |       |    |    |    |    |    |    |    |    |    |    |  |  |
|       | 12 |       |    |    |    |    |    |    |    |    |    |    |  |  |
|       | 14 |       |    |    |    |    |    |    |    |    |    |    |  |  |
|       | 15 |       |    |    |    |    |    |    |    |    |    |    |  |  |
|       | 16 |       |    |    |    |    |    |    |    |    |    |    |  |  |
|       | 17 |       |    |    |    |    |    |    |    |    |    |    |  |  |
|       | 18 |       |    |    |    |    |    |    |    |    |    |    |  |  |
|       | 19 |       |    |    |    |    |    |    |    |    |    |    |  |  |
|       | 20 |       |    |    |    |    |    |    |    |    |    |    |  |  |
|       | 22 |       |    |    |    |    |    |    |    |    |    |    |  |  |
|       | 24 |       |    |    |    |    |    |    |    |    |    |    |  |  |

Select your size, then complete the Part Number by adding the Width first then Depth.

## ADD SOLID BASE

### ABS FLAT RIGIDISER 22" WIDTH AND GREATER CUSHIONS INCLUDE A FLAT RIGIDISER AS STANDARD.

|                                   |     |                                      |     |                                      |     |
|-----------------------------------|-----|--------------------------------------|-----|--------------------------------------|-----|
| 10"W x 10"D. <i>DWR1010E188A.</i> | £39 | 16"W x 16"D. <i>DWR1010E191G.</i>    | £39 | 19"W x 19"D. <i>MEC1804152-1919.</i> | £39 |
| 10"W x 12"D. <i>DWR1010E188C.</i> | £39 | 16"W x 18"D. <i>DWR1010E191I.</i>    | £39 | 20"W x 16"D. <i>DWR1010E193G.</i>    | £39 |
| 12"W x 10"D. <i>DWR1010E189A.</i> | £39 | 16"W x 20"D. <i>DWR1010E191K.</i>    | £39 | 20"W x 18"D. <i>DWR1010E193I.</i>    | £39 |
| 12"W x 12"D. <i>DWR1010E189C.</i> | £39 | 17"W x 17"D. <i>DWR1010E234H.</i>    | £39 | 20"W x 20"D. <i>DWR1010E193K.</i>    | £39 |
| 12"W x 14"D. <i>DWR1010E189E.</i> | £39 | 17"W x 19"D. <i>DWR1010E234J.</i>    | £39 | 20"W x 22"D. <i>DWR1010E193M.</i>    | £39 |
| 14"W x 12"D. <i>DWR1010E190C.</i> | £39 | 18"W x 16"D. <i>DWR1010E192G.</i>    | £39 | 22"W x 18"D. <i>DWR1010E194I.</i>    | £39 |
| 14"W x 14"D. <i>DWR1010E190E.</i> | £39 | 18"W x 18"D. <i>DWR1010E192I.</i>    | £39 | 22"W x 20"D. <i>DWR1010E194K.</i>    | £39 |
| 14"W x 16"D. <i>DWR1010E190G.</i> | £39 | 18"W x 20"D. <i>DWR1010E192K.</i>    | £39 | 22"W x 22"D. <i>DWR1010E194M.</i>    | £39 |
| 14"W x 18"D. <i>DWR1010E190I.</i> | £39 | 18"W x 22"D. <i>DWR1010E192M.</i>    | £39 | 24"W x 20"D. <i>DWR1010E236K.</i>    | £39 |
| 15"W x 18"D. <i>DWR1010E232H.</i> | £39 | 19"W x 17"D. <i>MEC1804152-1917.</i> | £39 | 24"W x 22"D. <i>DWR1010E236M.</i>    | £39 |

### ACM CONTOURED RIGIDISER FOR USE WITH SLING SEATS ONLY.

|                                |     |                                |     |                                |     |
|--------------------------------|-----|--------------------------------|-----|--------------------------------|-----|
| 15"W x 18"D. <i>CUS142467.</i> | £72 | 18"W x 16"D. <i>CUS142466.</i> | £72 | 20"W x 18"D. <i>CUS142465.</i> | £72 |
| 16"W x 16"D. <i>CUS142467.</i> | £72 | 18"W x 18"D. <i>CUS142466.</i> | £72 | 20"W x 19"D. <i>CUS142465.</i> | £72 |
| 16"W x 17"D. <i>CUS142467.</i> | £72 | 18"W x 19"D. <i>CUS142466.</i> | £72 | 20"W x 20"D. <i>CUS142465.</i> | £72 |
| 16"W x 18"D. <i>CUS142467.</i> | £72 | 18"W x 20"D. <i>CUS142466.</i> | £72 |                                |     |

## NOTES & SPECIAL REQUESTS

Special Order Solution SPP Cushion.

£400

Please detail your special order cushion request.

All special order cushion requests must be approved with the Quantum Custom Department.

Please contact Quantum Inside Sales or Custom Department with questions. Extended lead time will apply.

Cushion part numbers involved will be the cushion part number listed in previous sections modified to meet special order request stated in details below. Drawings are available for sketch details.


## PRICE

TOTAL

FRONT VIEW



RIGHT SIDE VIEW



TOP VIEW

