



RETAIL

SOLUTION SPP CUSHION

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
	Postcode
Email	
Phone	Fax


DELIVERY DETAILS

Dealers Address		
Postcode		
PO No		
Order Type	Quote	Order

DETAILS

This cushion offers skin protection and positioning with a Coolcore, 4-way stretch, breathable cover and an urthane liner.

Standard viscoelastic foam top with twin cell aqueous viscous gel (heavier gel) insert in ishcial well, mid-level contour.



Weight capacities
79.37 kg: 10"-14" width
158.75 kg: 15"-21" width,
226.79 kg: 22" width+
Overall weight: 2.4kg
Overall height 3"
10x10" up to 24x24"
Odd sizes in 17" and 19"
Custom sizes available

CUSHION

		ADD WIDTH	ADD DEPTH	
Solution SPP. 10"-20" W.		ST-SOLSPP		£190
Solution SPP. 22"-24" W		ST-SOLSPP		£210
Cushion Cover.	Coolcore.	SOLSPP-NBLK		STD
	Mesh.	SOLSPP-MBLK		FOC
	Stealth Tek.	SOLSPP-TBLK		FOC
Fluid Proof Liner.	Urethane.	URL1709070		STD
	Polyurethane.	PUL1709030		£48
Extra Covers And Liners.	Coolcore Cover.	SOLSPP-NBLK		£70
	Mesh Cover.	SOLSPP-MBLK		£44
	Stealth Tek Cover.	SOLSPP-TBLK		£50
	Urethane Liner.	URL1709070		£30
	Polyurethane Liner.	PUL1709030		£48

		Width												
		10	12	14	15	16	17	18	19	20	22	24		
Depth	10													
	12													
	14													
	15													
	16													
	17													
	18													
	19													
	20													
	22													
	24													

Select your size, then complete the Part Number by adding the Width first then Depth.

ADD SOLID BASE

ABS FLAT RIGIDISER

10"W x 10"D. <i>DWR1010E188A.</i>	£34	16"W x 16"D. <i>DWR1010E191G.</i>	£34	19"W x 19"D. <i>MEC1804152-1919.</i>	£34
10"W x 12"D. <i>DWR1010E188C.</i>	£34	16"W x 18"D. <i>DWR1010E191I.</i>	£34	20"W x 16"D. <i>DWR1010E193G.</i>	£34
12"W x 10"D. <i>DWR1010E189A.</i>	£34	16"W x 20"D. <i>DWR1010E191K.</i>	£34	20"W x 18"D. <i>DWR1010E193I.</i>	£34
12"W x 12"D. <i>DWR1010E189C.</i>	£34	17"W x 17"D. <i>DWR1010E234H.</i>	£34	20"W x 20"D. <i>DWR1010E193K.</i>	£34
12"W x 14"D. <i>DWR1010E189E.</i>	£34	17"W x 19"D. <i>DWR1010E234J.</i>	£34	20"W x 22"D. <i>DWR1010E193M.</i>	£34
14"W x 12"D. <i>DWR1010E190C.</i>	£34	18"W x 16"D. <i>DWR1010E192G.</i>	£34	22"W x 18"D. <i>DWR1010E194I.</i>	£34
14"W x 14"D. <i>DWR1010E190E.</i>	£34	18"W x 18"D. <i>DWR1010E192I.</i>	£34	22"W x 20"D. <i>DWR1010E194K.</i>	£34
14"W x 16"D. <i>DWR1010E190G.</i>	£34	18"W x 20"D. <i>DWR1010E192K.</i>	£34	22"W x 22"D. <i>DWR1010E194M.</i>	£34
14"W x 18"D. <i>DWR1010E190I.</i>	£34	18"W x 22"D. <i>DWR1010E192M.</i>	£34	24"W x 20"D. <i>DWR1010E236K.</i>	£34
15"W x 18"D. <i>DWR1010E232H.</i>	£34	19"W x 17"D. <i>MEC1804152-1917.</i>	£34	24"W x 22"D. <i>DWR1010E236M.</i>	£34

22" width and greater cushions include a flat rigidizer as standard.

ACM CONTOURED RIGIDISER SUITABLE FOR SLING SEATS

15"W x 18"D. <i>CUS142467.</i>	£56	18"W x 16"D. <i>CUS142466.</i>	£56	20"W x 18"D. <i>CUS142465.</i>	£56
16"W x 16"D. <i>CUS142467.</i>	£56	18"W x 18"D. <i>CUS142466.</i>	£56	20"W x 19"D. <i>CUS142465.</i>	£56
16"W x 17"D. <i>CUS142467.</i>	£56	18"W x 19"D. <i>CUS142466.</i>	£56	20"W x 20"D. <i>CUS142465.</i>	£56
16"W x 18"D. <i>CUS142467.</i>	£56	18"W x 20"D. <i>CUS142466.</i>	£56		

For use on manual wheelchairs only.

NOTES & SPECIAL REQUESTS

Special Order Solution SPP Cushion.

£350

Please detail your special order cushion request.

All special order cushion requests must be approved with the Quantum Custom Department.

Please contact Quantum Inside Sales or Custom Department with questions. Extended lead time will apply.

Cushion part numbers involved will be the cushion part number listed in previous sections modified to meet special order request stated in details below. Drawings are available for sketch details.


FRONT VIEW



RIGHT SIDE VIEW



TOP VIEW



PRICE

TOTAL