

SPARKY

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
	Postcode
Phone	Email
PO No	Order Ref

DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Clients Weight	

BASE MODEL

£3,950

- 2 Colour Options
- Powered Seat Tilt
- Height Adjustable Removable Armrests
- Swing-away Leg Supports
- Tension Adjustable Fold-Forward Back Support
- 8" Head Support
- 4" x 14" Pneumatic Drive Tyres
- 2" x 8" Pneumatic Front Castors



- 10" - 14" Seat Width
- 10" - 14" Seat Depth
- VR2 Control System
- Rear Suspension
- 2 x 40 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight: 8 st 13lbs. (56.7 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

4" x 14" Pneumatic Black Tyres.	Standard
4" x 14" Flat Free Black Tyres.	£125

CASTOR WHEELS

2" x 8" Pneumatic Black Tyres.	Standard
2" x 8" Flat Free Black Tyres.	£55

COLOURS

White	Standard
Blue	FOC

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT WIDTH

10-14" Width

FOC

Set width at:

SEAT DEPTH

10-14" Depth

FOC

Set depth at:

STEALTH SEAT CUSHIONS

2" Seat Cushion Set. <i>Supplied with 10" x 10", 12" x 12" & 14" x 14" seat cushions.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	£83
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	£165
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£209

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability. The nearest size will be provided unless otherwise specified

SEAT POSITIONING

Powered Seat Tilt. <i>0° to 30°. 18" seat pan to floor height.</i>	Standard
--	----------

BACK SUPPORTS & POSITIONING

Tension Adjustable Back Support. <i>Supplied with a 10", 12" & 14" wide overlay set.</i>	Standard
Attendant Controlled Manual Recline. <i>Range of adjustments 95° to 125°.</i>	Standard
ADI Aluminium Back. <i>Skin protection, 4 way stretch, Coolcore cover.</i>	POA

Enter ADI Back Support part number here for quotation.

ARM SUPPORTS

Paediatric Single Post Height Adjustable. <i>With 10"L x 3.5"W arm support pads. Adjustment range 2'.</i>	Standard
---	----------

LEG SUPPORTS

70° Swing-away Leg Supports. <i>Adjustment range 9.5" to 12.5" from seat pan.</i>	Standard
---	----------

Manual Elevating Swing-Away Leg Supports.	LEFT	£60	£60	RIGHT
---	-------------	-----	-----	--------------

Adjustment ranges above are from seat pan to footplates. Individual legrests are supplied with heel loops.

STEALTH HEAD SUPPORTS

Sparky 8" Head Support.	Standard
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£230
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£92

Please visit stealthproducts.com for alternative head supports.

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.	Standard
--	----------

Stealth Positioning 2-Point Padded Pelvic Belt.	Extra Small	£84	Small	£84	} Single Pull	STD*	} Dual Pull	FOC		
	Medium	£84	Large	£84		Rear Pull		FOC	Side Release	STD*
						Push Button		FOC		

Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull.	Extra Small	£131	Small	£131	} Push Button	STD*
	Medium	£131	Large	£131		Side Release

Please visit stealthproducts.com for alternative lateral, pelvic & lower extremity solutions. *STD = Standard option.

TRUNK SOLUTIONS

Stealth 4-Point Harness	Extra Small	£143	Classic	STD*	Contoured	FOC
	Small	£143	Dynamic	STD*	Structured	FOC
	Medium	£143	Vinyl Face	FOC	Zipper Face**	FOC
Chest Strap.	Extra Small	£91	Side Release		Standard	
	Small	£91	Velcro Release		£20	

** Classic Only *STD = Standard option

NON-EXPANDABLE CONTROLS

VR2 Controls.		Standard
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	£100
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

HAND CONTROL MOUNTS

Swing-away.						Standard
Mounting Position.	Left	FOC	Right	Standard		
	Standard Carrot	STD	T-Bar	£35	Extended I-Handle	£35
Hand Control Joystick Handle.	Mushroom	£35	Medium Ball	£35	3" Goal Post	£50
	Chin Cup	£35	Soft Ball	£35	4" Goal Post	£50

ACCESSORIES

Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379
Quantum Back Pack.	£110
XLR USB Charger Port.	£17

NOTES & SPECIAL REQUESTS

PRICE

SPARKY



BASE PRICE

£3,950

COST OPTIONS

TOTAL

3



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
RETAIL UK - 01/04/25 V1

