

SPARKY

DEALERS DETAILS

Dealers Name	
<hr/>	
Order Date	
<hr/>	
Address	
<hr/>	
	Postcode
<hr/>	<hr/>
Phone	Email
<hr/>	<hr/>
PO No	Order Ref
<hr/>	<hr/>

DELIVERY DETAILS

Dealers Address	
<hr/>	
<hr/>	
<hr/>	
	Postcode
<hr/>	<hr/>
Quantum Specialist	
<hr/>	

BASE MODEL

£3,475

- 2 Colour Options
- Powered Tilt
- Height Adjustable Removable Armrests
- Swing-away Legrests
- Tension Adjustable Fold-Forward Back Support
- 8" Head Support
- 4" x 14" Pneumatic Drive Tyres
- 2" x 8" Pneumatic Front Castors



- 10" - 14" Seat Width
- 10" - 14" Seat Depth
- VR2 Control System
- Rear Suspension
- 2 x 40 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:
8 st 13lbs. (56.7 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

4" x 14" Pneumatic Black Tyres.	Standard
4" x 14" Flat Free Black Tyres.	£125

CASTOR WHEELS

2" x 8" Pneumatic Black Tyres.	Standard
2" x 8" Flat Free Black Tyres.	£55

COLOURS

White	Standard
Blue	FOC

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT WIDTH

10-14" Width

FOC

Set width at:

SEAT DEPTH

10-14" Depth

FOC

Set depth at:

STEALTH SEAT CUSHIONS

2" Seat Cushion Set. *Supplied with 10" x 10", 12" x 12" & 14" x 14" seat cushions.*

Standard

Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. *158.7 Kg weight limit.*

£75

Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. *160 Kg weight limit.*

£150

Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. *158.7 Kg weight limit.*

£190

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability. The nearest size will be provided unless otherwise specified

SEAT POSITIONING

Powered Seat Tilt. *0° to 30°. 18" seat pan to floor height.*

Standard

BACK SUPPORTS & POSITIONING

Tension Adjustable Back Support. *Supplied with a 10", 12" & 14" wide overlay set.*

Standard

Attendant Controlled Manual Recline. *Range of adjustments 95° to 125°.*

Standard

ADI Aluminium Back. *Skin protection, 4 way stretch, Coolcore cover.*

POA

Enter ADI Back Support part number here for quotation.

ARM SUPPORTS

Paediatric Single Post Height Adjustable. *With 10"L x 3.5"W arm support pads. Adjustment range 2".*

Standard

LEG SUPPORTS

70° Swing-away Legrests. *Adjustment range 9.5" to 12.5" from seat pan.*

Standard

Manual Elevating Swing-Away Legrests.

LEFT

£60

£60

RIGHT

Adjustment ranges above are from seat pan to footplates. Individual legrests are supplied with heel loops.

STEALTH HEAD SUPPORTS

Sparky 8" Head Support.

Standard

UniLink 10" Comfort Plus. *HMO475-17M / CP10CCFPM.*

£219

Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.

£88

Please visit stealthproducts.com for alternative head supports.

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.

Standard (STD)

Stealth Positioning
2-Point Padded
Pelvic Belt.

Extra Small £78

Small £78

Medium £78

Large £78

Single Pull STD*

Dual Pull FOC

Rear Pull FOC

Side Release STD*

Push Button FOC

Stealth Positioning 4-Point Linear
Padded Pelvic Belt, Single Pull.

Extra Small £120

Small £120

Medium £120

Large £120

Push Button STD*

Side Release £20

Please visit stealthproducts.com for alternative lateral, pelvic & lower extremity solutions. *STD = Standard option.

TRUNK SOLUTIONS

Stealth 4-Point Harness	Extra Small	£137	} Classic	STD*	Contoured	FOC	
	Small	£137		Dynamic	STD*	Structured	FOC
	Medium	£137		Vinyl Face	FOC	Zipper Face**	FOC
Chest Strap.	Extra Small	£87	} Side Release	Standard			
	Small	£87		Velcro Release	£19		

*STD = Standard option

NON-EXPANDABLE CONTROLS

VR2 Controls.				Standard
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>			Standard
	Attendant Control. <i>Hand control mounted at rear.</i>			£100
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>			£380

HAND CONTROL MOUNTS

Swing-away.						Standard
Mounting Position.	Left	FOC	Right	Standard		
	Standard Carrot	STD	T-Bar	£35	Extended I-Handle	£35
Hand Control Joystick Handle.	Mushroom	£35	Medium Ball	£35	3" Goal Post	£50
	Chin Cup	£35	Soft Ball	£35	4" Goal Post	£50

ACCESSORIES

Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379
Quantum Back Pack.	£110
Additional XLR USB Charger Port.	£17

NOTES & SPECIAL REQUESTS

PRICE

SPARKY



BASE PRICE **£3,475**

COST OPTIONS

TOTAL

3



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
RETAIL UK - 18/01/23 V1



Crash Tested



V1: 14/06/2022