



RETAIL

SPECTRUM FOAM SPP CUSHION

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
Postcode	
Email	
Phone	Fax

DELIVERY DETAILS

Dealers Address		
Postcode		
PO No		
Order Type	Quote	Order

DETAILS

This cushion offers skin protection and positioning with a Coolcore, 4-way stretch, breathable cover and an urthane liner.

SunMate FRG foam top (SunMate FRG XXSoft (blue); Medsoft (green)) with multi-layered foam insert or twin cell ultra-lightweight viscous gel insert in ischial well, sacral relief



Weight capacities
158.75 kg: 16"-20" width,
294.83 kg: 22" width+
Overall weight: 1.3 kg
Overall height 3"
14x18" up to 32x28"
Odd sizes in 17" up to 19"
Custom sizes available

CUSHION

		ADD WIDTH		ADD DEPTH		Width															
						14	16	17	18	19	20	22	24	26	28	30	32				
Spectrum Foam SPP. 16"-20"W.		ST-SPCFSP			£265	15															
Spectrum Foam SPP. 22"-24"W		ST-SPCFSP			£310	16															
Spectrum Foam SPP. 26"-32"W		ST-SPCFSP			£310	17															
Cushion Cover.	Coolcore.	SPCFSP-NBLK			STD	18															
	Mesh.	SPCFSP-MBLK			FOC	19															
	Stealth Tek.	SPCFSP-TBLK			FOC	20															
Fluid Proof Liner.	Urethane.	URL1709070			STD	21															
	Polyurethane.	PUL1709030			£48	22															
Extra Covers And Liners.	Coolcore Cover.	SPCFSP-NBLK			£70	24															
	Mesh Cover.	SPCFSP-MBLK			£44	26															
	Stealth Tek Cover.	SPCFSP-TBLK			£50	28															
	Urethane Liner.	URL1709070			£30																
	Polyurethane Liner.	PUL1709030			£48																

Select your size, then complete the Part Number by adding the Width first then Depth.

ADD SOLID BASE

ABS FLAT RIGIDISER

14"W x 18"D. DWR1010E1910I . £34	19"W x 17"D. MEC1804152-1917 . £34	24"W x 21"D. DWR1010E236K . £34
16"W x 15"D. MEC1804152-1615 . £34	19"W x 18"D. MEC1804152-1917 . £34	24"W x 22"D. DWR1010E236M . £34
16"W x 16"D. DWR1010E191G . £34	19"W x 19"D. MEC1804152-1919 . £34	24"W x 24"D. DWR1010E236O . £34
16"W x 17"D. DWR1010E191G . £34	19"W x 20"D. MEC1804152-1919 . £34	24"W x 26"D. DWR1010E236Q . £34
16"W x 18"D. DWR1010E191I . £34	20"W x 16"D. DWR1010E193G . £34	26"W x 18"D. DWR1010E240I . £34
16"W x 20"D. DWR1010E191I . £34	20"W x 18"D. DWR1010E193I . £34	26"W x 20"D. DWR1010E240K . £34
17"W x 16"D. MEC1804152-1715 . £34	20"W x 19"D. DWR1010E193I . £34	26"W x 22"D. DWR1010E240M . £34
17"W x 17"D. DWR1010E234H . £34	20"W x 20"D. DWR1010E193K . £34	26"W x 24"D. DWR1010E240O . £34
17"W x 18"D. DWR1010E234H . £34	20"W x 21"D. DWR1010E193K . £34	28"W x 18"D. DWR1010E244I . £34
17"W x 19"D. DWR1010E234J . £34	20"W x 22"D. DWR1010E193M . £34	28"W x 20"D. DWR1010E244K . £34
17"W x 20"D. DWR1010E234J . £34	20"W x 24"D. DWR1010E19O . £34	28"W x 22"D. DWR1010E244M . £34
18"W x 16"D. DWR1010E192G . £34	22"W x 16"D. DWR1010E194G . £34	28"W x 24"D. DWR1010E244O . £34
18"W x 17"D. DWR1010E192G . £34	22"W x 18"D. DWR1010E194I . £34	28"W x 28"D. DWR1010E244S . £34
18"W x 18"D. DWR1010E192I . £34	22"W x 19"D. DWR1010E194I . £34	30"W x 18"D. DWR1010E341II . £34
18"W x 19"D. DWR1010E192I . £34	22"W x 20"D. DWR1010E194K . £34	30"W x 20"D. DWR1010E341IK . £34
18"W x 20"D. DWR1010E192K . £34	22"W x 22"D. DWR1010E194M . £34	30"W x 22"D. DWR1010E341IM . £34
18"W x 21"D. DWR1010E192K . £34	22"W x 24"D. DWR1010E194O . £34	30"W x 24"D. DWR1010E341IO . £34
18"W x 22"D. DWR1010E192M . £34	24"W x 16"D. DWR1010E194G . £34	30"W x 28"D. DWR1010E341IS . £34
18"W x 24"D. DWR1010E192O . £34	24"W x 18"D. DWR1010E236I . £34	32"W x 26"D. MEC1804152-3226 . £34
19"W x 16"D. MEC1804152-1915 . £34	24"W x 20"D. DWR1010E236K . £34	32"W x 28"D. MEC1804152-3228 . £34

22" width and greater cushions include a flat rigidizer as standard.

ACM CONTOURED RIGIDISER SUITABLE FOR SLING SEATS

16"W x 16"D. CUS142467 . £56	18"W x 18"D. CUS142466 . £56	20"W x 19"D. CUS142465 . £56
16"W x 17"D. CUS142467 . £56	18"W x 19"D. CUS142466 . £56	20"W x 20"D. CUS142465 . £56
16"W x 18"D. CUS142467 . £56	18"W x 20"D. CUS142466 . £56	
16"W x 16"D. CUS142466 . £56	20"W x 18"D. CUS142465 . £56	

For use on manual wheelchairs only.


NOTES & SPECIAL REQUESTS

Special Order Spectrum Foam SPP Cushion.

£450

Please detail your special order cushion request. All special order cushion requests must be approved with the Quantum Custom Department. Please contact Quantum Inside Sales or Custom Department with questions. Extended lead time will apply.
Cushion part numbers involved will be the cushion part number listed in previous sections modified to meet special order request stated in details below. Drawings are available for sketch details.


TOP VIEW



FRONT VIEW



RIGHT SIDE VIEW



PRICE

TOTAL

2



Stealth Products

32 Wedgwood Road, Bicester,
Oxon, OX26 4UL

Phone: 01869 324600

Fax: 01869 323070

Email: sales@quantumrehab.co.uk

Web: www.quantumrehab.co.uk

SPECTRUM FOAM SPP - 02/11/23 V1



V1: 14/06/2022

