



RETAIL

ZEN SP CUSHION

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Email _____

Phone _____ Fax _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____


PO No _____

Order Type _____ Quote _____ Order _____

DETAILS

This cushion offers skin protection and a Coolcore 4-way stretch, breathable cover.

New SIM G base, SunMate FRG XXSoft (blue topper), low profile, lightweight



Weight capacities
 79.37 kg: 10"-14" width
 160.11 kg: 15"-21" width,
 181.43: 22" width+
 Overall weight: 1.5kg
 Overall height 2.5"
 10x10" up to 24x24"
 Odd & Custom sizes available

CUSHION

		ADD WIDTH		ADD DEPTH		Width											
						10	12	14	15	16	17	18	19	20	21	22	24
Zen SP.	10"-21" W.	ST-ZENSP	<input type="text"/>	<input type="text"/>	£150	10	12	14	15	16	17	18	19	20	21	22	24
Zen SP.	22"-24" W	ST-ZENSP	<input type="text"/>	<input type="text"/>	£170	10	12	14	15	16	17	18	19	20	21	22	24
Cushion Cover.	Coolcore.	ZENSP-NBLK	<input type="text"/>	<input type="text"/>	STD	10	12	14	15	16	17	18	19	20	21	22	24
	Mesh.	ZENSP-MBLK	<input type="text"/>	<input type="text"/>	FOC	10	12	14	15	16	17	18	19	20	21	22	24
	Stealth Tek.	ZENSP-TBLK	<input type="text"/>	<input type="text"/>	FOC	10	12	14	15	16	17	18	19	20	21	22	24
Fluid Proof Liner.	Urethane.	URL1709070	<input type="text"/>	<input type="text"/>	£40	10	12	14	15	16	17	18	19	20	21	22	24
	Polyurethane.	PUL1709030	<input type="text"/>	<input type="text"/>	£48	10	12	14	15	16	17	18	19	20	21	22	24
Extra Covers And Liners.	Coolcore Cover.	ZENSP-NBLK	<input type="text"/>	<input type="text"/>	£70	10	12	14	15	16	17	18	19	20	21	22	24
	Mesh Cover.	ZENSP-MBLK	<input type="text"/>	<input type="text"/>	£44	10	12	14	15	16	17	18	19	20	21	22	24
	Stealth Tek Cover.	ZENSP-TBLK	<input type="text"/>	<input type="text"/>	£50	10	12	14	15	16	17	18	19	20	21	22	24
	Urethane Liner.	URL1709070	<input type="text"/>	<input type="text"/>	£40	10	12	14	15	16	17	18	19	20	21	22	24
	Polyurethane Liner.	PUL1709030	<input type="text"/>	<input type="text"/>	£48	10	12	14	15	16	17	18	19	20	21	22	24

Select your size, then complete the Part Number by adding the Width first then Depth.

ADD SOLID BASE

ABS FLAT RIGIDISER

10"W x 10"D. <i>DWR1010E188A.</i>	£34	16"W x 20"D. <i>DWR1010E191K.</i>	£34	20"W x 22"D. <i>DWR1010E193M.</i>	£34
10"W x 12"D. <i>DWR1010E188C.</i>	£34	17"W x 17"D. <i>DWR1010E234H.</i>	£34	21"W x 19"D. <i>MEC1804152-2119.</i>	£34
10"W x 14"D. <i>DWR1010E188E.</i>	£34	17"W x 18"D. <i>DWR1010E234H.</i>	£34	21"W x 20"D. <i>MEC1804152-2119.</i>	£34
12"W x 10"D. <i>DWR1010E189A.</i>	£34	17"W x 19"D. <i>DWR1010E234J.</i>	£34	21"W x 21"D. <i>MEC1804152-2121.</i>	£34
12"W x 12"D. <i>DWR1010E189C.</i>	£34	17"W x 20"D. <i>DWR1010E234J.</i>	£34	22"W x 16"D. <i>DWR1010E193G.</i>	£34
12"W x 14"D. <i>DWR1010E189E.</i>	£34	18"W x 16"D. <i>DWR1010E192G.</i>	£34	22"W x 18"D. <i>DWR1010E194I.</i>	£34
14"W x 10"D. <i>MEC1804152-1410.</i>	£34	18"W x 18"D. <i>DWR1010E192I.</i>	£34	22"W x 20"D. <i>DWR1010E194K.</i>	£34
14"W x 12"D. <i>DWR1010E190C.</i>	£34	18"W x 20"D. <i>DWR1010E192K.</i>	£34	22"W x 22"D. <i>DWR1010E194M.</i>	£34
14"W x 14"D. <i>DWR1010E190E.</i>	£34	18"W x 22"D. <i>DWR1010E192M.</i>	£34	22"W x 24"D. <i>DWR1010E193O.</i>	£34
14"W x 16"D. <i>DWR1010E190G.</i>	£34	19"W x 17"D. <i>MEC1804152-1917.</i>	£34	24"W x 18"D. <i>DWR1010E236I.</i>	£34
14"W x 16"D. <i>DWR1010E190G.</i>	£34	19"W x 18"D. <i>MEC1804152-1917.</i>	£34	24"W x 19"D. <i>DWR1010E236I.</i>	£34
15"W x 15"D. <i>MEC1804152-1515.</i>	£34	19"W x 19"D. <i>MEC1804152-1919.</i>	£34	24"W x 20"D. <i>DWR1010E236K.</i>	£34
15"W x 16"D. <i>MEC1804152-1515.</i>	£34	19"W x 20"D. <i>MEC1804152-1919.</i>	£34	24"W x 21"D. <i>DWR1010E236K.</i>	£34
15"W x 17"D. <i>DWR1010E232H.</i>	£34	20"W x 16"D. <i>DWR1010E193G.</i>	£34	24"W x 22"D. <i>DWR1010E236M.</i>	£34
16"W x 16"D. <i>DWR1010E191G.</i>	£34	20"W x 18"D. <i>DWR1010E193I.</i>	£34	24"W x 23"D. <i>DWR1010E236M.</i>	£34
16"W x 18"D. <i>DWR1010E191I.</i>	£34	20"W x 20"D. <i>DWR1010E193K.</i>	£34	24"W x 24"D. <i>DWR1010E236O.</i>	£34

ACM CONTOURED RIGIDISER SUITABLE FOR SLING SEATS

16"W x 16"D. <i>CUS142467-1616.</i>	£56	17"W x 20"D. <i>CUS142466.</i>	£56	19"W x 19"D. <i>CUS142465-2018.</i>	£56
16"W x 18"D. <i>CUS142467.</i>	£56	18"W x 16"D. <i>CUS142466-1816.</i>	£56	19"W x 20"D. <i>CUS142465.</i>	£56
16"W x 20"D. <i>CUS142467.</i>	£56	18"W x 18"D. <i>CUS142466-1818.</i>	£56	20"W x 18"D. <i>CUS142465-2018.</i>	£56
17"W x 17"D. <i>CUS142466-1816.</i>	£56	18"W x 20"D. <i>CUS142466.</i>	£56	20"W x 20"D. <i>CUS142465.</i>	£56
17"W x 18"D. <i>CUS142466-1818.</i>	£56	18"W x 22"D. <i>CUS142466.</i>	£56	20"W x 22"D. <i>CUS142465.</i>	£56
17"W x 19"D. <i>CUS142466-1818.</i>	£56	19"W x 18"D. <i>CUS142465-2018.</i>	£56		


For use on manual wheelchairs only.

NOTES & SPECIAL REQUESTS

Special Order Zen SP Cushion.

£300

TOP VIEW



FRONT VIEW



RIGHT SIDE VIEW



PRICE

TOTAL

2



Stealth Products

32 Wedgwood Road, Bicester,
Oxon, OX26 4UL

Phone: 01869 324600

Fax: 01869 323070

Email: sales@quantumrehab.co.uk

Web: www.quantumrehab.co.uk

ZEN SP - 15/01/24 V1



V1: 14/06/2022