

ADI



Designed with you in mind

Stealth Products
Solutions For Life



Lightweight Positioning Backs

The lightest back supports on the market

CF SERIES

3 Heights / 2 Contours



Fits Chair Widths

12in (30.48 cm) to 19in (48.26 cm)

Support Heights**

Low (10in - 25.40cm)

Medium (13in - 33.02cm)

Standard (16in - 40.64cm)

Contour Depth

Active (2.5in - 6.35cm)

Deep (4.5in - 11.43cm)

Material

Carbon Fiber Material

Weight

.820lb (370 grams)

approx. weight of a 16in -40.64cm- standard height back

Note:

Custom options available

** Support Heights are measured on a base shell, cushion and cover will add approximately 2.5in (6.35cm) of overall height.

AL SERIES

5 Heights 3 Contours



Fits Chair Widths

10in (25.40cm) through 20in (50.80cm)

Support Heights**

Low (10in - 25.40cm)

Medium (13in - 33.02cm)

Standard (16in - 40.64cm)

Tall (18in - 45.72cm)

Extra Tall (20in - 50.80cm)

Contour Depth

Active (2.5in - 6.35cm)

Deep (4.5in - 11.43cm)

Extra Deep (6in - 15.24cm)

HPCPS Coding (USA only)

Active: E2613

Deep / Extra Deep: E2615

Material

Lightweight aircraft quality aluminum

Weight

1.5lb (680 grams)

approx. weight of a 16in -40.64cm- standard height back

Note:

Custom options available

** Support Heights are measured on a base shell, cushion and cover will add approximately 2.5in (6.35cm) of overall height.

AFT SERIES

5 Heights 3 Contours



Material

Lightweight aircraft quality aluminum

Weight

1.5lb (680 grams)

approx. weight of a 16in -40.64cm- standard height back

Note:

Custom options available

** Support Heights are measured on a base shell, cushion and cover will add approximately 2.5in (6.35cm) of overall height.

MOUNTING

REMOVABLE



Quick Release

All

1.3 lbs (589gm)

Removable mounting with adjustable positioning



Basic Release

AL/AFT Series Only

1.8 lbs (816gm)

ADJUSTABLE



2-Point PRO™

All

1.25 lbs (566gm)

Stationary mounting with adjustable positioning



4-Point PRO™

ALL

1.32 lbs (596gm)

FIXED



Fixed ELITE™

CF Series Only

.60 lbs (272gm)

Stationary mounting with fixed positioning

No Adjustments



Width Adjustable ELITE

ALL

.63 lbs (285gm)

Width: +/- 1in (2.54cm)

Angle: 20° | Depth: 1.25in (31.75cm) | Width: +/- 1in (2.54cm)

For folding frame chair users that require angle, depth and width adjustment.