



DYNO™



- Folds Compactly ◇
- 3 Wheel Configuration ◇
- Swiveling Front Wheel ◇
- ACTIVATOR™ Dynamic Seating ◇
- 30° Hip Angle Adjustment ◇
- Lightweight ◇
- Contoured Seating System ◇



Call: 01869 324600 • Web: www.quantumrehab.co.uk • Email: sales@quantumrehab.co.uk

© 2019 LEGGERO LLC • Prices and specifications are subject to change without notice

DYNO™



The only WC-19
transit tested
compact folding
3-wheeler!

Available with
Activator™
Dynamic Seating



SPECIFICATIONS & ADJUSTMENT RANGE

| | |
|---|---------------------------------------|
| Weight Limit | 80 lbs / 36 kg |
| Weight of Chair | 26.5 lbs / 12 kg |
| Chair Width | 26.5 in / 67 cm |
| Seat Depth | 6-14 in / 15-35 cm |
| Seat Width | 14 in / 35 cm |
| Back Height | 13-21 in / 33-53 cm |
| Lower Leg Length | 8-14 in / 20-35 cm |
| Chest Width Range | 7- 14 in / 18-35 cm |
| Seat to Floor | 19 in / 48 cm |
| Hip Angle Adjustment | 85°, 90°, 105° |
| Adjustable Seat Tilt | 25-30° |
| Knee Active Range (optional) | 10° |
| Seat to Back Active Range (optional) | 30° |
| Push Handle Height | 39-46 in / 99-117 cm |
| Front Wheel Diameter | 10 in / 25 cm |
| Rear Wheel Diameter | 16 in / 41 cm |
| Folded with Wheels | 42 x 17 x 26.5 in 107 x 43 x 67 cm |
| Folded without Wheels | 35 x 17 x 25 in 89 x 43 x 64 cm |

Freedom to go where you want and the
positioning options you need.



For more information visit www.quantumrehab.co.uk or call 01869 324600