



# Dyno SE UK Retail Order Form

**DYNO**

June 2021

Stealth UK Customer Service: 01869324600

Email: sales@quantumrehab.co.uk

Date: \_\_\_\_\_ ☐ QUOTE ☐ ORDER - PO# \_\_\_\_\_

Company Name: \_\_\_\_\_

Ship To: \_\_\_\_\_

Contact Name: \_\_\_\_\_

Attention: \_\_\_\_\_

Address: \_\_\_\_\_

Billing Address: \_\_\_\_\_

Post Code: \_\_\_\_\_

Post Code: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Email: \_\_\_\_\_

Mark for: \_\_\_\_\_



Model shown with accessories.

## DYNO SE™ INCLUDES:


- Frame Suspension
- Hip Angle Adjustment
- QR 16" Rear Wheels
- 10" Front Wheel with Swivel Lock
- Height Adjustable Push Handle
- Attendant Wheel Lock
- Under Seat Storage Bag
- Rear Pocket
- Sunshade

## Specifications

|                            |               |             |
|----------------------------|---------------|-------------|
| Max Occupant Weight        | 80 lbs.       | 36 kg       |
| Push Handle Height         | 30-43 in      | 76-109 cm   |
| Weight                     | 26.5 lbs.     | 12 kg       |
| Overall Width              | 26.5 in       | 67 cm       |
| Front Wheel Diameter       | 10 in         | 25 cm       |
| Rear Wheel Diameter        | 16 in         | 41 cm       |
| Folded without rear wheels | 32x17x25 in   | 82x43x64 cm |
| Folded with wheels         | 39x17x26.5 in | 99x43x67 cm |
| Reversing Width            | 41 in         | 104 cm      |
| Tilt Positions             | 25-30°        |             |
| Maximum Safe Slope         | 5°            |             |
| Knee Active Range          | 10°           |             |
| Seat-to-back Active Range  | 30°           |             |

## MEASUREMENT WORKSHEET

- A. Seat Depth \_\_\_\_\_
- B. Trunk Depth \_\_\_\_\_
- C. Trunk Width \_\_\_\_\_
- D. Hip Width \_\_\_\_\_
- E. Seat to Top of Shoulder \_\_\_\_\_
- F. Lower Leg Length \_\_\_\_\_
- G. Seat to Axilla \_\_\_\_\_



## Adjustability

|                      |                |          |
|----------------------|----------------|----------|
| Seat Depth           | 8.5-13.5 in    | 22-34 cm |
| Seat Width           | 14 in          | 36 cm    |
| Lower Leg Length     | 7.5-14 in      | 19-36 cm |
| Back Height          | 11.5-21 in     | 29-53 cm |
| Hip Angle Adjustment | 85°, 90°, 105° |          |
| Chest Width Range    | 7-14 in        | 18-35 cm |

|                    |
|--------------------|
| <b>DYNO™ Chair</b> |
|--------------------|

|                 |
|-----------------|
| <b>DYNO SE™</b> |
|-----------------|

|                          |           |       |          |           |
|--------------------------|-----------|-------|----------|-----------|
| <input type="checkbox"/> | 3001-1002 | E1238 | DYNO SE™ | £1,499.00 |
|--------------------------|-----------|-------|----------|-----------|

|                         |
|-------------------------|
| <b>Upholstery Color</b> |
|-------------------------|

|                          |           |                             |     |
|--------------------------|-----------|-----------------------------|-----|
| <input type="checkbox"/> | 3001-0101 | BIG SKY BLUE - Blue/Black   | N/C |
| <input type="checkbox"/> | 3001-0102 | ROSSO RED - Red/Black       | N/C |
| <input type="checkbox"/> | 3001-0103 | ICE SILVER - Silver/Black   | N/C |
| <input type="checkbox"/> | 3001-0104 | SNOWBERRY PINK - Pink/Black | N/C |

*Upholstery Color applies to back, seat cushions, and sunshade*

|                                                       |
|-------------------------------------------------------|
| <b>Transit Option &amp; Transit Harness Selection</b> |
|-------------------------------------------------------|

|                       |
|-----------------------|
| <b>Transit Option</b> |
|-----------------------|

*Requirements: Sure-Lok belt, 3001-4001 back rest and a head rest.*

|                          |           |                |         |
|--------------------------|-----------|----------------|---------|
| <input type="checkbox"/> | 3001-7007 | Transit Option | £187.00 |
|--------------------------|-----------|----------------|---------|

*This option is not retrofittable.*

*The transit option includes 4 installed transportation anchors and belt mount studs*

|                                |
|--------------------------------|
| <b>Transit Belt or Harness</b> |
|--------------------------------|

|                          |           |       |                           |        |
|--------------------------|-----------|-------|---------------------------|--------|
| <input type="checkbox"/> | 3001-7008 | K0108 | Transit Sure-Lok Lap Belt | £68.00 |
|--------------------------|-----------|-------|---------------------------|--------|

|                                   |
|-----------------------------------|
| <b>Features &amp; Accessories</b> |
|-----------------------------------|

|                                     |           |                              |     |
|-------------------------------------|-----------|------------------------------|-----|
| <input checked="" type="checkbox"/> | 3001-5001 | Height Adjustable Foot Plate | N/C |
| <input checked="" type="checkbox"/> | 3001-7009 | Attendant Wheel Lock         | N/C |
| <input checked="" type="checkbox"/> | 3001-7010 | Under Seat Storage           | N/C |
| <input checked="" type="checkbox"/> | 3001-7012 | Deluxe Sunshade              | N/C |

|                                     |           |                                            |         |
|-------------------------------------|-----------|--------------------------------------------|---------|
| <input type="checkbox"/>            | 3001-7015 | DYNO Grab Bar                              | £101.00 |
| <input type="checkbox"/>            | 3001-7020 | DYNO Rain Cover                            | £146.00 |
| <input checked="" type="checkbox"/> | 3001-2001 | 10 in (25 cm) Front Wheel with Pneu. Tire  | N/C     |
| <input checked="" type="checkbox"/> | 3001-2002 | 16 in (41 cm) Rear Wheels with Pneu. Tires | N/C     |

|                                     |
|-------------------------------------|
| <b>ACTIVATOR™ Dynamic Component</b> |
|-------------------------------------|

*If ACTIVATOR™ Dynamic Component is selected, select spring to right.*

|                                        |
|----------------------------------------|
| <b>A. ACTIVATOR™ Dynamic Component</b> |
|----------------------------------------|

|                          |           |       |                              |         |
|--------------------------|-----------|-------|------------------------------|---------|
| <input type="checkbox"/> | 3001-7018 | E2295 | ACTIVATOR™ Dynamic Component | £295.00 |
|--------------------------|-----------|-------|------------------------------|---------|

|                                                   |
|---------------------------------------------------|
| <b>ACTIVATOR™ Springs (Only available with A)</b> |
|---------------------------------------------------|

|                          |           |         |                               |        |
|--------------------------|-----------|---------|-------------------------------|--------|
| <input type="checkbox"/> | 3001-7003 |         | Light Springs for ACTIVATOR™  | N/C    |
| <input type="checkbox"/> | 3001-7004 | Default | Medium Springs for ACTIVATOR™ | N/C    |
| <input type="checkbox"/> | 3001-7005 |         | Stiff Springs for ACTIVATOR™  | N/C    |
| <input type="checkbox"/> | 3001-7006 |         | Full ACTIVATOR™ Spring Kit    | £38.00 |

|                                       |
|---------------------------------------|
| <b>Back Rest &amp; Trunk Supports</b> |
|---------------------------------------|

*Select back rest and review related options to the right.*

|                          |
|--------------------------|
| <b>A. Contoured Back</b> |
|--------------------------|

|                          |           |       |                             |     |
|--------------------------|-----------|-------|-----------------------------|-----|
| <input type="checkbox"/> | 3001-4001 | E2293 | Contoured Configurable Back | N/C |
|--------------------------|-----------|-------|-----------------------------|-----|

|                                                             |
|-------------------------------------------------------------|
| <b>Fixed Lateral Trunk Supports (Only available with A)</b> |
|-------------------------------------------------------------|

|                          |           |       |                                                    |     |
|--------------------------|-----------|-------|----------------------------------------------------|-----|
| <input type="checkbox"/> | 3003-6001 | E0956 | LEFT, Lateral Trunk Support<br>3x4 in (7.5x10 cm)  | N/C |
| <input type="checkbox"/> | 3003-6002 | E0956 | RIGHT, Lateral Trunk Support<br>3x4 in (7.5x10 cm) | N/C |

|                                                                  |
|------------------------------------------------------------------|
| <b>Swing Away Lateral Trunk Supports (Only available with A)</b> |
|------------------------------------------------------------------|

*Must select hardware & pad.*

|                          |            |           |       |                                     |         |
|--------------------------|------------|-----------|-------|-------------------------------------|---------|
| <input type="checkbox"/> | 3002-6111  |           | E1028 | LEFT, Leggero Swing In/Out Lateral  | £168.00 |
| <input type="checkbox"/> | 3002-6112  |           | E1028 | RIGHT, Leggero Swing In/Out Lateral | £168.00 |
| <input type="checkbox"/> | TWBTL-L-01 |           | E1028 | LEFT, Stealth Thin Lateral          | £152.00 |
| <input type="checkbox"/> | TWBTL-R-01 |           | E1028 | RIGHT, Stealth Thin Lateral         | £152.00 |
| L                        | R          | 3004-6035 | E0956 | Lateral Pad 3x5", Each              | £34.00  |
| L                        | R          | 3004-6045 | E0956 | Lateral Pad 4x5", Each              | £34.00  |
| L                        | R          | 3004-6046 | E0956 | Lateral Pad 4x6", Each              | £34.00  |

### Seat Cushion

|                          |           |       |                                       |     |
|--------------------------|-----------|-------|---------------------------------------|-----|
| <input type="checkbox"/> | 3001-3001 | E2292 | DYNO Cascading Contoured Seat Cushion | N/C |
|--------------------------|-----------|-------|---------------------------------------|-----|

### Additional Positioning

#### DYNO Lateral/Pelvic Thigh Supports

|                                                    |           |       |                                     |     |
|----------------------------------------------------|-----------|-------|-------------------------------------|-----|
| <input type="checkbox"/>                           | 3001-6001 | E0956 | LEFT, Lateral/Pelvic Thigh Support  | N/C |
| Select Pad Thickness:                              |           |       |                                     |     |
| <div><div>1"</div><div>2"</div><div>3"</div></div> |           |       |                                     |     |
| <input type="checkbox"/>                           | 3001-6002 | E0956 | RIGHT, Lateral/Pelvic Thigh Support | N/C |
| Select Pad Thickness:                              |           |       |                                     |     |
| <div><div>1"</div><div>2"</div><div>3"</div></div> |           |       |                                     |     |

Dyno hip width starts at 14". Adding 2 ea. x 2" pads reduces hip width to 14"-4"= 10".

#### Abductor/Medial Knee Block

Must select hardware & one size of abductor pad. Wedge sizes are width x depth x height.

|                          |           |       |                               |         |
|--------------------------|-----------|-------|-------------------------------|---------|
| <input type="checkbox"/> | 3005-7001 | E1028 | Leggero Flip Down Medial Hdw. | £180.00 |
| <input type="checkbox"/> | MTMD-12   | E0957 | Small Wedge Pad (3x5x3")      | £96.00  |
| <input type="checkbox"/> | MTMD-13   | E0957 | Medium Wedge Pad (5x5x3")     | £96.00  |

#### EVA Positioning Wedges for Seat and/or Back

|                          |           |       |                               |        |
|--------------------------|-----------|-------|-------------------------------|--------|
| <input type="checkbox"/> | 3001-7011 | K0108 | Wedges QTY 4 and Lumbar QTY 1 | £38.00 |
|--------------------------|-----------|-------|-------------------------------|--------|

### Positioning Tray

|                          |           |       |                               |         |
|--------------------------|-----------|-------|-------------------------------|---------|
| <input type="checkbox"/> | 3001-7014 | E1028 | Tray Hardware, Adj./Removable | £157.00 |
| <input type="checkbox"/> | 3001-7013 | E0950 | DYNO Clear Positioning Tray   | £112.00 |

3001-7014 Tray hardware is required

### Ankle Supports

|                          |           |       |                               |        |
|--------------------------|-----------|-------|-------------------------------|--------|
| <input type="checkbox"/> | 3001-6007 | E0951 | LEFT, Extra Small (Pad 2x8")  | £68.00 |
| <input type="checkbox"/> | 3001-6008 | E0951 | RIGHT, Extra Small (Pad 2x8") | £68.00 |
| <input type="checkbox"/> | 3001-6009 | E0951 | LEFT, Small (Pad 3x10")       | £68.00 |
| <input type="checkbox"/> | 3001-6010 | E0951 | RIGHT, Small (Pad 3x10")      | £68.00 |

### Stealth Products Trunk & Pelvic Positioning

#### SP™ Anterior Posture Supports

|                          |             |       |                               |         |
|--------------------------|-------------|-------|-------------------------------|---------|
| <input type="checkbox"/> | IFPS1101-01 | E0960 | SP™ Classic, Dynamic, XS      | £100.00 |
| <input type="checkbox"/> | IFPS1102-01 | E0960 | SP™ Classic, Dynamic, S       | £100.00 |
| <input type="checkbox"/> | IFPS1103-01 | E0960 | SP™ Classic Dynamic, M        | £100.00 |
| <input type="checkbox"/> | IFPS1201-01 | E0960 | SP™ Classic, Structured, XS   | £100.00 |
| <input type="checkbox"/> | IFPS1202-01 | E0960 | SP™ Classic, Structured, S    | £100.00 |
| <input type="checkbox"/> | IFPS1203-01 | E0960 | SP™ Classic, Structured, M    | £100.00 |
| <input type="checkbox"/> | IFPS2101-01 | E0960 | SP™ Contoured, Dynamic, XS    | £100.00 |
| <input type="checkbox"/> | IFPS2102-01 | E0960 | SP™ Contoured, Dynamic, S     | £100.00 |
| <input type="checkbox"/> | IFPS2103-01 | E0960 | SP™ Contoured, Dynamic, M     | £100.00 |
| <input type="checkbox"/> | IFPS2201-01 | E0960 | SP™ Contoured, Structured, XS | £100.00 |
| <input type="checkbox"/> | IFPS2202-01 | E0960 | SP™ Contoured, Structured, S  | £100.00 |
| <input type="checkbox"/> | IFPS2203-01 | E0960 | SP™ Contoured, Structured, M  | £100.00 |

#### SP™ Pelvic Positioning Belts

|                          |          |       |                               |        |
|--------------------------|----------|-------|-------------------------------|--------|
| <input type="checkbox"/> | IFPB5141 | E0978 | 4- Point Push Button 1", XS   | £90.00 |
| <input type="checkbox"/> | IFPB5242 | E0978 | 4- Point Push Button 1.5", SM | £90.00 |
| <input type="checkbox"/> | IFPB5111 | E0978 | 2-Point, Push Button 1", XS   | £76.00 |
| <input type="checkbox"/> | IFPB5112 | E0978 | 2-Point, Push Button 1", SM   | £76.00 |

Side release available on request.

#### SP™ Buckle Covers

|                          |          |  |                              |        |
|--------------------------|----------|--|------------------------------|--------|
| <input type="checkbox"/> | IFMDBC10 |  | Buckle Cover Kit 1" Buckle   | £21.00 |
| <input type="checkbox"/> | IFMDBC15 |  | Buckle Cover Kit 1.5" Buckle | £21.00 |

#### Posture Support Strap Risers

|                          |           |       |                    |        |
|--------------------------|-----------|-------|--------------------|--------|
| <input type="checkbox"/> | 3001-6011 | K0108 | Strap Risers, Pair | £33.00 |
|--------------------------|-----------|-------|--------------------|--------|

### Leggero Head Positioning

|                          |           |       |                                       |         |
|--------------------------|-----------|-------|---------------------------------------|---------|
| <input type="checkbox"/> | 3001-6005 | E0955 | 8" Oval Headrest and Hardware         | £146.00 |
| <input type="checkbox"/> | 3001-6006 | E0955 | RESPOND™ Cradle Headrest and Hardware | £146.00 |
| <input type="checkbox"/> | 3005-6099 |       | Omit Head Positioning                 | N/C     |

Hardware is included with 8" Multi-Adjustable & RESPOND™ Headrests

### Stealth Products Head Positioning

Select **A or B**. Then continue selections to right to complete the package.

#### A. Link Headrest Hardware 1" Ball

|                          |             |       |                       |         |
|--------------------------|-------------|-------|-----------------------|---------|
| <input type="checkbox"/> | TWB480-1-01 | E1028 | Link Hardware 1" Ball | £146.00 |
|--------------------------|-------------|-------|-----------------------|---------|

#### B. Link Headrest Hardware 17mm Ball/Flip Down

|                          |              |       |                               |         |
|--------------------------|--------------|-------|-------------------------------|---------|
| <input type="checkbox"/> | TWB480-17-01 | E1028 | Link Hdw. 17mm Ball/Flip Down | £201.00 |
|--------------------------|--------------|-------|-------------------------------|---------|

#### Comfort Plus Headrest (Only if A is selected)

|                          |          |       |                  |        |
|--------------------------|----------|-------|------------------|--------|
| <input type="checkbox"/> | CP250-SR | E0955 | Comfort Plus 10" | £96.00 |
| <input type="checkbox"/> | CP450-SR | E0955 | Comfort Plus 8"  | £96.00 |

#### i2i Head Positioning System (Only if B is selected)

|                          |        |       |                                |         |
|--------------------------|--------|-------|--------------------------------|---------|
| <input type="checkbox"/> | i2i250 | K0108 | i2i Low Pro Nino XS, Pad Only  | £213.00 |
| <input type="checkbox"/> | i2i350 | K0108 | i2i Low Pro Nino SM, Pad Only  | £213.00 |
| <input type="checkbox"/> | i2i450 | K0108 | i2i Low Pro Adult SM, Pad Only | £213.00 |

#### Chin Prompts (Only if i2i Head Positioning is selected)

|                          |        |       |                     |        |
|--------------------------|--------|-------|---------------------|--------|
| <input type="checkbox"/> | i2i460 | K0108 | i2i SM Chin Prompt  | £49.00 |
| <input type="checkbox"/> | i2i560 | K0108 | i2i MED Chin Prompt | £49.00 |

#### Combo Series Headrest

|                          |       |       |                      |         |
|--------------------------|-------|-------|----------------------|---------|
| <input type="checkbox"/> | CB315 | E0955 | Combo Series, Small  | £185.00 |
| <input type="checkbox"/> | CB305 | E0955 | Combo Series, Medium | £185.00 |
| <input type="checkbox"/> | CB325 | E0955 | Combo Series, Large  | £185.00 |