



REACH SE™ Retail Order Form *REACH*

July 2021

Stealth Customer Service: 01872 324600

Email: sales@quantumrehab.co.uk

Date: _____ ☐ QUOTE ☐ ORDER - PO# _____

Company Name: _____

Ship To: _____

Contact Name: _____

Attention: _____

ATP Name: _____

Address: _____

Billing Address: _____

Post Code: _____

Post Code: _____

Phone: _____ Fax: _____

Email: _____

Mark for: _____



REACH SE™ INCLUDES:


- Fixed Seat Tilt 10°, 20°, 30°
- Hip Angle Adjustment
- QR 12" Poly Rear Wheels
- 7" Poly Front Casters
- Height Adjustable Push Handle
- Height Adjustable Foot Plate
- Foot Activated Hub Lock
- Under Seat Storage Bag
- Sunshade

Specifications & Adjustability

	REACH™ 12	REACH™ 14	REACH™ 16
Max Occupant Weight	80 lbs/36 kg	110 lbs/50 kg	150 lbs/68 kg
Seat Width	8-12 in/20-30 cm	10-14 in/25-35 cm	12-16 in/30-40 cm
Seat Depth	8-13 in/20-33 cm	11-16 in/28-40 cm	13-18 in/33-46 cm
Back Height	11-18 in/28-46 cm	13-21 in/33-53 cm	16-24 in/40-61 cm
Seat to Footplate	8-13 in/20-33 cm	11-16 in/28-40 cm	13-18 in/33-46 cm
Seat to Floor	19.5 in/50 cm	19.5 in/50 cm	19.5 in/50 cm
Fixed Tilt Positions	10, 20, or 30°		
Seat to Back Angle (Range depends on seat angle)			
10° seat angle	85, 90, 95, 100, 105, 110, 115°		
20° seat angle	85, 90, 95, 100, 105°		
30° seat angle	85, 90, 95°		
Weight of Chair	27 lbs/12 kg	29 lbs/13 kg	31 lbs/14 kg
Chair Width	21.5 in/55cm	23.5 in/60 cm	25.5 in/65 cm
Folded Dimensions	31 x 21.5 x 18 in 79 x 55 x 46 cm	31 x 23.5 x 18 in 79 x 60 x 46 cm	33 x 25.5 x 18 in 84 x 65 x 46 cm

MEASUREMENT WORKSHEET

- A. Seat Depth
- B. Trunk Depth
- C. Trunk Width
- D. Hip Width
- E. Seat to Top of Shoulder
- F. Lower Leg Length
- G. Seat to Axilla



REACH™ Chair

REACH SE™

☐	3004-1004	E1238	REACH SE™ 12 - 12 in./30.5 cm width	£1,635.00
☐	3004-1005	E1238	REACH SE™ 14 - 14 in./35.5 cm width	£1,635.00
☐	3004-1006	E1238	REACH SE™ 16 - 16 in./40.5 cm width	£1,635.00

Upholstery Color

Upholstery color applies to back and seat cushions, sunshade.

☐	3004-0101	BIG SKY BLUE - Blue/Black	N/C
☐	3004-0102	ROSSO RED - Red/Black	N/C
☐	3004-0103	ICE SILVER - Silver/Black	N/C
☐	3004-0104	SNOWBERRY PINK - Pink/Black	N/C
☐	3004-0105	HUCKLEBERRY PURPLE - Purple/Black	N/C

Transit Option & Transit Harness Selection

Transit Option

Requires belt/harness, back rest and a head rest.

☐	3004-7001	Transit Option (RESNA WC-19)	£187.00
--------------------------	-----------	------------------------------	---------

This option is not retrofittable.

Transit Belt or Harness

Transit 5- Point Harness NOT COMPATIBLE with Reach 16

☐	3004-7002	K0108	Transit 5-Point Harness for 12/14	£124.00
		☐	Add pad set for 5-point harness (3002-7022)	£68.00
☐	3001-7008	K0108	Transit Sure-Lok Lap Belt for 12/14	£68.00
☐	3004-7004	K0108	Transit Sure-Lok Lap Belt for 16	£68.00

Select 5-PT Harness for occupants less than 50 lbs./23kg

Select Sure-Lok Lap Belt for occupants over 50 lbs./23kg

The transit option includes 4 installed transportation anchors & has been tested in compliance with WC-19 standards. For more information on RESNA WC-19 please visit www.wc-transportation-safety.com.

Seat Pan and Front Rigging

A. Standard Front Rigging and Solid Seat Pan

Select both items to order the standard version

☐	3004-3001	E2231	Solid Seat Pan	n/c
☐	3004-7032	E0990	Angle Adjustable Legrest	n/c

B. Swing-Away Leg Hangers and Solid Seat Pan

Seat Pan 3004-3002 Required for use with Swing-Away Leg Hangers

☐	3004-3002	E2231	Solid Seat Pan for Swing-Aways	n/c
--------------------------	-----------	-------	--------------------------------	-----

Select rigging style, then select length from table to the right.

☐	3004-5007		70° Swing Away Leg Hangers	n/c
☐	3004-5003		80° Swing Away Leg Hangers	n/c
☐	3004-5009	E0990	Swing Away ELR	£191.00

Option B is Not Compatible with Reach 12

Front Rigging Length, Foot Plate, Calf Pad (only if B is selected)

Append the rigging length combination's two digit part number extension to the rigging part number.

70°	Short Extension	Standard Extension
Short Upper*	10-14"/25.4-35.56 cm -11 ☐	11-16"/27.94-40.64 -12 ☐
Medium Upper	11.5-16"/29.21-40 cm -21 ☐	11.5-18"/19.21-45.72 -22 ☐
Long Upper	n/a	13-19"/33-48.26 cm -32 ☐

80°	Standard Extension
Short Upper*	10-14"/25-35.56 cm -12 ☐
Medium Upper	11.5-14"/29.21-35.56 -22 ☐
Long Upper	11-17"/28-43 cm -32 ☐

Certain measurements may cause caster interference

Foot Plate and Front Rigging Accessories

☐	3004-5199		Omit Foot Plate	n/c
☐	3004-5010	K0040	Flip Up, Height/Angle Adjustable, Alum.	n/c
☐	3004-5011	E0951	Heel Loop for Foot Plate	£30.00
☐	3004-5013	E0995	Calf Pad, Left	£30.00
☐	3004-5014	E0995	Calf Pad, Right	£30.00

**Calf Pad Not Compatible with Short Upper Extension*

Ankle Supports

☐	3001-6007	E0951	LEFT, Extra Small (Pad 2x8"), REACH 12/14	£68.00
☐	3001-6008	E0951	RIGHT, Extra Small (Pad 2x8"), REACH 12/14	£68.00
☐	3001-6009	E0951	LEFT, Small (Pad 3x10"), REACH 14/16	£68.00
☐	3001-6010	E0951	RIGHT, Small (Pad 3x10"), REACH 14/16	£68.00

Small ankle supports are compatible with REACH 14 and 16 ONLY

XS ankle supports are NOT compatible with swing away leg hangers

Seat Cushion

LT Seat Cushion

☐	3002-7099		Omit Seat Cushion	n/c
☐	3005-3001	E2292	REACH™ 12 Planar Seat Cushion	n/c
☐	3005-3002	E2292	REACH™ 14 Planar Seat Cushion	n/c
☐	3005-3003	E2292	REACH™ 16 Planar Seat Cushion	n/c

EVA Positioning Wedges for Seat and/or Back

☐	3001-7011	K0108	Wedges QTY 4 and Lumbar QTY 1	£38.00
--------------------------	-----------	-------	-------------------------------	--------

Lateral Pelvic/Thigh Supports

☐	3005-6010	E0956	Reach Adj. Thigh Support, LEFT	£90.00
☐	3005-6011	E0956	Reach Adj. Thigh Support, RIGHT	£90.00

Armrests & Positioning Tray

☐	3004-6009	E0973	Reach Height/Angle Adj. Armrests	n/c
--------------------------	-----------	-------	----------------------------------	-----

☐	3004-6113	E0950	REACH™ 12 Clear Positioning Tray	£112.00
☐	3004-6114	E0950	REACH™ 14 Clear Positioning Tray	£112.00
☐	3004-6115	E0950	REACH™ 16 Clear Positioning Tray	£112.00
☐	3004-6008	E1028	Tray Hardware, Adj./Removable	£157.00

3004-6009 Armrests & 3004-6008 Hardware are required when ordering Tray Platform

Anterior Posture Supports & Pelvic Positioning

SP™ Pelvic Positioning Belts

☐	IFPB5141	E0978	4-Point, Push Button 1", XS	£90.00
☐	IFPB5242	E0978	4-Point, Push Button 1.5", SM	£90.00
☐	IFPB5243	E0978	4-Point, Push Button 1.5", M	£90.00
☐	IFPB5111	E0978	2-Point, Push Button 1", XS	£76.00
☐	IFPB5112	E0978	2-Point, Push Button 1", SM	£76.00

Abductor/Medial Knee Block

Must select hardware & one size of abductor pad. Wedge sizes are width x depth x height.

☐	3005-7001	E1028	Leggero Flip Down Medial Hdw.	£180.00
☐	MTMD-12	E0957	Small Wedge Pad (3x5x3")	£96.00
☐	MTMD-13	E0957	Medium Wedge Pad (5x5x3")	£96.00

Harness/Belt Accessories

☐	3001-6011	K0108	Strap Risers, Pair for 1.25" Webbing	£33.00
☐	IFMDBC10		SP™ Buckle Cover Kit 1" Buckle	£21.00
☐	IFMDBC15		SP™ Buckle Cover Kit 1.5" Buckle	£21.00

Anterior Posture Supports

☐	IFPS1101-02	E0960	SP™ Classic, Dynamic, XS	£100.00
☐	IFPS1102-02	E0960	SP™ Classic, Dynamic, S	£100.00
☐	IFPS1103-02	E0960	SP™ Classic Dynamic, M	£100.00
☐	IFPS1201-02	E0960	SP™ Classic, Structured, XS	£100.00
☐	IFPS1202-02	E0960	SP™ Classic, Structured, S	£100.00
☐	IFPS1203-02	E0960	SP™ Classic, Structured, M	£100.00
☐	IFPS2101-02	E0960	SP™ Contoured, Dynamic, XS	£100.00
☐	IFPS2102-02	E0960	SP™ Contoured, Dynamic, S	£100.00
☐	IFPS2103-02	E0960	SP™ Contoured, Dynamic, M	£100.00
☐	IFPS2201-02	E0960	SP™ Contoured, Structured, XS	£100.00
☐	IFPS2202-02	E0960	SP™ Contoured, Structured, S	£100.00
☐	IFPS2203-02	E0960	SP™ Contoured, Structured, M	£100.00
☐	3005-6016	E0960	Leggero 'H' Harness, S	£124.00
☐	3005-6017	E0960	Leggero 'H' Harness, M	£124.00
☐	3005-6018	E0960	Leggero 'H' Harness, L	£124.00
☐	3005-6014-12	E0960	5-Point Positioning Harness	£124.00
☐	3005-6014-14	E0960	5-Point Positioning Harness	£124.00
☐	3005-6014-16	E0960	5-Point Positioning Harness	£124.00

☐ Add pad set for 5-point harness (3002-7022) £68.00

Back Rest & Accessories

Select a system from options **A, B, or C**. Then continue your selections to the right to complete the package.

A. Reach Planar Back

☐	3005-4004	E2291	REACH™ 12 Planar Back	n/c
☐	3005-4005	E2291	REACH™ 14 Planar Back	n/c
☐	3005-4006	E2291	REACH™ 16 Planar Back	n/c
☐	3003-4099		Omit Back	n/c

When omitting back, crossbars between canes will remain in place.

B. LT Contoured Back

☐	3005-4001	E2293	REACH™ 12 Contoured Adj. Back	n/c
☐	3005-4002	E2293	REACH™ 14 Contoured Adj. Back	n/c
☐	3005-4003	E2293	REACH™ 16 Contoured Adj. Back	n/c

Soft Lateral w/ LT 'H' Harness Options (Only available with B)

☐	3005-6005	E0960	Soft Laterals w/ SM 'H' Harness	n/c
☐	3005-6006	E0960	Soft Laterals w/ MED 'H' Harness	n/c
☐	3005-6007	E0960	Soft Laterals w/ LRG 'H' Harness	n/c

Soft Laterals' are pads that integrate with back and harness

LT Lumbar Support (Only available with B)

☐	3005-6002	K0108	REACH™ 12 Adj. Lumbar Support	n/c
☐	3005-6003	K0108	REACH™ 14 Adj. Lumbar Support	n/c
☐	3005-6004	K0108	REACH™ 16 Adj. Lumbar Support	n/c

C. REACH RCB™ Rigid Contoured Back

REACH RCB™ 12 Contoured Back Options (includes hardware, plate, and cushion)

☐	3004-4212	E2293	12" W x 12" H (Adj. 12-14")	£392.00
☐	3004-4214	E2293	12" W x 14" H (Adj. 14-16")	£392.00
☐	3004-4216	E2293	12" W x 16" H (Adj. 16-18")	£392.00
☐	3004-4218	E2293	12" W x 18" H (Fixed Height)	£392.00

REACH RCB™ 14 Contoured Back Options (includes hardware, plate, and cushion)

☐	3004-4414	E2293	14" W x 14" H (Adj. 14-16")	£392.00
☐	3004-4416	E2293	14" W x 16" H (Adj. 16-18")	£392.00
☐	3004-4418	E2293	14" W x 18" H (Adj. 18-20")	£392.00
☐	3004-4420	E2293	14" W x 20" H (Adj. 20-21")	£392.00

REACH RCB™ 16 Contoured Back Options (includes hardware, plate, and cushion)

☐	3004-4616	E2293	16" W x 16" H (Adj. 16-18")	£392.00
☐	3004-4618	E2293	16" W x 18" H (Adj. 18-20")	£392.00
☐	3004-4620	E2293	16" W x 20" H (Adj. 20-22")	£392.00
☐	3004-4622	E2293	16" W x 22" H (Adj. 22-24")	£392.00
☐	3004-4624	E2293	16" W x 24" H (Fixed Height)	£392.00

Swing Away Lateral Hardware & Pads (Only available with C)

Must select hardware & pad. Leggero Swing In/Out Lateral includes rail and track slider.

☐	3004-6111		E1028	LEFT, Leggero Swing In/Out Lateral	£168.00
☐	3004-6112		E1028	RIGHT, Leggero Swing In/Out Lateral	£168.00
L	R	3004-6035	E0956	Lateral Pad 3x5", Each	£34.00
L	R	3004-6045	E0956	Lateral Pad 4x5", Each	£34.00
L	R	3004-6046	E0956	Lateral Pad 4x6", Each	£34.00

Leggero and Stealth Products Head Positioning

Leggero Head Positioning

☐	3005-6008	E0966	LT Planar Rectangle Headrest	£122.00
☐	3001-6005	E0955	8" Oval Headrest	£146.00
☐	3002-6003	E0955	RESPOND™ Cradle Headrest and Hardware	£146.00
☐	3005-6099		Omit Head Positioning	n/c

Hardware is included with LT, 8" Multi-Adjustable & RESPOND™ Headrests

A. Link Headrest Hardware 1" Ball

☐	TWB480-1-01	E1028	Link Hardware 1" Ball	£146.00
--------------------------	-------------	-------	-----------------------	---------

B. Link Headrest Hardware 17mm Ball/Flip Down

☐	TWB480-17-01	E1028	Link Hdw. 17mm Ball/Flip Down	£201.00
--------------------------	--------------	-------	-------------------------------	---------

Comfort Plus Headrest (Only if A is selected)

☐	CP250-SR	E0955	Comfort Plus 10"	£96.00
☐	CP450-SR	E0955	Comfort Plus 8"	£96.00

i2i Head Positioning System (Only if B is selected)

☐	i2i250	K0108	i2i Low Pro Nino XS, Pad Only	£213.00
☐	i2i350	K0108	i2i Low Pro Nino SM, Pad Only	£213.00
☐	i2i450	K0108	i2i Low Pro Adult SM, Pad Only	£213.00

Chin Prompts (Only if i2i Head Positioning is selected)

☐	i2i460	K0108	i2i SM Chin Prompt	£49.00
☐	i2i560	K0108	i2i MED Chin Prompt	£49.00

Combo Series Headrest

☐	CB315	E0955	Combo Series, Small	£185.00
☐	CB305	E0955	Combo Series, Medium	£185.00
☐	CB325	E0955	Combo Series, Large	£185.00

Accessories

☐	3004-7020		Rain Cover, REACH™ 12	£146.00
☐	3004-7021		Rain Cover, REACH™ 14	£146.00
☐	3004-7022		Rain Cover, REACH™ 16	£146.00

☐	3004-7011	K0108	Calf Support, REACH™ 12	£56.00
☐	3004-7027	K0108	Calf Support, REACH™ 14	£56.00
☐	3004-7028	K0108	Calf Support, REACH™ 16	£56.00

May interfere with folding

☐	3004-5002		Custom set lower leg length to (F) dimension	n/c
--------------------------	-----------	--	--	-----

When 3004-5002 is selected, most leg lengths can be accommodated.

☐	3004-7012		Grab Bar, REACH™ 12	£101.00
☐	3004-7014		Grab Bar, REACH™ 14	£101.00
☐	3004-7016		Grab Bar, REACH™ 16	£101.00

☒	3004-7009	E0971	Rear Anti-Tips (required)	n/c
☐	3004-7025	K0105	IV Pole	£112.00
☐	3004-7026	E2208	Oxygen Tank Holder (D size)	£131.00