

RETAIL UK

TRAK SE™

ORGANISATION DETAILS

Organisation Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____


DELIVERY DETAILS

Organisation Address _____

Postcode _____

Stealth Specialist _____

① Clients Weight _____



Adjustability	TRAK™ 12	TRAK™ 14
Seat Depth	6.5-15 in / 17-38 cm	6.5-15 in / 17-38 cm
Seat Width	6-12 in / 15-30 cm	8-14 in / 20-36 cm
Lower Leg Width	4-17 in / 10-43 cm	4-17 in / 10-43 cm
Back Height	12-19 in / 30-48 cm	13-19 in / 33-48 cm
Hip Angles	90°, 105°, 120°	90°, 105°, 120°

BASE MODEL

- Integrated Suspension
- 45° Tilt -in -Space Frame
- 30 °Hip Angle Adjustment
- QR 16" Poly Rear Wheels
- 7" Poly Front Casters
- Height Adjustable Push Handle
- Height Adjustable Foot Plate
- Foot Activated Wheel Lock
- Caster Swivel Locks
- Under Seat Storage Bag
- Sunshade



Specifications

Max Chair Capacity	TRAK 12 & 14 - 80lbs / 36kg	
Push Handle Height	40-50 in / 101-127 cm	
Weight	29 lbs / 13 kg	
Rear Wheel Diameter	16 in / 41 cm	
Front Caster Diameter	7 in / 25 cm	
Reversing Width	41 in / 104 cm	
Maximum Safe Slope	5°	

	TRAK™ 12	TRAK™ 14
Overall Width	22.5in / 57cm	24.5in / 62cm
Folded no rear wheels	89x43x57cm	89x43x62cm
Folded with wheels	92x51x57cm	92x51x62cm

TRAK SE™ CHAIR

TRAK SE™ 12.
3002-1004. £2,179

TRAK SE™ 14.
3002-1002. £2,179

UPHOLSTERY COLOR

BIG SKY BLUE - Blue/Black. 3002-0101.	FOC
ROSSO RED - Red/Black. 3002-0102.	FOC
ICE SILVER - Silver/Black. 3002-0103.	FOC
SNOWBERRY PINK - Pink/Black. 3002-0104.	FOC
HUCKLEBERRY PURPLE - Purple/Black. 3002-0105.	FOC

Upholstery Color applies to RESPOND™ back and seat cushions, sunshade.

TRANSIT OPTION

Transit Option (RESNA WC-19). *This option is not retrofittable. 3002-7001.*

£187

The transit option includes 4 installed anchors & has been tested in compliance with WC-19 standards. Tested with Stealth WC-20 Transit Knob Hardware. For more information on please visit www.wc-transportation-safety.com.

SEAT TILT RANGE

Transit - 5-Point Harness. *3002-7002.*

£124

Add pad set for 5-point harness. *3002-7022.*

£68

Transit - Sure-Lok Lap Belt. *3002-7008.*

£68

Transit Belt is Required. Select 5-PT Harness for occupants less than 50 lbs./23kg. Select Sure-Lok Lap Belt for occupants over 50 lbs./23kg

FRONT RIGGING

90° Fixed Leg Rest. *3002-7015.*

FOC

Angle Adjustable Leg Rest. *3002-7006.*

FOC

Foot
Plate.

Standard Height Adjustable. *3002-5001.*

FOC

Height & Angle Adjustable. *3002-5002.*

£124

RESPOND™ SEAT CUSHIONS

Solid Seat Pan. *3002-3006.*

FOC

Omit Seat Cushion. *3002-7099.*

£168

RESPOND™
Planar
Seat
Cushion

TRAK 12 Cascading Planar Cushion. *3002-5001.*

£221

TRAK 12 Extra Planar Cushion Cover. *3002-5002.*

£53

TRAK 14 Cascading Planar Cushion. *3002-5001.*

£221

TRAK 14 Extra Planar Cushion Cover. *3002-5002.*

£53

RESPOND™
Contoured
Seat
Cushion

TRAK 12 Cascading Contoured
Cushion. *3002-3005.*

FOC

Contoured Seat Ischial Insert Options

TRAK 12 Extra Cushion Cover.
3002-7025.

£53

TRAK 12 PU Foam Ischial Insert.
3003-3010.

FOC

TRAK 14 Cascading Contoured
Cushion. *3002-3002.*

FOC

TRAK 14 PU Foam Ischial Insert.
3003-3006.

FOC

TRAK 14 Extra Cushion Cover.
3002-7011.

£30

TRAK 14 Soft Mem.
Foam Ischial Insert. *3003-3007.*

£38

Lumbar EVA & 4 EVA Wedges.
3001-7011.

£53

TRAK 14 Med. Mem.
Foam Ischial Insert. *3003-3008.*

£38

RESPOND™
Back

TRAK 12 Contoured Config. Back. *3003-4002.*

FOC

TRAK 12 Extra Back Cushion Cover. *3003-3013.*

£53

TRAK 14 Contoured Config. Back. *3003-4001.*

FOC

TRAK 14 Extra Back Cushion Cover. *3003-3009.*

£53

Lumbar EVA & 4 EVA Wedges. *3001-7011.*

£38

Omit Back. *3003-4099.*

FOC

When omitting back, crossbars between canes will remain in place.

LATERAL/PELVIC THIGH SUPPORTS

TRAK Adj. Thigh Support, LEFT. *3003-6006.*

FOC

TRAK Adj. Thigh Support, RIGHT. *3003-6007.*

FOC

ARMRESTS & POSITIONING TRAY

Armrest.	Height/Angle Adj. Armrests. <i>3002-6017.</i>	FOC
Positioning Tray.	TRAK 12 Clear Positioning Tray Platform. <i>3002-6021.</i>	£112
	TRAK 14 Clear Positioning Tray Platform. <i>3002-6015.</i>	£112
	Tray Hdw., Height & Angle Adjustable. <i>3002-6016.</i>	£157

3002-6017 Armrests & 3002-6016 Hardware are required when ordering Tray Platform

ANKLE SUPPORTS

Ankle Supports.	LEFT	£68	<i>3001-6007 < Extra Small (Pad 2x8"). > 3001-6008</i>	£68	RIGHT
		£68	<i>3001-6009 < Small (Pad 3x10"). > 3001-6010</i>	£68	

TRUNK & PELVIC POSITIONING

SP™ Pelvic Positioning Belts.	4-Point, Push Button 1", XS. <i>IFPB5141.</i>	£125			
	4-Point, Push Button 1.5", SM. <i>IFPB5242.</i>	£125			
	4-Point, Push Button 1.5", M. <i>IFPB5243.</i>	£125			
	2-Point, Push Button 1", XS. <i>IFPB5111.</i>	£83			
	2-Point, Push Button 1", SM. <i>IFPB5112.</i>	£83			
Abductor/ Medial Knee Block.	Leggero Flip Down Medial Hdw. <i>3005-7001.</i>	£180			
	Small Wedge Pad (3x5x3"). <i>MTMD-12.</i>	£103			
	Medium Wedge Pad (5x5x3"). <i>MTMD-13.</i>	£103			
Swing Away Lateral Hardware & Pads.	£168	<i>3002-6111 < Leggero Swing In/Out Lateral. > 3002-6112</i>	£168		
	£179	<i>TWBTL-L-01 < Stealth Thin Lateral. > TWBTL-R-01</i>	£179		
	£40	<i>3004-6035 < Lateral Pad 3x5", Each. > 3004-6035</i>	£40		
	£40	<i>3004-6045 < SLateral Pad 4x5", Each > 3004-6045</i>	£40		
	£40	<i>3004-6046 < Lateral Pad 3x6", Each. > 3004-6046</i>	£40		
Fixed Trunk Supports.	LEFT	FOC	<i>3003-6001 < 3x4" (7.5x10cm). > 3003-6002</i> <i>Requires RESPOND™ Back 3003-4001</i>	FOC	RIGHT
Anterior Posture Supports.	SP™ Classic, Dynamic.	XS. <i>IFPS1101-02.</i> £137	S. <i>IFPS1102-02.</i> £137	M. <i>IFPS1103-02.</i> £137	
	SP™ Classic, Structured.	XS. <i>IFPS1201-02.</i> £137	S. <i>IFPS1202-02.</i> £137	M. <i>IFPS1203-02.</i> £137	
	SP™ Contoured, Dynamic.	XS. <i>IFPS2101-02.</i> £137	S. <i>IFPS2102-02.</i> £137	M. <i>IFPS2103-02.</i> £137	
	SP™ Contoured, Structured.	XS. <i>IFPS2201-02.</i> £137	S. <i>IFPS2202-02.</i> £137	M. <i>IFPS2203-02.</i> £137	
	Leggero 'H' Harness.	S. <i>3005-6016.</i> £124	M. <i>3005-6017.</i> £124	L. <i>3005-6018.</i> £124	
	5-Point Positioning Harness.	S <i>3005-6014-12.</i> £124	M <i>3005-6014-14.</i> £124	L <i>3005-6014-16.</i> £124	
	Pad set for 5-point harness. <i>(3002-7022).</i>				£68
Harness/Belt Accessories.	Strap Risers, Pair for 1.25" Webbing. <i>3001-6011.</i>				£35
	SP™ Buckle Cover Kit 1" Buckle. <i>IFMDBC10.</i>				£24
	SP™ Buckle Cover Kit 1.5" Buckle. <i>IFMDBC15.</i>				£24

LEGGERO HEAD POSITIONING

8" Oval Headrest and Hardware. <i>3001-6005.</i>	£146
RESPOND™ Cradle Headrest and Hardware. <i>3001-6006.</i>	£146
Omit Head Positioning. <i>3005-6099.</i>	FOC

STEALTH PRODUCTS HEAD POSITIONING

Link Headrest.	Link Hardware 1" Ball. <i>TWB480-1-01.</i>	£167	▶	Comfort Plus Headrest 10". <i>CP250-SR.</i>	£108	
				Comfort Plus Headrest 8". <i>CP450-SR.</i>	£108	
Link Headrest Hardware.	17mm Ball/Flip Down. <i>TWB480-1-01.</i>	£201	▶	i2i Head Pos System Low Pro Nino XS, Pad Only. <i>i2i250.</i>	£228	
				i2i Head Pos System Low Pro Nino SM, Pad Only. <i>i2i350.</i>	£228	
				i2i Head Pos System Low Pro Adult SM, Pad Only. <i>i2i450.</i>	£228	
				i2i SM Chin Prompt. <i>i2i460.</i>	£53	
				i2i MED Chin Prompt. <i>i2i560.</i>	£53	
Combo Series Headrest	S. <i>CB315.</i>	£193	M. <i>CB305.</i>	£193	L. <i>CB325.</i>	£193

REAR ANTI-TIPS

Rear Anti-Tips. <i>3002-7007.</i>	FOC
-----------------------------------	-----

MEDICAL MOUNTING & ACCESSORIES

IV Pole. <i>3002-7016. Must be factory installed.</i>	£112
Oxygen Tank Holder (D size). <i>3002-7017. Must be factory installed.</i>	£131
LTV Bracket for LTV950 ¹ . <i>3002-7019. Must be factory installed.</i>	£187
TRAK 12 Equipment Tray. <i>3002-7027.</i>	£262
TRAK 12 LTV Bracket ¹ and Equip. Tray. <i>3002-7028.</i>	£370
TRAK 14 Equipment Tray. <i>3002-7018.</i>	£262
TRAK 14 LTV Bracket ¹ and Equip. Tray. <i>3002-7021.</i>	£370
Under Seat Storage. <i>3002-7004.</i>	FOC
Deluxe Sunshade. <i>3002-7012.</i>	FOC
TRAK 12 Grab Bar. <i>3002-6113.</i>	£101
TRAK 12 Rain Cover. <i>3002-7026.</i>	£146
TRAK 14 Grab Bar. <i>3002-6019.</i>	£101
TRAK 14 Rain Cover. <i>3002-7014.</i>	£146

Combined weight of accessories & occupant must not exceed chair capacity. Ventilator, oxygen tank, and any other accessories must be removed and secured off chair during transit in a motor vehicle.

PRICE

TRAK SE™	BASE PRICE	£2,179
	COST OPTIONS	
	TOTAL	



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
RETAIL UK - 30/07/24 V1



V1: 14/06/2022



TRAK SE™ GALLERY



5-POINT HARNESS ADJUSTMENT



CASCADING SEAT CUSHION



PLANAR SEAT CUSHION



LT HARNESS WITH PADS



RESPOND SEATING



8 INCH OVAL HEAD SUPPORT



TILT LEVER




2 POINT PADDED BELTS (DUAL PULL)
4 POINT BELTS (V-CONTOUR STYLE)

BELT



CLASSIC STYLE CONTOURED STYLE


HARNESS



ANTI-TIPS



COMBO SERIES HEAD SUPPORT



COMFORT PLUS HEAD SUPPORT



I2I HEAD SUPPORT



OXYGEN BOTTLE HOLDER



SUNSHADE