



2 SIZES
AVAILABLE



Facebook.com/
QuantumRehabUK



14



12

TRAK™

- ✔ Free Sunshade & Storage Bag
- ✔ Lightweight (29 lbs/13 kg)
- ✔ 30° Hip Angle Adjustment
- ✔ 45° Tilting Frame
- ✔ Configurable Seating
- ✔ Folds Compactly
- ✔ Rear Suspension



Call: 01869 324600 • Web: www.quantumrehab.co.uk • Email: sales@quantumrehab.co.uk

© 2019 LEGGERO LLC • Prices and specifications are subject to change without notice

SPECIFICATIONS & ADJUSTMENT RANGE

	TRAK™ 12	TRAK™ 14
Weight Limit	80 lbs / 36 kg	80 lbs / 36 kg
Seat Depth	6.5-15 in / 17-38 cm	6.5-15 in / 17-38 cm
Seat Width	6-12 in / 15-30 cm	8-14 in / 20-36 cm
Back Height	12-19 in / 30-48 cm	13-19 in / 33-48 cm
Seat to Footplate	4-17 in / 10-43 cm	4-17 in / 10-43 cm
Tilt Range	45°	45°
Hip Angle Positions	90°, 105°, 120°	90°, 105°, 120°
Maximum Safe Slope	5°	5°
Overall Width	22.5 in / 57 cm	24.5 in / 62 cm
Push Handle Height	40-50 in / 101-127 cm	40-50 in / 101-127 cm
Overall Length	35 in / 89 cm	35 in / 89 cm
Caster Diameter	7 in / 18 cm	7 in / 18 cm
Rear Wheel Diameter	16 in / 41 cm	16 in / 41 cm
Weight of Chair	29 lbs / 13 kg	29 lbs / 13 kg
Folded without Rear Wheels	35 x 17 x 22.5 in 89 x 43 x 57 cm	35 x 17 x 24.5 in 89 x 43 x 62 cm
Folded with Wheels	36 x 21 x 22.5 in 92 x 51 x 57 cm	36 x 21 x 24.5 in 92 x 51 x 62 cm

30°
Hip
Angle
Range

45° Tilt Range



Easy folding with
seating components
in place.

FEATURES & OPTIONS

- ✔ Easy to operate foot tilt lever
- ✔ Accessories include Oxygen Tank Holder, IV Pole, LTV Vent Bracket & Ventilator/Battery Tray
- ✔ Numerous Positioning Components
- ✔ Height Adjustable Push Handle
- ✔ Foot Activated Wheel Locks
- ✔ Polyurethane Wheels & Casters
- ✔ Height Adjustable Footplate
- ✔ Rear Anti-tips

RESPOND™ SEATING

- ✔ Contoured Back Cushion
- ✔ Contoured Seat Cushion
- ✔ Adjustable Lateral Thoracic Supports
- ✔ Adjustable Hip Guides

