



TRAK SE™ UK Retail Order Form



June 2021

Stealth UK Customer Service: 01869 324600

Email: sales@quantumrehab.co.uk

Date: _____ QUOTE ORDER - PO# _____

Service Name: _____

Ship To: _____

Contact Name: _____

Attention: _____

Billing Address: _____

Address: _____

Postcode: _____

Postcode: _____

Phone: _____ Fax: _____

Email: _____

Mark for: _____

TRAK SE™ INCLUDES:



Model shown with some accessories


- Integrated Suspension
- 45° Tilt-in-Space Frame
- 30° Hip Angle Adjustment
- QR 16" Poly Rear Wheels
- 7" Poly Front Casters
- Height Adjustable Push Handle
- Height Adjustable Foot Plate
- Foot Activated Wheel Lock
- Caster Swivel Locks
- Under Seat Storage Bag
- Sunshade

Specifications

Maximum Chair Capacity	(TRAK 12 & TRAK 14) 80 lbs / 36 kg	
Push Handle Height	40-50 in / 101-127 cm	
Weight	29 lbs / 13 kg	
Rear Wheel Diameter	16 in / 41 cm	
Front Caster Diameter	7 in / 25 cm	
Reversing Width	41 in / 104 cm	
Maximum Safe Slope	5°	
	TRAK™ 12	TRAK™ 14
Overall Width	22.5 in / 57 cm	24.5 in / 62 cm
Folded without rear wheels	35x17x22.5 in 89x43x57 cm	35x17x24.5 in 89x43x62 cm
Folded with wheels	36x21x22.5 in 92x51x57 cm	36x21x24.5 in 92x51x62 cm

MEASUREMENT WORKSHEET

- A. Seat Depth _____
- B. Trunk Depth _____
- C. Trunk Width _____
- D. Hip Width _____
- E. Seat to Top of Shoulder _____
- F. Lower Leg Length _____
- G. Seat to Axilla _____



Adjustability

	TRAK™ 12	TRAK™ 14
Seat Depth	6.5-15 in / 17-38 cm	6.5-15 in / 17-38 cm
Seat Width	6-12 in / 15-30 cm	8-14 in / 20-36 cm
Lower Leg Length	4-17 in / 10-43 cm	4-17 in / 10-43 cm
Back Height	12-19 in / 30-48 cm	13-19 in / 33-48 cm
Hip Angles	90°, 105°, 120°	90°, 105°, 120°

TRAK™ Chair

TRAK SE™

<input type="checkbox"/>	3002-1004	E1234	TRAK SE™ 12	£2,179.00
<input type="checkbox"/>	3002-1002	E1234	TRAK SE™ 14	£2,179.00

Upholstery Color

Upholstery Color applies to RESPOND™ back and seat cushions, sunshade.

<input type="checkbox"/>	3002-0101	BIG SKY BLUE - Blue/Black	N/C
<input type="checkbox"/>	3002-0102	ROSSO RED - Red/Black	N/C
<input type="checkbox"/>	3002-0103	ICE SILVER - Silver/Black	N/C
<input type="checkbox"/>	3002-0104	SNOWBERRY PINK - Pink/Black	N/C
<input type="checkbox"/>	3002-0105	HUCKLEBERRY PURPLE - Purple/Black	N/C

Transit Option & Transit Harness Selection

Transit Option

Select belt or harness to right, back rest and a head rest.

<input type="checkbox"/>	3002-7001	Transit Option (RESNA WC-19)	£187.00
--------------------------	-----------	------------------------------	---------

This option is not retrofittable.

Transit Belt or Harness

Transit Belt is Required

<input type="checkbox"/>	3002-7002	K0108	Transit - 5-Point Harness	£124.00
		<input type="checkbox"/>	Add pad set for 5-point harness (3002-7022)	£68.00
<input type="checkbox"/>	3001-7008	K0108	Transit - Sure-Lok Lap Belt	£68.00

Select 5-PT Harness for occupants less than 50 lbs./23kg

Select Sure-Lok Lap Belt for occupants over 50 lbs./23kg

The transit option includes 4 installed transportation anchors & has been tested in compliance with WC-19 standards. For more information on RESNA WC-19 please visit www.wc-transportation-safety.com.

Front Rigging & Foot Plate

Front Rigging

<input type="checkbox"/>	3002-7015		90° Fixed Leg Rest	N/C
<input type="checkbox"/>	3002-7006	E0990	Angle Adjustable Leg Rest	£168.00

Foot Plate

<input type="checkbox"/>	3002-5001		Standard Height Adjustable	N/C
<input type="checkbox"/>	3002-5002	K0040	Height & Angle Adjustable	£124.00

RESPOND™ Seat Cushions

<input type="checkbox"/>	3002-3006	E2231	Solid Seat Pan	N/C
<input type="checkbox"/>	3002-7099		Omit Seat Cushion	N/C

Select from A or B. If B is chosen, continue selection to right.

A. RESPOND™ Planar Seat Cushion

<input type="checkbox"/>	3002-3004	E2292	TRAK 12 Cascading Planar Cushion	£221.00
<input type="checkbox"/>	3002-7024	E0981	TRAK 12 Extra Planar Cushion Cover	£53.00
<input type="checkbox"/>	3002-3001	E2292	TRAK 14 Cascading Planar Cushion	£221.00
<input type="checkbox"/>	3002-7010	E0981	TRAK 14 Extra Planar Cushion Cover	£53.00

B. RESPOND™ Contoured Seat Cushions

<input type="checkbox"/>	3002-3005	E2292	TRAK 12 Cascading Contoured Cushion	N/C
<input type="checkbox"/>	3002-7025	E0981	TRAK 12 Extra Cushion Cover	£53.00
<input type="checkbox"/>	3002-3002	E2292	TRAK 14 Cascading Contoured Cushion	N/C
<input type="checkbox"/>	3002-7011	E0981	TRAK 14 Extra Cushion Cover	£30.00
<input type="checkbox"/>	3001-7011	K0108	Lumbar EVA & 4 EVA Wedges	£53.00

Ischial Insert Options (Only if B is selected)

<input type="checkbox"/>	3003-3010		TRAK 12 PU Foam Ischial Insert	N/C
<input type="checkbox"/>	3003-3006		TRAK 14 PU Foam Ischial Insert	N/C
<input type="checkbox"/>	3003-3007	K0108	TRAK 14 Soft Mem. Foam Ischial Insert	£38.00
<input type="checkbox"/>	3003-3008	K0108	TRAK 14 Med. Mem. Foam Ischial Insert	£38.00

RESPOND™ Back

<input type="checkbox"/>	3003-4002	E2293	TRAK 12 Contoured Config. Back	N/C
<input type="checkbox"/>	3003-3013	E0982	TRAK 12 Extra Back Cushion Cover	£53.00
<input type="checkbox"/>	3003-4001	E2293	TRAK 14 Contoured Config. Back	N/C
<input type="checkbox"/>	3003-3009	E0982	TRAK 14 Extra Back Cushion Cover	£53.00
<input type="checkbox"/>	3001-7011	K0108	Lumbar EVA & 4 EVA Wedges	£38.00

<input type="checkbox"/>	3003-4099		Omit Back	N/C
--------------------------	-----------	--	-----------	-----

When omitting back, crossbars between canes will remain in place.

Lateral/Pelvic Thigh Supports

<input type="checkbox"/>	3003-6006	E0956	TRAK Adj. Thigh Support, LEFT	N/C
<input type="checkbox"/>	3003-6007	E0956	TRAK Adj. Thigh Support, RIGHT	N/C

Armrests & Positioning Tray

3002-6017 Armrests & 3002-6016 Hardware are required when ordering Tray Platform

<input type="checkbox"/>	3002-6017	E0973	Height/Angle Adj. Armrests	N/C
--------------------------	-----------	-------	----------------------------	-----

<input type="checkbox"/>	3002-6021	E0950	TRAK 12 Clear Positioning Tray Platform	£112.00
<input type="checkbox"/>	3002-6015	E0950	TRAK 14 Clear Positioning Tray Platform	£112.00
<input type="checkbox"/>	3002-6016	E1028	Tray Hdw., Height & Angle Adjustable	£157.00

Ankle Supports

<input type="checkbox"/>	3001-6007	E0951	LEFT, Extra Small (Pad 2x8")	£68.00
<input type="checkbox"/>	3001-6008	E0951	RIGHT, Extra Small (Pad 2x8")	£68.00
<input type="checkbox"/>	3001-6009	E0951	LEFT, Small (Pad 3x10")	£68.00
<input type="checkbox"/>	3001-6010	E0951	RIGHT, Small (Pad 3x10")	£68.00

Trunk & Pelvic Positioning**SP™ Pelvic Positioning Belts**

<input type="checkbox"/>	IFPB5141	E0978	4-Point, Push Button 1", XS	£90.00
<input type="checkbox"/>	IFPB5242	E0978	4-Point, Push Button 1.5", SM	£90.00
<input type="checkbox"/>	IFPB5243	E0978	4-Point, Push Button 1.5", M	£90.00
<input type="checkbox"/>	IFPB5111	E0978	2-Point, Push Button 1", XS	£76.00
<input type="checkbox"/>	IFPB5112	E0978	2-Point, Push Button 1", SM	£76.00

Side release available on request.

Anterior Posture Supports

<input type="checkbox"/>	IFPS1101-02	E0960	SP™ Classic, Dynamic, XS	£100.00
<input type="checkbox"/>	IFPS1102-02	E0960	SP™ Classic, Dynamic, S	£100.00
<input type="checkbox"/>	IFPS1103-02	E0960	SP™ Classic Dynamic, M	£100.00
<input type="checkbox"/>	IFPS1201-02	E0960	SP™ Classic, Structured, XS	£100.00
<input type="checkbox"/>	IFPS1202-02	E0960	SP™ Classic, Structured, S	£100.00
<input type="checkbox"/>	IFPS1203-02	E0960	SP™ Classic, Structured, M	£100.00
<input type="checkbox"/>	IFPS2101-02	E0960	SP™ Contoured, Dynamic, XS	£100.00
<input type="checkbox"/>	IFPS2102-02	E0960	SP™ Contoured, Dynamic, S	£100.00
<input type="checkbox"/>	IFPS2103-02	E0960	SP™ Contoured, Dynamic, M	£100.00
<input type="checkbox"/>	IFPS2201-02	E0960	SP™ Contoured, Structured, XS	£100.00
<input type="checkbox"/>	IFPS2202-02	E0960	SP™ Contoured, Structured, S	£100.00
<input type="checkbox"/>	IFPS2203-02	E0960	SP™ Contoured, Structured, M	£100.00
<input type="checkbox"/>	3005-6016	E0960	Leggero 'H' Harness, S	£124.00
<input type="checkbox"/>	3005-6017	E0960	Leggero 'H' Harness, M	£124.00
<input type="checkbox"/>	3005-6018	E0960	Leggero 'H' Harness, L	£124.00
<input type="checkbox"/>	3005-6014-12	E0960	5-Point Positioning Harness	£124.00
<input type="checkbox"/>	3005-6014-14	E0960	5-Point Positioning Harness	£124.00
<input type="checkbox"/>	3005-6014-16	E0960	5-Point Positioning Harness	£124.00

Add pad set for 5-point harness (3002-7022) £68.00

Abductor/Medial Knee Block

<input type="checkbox"/>	3005-7001	E1028	Leggero Flip Down Medial Hdw.	£180.00
<input type="checkbox"/>	MTMD-12	E0957	Small Wedge Pad (3x5x3")	£96.00
<input type="checkbox"/>	MTMD-13	E0957	Medium Wedge Pad (5x5x3")	£96.00

Swing Away Lateral Hardware & Pads

<input type="checkbox"/>	3002-6111	E1028	LEFT, Leggero Swing In/Out Lateral	£168.00	
<input type="checkbox"/>	3002-6112	E1028	RIGHT, Leggero Swing In/Out Lateral	£168.00	
<input type="checkbox"/>	TWBTL-L-01	E1028	LEFT, Stealth Thin Lateral	£152.00	
<input type="checkbox"/>	TWBTL-R-01	E1028	RIGHT, Stealth Thin Lateral	£152.00	
L	R	3004-6035	E0956	Lateral Pad 3x5", Each	£34.00
L	R	3004-6045	E0956	Lateral Pad 4x5", Each	£34.00
L	R	3004-6046	E0956	Lateral Pad 4x6", Each	£34.00

Fixed Trunk Supports

Requires RESPOND™ Back 3003-4001

<input type="checkbox"/>	3003-6001	E0956	LEFT, 3x4" (7.5x10cm)	N/C
<input type="checkbox"/>	3003-6002	E0956	RIGHT, 3x4" (7.5x10cm)	N/C

Harness/Belt Accessories

<input type="checkbox"/>	3001-6011	K0108	Strap Risers, Pair for 1.25" Webbing	£33.00
<input type="checkbox"/>	IFMDBC10		SP™ Buckle Cover Kit 1" Buckle	£21.00
<input type="checkbox"/>	IFMDBC15		SP™ Buckle Cover Kit 1.5" Buckle	£21.00

Leggero Head Positioning

<input type="checkbox"/>	3001-6005	E0955	8" Oval Headrest and Hardware	£146.00
<input type="checkbox"/>	3001-6006	E0955	RESPOND™ Cradle Headrest and Hardware	£146.00
<input type="checkbox"/>	3005-6099		Omit Head Positioning	N/C

Stealth Products Head Positioning

A. Link Headrest Hardware 1" Ball

<input type="checkbox"/>	TWB480-1-01	E1028	Link Hardware 1" Ball	£146.00
--------------------------	-------------	-------	-----------------------	---------

B. Link Headrest Hardware 17mm Ball/Flip Down

<input type="checkbox"/>	TWB480-17-01	E1028	Link Hdw. 17mm Ball/Flip Down	£201.00
--------------------------	--------------	-------	-------------------------------	---------

Comfort Plus Headrest (Only if A is selected)

<input type="checkbox"/>	CP250-SR	E0955	Comfort Plus 10"	£96.00
<input type="checkbox"/>	CP450-SR	E0955	Comfort Plus 8"	£96.00

i2i Head Positioning System (Only if B is selected)

<input type="checkbox"/>	i2i250	K0108	i2i Low Pro Nino XS, Pad Only	£213.00
<input type="checkbox"/>	i2i350	K0108	i2i Low Pro Nino SM, Pad Only	£213.00
<input type="checkbox"/>	i2i450	K0108	i2i Low Pro Adult SM, Pad Only	£213.00

Chin Prompts (Only if i2i Head Positioning is selected)

<input type="checkbox"/>	i2i460	K0108	i2i SM Chin Prompt	£49.00
<input type="checkbox"/>	i2i560	K0108	i2i MED Chin Prompt	£49.00

Combo Series Headrest

<input type="checkbox"/>	CB315	E0955	Combo Series, Small	£185.00
<input type="checkbox"/>	CB305	E0955	Combo Series, Medium	£185.00
<input type="checkbox"/>	CB325	E0955	Combo Series, Large	£185.00

Medical Mounting Accessories

Must be factory installed.

<input type="checkbox"/>	3002-7016	K0105	IV Pole	£112.00
<input type="checkbox"/>	3002-7017	E2208	Oxygen Tank Holder (D size)	£131.00
<input type="checkbox"/>	3002-7019	E1029	LTV Bracket for LTV950 ¹	£187.00

<input type="checkbox"/>	3002-7027	E1029	TRAK 12 Equipment Tray	£262.00
<input type="checkbox"/>	3002-7028	E1029	TRAK 12 LTV Bracket ¹ and Equip. Tray	£370.00
<input type="checkbox"/>	3002-7018	E1029	TRAK 14 Equipment Tray	£262.00
<input type="checkbox"/>	3002-7021	E1029	TRAK 14 LTV Bracket ¹ and Equip. Tray	£370.00

Combined weight of accessories & occupant must not exceed chair capacity.

Ventilator, oxygen tank, and any other accessories must be removed and secured off chair during transit in a motor vehicle.

Anti-Tips and Accessories

<input checked="" type="checkbox"/>	3002-7007	E0971	Rear Anti-Tips	N/C
<input checked="" type="checkbox"/>	3002-7004		Under Seat Storage	N/C
<input checked="" type="checkbox"/>	3002-7012		Deluxe Sunshade	N/C
<input type="checkbox"/>	3002-6113		TRAK 12 Grab Bar	£101.00
<input type="checkbox"/>	3002-7026		TRAK 12 Rain Cover	£146.00
<input type="checkbox"/>	3002-6019		TRAK 14 Grab Bar	£101.00
<input type="checkbox"/>	3002-7014		TRAK 14 Rain Cover	£146.00