

FUSION E

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Ⓢ Clients Weight _____

BASE MODEL

OPTION 1

£2,519

4 mph Motor Package
 Height Adjustable Removable Armrests
 70° Swing-Away Legrests
 Adjustable Seat Width & Depth Package
 Fold-forward Back Support with Tension
 Adjustable Upholstery
 14" Solid Drive Wheels
 10" Solid Castors
 Rear Suspension System
 Black Frame
 Yellow Wheel Accent Inserts
 VR2 Control System
 2 x 50 ah Batteries
 Max User Weight: 23 st 6lbs. (150 kg)



OPTION 2

£2,749

6 mph Motor Package
 Lights & Indicators
 Height Adjustable Removable Armrests
 70° Swing-Away Legrests
 Adjustable Seat Width & Depth Package
 Fold-forward Back Support with Tension
 Adjustable Upholstery
 14" Pneumatic Drive Wheels
 10" Pneumatic Castors
 Rear Suspension System
 Black Frame
 Yellow Wheel Accent Inserts
 VR2 Control System
 2 x 50 ah Batteries
 Max User Weight: 23 st 6lbs. (150 kg)

SEAT PACKAGE

18"W x 20"D. *Seat width can be adjusted between 16-23". Seat depth can be adjusted between 16-21".*

Standard

SEAT CUSHION

2" Foam Seat Cushions. *Supplied with 1 x 16"W x 18"D & 1 x 20" x 22".*

Standard

BACK SUPPORT POSITIONING

Angle Adjustable, Fold-Forward. *3° - 33° angle adjustment in 4 steps of 10°.*

Standard

BACK SUPPORTS

Tension Adjustable Upholstery. *23.5" tall from seat pan. Supplied with 1 x 16.5"W & 1 x 18.5"W overlays.* Standard

SEAT POSITIONING

Static Seat. *18.5" front seat to floor and 17" rear seat to floor height from seat pan.* Standard

ARM SUPPORTS

Height Adjustable, Removable. *Includes 15"L x 3"W arm pads. Adjustment range 10.5" to 14.5" from seat pan.* Standard

LEG SUPPORTS

70° Swing-away Legrests. *Adjustment range 13" to 16" from seat pan.* Standard

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release. Standard

NON-EXPANDABLE CONTROLS

VR2 Hand Control. Standard

HAND CONTROL MOUNTS

Right Side with Fixed Position Mount. *Can be moved to left side.* Standard

ADDITIONAL OPTIONS

Please note that all additional options below will be supplied uninstalled.

SEAT CUSHIONS

Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. *158.7 Kg weight limit.* £75

Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. *160 Kg weight limit.* £150

TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. *158.7 Kg weight limit.* £165

Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. *158.7 Kg weight limit.* £190

Please enter cushion size here

WIDTH

DEPTH

Refer to cushion order forms for size availability.

STEALTH PRO BACK SUPPORTS

Standard, Coolcore Cover. *Solid back, 9 cm contour.* £227

Active, Coolcore Cover. *Solid back, 12 cm contour.* £341

Extra Deep Active, Coolcore Cover. *Solid back, 17 cm contour.* £448

Active Flex, Coolcore Cover. *Tension adjustable straps. 12 cm contour.* £559

Back Height 50 cm. *Actual height.* Standard 57 cm. *Actual height.* FOC

STEALTH PRO BACK SUPPORT WIDTH

16"W.	16.5"W.	17"W.
18"W.	19"W.	20"W.
21"W.	22"W.	

STANDARD HEAD SUPPORTS

Fusion 14" Adjustable Comfort. <i>FRMASMB16245..</i>	£65
--	-----

STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM. Not compatible with TAB or Active Flex Pro Back Support.</i>	£219
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM. Not compatible with TAB or Active Flex Pro Back Support.</i>	£219
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M/ ACP150M. Not compatible with TAB or Active Flex Pro Back Support.</i>	£327

CONTROL OPTIONS

Attendant Control. <i>Bracket to mount hand control to backrest.</i>	£100
Dual Controls. <i>Hand control and secondary attendant control mounted to back bar. Only available with option 2 (6 mph).</i>	£430

NOTES & SPECIAL REQUESTS

PRICE

FUSION E



OPTION 1

BASE PRICE	£2,519
COST OPTIONS	
TOTAL	

OPTION 2

BASE PRICE	£2,749
COST OPTIONS	
TOTAL	



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
RETAIL UK - 05/12/24 V1

