

PRODUCT SPECIFICATIONS SHEET

Class of Use	A (Group 2)
Weight Capacity	200 lbs. (91 kg/14 stone 4 lbs.) 250 lbs. (113 kg/17 stone 12 lbs.) United Kingdom with Swing-away Leg Rests only
Maximum Speed ¹	Up to 3.6 mph (5.8 km/h)
Ground Clearance ²	1.6 in. (42 mm) at anti-tip bracket 1.9 in. (45 mm) at motor
Turning Radius ²	25.45 in. (645 mm) with foot platform
Overall Size ^{2, 6}	Length: 33.09 in. (838 mm) with front riggings Width: 22.6 in. (574 mm)
Drive Wheels	9 in. (229 mm), solid
Front Caster Wheels	6 in. (152 mm), solid
Rear Anti-tip Wheels	3 in. (76 mm), solid
Suspension	Limited
Range ^{1, 4}	Up to 13.8 miles (22.21 km)
Component Weights ^{2, 5}	Base: 67 lbs. (30.4 kg) Seat: 48.5 lbs. (22 kg) Batteries: 18 Ah: 31 lbs. (14.06 kg) (battery box)
Drivetrain	Two motor, rear-wheel drive
Brakes	“Intelligent Braking” electronic regenerative, disc park brake
Batteries ³	Two 12 volt, deep-cycle: 18 Ah (standard)
Battery Charger	3.5 amp, off-board
Freewheel Operational Force	Less than 60 N
Drive Control Operational Force	Less than 5 N
Maximum Rated Slope	6° (10.5%)
Maximum Obstacle Climbing Ability (Ascending/Descending)	1.57 in. (40 mm)
Measured Sound Level (A-weighted)	Less than 65 dB

- ¹ Varies with base model, user weight, terrain type, battery amp hour rating (Ah), battery charge, battery condition, motors, controller type, tire type, and tire condition. This specification is subject to a variance of +10%, -5%.
- ² Due to manufacturing tolerances and continued product improvement, this specification can be subject to a variance of (+ or -) 3%.
- ³ AGM, Gel-Cell, or Hybrid Gel (Hybrid Gel is defined as Nano Silica Gel Electrolyte with High Quality Glass Mat Separators) type required. See “Batteries and Charging.”
- ⁴ Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- ⁵ Battery weight may vary based on manufacturer.
- ⁶ Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

NOTE: This product conforms to all applicable ANSI/RESNA, ISO 7176 Series, AS/NZS 3695.2, and EN12184 standards. All specifications subject to change without notice.



Rev E/March 2020






KOZMO / KOZMO PORTABLE


PRODUCT SPECIFICATIONS SHEET

LENGTH (WITH FRONT RIGGINGS)

33.09 in. (838 mm)




TURNING RADIUS



WIDTH

22.6 in. (574 mm)



1.9 in. (45 mm)

**GROUND CLEARANCE
AT MOTOR**

1.6 in. (42 mm)

**GROUND CLEARANCE
AT ANTI-TIP BRACKET**

Rev E/March 2020



I N F I N F B 3 8 7 8