

# Quantum®

32 Wedgwood Road, Bicester,  
Oxon, OX26 4UL  
Phone: 01869 324600  
Fax: 01869 323070  
Email: sales@quantumrehab.co.uk  
Web: www.quantumrehab.co.uk



## QUANTUM®

QUANTUM  
**KOZMO**  
PORTABLE



April 2022 (V1.4)



### Kozmo Portable

Dealer Name: \_\_\_\_\_ Order Date: \_\_\_\_\_

Address: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_ Post Code: \_\_\_\_\_

Phone \_\_\_\_\_ Email: \_\_\_\_\_

PO No: \_\_\_\_\_ Order Ref: \_\_\_\_\_

Delivery Address: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_ Post Code: \_\_\_\_\_

Quantum Specialist: \_\_\_\_\_

### Standard Features

- Available Seat Widths = 12" - 18"
- Available Seat Depths = 12" - 18"
- 2 x 18Ah batteries
- GC3 Joystick
- Max User Weight 17 Stone 8lbs (113 kg)
- 6" Flat Free Front Castors
- 9" Flat Free Drive Tyres
- 25.4" Turning Radius
- 3.6 Mph (5.79 Km/h) Motors
- Dismantles into 4 pieces
- Quick Height Adjustable Armrests
- Flip-up foot Platform
- 6 Base Colours
- Padded Backrest
- Auto Buckle Lap Belt, Push Button Release

**Standard = Standard Factory Feature | FOC = Select instead of Standard Factory Feature**

### Base Model

**Kozmo - £2,250.00**

### Drive Wheels

9" Flat Free, Black Tyres Standard

### Castor Wheels

6" Flat Free, Black Tyres Standard

### Colours



Ruby Red



Pearl White



Rose Quartz



Citrine Yellow



Amber Orange



Sapphire Blue

### Seat Dimensions

Select width and depth package and enter desired seat width and depth in that package range. Available in 1" increments.

#### Seat Width Packages

12"W-16"W       16"W-18"W

Set Width at: \_\_\_\_\_

#### Seat Depth Packages

12"D-16"D       16"D-18"D

Set Depth at: \_\_\_\_\_

### Seat Heights

Standard Seat Post For 19" or 20" seat pan to floor heights

Short Seat Post For 16" or 17" seat pan to floor heights

#### Seat to Floor Heights

16"       17"       19"       20"

