

## Q6 EDGE 2.0 WITH TB FLEX SEATING

### DEALERS DETAILS

Dealers Name \_\_\_\_\_

Order Date \_\_\_\_\_

Address \_\_\_\_\_

Postcode \_\_\_\_\_

Phone \_\_\_\_\_ Email \_\_\_\_\_

PO No \_\_\_\_\_ Order Ref \_\_\_\_\_

### DELIVERY DETAILS

Dealers Address \_\_\_\_\_

Postcode \_\_\_\_\_

Quantum Specialist \_\_\_\_\_

Clients Weight \_\_\_\_\_

### BASE MODEL

**£4,450**

- 13 Colour Options
- Height Adjustable Removable Armrests
- 70° Swing-Away Legrests
- Angle Adjustable Fold-Forward Backrest
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- ATX Suspension System
- QLNES Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature  
FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 14" Flat Free Black Tyres.	Standard
3" x 14" Pneumatic Black Tyres.	FOC

### COLOURS

Back in Black	Standard	Lemon Crush*	FOC
Green Machine	FOC	Night Train*	FOC
Midnight Blue	FOC	Evolution Orange*	FOC
Electric Blue	FOC	Military Green*	FOC
Pretty in Pink	FOC	Grey Street*	FOC
Deep Purple	FOC	Custom Hydro Dip*	£600
Red Corvette	FOC		

\* Extended lead-time

Please download our free PQ-365 app to view colours and other product features in more detail.



### SEAT WIDTH

15" Width	FOC	16" Width	FOC
17" Width	FOC	18" Width	Standard
19" Width	FOC	20" Width	FOC

### SEAT DEPTH

15" Depth	FOC	16" Depth	FOC
17" Depth	FOC	18" Depth	Standard
19" Depth	FOC	20" Depth	FOC

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	£100
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£115
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£115
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£135
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i>	£135
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability. The nearest size will be provided unless otherwise specified

## SEAT POSITIONING

Static Seat.	Standard
Powered Seat Tilt. <i>0° to 30°. Uses 1 power option.</i>	£300
Static Seat to Floor Height. <i>From seat pan.</i>	16"      17"      18"
Powered Tilt Seat to Floor Height. <i>From seat pan.</i>	16"      17"      18"

## STANDARD BACK SUPPORT POSITIONING

Manual Angle Settable, Fold-Forward. <i>Adjustment range 95° to 125°.</i>	Standard
Attendant Controlled Manual Recline. <i>Adjustment range 95° to 125°. Not suitable for occupied transportation.</i>	£150
Powered Back Recline. <i>Adjustment range 95° to 125°. Uses 1 power option.</i>	£200

## STANDARD BACK SUPPORTS

TRU-Comfort 2. <i>Plastic shell, 4-way stretch cover. Adjustable range 24.5" to 28" from seat board.</i>	Standard
Tension Adjustable. <i>Air-mesh cover. Adjustment range 19" to 23" from seat board. OBAS** required to mount headrest.</i>	FOC
ADI Aluminium Back. <i>Skin protection, 4 way stretch, Coolcore cover.</i>	POA
Enter ADI Back Support part number here for quotation.	

## STEALTH PRO BACK SUPPORTS

Standard, Coolcore Cover. <i>Solid back, 9 cm contour.</i>	£227
Active, Coolcore Cover. <i>Solid back, 12 cm contour.</i>	£341
Extra Deep Active, Coolcore Cover. <i>Solid back, 17 cm contour.</i>	£448
Active Flex, Coolcore Cover. <i>Tension adjustable straps. 12 cm contour.</i>	£559
Back Height      50 cm. <i>Actual height.</i>	Standard      57 cm. <i>Actual height.</i>
The back support width will match the seat width selected in seat dimensions.	

## ARM SUPPORTS

TB Flex Height Adjustable, Removable, Drop-in T Arm. <i>Adjustment range 9" to 14" from seat pan.</i>	Standard
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>	£279
Arm Support Pads      10" Standard	Standard      10" Waterfall.
14" Standard	FOC      14" Waterfall.
	FOC

## LEG SUPPORTS

70° Swing-away Legrests. <i>Adjustment range 12" to 15".</i>	Standard
Base Mounted Foot Platform. <i>Mounted to chassis. Will not move with power tilt.</i>	FOC
60° to 90° Tool Adjustable Swing-away Legrests. <i>Adjustment range 11" to 16.5".</i>	£330
Manual Elevating Swing-Away Legrests. <i>Adjustment range 15" to 20".</i>	£250
Manually Adjustable Centre Mount Foot Platform. <i>Adjustment range range 10" to 16". 5"W x 11"L Footplates.</i>	£340
Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away legrest hanger for mounting.</i>	£300

Limb Support	<b>LEFT</b>	<b>RIGHT</b>
Seat pan to heel measurement		(Leave blank for standard factory)

Adjustment ranges above are from seat pan to footplates. Individual legrests supplied with heel loops. Elevating legrests supplied with calf supports.

## STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£219
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	£219
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i>	£327
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£88
Push Handles Head Mount Adapter OBAS**. <i>Required to mount head support with TA Back or Active Flex Pro Back.</i>	£266

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative head supports.

## TRUNK SOLUTIONS

Thoracic Lateral Supports.	<b>LEFT</b>	£150	Lateral Thoracic Swing-away Hardware.	£150	<b>RIGHT</b>			
		£35	3" H x 5" L Curved Lateral Pad.	£35				
		£35	4" H x 6" L Curved Lateral Pad.	£35				
Stealth 4-Point Harness	Extra Small	£137	Small	£137	Classic	STD*	Contoured	FOC
	Medium	£137	Large	£137	Dynamic	STD*	Structured	FOC
	Extra Large	£137			Vinyl Face	FOC	Zipper Face**	FOC
Chest Strap.	Extra Small	£87	Small	£87	Side Release			STD*
	Medium	£87	Large	£87	Velcro Release			£19
	Extra Large	£87						

\*STD = Standard option \*\*Classic Only

## PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.	Standard (STD)							
Retractable Lap Belt, Push Button Release.	£156							
Stealth Positioning 2-Point Padded Pelvic Belt.	Extra Small	£78	Small	£78	Single Pull	STD*	Dual Pull	FOC
	Medium	£78	Large	£78	Rear Pull	FOC	Side Release	STD*
	Extra Large	£78			Push Button	FOC		
Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull.	Extra Small	£120	Small	£120	Push Button	STD*		
	Medium	£120	Large	£120	Side Release	£20		

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative lateral, pelvic & lower extremity solutions \*STD = Standard option

## NON-EXPANDABLE CONTROLS

QLNES Controls.	FOC	
QLNES Controls. <i>2 power option ready.</i>	£125	
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	£100
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

## HAND CONTROL MOUNTS

Fixed.	Standard					
Basic swing away.	£75					
Advanced swing away. <i>Chain driven. Swings in &amp; out with height, tension and pivot adjustment.</i>	£100					
Mounting Position.	Left	FOC	Right	Standard		
Hand Control Joystick Handle.	Standard Carrot	STD	T-Bar	£35	Extended I-Handle	£35
	Mushroom	£35	Medium Ball	£35	3" Goal Post	£50
	Chin Cup	£35	Soft Ball	£35	4" Goal Post	£50

## SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
4 Mph Program with Lights and indicators.	£450
6 Mph Program with Lights and indicators.	£450
Lights & Indicator Surcharge. <i>Required if no power options have been selected.</i>	£260

## ACCESSORIES

Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379
Quantum Back Pack.	£110
XLR USB Charger Port.	£17

## NOTES & SPECIAL REQUESTS

## PRICE

**Q6 EDGE 2.0 WITH  
TB FLEX SEATING**



**BASE PRICE**

**£4,450**

**COST OPTIONS**

**TOTAL**

**4**



**Quantum®**  
32 Wedgwood Road, Bicester,  
Oxon, OX26 4UL  
Phone: 01869 324600  
Fax: 01869 323070  
Email: sales@quantumrehab.co.uk  
Web: www.quantumrehab.co.uk  
**RETAIL UK - 02/12/24**

