

Q6 EDGE 2.0 WITH TB FLEX SEATING

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Phone _____ Email _____

PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

Quantum Specialist _____

Clients Weight _____

BASE MODEL

£4,750

- 13 Colour Options
- Height Adjustable Removable Armrests
- 70° Swing-Away Legrests
- Angle Adjustable Fold-Forward Backrest
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- ATX Suspension System
- QLNES Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres.	Standard
3" x 14" Pneumatic Black Tyres.	FOC

COLOURS

Back in Black	Standard	Lemon Crush*	FOC
Green Machine	FOC	Night Train*	FOC
Midnight Blue	FOC	Evolution Orange*	FOC
Electric Blue	FOC	Military Green*	FOC
Pretty in Pink	FOC	Grey Street*	FOC
Deep Purple	FOC	Custom Hydro Dip*	£600
Red Corvette	FOC		

* Extended lead-time

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT WIDTH

15" Width	FOC	16" Width	FOC
17" Width	FOC	18" Width	Standard
19" Width	FOC	20" Width	FOC

SEAT DEPTH

15" Depth	FOC	16" Depth	FOC
17" Depth	FOC	18" Depth	Standard
19" Depth	FOC	20" Depth	FOC

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	£99
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£126
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£126
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£148
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i>	£148
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£209
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£209

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability. The nearest size will be provided unless otherwise specified

SEAT POSITIONING

Static Seat.	Standard
Powered Seat Tilt. <i>0° to 30°. Uses 1 power option.</i>	£500
Static Seat to Floor Height. <i>From seat pan.</i>	16" 17" 18"
Powered Tilt Seat to Floor Height. <i>From seat pan.</i>	16" 17" 18"

STANDARD BACK SUPPORT POSITIONING

Manual Angle Settable, Fold-Forward. <i>Adjustment range 95° to 125°.</i>	Standard
Attendant Controlled Manual Recline. <i>Adjustment range 95° to 125°. Not suitable for occupied transportation.</i>	£150
Powered Back Recline. <i>Adjustment range 95° to 125°. Uses 1 power option.</i>	£400

STANDARD BACK SUPPORTS

TRU-Comfort 2. <i>Plastic shell, 4-way stretch cover. Adjustable range 24.5" to 28" from seat board.</i>	Standard
Tension Adjustable. <i>Air-mesh cover. Adjustment range 19" to 23" from seat board. OBAS** required to mount headrest.</i>	FOC
ADI Aluminium Back. <i>Skin protection, 4 way stretch, Coolcore cover.</i>	POA
Enter ADI Back Support part number here for quotation.	

STEALTH PRO BACK SUPPORTS

Standard, Coolcore Cover. <i>Solid back, 9 cm contour.</i>	£227
Active, Coolcore Cover. <i>Solid back, 12 cm contour.</i>	£341
Extra Deep Active, Coolcore Cover. <i>Solid back, 17 cm contour.</i>	£448
Active Flex, Coolcore Cover. <i>Tension adjustable straps. 12 cm contour.</i>	£559
Back Height 50 cm. <i>Actual height.</i>	Standard 57 cm. <i>Actual height.</i>
The back support width will match the seat width selected in seat dimensions.	

ARM SUPPORTS

TB Flex Height Adjustable, Removable, Drop-in T Arm. <i>Adjustment range 9" to 14" from seat pan.</i>	Standard
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>	£279
Arm Support Pads	10" Standard Standard 10" Waterfall. FOC
	14" Standard FOC 14" Waterfall. FOC

LEG SUPPORTS

Base Mounted Foot Platform. <i>Mounted to chassis. Will not move with powered seat.</i>	Standard
70° Swing-Away. <i>Extension range 12" to 15". Minimum seat width 15".</i>	FOC
60° to 90° Tool Adjustable Swing-Away. <i>6"W x 8"D foot supports. Extension range 11" to 19.5". Min. seat width 15".</i>	£330
Manually Adjustable Centre Mount Foot Platform. <i>5"W x 11"L Footplates. Extension range 10" to 19". Min. seat width 15".</i>	£340
Manual Elevating Swing-Away. <i>Extension range 15" to 20". Min. seat width 15"</i>	£425
Power Elevating Articulating Swing-Away. <i>5"W x 7"D foot supports. Ext. range 9" to 21". Min. seat width 15". Uses 2 power options</i>	£850
Powered Leg Support Operation Method. Independent. <i>Uses 2 power options.</i> Combined. <i>Uses 1 power option.</i>	
Power Elevating Articulating Centre Mount Foot Platform. <i>5"W x 11"L Footplates. Extension range 10" to 19". Min. seat width 15"</i>	£950

① Enter seat pan to heel measurement.

If NO measurement is supplied, standard factory length option will be fitted.

Stealth Residual Limb Support (each). *Includes a 60-90° Swing-away legrest hanger for mounting.*

LEFT

£300

Limb Support

£300

RIGHT

Extension ranges above are from seat pan to foot supports. Individual leg supports supplied with heel loops. Elevating leg supports supplied with standard calf supports.

STEALTH HEAD SUPPORTS

UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£230
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	£230
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i>	£342
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£92
Push Handles Head Mount Adapter OBAS**. <i>Required to mount head support with TA Back or Active Flex Pro Back.</i>	£266

Please visit stealthproducts.com for alternative head supports.

TRUNK SOLUTIONS

Thoracic Lateral Supports.	LEFT	£175	Lateral Thoracic Swing-away Hardware.	£175	RIGHT	
		£42	3" H x 5" L Curved Lateral Pad.	£42		
		£42	4" H x 6" L Curved Lateral Pad.	£42		
Stealth 4-Point Harness	Extra Small	£143	Small	£143	} Classic STD* Contoured FOC	
	Medium	£143	Large	£143		} Dynamic STD* Structured FOC
	Extra Large	£143				
Chest Strap.	Extra Small	£91	Small	£91	} Side Release STD*	
	Medium	£91	Large	£91		} Velcro Release £20
	Extra Large	£91				

*STD = Standard option **Classic Only

PELVIC SOLUTIONS

Auto Buckle Lap Belt, Push Button Release.	Standard					
Retractable Lap Belt, Push Button Release.	£169					
Stealth Positioning 2-Point Padded Pelvic Belt.	Extra Small	£84	Small	£84	} Single Pull STD* Dual Pull FOC	
	Medium	£84	Large	£84		} Rear Pull FOC Side Release STD*
	Extra Large	£84				
Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull.	Extra Small	£143	Small	£143	} Push Button STD*	
	Medium	£143	Large	£143		} Side Release £21

Please visit stealthproducts.com for alternative lateral, pelvic & lower extremity solutions *STD = Standard option

NON-EXPANDABLE CONTROLS

QLNES Controls.	FOC	
QLNES Controls. <i>2 power option ready.</i>	£125	
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	£100
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

HAND CONTROL MOUNTS

Fixed.	Standard					
Basic swing away.	£75					
Advanced swing away. <i>Chain driven. Swings in & out with height, tension and pivot adjustment.</i>	£100					
Mounting Position.	Left	FOC	Right	Standard		
Hand Control Joystick Handle.	Standard Carrot	STD	T-Bar	£35	Extended I-Handle	£35
	Mushroom	£35	Medium Ball	£35	3" Goal Post	£50
	Chin Cup	£35	Soft Ball	£35	4" Goal Post	£50

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
4 Mph Program with Lights and indicators.	£450
6 Mph Program with Lights and indicators.	£450
Lights & Indicator Surcharge. <i>Required if no power options have been selected.</i>	£260

ACCESSORIES

Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379
Quantum Back Pack.	£110
XLR USB Charger Port.	£17

NOTES & SPECIAL REQUESTS

PRICE

**Q6 EDGE 2.0 WITH
TB FLEX SEATING**



BASE PRICE **£4,750**

COST OPTIONS

TOTAL

4



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
RETAIL UK - 01/04/25

