

Quantum®

32 Wedgwood Road, Bicester,
Oxon, OX26 4UL

Phone: 01869 324600

Fax: 01869 323070

Email: sales@quantumrehab.co.uk

Web: www.quantumrehab.co.uk



Dealer April 2021



**Crash
Tested**



Q6 Edge 2.0 with Synergy® Seating & VR2 or R-Net Controls

Dealer Name:

Order Date:

Address:

Post Code:

Phone

Email:

PO No:

Order Ref:

Delivery Address:

Post Code:

Quantum Specialist:

Standard Features

- Synergy Seating System
- 4" of Seat Width Adjustment
- 4" of Seat Depth Adjustment
- 2 x 55ah batteries
- VR2 Control System
- 2" x 6" Front & Rear Castors, Flat Free
- 3" x 14" Pneumatic Drive Tyres
- 12 Colour Options
- 4 Mph (6 Km/h) Motors
- ATX Suspension System
- Single Post Height Adj. Removable Armrests
- Centre Mount Foot Platform
- Maximum User Weight 21 Stone 6lbs
- Synergy Simplicity Seat Cushion
- Tie-Down System for Transportation

Base Model

Q6 Edge 2.0 with Synergy® Seating and VR2 Controls - £4,595.00

Drive Wheels

- | | |
|---|----------|
| ☐ 3" x 14" Pneumatic, Black Tyres | Standard |
| ☐ 3" x 14" Flat Free, Black Tyres | £150.00 |
| ☐ 4" x 14" Pneumatic, Black Tyres, Knobby, Extended Traction Tyres <i>Increases overall base width by 1.5"</i> | £220.00 |
| ☐ 4" x 14" Flat Free, Black Tyres, Knobby, Extended Traction Tyres <i>Increases overall base width by 1.5"</i> | £320.00 |

Castor Wheels

- | | |
|--|----------|
| ☐ 2" x 6" Flat Free, Black Tyres, Front and Rear | Standard |
| ☐ 2" x 8" Flat Free, Black Tyres, Front and Rear <i>Includes extended castor arms</i> | £364.00 |

Colours



Matching back shroud included where applicable.

Q6 EDGE 2.0 SYNERGY SEATING - DEALER

Seat Dimensions

Select width and depth package and enter desired seat width and depth in that package range. Available in 1" increments. Maximum seat size is 22"W x 22"D.

Seat Width Packages

☐ 14"W-18"W ☐ 16"W-20"W ☐ 20"W-22"W

Set Width at:_____

Seat Depth Packages

☐ 14"D-18"D ☐ 16"D-20"D ☐ 20"D-22"D

Set Depth at:_____

Seat Cushions

Please refer to our Stealth Cushion Order form for available sizes and enter your preferred size here.....W xD

- ☐ No Seat Cushion Required FOC
- ☐ Stealth Simplicity G: General use, Fluid resistant Stealth-Tek cover - 25 stone weight limit FOC
- ☐ Stealth TRU-Comfort 2 SPP: Skin protection, positioning, 4-way stretch Coolcore cover - 32 stone weight limit Standard
- ☐ Stealth Solution® SPP: Skin protection, positioning, 4-way stretch, Coolcore cover and urethane liner - 25 stone weight limit £100.00
- ☐ Stealth Essence® SPP: Skin protection, positioning, 4-way stretch and Coolcore cover - 35 stone 7lbs weight limit £120.00
- ☐ Stealth Spectrum® Gel SPP: Skin protection, positioning, 4-way stretch, Coolcore cover and urethane liner - 25 stone weight limit £160.00
- ☐ Stealth Spectrum® Foam SPP: Skin protection, positioning, 4-way stretch, Coolcore cover and urethane liner - 42 stone weight limit £160.00

Synergy® Backrest Specifications

- ☐ Synergy® Backrest Assembly Including Visco padded backrest Standard
- ☐ Tension Adjustable Backrest Upholstery This option may increase the delivery time £120.00
- ☐ Back Canes Only Suitable for third party backrests. A backrest will not be supplied FOC

Backrest/Back Cane Height:

Set Backrest/Back Cane Height at:_____ Available in 1" increments. Please refer to the Synergy Backrest Matrix below for seat width and backrest height compatibility. Backrest Width must match Seat Width selected unless using 1" or 2" Endo brackets in Backrest Positioning Adjustments section, in which case Backrest Width can be 2" or 4" narrower than Seat Width.

Back Height = Seat Pan to Top of Back | Odd heights = 3" gap at bottom, Even heights = 4" gap at bottom.

		Synergy Static Back Width																				
		10"	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"	25"	26"	27"	28"	29"	30"	32"
Synergy Back Height	12"																					
	13"	✓																				
	14"	✓																				
	15"		✓		✓																	
	16"		✓		✓																	
	17"			✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓							
	18"			✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓							
	19"				✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓
	20"				✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓
	21"				✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	
	22"				✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	
	23"				✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	
24"				✓		✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓		✓		
25"									✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			✓	

ADI Back

- ☐ ADI Backrest Refer to ADI Backrest configurator on page 7 and enter part number and price here. Part No:..... £_____

TRU-Comfort Plus® Backrest Specifications

- ☐ TRU-Comfort Plus® Backrest Assembly 4-way stretch fabric (Positioning). Only available with seat widths 16", 18", 20" and 22" £220.00

Backrest Height:

- ☐ 27" Tall Standard
- ☐ 24" Tall FOC

Back Height = Seat Pan to Top of Back | 3" gap at bottom.

Seat Positioning Adjustments

- ☐ Static Seat Standard
- ☐ TRU-Balance®2 Power Tilt 0-55 degree centre of gravity (CG) power tilt seating system. Uses 1 power option £1100.00
- ☐ Tru-Balance®2 10" Power Adjustable Seat Lift and Tilt Requires R-Net Controls. Uses 2 power options. £1400.00

Seat to Floor Height

- Static Seat ☐ 15.875" Seat to Floor ☐ 16.875" Seat to Floor ☐ 17.875" Seat to Floor
- TB2 Power Tilt ☐ 17.125" Seat to Floor ☐ 18.125" Seat to Floor ☐ 19.125" Seat to Floor
- 10" Seat Lift and Tilt ☐ 17.875" Seat to Floor ☐ 18.875" Seat to Floor ☐ 19.875" Seat to Floor

Backrest Positioning Adjustments

☐ Static Backrest Tool angle adjustable from 85° to 120° in 5° increments	Standard
Additional Static Back Options:	
☐ Flip Forward Backrest <i>Not suitable for occupied transportation</i>	£200.00
☐ 1" Endo Brackets <i>Reduces backrest width by 2". Not available with TRU-Comfort Plus backrest on seat widths 14"-17"</i>	£40.00
☐ 2" Endo Brackets <i>Reduces backrest width by 4". Not available with TRU-Comfort Plus backrest on seat widths 14"-19"</i>	£40.00

Armrests

☐ Single Post, Quick Height Adjustable Removable Armrests <i>Height adjustment from 11" - 17" in 1/2" increments</i>	Standard
☐ Paediatric, Removable, Height Adjustable <i>Height adjustment from 7.75" to 9.75" in 1/2" increments. Not available with full length armpads</i>	FOC
☐ Single Post, Tool Height Adjustable Armrests, Heavy Duty, Removable <i>Height adjustment from 10" - 14"</i>	£220.00
☐ 2-Post, Flip Back, Height Adjustable, Removable <i>Height range 9"-13". Not available with a seat depth less than 16"</i>	£145.00
☐ Cantilever Flip-Up Armrests <i>Not available with a tension adjustable backrest. Increases overall width by 7". Height adj. 11"-17" depending on backrest height. Not available with Flip Forward Backrest.</i>	£280.00

Armrests Pad Options

☐ 14"L x 2"W Straight, Pair	☐ 10"L x 2"W Straight, Pair	Standard
☐ 14"L x 3"W Waterfall, Pair	☐ 10"L x 3"W Waterfall, Pair	FOC
☐ 14"L x 3.5"W Gel Armpad, Pair		£160.00
☐ 10"L x 3.5"W Gel Armpad, Pair		£130.00

Headrests

☐ Stealth 10" Comfort Plus CP250	£140.00
☐ Stealth 14" Comfort Plus CP150	£140.00
☐ Stealth 14" Adjustable Comfort Plus ACP150	£250.00
☐ Stealth Small Contoured STL-OTSM250	£200.00
☐ Stealth Large Contoured STL-OTLG150	£200.00

Lateral and Trunk Supports

Stealth Thoracic Lateral Supports

Left Lateral Package

☐ Thin Hardware, Short Link TWBTL-L	£175.00
☐ Thin Hardware, Standard Link TWBSTL-L	£175.00

Left Lateral Pads

☐ 3" H x 5" L Curved Lateral Pad TWBL-3X5	£35.00
☐ 4" H x 6" L Curved Lateral Pad TWBL-4X6	£35.00

Right Lateral Package

☐ Thin Hardware, Short Link TWBTL-R	£175.00
☐ Thin Hardware, Standard Link TWBSTL-R	£175.00

Right Lateral Pads

☐ 3" H x 5" L Curved Lateral Pad TWBL-3X5	£35.00
☐ 4" H x 6" L Curved Lateral Pad TWBL-4X6	£35.00

i-Fit Chest Harness *Supplied uninstalled*

☐ Classic Style <i>Offers a large padded surface for the thoracic area</i>	£100.00
☐ Contoured Style <i>Narrowed sectioned over the sternum that allows a better fit for females, G-tube users, or those requiring a smaller support surface</i>	£100.00

Size

☐ Extra Small	☐ Small	☐ Medium	☐ Large	☐ Extra Large
--------------------------------------	--------------------------------	---------------------------------	--------------------------------	--------------------------------------

Construction

☐ Dynamic <i>Allows for a degree of trunk movement and softer fit</i>	☐ Structured <i>More resistant to high tone, aiding maintenance of correct posture</i>
--	---

Pelvic Supports

Stealth Pelvic Supports *May conflict with armrest option. Contact Quantum for more information*

Left Pelvic Package

☐ Quick Release Removable Hardware <i>INDPART2805</i>	£100.00
☐ Flip-Down Hardware TWBL-L	£150.00
☐ Additional link for Flip-down Hardware LHW-129	£35.00
☐ Full Surface Contact Feature for Flip-Down Hardware <i>LTHW511</i>	£50.00

Right Pelvic Package

☐ Quick Release Removable Hardware <i>INDPART2805</i>	£100.00
☐ Flip-Down Hardware TWBL-R	£150.00
☐ Additional link for Flip-down Hardware LHW-129	£35.00
☐ Full Surface Contact Feature for Flip-Down Hardware <i>LTHW511</i>	£50.00

Pelvic Supports (continued)

Left Pelvic Pads		Right Pelvic Pads	
☐ 4" H x 4" L Pelvic Pad TWBADD-4x4	£35.00	☐ 4" H x 4" L Pelvic Pad TWBADD-4x4	£35.00
☐ 4" H x 8" L Pelvic Pad TWBADD-4x8	£35.00	☐ 4" H x 8" L Pelvic Pad TWBADD-4x8	£35.00
☐ 4" H x 12" L Pelvic Pad TWBADD-4x12	£35.00	☐ 4" H x 12" L Pelvic Pad TWBADD-4x12	£35.00

Stealth Knee Adductors		Right Knee Support Package	
Left Knee Support Package			
☐ Cane Mounted Hardware <i>Only available with individual legrests</i> KN1ADD	£145.00	☐ Cane Mounted Hardware <i>Only available with individual legrests</i> KN1ADD	£145.00
☐ Swing-Away Track Mount Hardware <i>Mounts to seat rail</i> KN1ADD-TRK <i>May conflict with tie-down bracket</i>	£160.00	☐ Swing-Away Track Mount Hardware <i>Mounts to seat rail</i> KN1ADD-TRK <i>May conflict with tie-down bracket</i>	£160.00
Left Knee Pads <i>With 1/2" medium sunmate</i>		Right Knee Pads <i>With 1/2" medium sunmate</i>	
☐ 3" Round Knee Adductor Pad KN1ADD3	£17.00	☐ 3" Round Knee Adductor Pad KN1ADD3	£17.00
☐ 4" Round Knee Adductor Pad KN1ADD4	£17.00	☐ 4" Round Knee Adductor Pad KN1ADD4	£17.00

i-Fit Pelvic Belts

☐ Padded <i>Offers a large padded surface area</i>	£70.00
☐ Non-Padded <i>Narrowed sectioned for smaller support surface</i>	£65.00
Size	
☐ Extra Small <i>Not available with non-padded</i>	☐ Small
☐ Medium	☐ Large
☐ Extra Large <i>Not available with non-padded or 4-point.</i>	
Attachment Method	
☐ 2-Point	☐ 4-Point <i>Single pull only</i>
Release Mechanism	
☐ Push Button	☐ Side Release
Adjustment	
☐ Single Pull	☐ Dual Pull

Medial (Pommel) Supports

Stealth Pommel Mounting Hardware <i>May conflict with centre mount foot platform. Please contact Quantum for more information</i>	
☐ Fixed Modular Removable Medial Hardware <i>Push button release. MTHW-650</i>	£225.00
☐ Flip-Down Link Style Medial Hardware. MTHW-620	£170.00
☐ Fixed Flip-Down Medial Hardware. MTHW-640	£135.00
Stealth Pommel Pads	
☐ Small Wedge Pad 3" W x 4" L x 3" H x 2" W. MTMD-12	£95.00
☐ Medium Wedge Pad 4.5" W x 5" L x 3" H x 2" W. MTMD-13	£95.00
☐ Large Wedge Pad 5" W x 5" L x 4" H x 3" W. MTMD-14	£95.00
☐ Adjustable Pommel with Dual 3" x 5" pads <i>Includes slip cover. MTMD-42</i>	£185.00

Centre Mount Foot Platforms

☐ Base Mounted Foot Platform <i>Mounted to chassis, will not move with Power Tilt or Seat Lift</i>	FOC
☐ Adjustable Centre Mount Foot Platform <i>Mounted to Seat.</i>	Standard
☐ Power Articulating Centre Mount Foot Platform <i>Mounted to Seat. Uses 1 power option.</i>	£642.00

A 16"D-20"D Seat Depth Package will be supplied if TB2 Power Tilt has been selected in combination with a Centre Mount Foot Platform. This is to accommodate the Power Tilt Module and the Centre Mount Foot Platform attachment. Seat Depths of 14"D-18"D can still be accommodated using a 16"D-20"D Seat Depth Package.

Legrest Options	
Lower Extensions <i>Seat pan to Heel</i>	
☐ 10-16"L <i>Seat pan to heel</i>	☐ 13-19"L <i>Seat pan to heel</i>
Footplates	Standard
☐ 10"W x 11"D Split Footplates	☐ 10"W x 8"D Split Footplates
Calf Pads	Standard
☐ Small 6"H x 6"W	☐ Large 8"H x 8"W <i>tapered to 6"W</i>
Calf Pad Covers (Pair) <i>1/2" Soft Visco Foam</i>	
☐ Air Mesh	☐ Reversed Dartex
	£42.00

Centre Mount Foot Platforms (continued)

- ☐ Calf Panel 8"H x 18"W. Not available with Calf Pads £42.00
- ☐ Calf Panel Extension (each) £30.00
- Calf Panel Cover 1/2" Soft Visco Foam
- ☐ Air Mesh or ☐ Reversed Dartex £72.00

Footplate Extension Kits



- ☐ 4" Width Extension Kit. Not available with 2" Bell Shaped, Length Extension Kit £16.00



- ☐ 2" Bell Shaped, Length Extension Kit. Not available with 4" Width Extension Kit £16.00

Additional Footplate Accessories:

- ☐ Footplate Outside Edge, pair £10.00
- ☐ Wheel Bumpers, pair £10.00
- ☐ Footplate Connector £10.00

Individual Legrests

Choose instead of centre mount foot platform options

- ☐ 60-90° Swing-Away Legrests Tool adjustable in 5° Increments. Factory set at 70°. Not available with a seat width less than 16" £110.00

Footplate Selection Includes heel loops

- ☐ Flip-up, Extruded Footplates, 8.5"W x 6"D Standard
- ☐ Angle Adjustable Footplates, 6"W x 8"D Other sizes available £115.00

Please provide the seat pan to heel measurement required, minimum 11", maximum 19.5": _____

- ☐ 70° Heavy duty, drop in, swing-away legrests with angle adjustable footplates. Not available with a seat width less than 16" £130.00

Footplate Selection Includes heel loops

- ☐ 6"W x 8"D ☐ 8"W x 11"D Other sizes available

Please provide the seat pan to heel measurement required, minimum 12", maximum 18": _____

- ☐ 90° Paediatric Swing-Away Legrests Seat Pan to Heel range of adjustment 7"-10"L. Includes 4" x 6"D angle adjustable footplates and heel loops. Other sizes available £110.00

Elevating Legrests

Choose instead of centre mount foot platform options

- ☐ Manual Elevating Swing-Away Legrests Not suitable for occupied transportation. Not available with a seat width less than 16" £130.00

- ☐ Supply base mounted foot platform to allow power chair to remain crash tested Manual elevating legrests must be removed during transport FOC

Footplate Selection Includes heel loops

- ☐ Flip-up, Extruded Footplates, 8.5"W x 6"D & 14"-17" Lower leg extensions Standard
- ☐ Angle Adjustable Footplates, 6" W x 8"D & 14.5"-18" Lower leg extensions £115.00

- ☐ 70° Power Articulating Elevating Swing-Away Legrests 4" articulation range. Not available with a seat width less than 16". Uses 2 power options £1320.00

Footplate Selection Includes heel loops

- ☐ Angle Adjustable Footplates, 6"W x 8"D Other sizes available Standard

Please provide the seat pan to heel measurement required, minimum 11", maximum 18": _____

Accessories:

- ☐ Small Calf Pads, Pair 5.35"W x 6"H or ☐ Large Calf Pads, Pair 5.86"W x 8"H £88.00
- ☐ Gel Padded Adductor Buttons, Pair £30.00

Operation Method: ☐ Independent ☐ Combined

Additional Legrest Options

- ☐ Stump Board support (each) Includes a 60-90° Swing-Away Legrest hanger ☐ Left ☐ Right £300.00
- ☐ Residual Limb support (each) ☐ Left ☐ Right £180.00
- If both a left and right support are required, please amend price accordingly

VR2 Non-Expandable Controls

☐ VR2 - No Power Options	Standard
☐ VR2 - 2 Power Option Ready <i>Required when 1 or 2 power options are selected</i>	£125.00
☐ Mount Joystick as Attendant Control Only <i>Please select VR2 power option above</i>	£100.00
☐ Dual Control <i>Please select VR2 power option above</i>	£380.00

R-Net Expandable Controls

☐ R-Net - 1 Power Option Ready <i>Required when 0 or 1 power option is selected</i>	£350.00
☐ R-Net - 5 Power Option Ready <i>Required when 2 or more power options have been selected</i>	£550.00
☐ User Joystick not required <i>An R-Net hand control will not be supplied. An Omni display must be selected</i>	-£100.00
☐ Mount Joystick as Attendant Control Only <i>Please select R-Net power option above</i>	£100.00
☐ Dual Control <i>Please select R-Net power option above</i>	£430.00
☐ R-Net OMNI2 Display <i>Includes 2 egg switches and swing-away display mount. Does not include R-Net Joystick. Select R-Net power option above</i>	£763.00

Joystick Handles

☐ Standard Carrot Handle	Standard
☐ Mushroom Handle	£35.00 (each)
☐ Extended I-Handle	
☐ Bodypoint 3 Inch Goal Post Handle	
☐ T-Bar	
☐ Soft Ball Handle	£35.00 (each)
☐ Bodypoint 4 Inch Goal Post Handle	
☐ Chin Cup	
☐ Medium Ball Handle	

Joystick Mounting Brackets

☐ Fixed Inline	Standard
☐ Swing-Away	£75.00
☐ Multi-Axis	£125.00
☐ Retract 4 Swing-Away Joystick Arm Kit <i>Full 180° Movement</i>	£150.00

Mounting Position

☐ Left	or	☐ Right
☐ Inline	or	☐ Offset <i>Not available with Fixed Joystick Bracket</i>

Speed Program and Lights Package

☐ 4 Mph program <i>Available with VR2 or R-Net. Lights and indicators not included</i>	Standard
☐ 4 Mph program, with Lights and indicators	£450.00
☐ 6 Mph program, with Lights and indicators	£450.00
☐ Lights and indicator kit surcharge <i>If no power options have been selected, you must tick this box in addition to the lights and indicator package to receive the correct lighting module.</i>	£260.00

Quantum Accessories

Some accessories will be shipped with the unit, not installed. Multiple accessories can be mounted along the various tracks, some interference may occur.

☐		Glove Box unit <i>uninstalled</i>	£30.00
☐		Clothing Guard (pair) <i>Only available with Flip Up armrests</i>	£100.00
☐		Quantum Backpack	£40.00

Price

	Sub Total:	£4,595.00
	Options:	
	Total:	

ALBR10LAF

Shell Type

Back Rest

Back Width

Height

Contour

Mounting

AL - AL Series	<u>Shell Assembly Weight</u>	10" / 25cm	L	A	F
AFT - AFT Series	AL Series - 1.5lb / 680gm	11" / 28cm	low (10"/25cm)	active (2.5"/63.5mm)	fixed elite (.60lb / 272gm)
CF - CF Series	CF Series - .820lb / 370gm	12" / 30.5cm	M	D	2
	above weights are of a 16" STD	13" / 33cm	medium (13"/33cm)	deep (4"/100mm)	2-point pro (1.25lb/566gm)
	- all weights approx.	14" / 35.5cm	S	6	4
		15" / 38cm	standard (16"/40.5cm)	extra deep (6"/150mm)	4-point pro (2.5lb/1132gm)
		16" / 40.5cm	T		R
		17" / 43cm	tall (18"/45.7cm)		quick release (1.30lb / 589gm)
		18" / 45.5cm	X		W
		19" / 48cm	x-tall (20"/51cm)		width adjustable elite (.63lb / 285gm)
		20" / 51cm			B
					Basic release (1.50lb / 680gm)

Notes

For additional seating options and alternative drive controls please refer to your Quantum Specialist.

Quantum Notes:

1. This product is suitable for occupied transportation in accordance with ISO 7176-19. We strongly recommend that all heavy accessories such as manual ELRs are removed and secured during occupied transportation. We also advise that the use of third party accessories including life Support Systems, Chin Controls and Trays be risk assessed prior to the wheelchair being transported. Full details of suitable restraints for this wheelchair can be obtained from Quantum Technical on 01869 324600 or Email us on: technical@quantumrehab.co.uk
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or – 3 %. Overall measurements will vary based on seating and accessory selections.