

The Fusion E is a cost effective rear wheel drive power chair with a modular state of the art design that provides improved obstacle climbing performance and a smooth ride. The Fusion E is a versatile chair that will meet the demands of even the most active user.



| Specifications | | |
|--|--|--|
| Power Chairs | | |
| Weight Capacity | 23 stone 6 lbs (150 kg) | |
| Overall Product Weight ⁵ | 211.64 lbs (96 kg) | |
| Maximum Speed ¹ | Up to 4 mph & 6 mph Speed Options | |
| Ground Clearance ² | 40 mm (1.57") at Lowest Point of Frame | |
| Turning Radius ² | 641.5 mm (25.25") | |
| Overall Size ² | Length: Width: | 1050 mm (41.33") with Front Riggings 630 mm (24.8") |
| Maximum Obstacle Climbing Ability | 100 mm (3.94") | |
| Seat Sizes | 15"W x 15"D to 22"W x 23"D | |
| Controls | VR2 Controls | |
| Angle Adjustable Backrest | - 5°, 15°, 25°, 35° | |
| Wheels | | |
| Drive Wheels | 355 mm (14") Solid/Pneumatic/ Rear Wheel Drive | |
| Caster Wheels | 254 mm (10") Pneumatic | |
| Anti-Tip Wheels | 50.8 mm (2") | |
| Suspension | Independent Rear Suspension | |
| Drivetrain | Two 350w Motor, Rear-Wheel Drive | |
| Brakes | Auto Magnetic Brakes | |
| Batteries ³ | 2 x 40 Ah Deep Cycle | |
| Battery Charger | 5-Amp, Off-Board | |
| Other Features | | |
| Versatile Design, Easily Adjustable Seat to Allow for Growth | | |
| Swing-Away Legrests | | |

Height Adjustable Armrests

- **1)** Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings.
- **3)** Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 21 stone 6 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- **6)** Maximum dimensions of the Manual Recline vary by powerbase and seating system. Call for details.
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.
- **8)** All Quantum Rehab® units can be ordered with specialty controls.
- 9) Calls for availability.

Colours

Please contact your representative for the full range of options and colours.



Black



REV2: 08/11/22

