

## 4FRONT 2 HD WITH TB3 SEATING

### DEALERS DETAILS

Dealers Name \_\_\_\_\_

Order Date \_\_\_\_\_

Address \_\_\_\_\_

Postcode \_\_\_\_\_

Phone \_\_\_\_\_ Email \_\_\_\_\_

PO No \_\_\_\_\_ Order Ref \_\_\_\_\_

### DELIVERY DETAILS

Dealers Address \_\_\_\_\_

Postcode \_\_\_\_\_

Quantum Specialist \_\_\_\_\_

Ⓢ Clients Weight (required for suspension settings) \_\_\_\_\_

### BASE MODEL

**£8,995**

- 14 Colour Options
- Height Adjustable Removable Armrests
- Centre Mount Foot Platform
- TB3 Seating
- 14" Flat Free Drive Tyres
- 2" x 9" Flat Free Rear Castors



- SRS Suspension System
- Smart Traction Control
- Q-Logic 3e Control System
- 2 x 70 ah Batteries
- 4 Mph (6 Km/h) 4-Pole Motors
- Max User Weight: 32 st 14lbs. (204 kg)

Standard = Standard Factory Feature  
FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	Standard
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC
4" x 14" Pneumatic Black Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£220
4" x 14" Flat Free Black Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£320

### COLOURS

Back in Black	Standard	Lemon Crush*	FOC
Green Machine	FOC	Night Train*	FOC
Midnight Blue	FOC	Evolution Orange*	FOC
Electric Blue	FOC	Military Green*	FOC
Blue Jean	FOC	Grey Street*	FOC
Deep Purple	FOC	Raspberry Beret*	FOC
Red Corvette	FOC	Custom Hydro Dip*	£600
Pretty in Pink	FOC		

\* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



## SEAT PACKAGES - SEAT WIDTH

20"-24"W.

Set width at:

## SEAT PACKAGES - SEAT DEPTH

20"-24"D.

Set depth at:

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>294 Kg weight limit.</i>	Standard
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>294 Kg weight limit.</i>	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	FOC
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	FOC
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i>	FOC

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

## SEAT POSITIONING

Static Seat.	Standard		
TRU-Balance 3 HD Powered Seat Tilt. <i>0° to 50° COG tilt. Uses 1 power option.</i>	£1200		
Powered Seat Lift & Tilt. <i>10" lift &amp; 50° tilt. Up to 3.5 Mph drive speed while elevated. Uses 2 power options.</i>	£2170		
Static Seat to Floor Height. <i>From seat pan.</i>	15"	16"	17"
Powered Tilt Seat to Floor Height. <i>From seat pan.</i>	17.75"	18.75"	19.75"
Powered iLevel Lift & Tilt Seat to Floor Height. <i>From seat pan.</i>	17.75"	18.75"	19.75"

## BACKREST SPECIFICATIONS

TRU-Balance 3 Back Support System. <i>Select TB3 back options below. Ignore all 'Back Cane' options.</i>	Standard
Back Canes Only System. <i>Not available with Power recline or TB3 flip-back arm supports. Move directly to 'Back Cane Only Support Options' if selected.</i>	FOC

Select one option above. Both options above include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a set of wedge inserts.

## TB3 BACK SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i>	Standard
Privacy Flap. <i>12.75"W. Closes gap between back support and seat pan.</i>	£50
Stealth Adapter Interface Plate. <i>SSAPT4. For mounting third party back supports to TB3 back system. Not required for Q-Back 2.</i>	£341
ADI Q-Back 2. <i>Skin protection, 4 way stretch, Coolcore cover. Mounts directly to TB3 system.</i>	POA

Enter ADI Q-Back 2 part number here for quotation.

## TB3 BACK SUPPORT POSTIONING

TRU-Balance 3 Static Back. <i>Angle can be preset manually.</i>	Standard
TRU-Balance 3 Flip-forward Back. <i>Not available with TB3 flip-back arm supports, powered recline, or 12"-14"W seats.</i>	£150
TRU-Balance 3 HD Powered Recline. <i>85° to 175° with sheer reduction. Uses 1 power option.</i>	£800

## BACK CANE ONLY SUPPORT

Active Contour Back Support with TRU-Comfort 2 Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions.</i>	FOC
Privacy Flap. <i>12.75"W. Closes gap between back support and seat pan.</i>	£50
ADI Back Support. <i>Skin protection, 4 way stretch, Coolcore cover. Mounts to canes with cane clamps.</i>	Standard

Enter ADI Back Support part number here for quotation

Back Cane Only Support Positioning.	Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i>					Standard
	Flip-forward Back Canes. <i>Not suitable for occupied transport or with Endo mounts.</i>					£200
	2" Endo Cane Mounts. <i>1" each side. Reduces back width by 2".</i>					£40
	4" Endo Cane Mounts. <i>2" each side. Reduces back width by 4".</i>					£40
Back Cane Height. <i>From seat pan.</i>	16"H	17"H	20"H	21"H	24"H	
	25"H					

## BACK SUPPORT DIMENSIONS

Compatible with TB3 and Back Canes options above. Back Support can be +/- 4" of seat width with TB3 Back Support System, but must be same as seat width with Back Canes Only System unless -2" or -4" Endo Mounts have been selected.

Backrest Support Width Package.	16"-20"W		20"-24"W		24"-28"W	
	28"-32"W					
Back Support Width.	16"W.	17"W.	18"W	19"W		
	20"W	21"W	22"W	24"W		
	26"W	28"W				
Back Support Height. <i>From seat pan.</i>	16"H. <i>Not available with power recline.</i>	17"H. <i>Not available with power recline.</i>	18"H. <i>Not available with power recline.</i>	19"H.		
	20"H.	21"H.	22"H.	23"H.		
	24"H.	25"H.	26"H.	27"H.		
	28"H.					

19" minimum backrest height when powered recline is selected.

## ARM SUPPORTS

Height Adjustable, Removable, Heavy Duty.				Standard
HD Removable Height Package. <i>From seat pan.</i>	8"-12"	FOC	10"-14".	Standard
TB3 Flip-back with Down Posts. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>				£450
Arm Support Pads.	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	17" Standard. <i>TB3 only. Minimum seat depth 16".</i>	£50	Falcon Gel Mounting Plates. <i>Required for 2-post, Single Post &amp; Paediatric arm supports</i>	£40
	10" Falcon Gel	£40		
	14" Falcon Gel	£40		

## CENTRE MOUNTED LEG SUPPORTS

Manually Adjustable Centre Mount Foot Platform.				Standard
Power Elevating Articulating Centre Mount Foot Platform. <i>Uses 1 power option.</i>				£700
Power Elevating Articulating & Lower To Floor Centre Mount Foot Platform. <i>Minimum seat width 14". Only with Q-Logic 3. Lowers to the floor to aid standing transfers. Uses 2 power options.</i>				£1112
Lower Leg Extensions.	10"-16". <i>From seat pan.</i>	Standard	13"-19".	FOC
Footplates.	10"W x 8"D Split Footplates.			FOC
	10"W x 11"D Split Footplates.			Standard
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>	Standard	Calf Panel. <i>18"W x 8"H.</i>	£42
	Reverse Dartex Calf Pad Covers.	£42	Reverse Dartex Calf Panel Cover.	£62
Footplate Accessories.	4" Width Extension Kit. <i>Not available with 2" length kit.</i>			£16
	2" Length Extension Kit. <i>Not available with 4" length kit.</i>			£16
	Footplate Connector.	£10	Outer Edge Supports.	£10
	Wheel Bumpers.	£10		

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.

## INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

70° Heavy Duty, Drop-in, Swing-away Leg Supports. <i>Select lowers and footplates below.</i>							£360
70° HD	Lower Leg Extensions.	8.5"-12"	FOC	13"-18"	Standard	16"-21"	FOC
	Footplate Type.	Angle & Depth Adjustable. <i>Select size below.</i>					Standard
	Footplate Size.	6"W x 8"D	Standard	6"W x 11"D	FOC	8"W x 11"D	FOC
	Accessories.	Wheel Bumpers					£10
Heavy Duty, Manual Elevating, Swing-away Leg Supports.							£360
ELEVATING HD	Lower Leg Extensions.	11"-14"	FOC	14"-17"	Standard	12.5"-16.5"	FOC
		15.5"-19.5"	FOC				
	Footplate Type.	Flip-up Extruded. <i>8.5"W x 6"D.</i>					Standard
		Angle & Depth Adjustable. <i>Select size below.</i>					£115
	Footplate Size.	6"W x 8"D					Standard
		6"W x 11"D					FOC
	8"W x 11"D					FOC	
Accessories.	Wheel Bumpers					£10	

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information. Adjustment ranges above are from seat pan to footplates. Individual leg supports supplied with heel loops. Elevating leg supports supplied with calf supports.

Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away leg support hanger for mounting.</i>	LEFT	£300	Limb Support	£300	RIGHT
--	------	------	--------------	------	-------

## STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>	£219
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£219
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	£219
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i>	£327
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£88

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative head supports.

## STEALTH TRUNK SOLUTIONS

UniLink Thoracic Lateral Supports. <i>Supplied as a pair.</i>	Lateral Thoracic Swing-away Hardware.	£425	Additional 1.5" Link (each).	£52		
	3" H x 5" L Curved Lateral Pad.	£72	Additional 2" Link (each).	£52		
	4" H x 6" L Curved Lateral Pad.	£72	Additional 3" Link (each).	£52		
4-Point Harness.	Medium	£137	Classic STD*	Dynamic STD*	Vinyl Face	FOC
	Large	£137	Contoured FOC	Structured FOC	Zipper Face**	FOC
	Extra Large	£137				
Chest Strap.	Medium	£87	Side Release	Standard		
	Large	£87	Velcro Release		£19	
	Extra Large	£87				

\*STD = Standard option \*\*Classic Only

## STEALTH PELVIC & THIGH SOLUTIONS

UniLink Thigh Supports.	LEFT	£176	UniLink Hip Guide Long Removable.	£176	RIGHT	
		£176	Unilink Hip Guide Short Removable.	£176		
		£48	4x4 Pelvic Thigh Pad.	£48		
		£48	4x8 Pelvic Thigh Pad.	£48		
		£48	4x12 Pelvic Thigh Pad.	£48		
Knee Adductors.	LEFT	£156	Fixed Cane Mounted.	£156	RIGHT	
		£170	Swing-away Track Mounted.	£170		
		£36	3" Round Knee Pad.	£36		
		£36	4" Round Knee Pad.	£36		
Pelvic Belts.	Auto Buckle Lap Belt, Push Button Release.			Standard		
	Retractable Lap Belt, Push Button Release.			£156		
Positioning 2-Point Padded Pelvic Belt.	Medium	£78	Single Pull	Standard	Push Button	Standard
	Large	£78	Rear Pull	FOC		
	Extra Large	£78	Dual Pull	FOC		
Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull.	Medium	£120	Push Button	Standard		
	Large	£120	Side Release		£20	
Medial Supports (Pommels).	Fixed Flip-down Hardware. <i>MTHW-660.</i>	£180	Small Wedge. <i>MTMD-12.</i>	£103		
	Medium Wedge. <i>MTMD-13.</i>	£103	Large Wedge. <i>MTMD-14.</i>	£103		
	Adjustable Wedge. <i>MTMD-42.</i>	£199				

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative lateral, pelvic & lower extremity solutions

## EXPANDABLE CONTROLS

Q-Logic 3e Controls. <i>No power options.</i>	Standard
Q-Logic 3e Controls. <i>2 power option ready.</i>	£300
Q-Logic 3e Controls. <i>5 power option ready.</i>	£500
Q-Logic 3 EX Controls. <i>No power options.</i>	£275
Q-Logic 3 EX Controls. <i>2 power option ready.</i>	£350
Q-Logic 3 EX Controls. <i>5 power option ready.</i>	£550
Q-Logic 3 Enhanced Display. <i>Includes 2 egg switches &amp; flexible display mount.</i>	£1074
SCIM Module Kit. <i>Allows Q-Logic 3 EX joystick display to be used for speciality controls. Not available with 3e controls.</i>	£316

iAccess OLED 4 Button Switch Module.  
*Mounted to opposite armrest unless specified.*

£300

Push Button iAccess  
Toggle Lever iAccess

If iAccess is selected choose either Push Button or Toggle lever option.

Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Omit Hand Control. <i>Often selected when ED is used as master remote with special device.</i>	-£50
	Attendant Control. <i>Hand control mounted at rear.</i>	£100
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative control options such as mini joysticks, head arrays & switches.

## HAND CONTROL MOUNTS

Fixed.	Standard				
Swing-away. <i>Swings in &amp; out with height, tension and pivot adjustment.</i>	£75				
Retract 4 Swing-away. <i>Full 180° Movement.</i>	£150				
Mounting Position.	Left	FOC	Right	Standard	
Hand Control Joystick Handle.	Standard Carrot	Standard	T-Bar	£35	Extended I-Handle £35
	Mushroom	£35	Medium Ball	£35	3" Goal Post £50
	Chin Cup	£35	Soft Ball	£35	4" Goal Post £50

## SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
4 Mph Program with Lights and indicators.	£450
Lights & Indicator Surcharge. <i>Required if no power options have been selected.</i>	£260
LED Fender Lights. <i>Includes single switch button mounted to side rail. Not road legal.</i>	Standard
Omit LED Fender Lights.	FOC

## ACCESSORIES

USB Mobile Device Charger. <i>Please specify mounting position. i.e arm support, backrest, seat rail.</i>	£110
Cup Holder. <i>Only with TB3 flip-back arm supports.</i>	£20
Personal Item Hooks. <i>Not with back cane option.</i>	£15
Mobile Phone Holder. <i>Only with TB3 arm supports.</i>	£40
Back Pack Holders. <i>Not with back cane option.</i>	£30
Clothes Guards. <i>Only with TB3 arm supports.</i>	£100

## ACCESSORIES CONTINUED

Glove Box.	£30
Hydration System.	£125
Oxygen Holder.	£50
Quick Release Transfer Bars.	£125
Stealth Tablet Holder.	£342
Push Handles. <i>Not with back cane option.</i>	£50
Accessory Bar.	£50
Back-up Camera.	£300
Quantum Back Pack.	£110
Additional XLR USB Charger Port.	£17
Magnetic to XLR Charger Adapter. <i>Only with Q-Logic 3e.</i>	£49
Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379
Wooden Trim Kit for Tray. <i>Only with swing-away tray table.</i>	£100

## NOTES & SPECIAL REQUESTS

### PRICE

**4FRONT 2 HD  
WITH TB3 SEATING**



**BASE PRICE** **£8,995**

**COST OPTIONS**

**TOTAL**

7



**Quantum®**  
32 Wedgwood Road, Bicester,  
Oxon, OX26 4UL  
Phone: 01869 324600  
Fax: 01869 323070  
Email: sales@quantumrehab.co.uk  
Web: www.quantumrehab.co.uk  
**RETAIL UK - 06/10/23 V1**



**Crash Tested**



V1: 14/06/2022

