





Weight capacity of up to 32 st 2 lbs (204 kg)



Two 73 Ah Batteries (as Standard)



Available 10" Lift for weights up to 28 st 8 lbs



14" (356 mm) flat free or pneumatic

4FRONT Ho

The 4Front® 2 HD is the next generation of the 4Front® powerchair with added heavy duty capabilities. Featuring new Smart Traction Control (STC), along with Smooth Ride Suspension (SRS), to provide stability for smooth obstacle transitions and unmatched driving performance.







Specifications		
Power Chair		
Weight Capacity	32 st 2 lbs (204 kg) 28 st 8 lbs (181 kg) with 10" lift	
Maximum Speed ¹	Up to 4 mph (6.4 kh/h) Up to 2.7 mph (4.3 km/h) with 10" lift	
Ground Clearence ²	80 mm (3.15") at battery tray	
Turning Radius ²	623 mm (24.5")	
Overall Size ^{2,6}	Length: Width:	1168 mm (45.9") 618 mm (24.3")
Component Weights ⁵	Base: Batteries:	63.5 kg (140 lbs) 23.36 kg (51.5 lbs) each (standard)
Wheels		
Drive Wheels	Standard: 356 mm (14") flat free	
Caster Wheels	Rear: 229 mm (9") flat free	
Anti-Tip Wheels	Front: 102 mm (4") flat free	
Suspension	SRS (smooth ride suspension)	
Batteries		
Batteries ³	2 x 73 Ah batteries (as standard)	
Battery Charger	8-amp, off-board	
Seating Sizes		
Width	20" - 24"	
Depth	20" - 24"	
Available Options		
10" of power adjustable safe seat elevation at speeds of up to 3.5 mph		
USB Mobile Device Charger Please specify fitting location i.e. left armrest or seat rail		
Cup Holder Only available with TB3 Flip-up armrests		
Personal Item Hook Not available with Back Cane Option		
Reflector Kit		
Mobile Phone/Media Holder Only available with TB3 Flip-up armrests		
Backpack Holder (pair) Not available with Back Cane Option		
Glove Box		
Hydration System		
Oxygen Holder		

Personal Item Hook Not available with Back Cane Option

- 1) Varies with base model, user weight, terrain type, battery amp hour rating (Ah), battery charge, battery condition, motors, controller type, tire type, and tire condition.

 This specification is subject to a variance of +10%, -5%.
- 2) Due to manufacturing tolerances and continued product improvement, this specification can be subject to a variance of (+ or -) 3%.
- **3)** AGM, Gel-Cell, or Hybrid Gel (Hybrid Gel is defined as Nano Silica Gel Electrolyte with High Quality Glass Mat Separators) type required. See "Batteries and Charging."
- 4) Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 5) Battery weight may vary based on manufacturer.
- **6)** Dimensions listed are for power base only.

 Overall measurements will vary based on seating and accessory selections.

Colours

Please contact your representative for the full range of options and colours.





VER1: 25/07/22

