



Bold Aesthetics



iLevel is Retrofittable on the Q6 Edge® 2.0 HD



iLevel® Option Allows 12" of Seat Elevation While at Walking Speed (3.5 mph)



Standard 4 mph 4-Pole Motors Offer Advanced Performance, Precise Control and Responsiveness via Innovative Single-Stage Gearbox

Q6 edge® 2.0 HD

With features like standard 4 mph 4-pole motors and advanced responsiveness, the Q6 Edge® 2.0 HD is engineered to meet the performance needs of the most active user. The Q6 Edge 2.0 HD accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value, and is available with our iLevel® technology.

Specifications

Power Chair

Weight Capacity	28 st 5 lbs. (181 kg)
Maximum Speed ¹	Up to 4 mph / Up to 3.5 mph w/iLevel Activated
Ground Clearance ²	69 mm (2.7")
Turning Radius ²	521 mm (20.5")
Overall Size ^{2,6}	Length: 902 mm (35.5") Without Front Rigging Width: 622 mm (24.4") 735 mm (28.9") With Armrests
Component Weights ^{5,7}	Base: 138.6 lbs (63 kg) Batteries: 38 lbs (17.24 kg) Each
Drivetrain	Two Motor, Mid-Wheel 6 [®] Drive
Maximum Obstacle Climbing Ability	90 mm (3.5")

Wheels

Drive Wheels	356 mm (14") Pneumatic or Solid Available (14" Knobby Option Available)
Caster Wheels	Front: 152 mm (6") Solid Rear: 152 mm (6") Solid (8" Solid Front & Rear Optional)
Suspension	ATX Suspension (Active-Trac [®] With Extra Stability)
Brakes	"Intelligent Braking" Electronic Regenerative, Disc Park Brake

Batteries

Batteries ³	Two 12v, Deep-Cycle 62 Ah Batteries
Battery Charger	8-Amp, Off-Board

Seating Sizes

TRU-Balance [®] 3 Static Seat	Width: 20" - 24" Depth: 20" - 24"
TRU-Balance [®] 3 Power Tilt	Width: 20" - 24" Depth: 20" - 24"
TRU-Balance [®] 3 Tilt & iLevel [®]	Width: 20" - 24" Depth: 20" - 24"

Seat to Floor Heights

TRU-Balance [®] 3 Static Seat ⁵	15.875-17.875"
TRU-Balance [®] 3 Power Tilt ⁵	17.5-19.5"
TRU-Balance [®] 3 Tilt & iLevel [®]	17.5-19.5"
TRU-Balance [®]	15.875-19.875"

Available Options

TRU-Balance [®] 3 Power Recline
TRU-Balance [®] 3 Power Tilt and Power Recline
LED Fender Lights (Optional)
8" Solid Castors / 14" Knobby Drive Wheels
Powered Elevating Leg Rests or Footboard
NE*
Q-Logic 3 EX
VR2
R-Net
iLevel [®] 12" Power Adjustable Seat Height

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings.
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 21 stone 6 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Maximum dimensions of the Manual Recline vary by power base and seating system call for details.
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 8) All Quantum Rehab[®] units can be ordered with speciality controls.
- 9) Call for availability.

Colours

Please contact your representative for the full range of options and colours.



Electric Blue



Red Corvette



Yellow Blaze



Lime Light



Back in Black



Midnight Blue



VER 1 : 25/07/22



quantumrehab.co.uk
01869 324600