



Sleek new aesthetic without compromising power



Impressive new SRS suspension system



LED fender lights package

QUANTUM[®] R-TRAK

Hold the smoothest ride at your fingertips. Designed housing an impressive SRS suspension system while promising an all over exceptional performance. A new sleek design is further enhanced by the clean aesthetic detailing combining elegance and power. R-Trak stands as the best answer for all users with active lifestyles.

01869 324600 | quantumrehab.co.uk



Specifications

Power Chair	
Weight Capacity	300 lbs. (136 kg/21 stone 6 lbs.)
Maximum Speed	Up to 6 mph (11 km/h) Up to 3.5 mph (5.6 km/h) with 10-in. lift
Ground Clearance	3.97 in. (101 mm) at battery tray
Turning Radius	25.7 in. (653 mm)
Overall Size	Length: 34.8 in (884 mm) without front rigging 44.5 in. (1130 mm) with front riggings Width: Width: 24.2 in. (616 mm)
Component Weights	Base: 126.4 lbs (57.25 kg)
Wheels	
Drive Wheels	14 in. (356 mm), solid or pneumatic
Caster Wheels	9 in. (102 mm), flat free
Suspension	SRS (Smooth Ride Suspension)
Drivetrain	Two motor, rear-wheel drive
Brakes	"Electronic regenerative, disc park brake
Batteries	
Batteries	Two 12 volt, Flat top deep-cycle C20 70 Ah, C5 62 Ah
Battery Charger	8-Amp, Off-Board
Optional Extras	
Control Options	
	NE+ Drive Controls
	Q-Logic 3e
	Q-Logic 3 EX Drive Controls


Colours



Blue
Jean



Red
Corvette



Lemon
Crush



Military
Green



Raspberry
Beret



Evolution
Orange



Sugar
Plum



Green
Machine



Deep
Purple



Electric
Blue



Black in
Black



Grey
Street



White
Satin



Cherry
Bomb



Night
Train