

## R-TRAK WITH TB FLEX SEATING

### DEALERS DETAILS

Dealers Name \_\_\_\_\_  
 Order Date \_\_\_\_\_  
 Address \_\_\_\_\_  
 \_\_\_\_\_ Postcode \_\_\_\_\_  
 Phone \_\_\_\_\_ Email \_\_\_\_\_  
 PO No \_\_\_\_\_ Order Ref \_\_\_\_\_

### DELIVERY DETAILS

Dealers Address \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_ Postcode \_\_\_\_\_  
 Quantum Specialist \_\_\_\_\_  
 Ⓢ Clients Weight \_\_\_\_\_

### BASE MODEL

**£5,495**

- 14 Colour Options
- Height Adjustable Removable Armrests
- Centre Mount Foot Platform
- Angle Adjustable Fold-Forward Backrest
- 14" Flat Free Drive Tyres
- 3" x 9" Flat Free Front Castors



- SRS Suspension System
- NE+ Control System
- 2 x 70 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:  
21 st 6lbs. (136 kg)

Standard = Standard Factory Feature  
 FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	Standard
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC
4" x 14" Pneumatic Black Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£220
4" x 14" Flat Free Black Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£320

### COLOURS

Back in Black	Standard	Lemon Crush*	FOC
Green Machine	FOC	Night Train*	FOC
Midnight Blue	FOC	Evolution Orange*	FOC
Electric Blue	FOC	Military Green*	FOC
Blue Jean	FOC	Grey Street*	FOC
Deep Purple	FOC	Raspberry Beret*	FOC
Red Corvette	FOC	Custom Hydro Dip*	£600
Pretty in Pink	FOC		

\* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



## SEAT PACKAGES - SEAT WIDTH

15"W.	FOC	19"W.	FOC
16"W.	FOC	20"W.	FOC
17"W.	FOC	21"W.	FOC
18"W.	Standard	22"W.	FOC

## SEAT PACKAGES - SEAT DEPTH

15"D.	FOC	19"D.	FOC
16"D.	FOC	20"D.	FOC
17"D.	FOC		
18"D.	Standard		

## STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	£100
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£115
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£115
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£135
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i>	£135
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

## SEAT POSITIONING

Manual Settable Tilt. <i>0° to 30° tool adjustable with a telescopic tube.</i>	Standard		
Powered Seat Tilt. <i>0° to 30°. Uses 1 power option.</i>	£450		
Powered Seat Lift & Tilt. <i>110" lift &amp; 50° tilt. Uses 2 power options.</i>	£1895		
Manual Seat to Floor Height. <i>From seat pan.</i>	16.5"	17.5"	
Powered Tilt Seat to Floor Height. <i>From seat pan.</i>	17.5"	18" <small>Changes Tilt range to -5° to 25°</small>	18.5" <small>Changes Tilt range to -10° to 20°</small>
Powered Lift & Tilt Seat to Floor Height. <i>From seat pan.</i>	17.5"	18"	18.5"

## STANDARD BACK SUPPORT POSITIONING

Manual Angle Settable, Fold-Forward. <i>Adj. range 95° to 125°. Will not fold with 10" lift.</i>	Standard
Attendant Controlled Manual Recline. <i>Adj. range 95° to 125°. Not suitable for occupied transportation.</i>	£150
Powered Back Recline. <i>Adj. range 95° to 125°. Uses 1 power option.</i>	£345

## STANDARD BACK SUPPORT

Tension Adjustable. <i>Air-mesh cover. Adjustment range 19" to 23" from seat board. OBAS** required to mount headrest.</i>	Standard
TRU-Comfort 2. <i>4-way stretch cover. Adjustable range 24.5" to 28" from seat board.</i>	£135
ADI Aluminium Back. <i>Skin protection, 4 way stretch, Coolcore cover.</i>	POA
Enter ADI Back Support part number here for quotation.	

## STEALTH PRO BACK SUPPORTS

Standard, Coolcore Cover. <i>Solid back, 9 cm contour</i>	£227
Active, Coolcore Cover. <i>Solid back, 12 cm contour seat pan.</i>	£341
Extra Deep Active, Coolcore Cover. <i>Solid back, 17 cm contour.</i>	£448
Active Flex, Coolcore Cover. <i>Tension adjustable straps. 12 cm contour.</i>	£559
Back Height	50cm. <i>Actual height.</i> Standard 57cm. <i>Actual height.</i> FOC

The back support width will match seat width selected in seat dimensions.

## ARM SUPPORTS

Std Flip-up Height Adjustable. <i>Includes contour arm pads. Do not select arm pads below.</i>	Standard
TB Flex Removable, Quick Release. <i>Adj. range 9" to 14" from seat pan.</i>	£85
Single Post Quick Height Adjustable Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14"</i>	£100
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16"</i>	£279
Arm Support Pads.	10" Standard Standard 10" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC 14" Standard FOC 14" Waterfall. <i>Not with Std flip-up or TB Flex.</i> FOC

## CENTRE MOUNTED LEG SUPPORTS

Manually Adjustable Centre Mount Foot Platform. <i>Minimum seat width 14"</i>	Standard
Paediatric Adjustable Centre Mount Foot Platform. <i>12"-14" seat widths only.</i>	FOC
Lower Leg Extensions. <i>From seat pan.</i>	8"-12". <i>Paediatric centre mount only.</i> FOC 10"-16". <i>Manually Adj. &amp; powered elevating only</i> FOC 13"-19". <i>Manually Adj. &amp; powered elevating only</i> Standard
Footplates.	10"W x 6"D Split Footplates. <i>Paediatric centre mount only.</i> FOC 10"W x 8"D Split Footplates. FOC 10"W x 11"D Split Footplates. Standard
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i> Standard } Calf Panel. <i>18"W x 8"H.</i> £42 Reverse Dartex Calf Pad Covers. £42 } Reverse Dartex Calf Panel Cover. £62
Footplate Accessories.	4" Width Extension Kit. <i>Not available with 2" length kit or 10"W x 6"D footplates.</i> £16 2" Length Extension Kit. <i>Not available with 4" width kit or 10"W x 6"D footplates.</i> £16 Footplate Connector. £10 Outer Edge Supports. £10 Wheel Bumpers. £10

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.

## INDIVIDUALLY MOUNTED LEG SUPPORTS

60° to 90° Tool Adjustable Swing-away Legrests. <i>Adjustment range 11" to 16.5"</i>	FOC
Manual Elevating Swing-Away Legrests. <i>Adjustment range 15" to 20"</i>	£250
Power Elevating Swing-away Legrests. <i>Adj. range 14.5" to 19.5". Min. seat width 16". Uses 2 power options.</i>	£650
Enter seat pan to heel measurement. <i>If no measurement is supplied the standard factory length option will be fitted.</i>	
Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away leg support hanger for mounting.</i>	LEFT £300 Limb Support £300 RIGHT

Adjustment ranges above are from seat pan to footplates. Individual legrests supplied with heel loops. Elevating legrests supplied with calf supports.

## STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>	£219
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£219
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	£219
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i>	£327
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£88
Push Handles Head Mount Adapter OBAS**. <i>Required to mount head support with TA Back or Active Flex Pro Back.</i>	£266

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative head supports.

## TRUNK SOLUTIONS

Thoracic Lateral Supports. <i>Supplied as a pair.</i>	<b>LEFT</b>	£150	Lateral Thoracic Swing-away Hardware.	£150	<b>RIGHT</b>	
		£35	3" H x 5" L Curved Lateral Pad.	£35		
		£35	4" H x 6" L Curved Lateral Pad.	£35		
Stealth 4-Point Harness.	Extra Small	£137	Classic STD*	Dynamic STD*	Vinyl Face	FOC
	Small	£137	Contoured FOC	Structured FOC	Zipper Face**	FOC
	Medium	£137				
	Large	£137				
	Extra Large	£137				
Chest Strap.	Extra Small	£87	Small	£87	Side Release	Standard
	Medium	£87	Large	£87	Velcro Release	£19
	Extra Large	£87				

\*STD = Standard option \*\*Classic Only

## PELVIC SOLUTIONS

Pelvic Belts.	Auto Buckle Lap Belt, Push Button Release.				Standard	
	Retractable Lap Belt, Push Button Release.				£156	
Stealth Positioning 2-Point Padded Pelvic Belt.	Extra Small	£78	Single Pull	Standard	Side Release	Standard
	Small	£78	Rear Pull	FOC	Push Button	FOC
	Medium	£78	Dual Pull	FOC		
	Large	£78				
	Extra Large	£78				
Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull.	Extra Small	£120	Small	£120	Push Button	Standard
	Medium	£120	Large	£120	Side Release	£20

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative lateral, pelvic & lower extremity solutions.

## NON-EXPANDABLE CONTROLS

NE + Controls. <i>With no power options.</i>	Standard	
NE + Controls. <i>2 power option ready.</i>	£125	
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	£100
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

## EXPANDABLE CONTROLS

Q-Logic 3e and 3 EX hand controls are supplied with built-in Bluetooth. Q-Logic 3 Enhanced Display is supplied with built-in Bluetooth and Infrared.

Q-Logic 3e Controls. <i>No power options.</i>	£225
Q-Logic 3e Controls. <i>2 power option ready.</i>	£300
Q-Logic 3e Controls. <i>5 power option ready.</i>	£500
Q-Logic 3 EX Controls. <i>No power options.</i>	£275
Q-Logic 3 EX Controls. <i>2 power option ready.</i>	£350
Q-Logic 3 EX Controls. <i>5 power option ready.</i>	£550
Q-Logic 3 Enhanced Display. <i>Includes 2 egg switches &amp; flexible display mount.</i>	£1074
IR Blaster. <i>Infrared module for Q-Logic 3 EX hand control.</i>	£50
SCIM Module Kit. <i>Allows Q-Logic 3 EX joystick display to be used for speciality controls. Not available with 3e controls.</i>	£316

iAccess OLED 4 Button Switch Module.  
*Mounted to opposite armrest unless specified.*

£300

Push Button iAccess  
Toggle Lever iAccess

If iAccess is selected choose either  
Push Button or Toggle lever option.

Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Omit Hand Control. <i>Often selected when ED is used as master remote with special device.</i>	£-100
	Attendant Control. <i>Hand control mounted at rear.</i>	£125
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative control options such as mini joysticks, head arrays & switches.

## HAND CONTROL MOUNTS

Fixed.						Standard
Swing-away. <i>Swings in &amp; out with height, tension and pivot adjustment.</i>						£75
Retract 4 Swing-away. <i>Full 180° Movement.</i>						£150
Mounting Position.	Left	FOC	Right	Standard		
Hand Control Joystick Handle.	Standard Carrot	Standard	T-Bar	£35	Extended I-Handle	£35
	Mushroom	£35	Medium Ball	£35	3" Goal Post	£50
	Chin Cup	£35	Soft Ball	£35	4" Goal Post	£50

## SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
4 Mph Program with Lights and indicators.	£450
6 Mph Program with Lights and indicators.	£450
Lights & Indicator Surcharge. <i>Required if no power options have been selected.</i>	
LED Front Castor and Rear Fender Lights. <i>Includes single switch button mounted to side rail. Not road legal.</i>	Standard
Omit LED Front Castor and Rear Fender Lights.	FOC

## ACCESSORIES

Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379
Quantum Back Pack.	£110
XLR USB Charger Port.	£17
Magnetic to XLR Charger Adapter. <i>Only with Q-Logic 3e.</i>	£49
Back-up Camera.	£300

## NOTES & SPECIAL REQUESTS

## PRICE

**R-TRAK WITH  
TB FLEX SEATING**



BASE PRICE	<b>£5,495</b>
COST OPTIONS	
TOTAL	

6



**Quantum®**  
32 Wedgwood Road, Bicester,  
Oxon, OX26 4UL  
Phone: 01869 324600  
Fax: 01869 323070  
Email: sales@quantumrehab.co.uk  
Web: www.quantumrehab.co.uk  
**RETAIL UK - 21/02/24 V1**



**Crash  
Tested**



V1: 14/06/2022