

R-TRAK WITH TB FLEX SEATING

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
	Postcode
Phone	Email
PO No	Order Ref

DELIVERY DETAILS

Dealers Address	
Postcode	
Quantum Specialist	
Ⓢ Clients Weight	

BASE MODEL

£5,650

14 Colour Options
Height Adj. Removable Arm Supports
Centre Mount Foot Platform
Angle Adj. Fold-Forward Back Support
14" Flat Free Drive Tyres
3" x 9" Flat Free Front Castors



SRS Suspension System
QLNES Control System
2 x 70 ah Batteries
4 Mph (6 Km/h) Motors
Max User Weight:
21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres. <i>With matching colour wheel accents.</i>	Standard
3" x 14" Pneumatic Black Tyres. <i>With matching colour wheel accents.</i>	FOC
4" x 14" Pneumatic Black Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£220
4" x 14" Flat Free Black Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£320

COLOURS

Back in Black	Standard	Lemon Crush*	FOC
Green Machine	FOC	Night Train*	FOC
Midnight Blue	FOC	Evolution Orange*	FOC
Electric Blue	FOC	Military Green*	FOC
Blue Jean	FOC	Grey Street*	FOC
Deep Purple	FOC	Raspberry Beret*	FOC
Red Corvette	FOC	Custom Hydro Dip*	£600
Pretty in Pink	FOC		

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

15"W.	Not available with 10" Seat Lift.	FOC	19"W.	FOC
16"W.		FOC	20"W.	FOC
17"W.		FOC	21"W.	Not available with 10" Seat Lift.
18"W.	Standard		22"W.	Not available with 10" Seat Lift.

SEAT PACKAGES - SEAT DEPTH

15"D.	FOC	19"D.	FOC
16"D.	FOC	20"D.	FOC
17"D.	FOC		
18"D.	Standard		

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	Standard
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	FOC
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	£99
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£126
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£126
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£148
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i>	£148
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£209
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£209

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

Static Seat.			Standard
Powered Seat Tilt. <i>0° to 30°. Uses 1 power option.</i>			£500
Powered Seat Tilt. <i>0° to 45°. Uses 1 power option.</i>			£650
Powered Seat Lift & Tilt. <i>10" lift & 50° tilt. Maximum seat width 20". Uses 2 power options.</i>			£1950
Static Seat. <i>From seat pan.</i>		16.5"	17.5"
Powered 30° Tilt Seat to Floor Height. <i>From seat pan.</i>	17.5"	18" <small>Changes Tilt range to -5° to 25°</small>	18.5" <small>Changes Tilt range to -10° to 20°</small>
Powered 45° Tilt Seat to Floor Height. <i>From seat pan.</i>	17.5"	18" <small>Changes Tilt range to -5° to 40°</small>	18.5" <small>Changes Tilt range to -10° to 40°</small>
Powered Lift & Tilt Seat to Floor Height. <i>From seat pan.</i>		18.5"	19.5"

STANDARD BACK SUPPORT POSTIONING

Manual Angle Settable, Fold-Forward. <i>Adj. range 95° to 125°. Will not fold with 10" lift.</i>	Standard
Attendant Controlled Manual Recline. <i>Adj. range 95° to 125°. Not suitable for occupied transportation.</i>	£150
Powered Back Recline. <i>Adj. range 95° to 125°. Uses 1 power option.</i>	£400

STANDARD BACK SUPPORT

Tension Adjustable. <i>Air-mesh cover. Adjustment range 19" to 23" from seat board. OBAS** required to mount headrest.</i>	Standard
TRU-Comfort 2. <i>4-way stretch cover. Adjustable range 24.5" to 28" from seat board.</i>	£135
ADI Aluminium Back. <i>Skin protection, 4 way stretch, Coolcore cover.</i>	POA

Enter ADI Back Support part number here for quotation.

STEALTH PRO BACK SUPPORTS

Standard, Coolcore Cover. <i>Solid back, 9 cm contour.</i>				£227
Active, Coolcore Cover. <i>Solid back, 12 cm contour.</i>				£341
Extra Deep Active, Coolcore Cover. <i>Solid back, 17 cm contour.</i>				£448
Active Flex, Coolcore Cover. <i>Tension adjustable straps. 12 cm contour.</i>				£559
Back Height	50cm. <i>Actual height.</i>	Standard	57cm. <i>Actual height.</i>	FOC

The back support width will match seat width selected in seat dimensions.

ARM SUPPORTS

TB Flex Removable, Quick Release. <i>Adj. range 9" to 14" from seat pan.</i>				Standard
Single Post Quick Height Adjustable Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14"</i>				£100
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16"</i>				£279
Arm Support Pads.	10" Standard	Standard	10" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC
	14" Standard	FOC	14" Waterfall. <i>Not with Std flip-up or TB Flex.</i>	FOC

LEG SUPPORTS

70° Swing-Away. <i>Extension range 12" to 15". Minimum seat width 15".</i>	Standard				
60° to 90° Tool Adjustable Swing-Away. <i>6"W x 8"D foot supports. Extension range 11" to 19.5". Min. seat width 15".</i>	£330				
Manually Adjustable Centre Mount Foot Platform. <i>5"W x 11"L Footplates. Extension range 10" to 19". Min. seat width 15".</i>	£340				
Manual Elevating Swing-Away. <i>Extension range 15" to 20". Min. seat width 15"</i>	£425				
Power Elevating Articulating Swing-Away. <i>5"W x 7"D foot supports. Ext. range 9" to 21". Min. seat width 15". Uses 2 power options</i>	£850				
Powered Leg Support Operation Method.	Independent. <i>Uses 2 power options.</i>	Combined. <i>Uses 1 power option.</i>			
Power Elevating Articulating Centre Mount Foot Platform. <i>5"W x 11"L Footplates. Extension range 10" to 19". Min. seat width 15"</i>	£950				
❶ Enter seat pan to heel measurement. <i>If NO measurement is supplied, standard factory length option will be fitted.</i>					
Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away legrest hanger for mounting.</i>	LEFT	£300	Limb Support	£300	RIGHT

Extension ranges above are from seat pan to foot supports. Individual leg supports supplied with heel loops. Elevating leg supports supplied with standard calf supports.

STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>	£230
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£230
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	£230
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i>	£342
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£92
Push Handles Head Mount Adapter OBAS**. <i>Required to mount head support with TA Back or Active Flex Pro Back.</i>	£266

Please visit stealthproducts.com for alternative head supports.

TRUNK SOLUTIONS

Thoracic Lateral Supports. <i>Supplied as a pair.</i>	LEFT	£175	Lateral Thoracic Swing-away Hardware.				£175	RIGHT			
		£42	3" H x 5" L Curved Lateral Pad.				£42				
		£42	4" H x 6" L Curved Lateral Pad.				£42				
Stealth 4-Point Harness.	Extra Small	£143	}	Classic	STD*	}	Dynamic	STD*	}	Vinyl Face	FOC
	Small	£143		Contoured	FOC		Structured	FOC		Zipper Face**	FOC
	Medium	£143									
	Large	£143									
	Extra Large	£143									
Chest Strap.	Extra Small	£91	Small	£91	}	Side Release		Standard			
	Medium	£91	Large	£91		Velcro Release		£20			
	Extra Large	£91									

*STD = Standard option

**Classic Only

*STD = Standard option **Classic Only

PELVIC SOLUTIONS

Pelvic Belts.	Auto Buckle Lap Belt, Push Button Release.				Standard			
	Retractable Lap Belt, Push Button Release.				£169			
Stealth Positioning 2-Point Padded Pelvic Belt.	Extra Small	£84	Single Pull	Standard	Side Release	Standard		
	Small	£84	Rear Pull	FOC	Push Button	FOC		
	Medium	£84	Dual Pull	FOC				
	Large	£84						
	Extra Large	£84						
Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull.	Extra Small	£143	Small	£143	Push Button	Standard		
	Medium	£143	Large	£143	Side Release	£21		

Please visit stealthproducts.com for alternative lateral, pelvic & lower extremity solutions.

NON-EXPANDABLE CONTROLS

QLNES Controls. <i>With no power options.</i>		Standard
QLNES Controls. <i>2 power option ready.</i>		£125
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	£100
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

EXPANDABLE CONTROLS

Q-Logic 3e and 3 EX hand controls are supplied with built-in Bluetooth. Q-Logic 3 Enhanced Display is supplied with built-in Bluetooth and Infrared.

Q-Logic 3e Controls. <i>No power options.</i>	£300
Q-Logic 3e Controls. <i>3 power option ready.</i>	£350
Q-Logic 3e Controls. <i>5 power option ready.</i>	£500
Q-Logic 3 EX Controls. <i>No power options.</i>	£350
Q-Logic 3 EX Controls. <i>3 power option ready.</i>	£375
Q-Logic 3 EX Controls. <i>5 power option ready.</i>	£550
Q-Logic 3 Enhanced Display. <i>Includes 2 egg switches & flexible display mount. Not available with 3e controls.</i>	£1074
IR Blaster. <i>Infrared module for Q-Logic 3 EX hand control.</i>	£50
SCIM Module Kit. <i>Allows Q-Logic 3 EX joystick display to be used for speciality controls. Not available with 3e controls.</i>	£316

iAccess OLED 4 Button Switch Module.
Mounted to opposite armrest unless specified.

£300

Push Button iAccess
Toggle Lever iAccess

If iAccess is selected choose either
Push Button or Toggle lever option.

Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Omit Hand Control. <i>Often selected when ED is used as master remote with special device.</i>	-£100
	Attendant Control. <i>Hand control mounted at rear.</i>	£125
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

Please visit stealthproducts.com for alternative control options such as mini joysticks, head arrays & switches.

HAND CONTROL MOUNTS

Fixed.	Standard					
Basic swing-away.	£75					
Advanced swing-away. <i>Chain driven . Swings in & out with height, tension and pivot adjustment.</i>	£100					
Mounting Position.	Left	FOC	Right	Standard		
Hand Control Joystick Handle.	Standard Carrot	Standard	T-Bar	£35	Extended I-Handle	£35
	Mushroom	£35	Medium Ball	£35	3" Goal Post	£50
	Chin Cup	£35	Soft Ball	£35	4" Goal Post	£50

SPEED PROGRAM & LIGHTS PACKAGE

2-Pole Motors & 4 Mph Program. <i>Lights and indicators not included.</i>	Standard
2-Pole Motors & 4 Mph Program with Lights and indicators.	£450
2-Pole Motors & 6 Mph Program with Lights and indicators.	£450
4-Pole Motors & 7 Mph Program with Lights and indicators.	£900
Lights & Indicator Surcharge. <i>Required if no power options have been selected.</i>	£260
LED Front Castor and Rear Fender Lights. <i>Includes single switch button mounted to side rail. Not road legal.</i>	Standard
Omit LED Front Castor and Rear Fender Lights.	FOC

ACCESSORIES

Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379
Quantum Back Pack.	£110
XLR USB Charger Port.	£17
Magnetic to XLR Charger Adapter. <i>Only with Q-Logic 3e.</i>	£49
Back-up Camera.	£300

NOTES & SPECIAL REQUESTS

PRICE

**R-TRAK WITH
TB FLEX SEATING**



BASE PRICE **£5,650**

COST OPTIONS

TOTAL

6



Quantum®
32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk
RETAIL UK - 01/04/25 V1

